



Now with Dual Plasma Technology

Unsurpassed

Stability,

Selectivity,

and

Sensitivity

for your

GC Analysis

SIEVERS® Dual Plasma Sulfur and Nitrogen Chemiluminescence Detectors

Overview

The Sievers® Sulfur Chemiluminescence Detector (355 SCD) and Nitrogen Chemiluminescence Detector (255 NCD) are the world's most sensitive and selective chromatographic detectors for sulfur and nitrogen-containing compounds. Ionics Instruments has further enhanced 355 SCD and 255 NCD performance and ease of use by developing the Sievers Dual Plasma™ Technology with its patented detection method, now available in the new Sievers Dual Plasma Burner and Controller.


Sievers Dual Plasma Burner and Controller

The compact Sievers Dual Plasma Burner advances our established technology by improving performance, reducing maintenance, and adding a low temperature safety shroud. The Burner easily mounts on major GC brands and incorporates features for easier and less frequent maintenance. The Dual Plasma Controller has in-line electronic flow sensors, a digital display with added controls and functionality, and a smaller footprint. The new Dual Plasma system delivers the easiest and fastest start-up times of all chemiluminescence detector systems.

Dual Plasma Technology

The Sievers Dual Plasma Technology harnesses the power of two-flame plasma combustion, optimizing combustion of the sample matrix and formation of either sulfur oxide (SO) or nitric oxide (NO). This results in unsurpassed:

- Stability
- Selectivity over carbon
- Sensitivity
- Equimolar and linear response
- Absence of quenching



Ionics Instruments

Sievers



SIEVERS 355 SCD Sulfur Chemiluminescence Detector




Overview

The Sievers 355 Sulfur Chemiluminescence Detector (355 SCD) is the world's most sensitive and selective chromatographic detector for sulfur containing compounds. The 355 SCD is based on Sievers patented technology, and it provides a linear and equimolar response to sulfur compounds without interference from most sample matrices. The exceptional performance of the 355 SCD has resulted in its widespread use and acceptance for the analysis of sulfur compounds in a wide range of applications. The new Sievers Dual Plasma Burner and Controller significantly enhance the performance of the 355 SCD, increase ease of use, and reduce maintenance.

The Sievers Patented Technology

The 355 SCD utilizes the combustion of sulfur compounds to form sulfur monoxide (SO) and the subsequent chemiluminescence reaction of SO with ozone (O₃). The unique combustion process achieves high temperatures (> 1,800 °C) which are unattainable by standard pyrolysis methods. This patented technology allows the 355 SCD to make ultra-sensitive measurements of any sulfur containing compound that can be analyzed by gas chromatography (GC) or supercritical fluid chromatography (SFC).




The light (hv) passes through an optical filter and is detected by a photomultiplier tube (see 355 SCD Fluidics diagram). This mechanism provides selective sulfur detection which is described in the following US and foreign patents: 5,330,714; 5,227,135; 5,310,683; 5,501,981; 5,424,217; 5,661,032; 6,130,095; WO 95/22049 and patents pending.


Method Approval

The 355 SCD is the detector of choice for ASTM Standard Test Method D 5504: *Determination of Sulfur Compounds in Natural Gas and Gaseous Fuels by Gas Chromatography*

Chrom 1. Sulfur Compounds in Gasoline



and Chemiluminescence, and ASTM D 5623: *Sulfur Compounds in Light Petroleum Liquids by Gas Chromatography and Sulfur Selection Detection*, and is the original detector of choice for ASTM Standard Test Method D 7011: *Standard Test Method for Determination of Trace Thiophene in Refined Benzene by Gas Chromatography and Sulfur Selective Detection*. The SCD is the only detector tested for ASTM D 5623-95 in which data was sufficient for determining method precision. (ASTM Research Report: RR:D02-1335).



Dual Plasma™



Unsurpassed
Stability • Selectivity • Sensitivity
for a wide range of applications


Applications

Natural Gas, Petroleum, and Petrochemical Products

Gas chromatography with sulfur chemiluminescence detection provides a rapid means of identification and quantification of sulfur compounds in petroleum feeds and products. Examples include the analysis of sulfur compounds in monomers such as ethylene and propylene; solvents such as paraffins, benzene, toluene and xylenes; and fuels such as natural gas, LPG, gasoline, kerosene, jet, and diesel fuels (see *Chrom 1, Analysis of Sulfur Compounds in Gasoline*).

Analysis of trace levels of thiophene (ASTM D 7011) and carbon disulfide in refined benzene serves as a good example of the superior performance of the Dual Plasma 355 SCD. Benzene is an important aromatic solvent and chemical intermediate and component that is widely produced and used in the petrochemical and other industries.


Sulfur containing compounds are notorious for their detrimental effects as catalyst poisons, and the use of more selective catalysts generally require higher purity starting materials. Thiophene is of particular concern as a source of sulfur because its boiling point is close to that of benzene and it is a known potential contaminant.



Chrom 2. CP Wax 52 GC Column as per ASTM D 7011, 1 µL Injection Split 1:10

Short-term Precision: 1.4% RSD for 1 ppm S Thiophene, 2.6% RSD for CS₂ at 90 ppb S, N=10
Long-term Precision: 3.6% RSD for 1 ppm S Thiophene, 10.4% RSD for CS₂ at 90 ppb S, N=42 over 72 hours

Chrom 3. Volatile Sulfur Compounds in Beer



Chrom 2 below shows an isometric plot of ten repetitive injections of a refined benzene sample that contained about 90 ppb carbon disulfide and 1 ppm thiophene (each as sulfur). The chromatograms demonstrate excellent sensitivity and selectivity, and these ten injections yielded 2.6% RSD for carbon disulfide and 1.4% RSD for thiophene. Repetitive analyses made over 72 hours yielded 10.4 and 3.6% RSD respectively.

Food and Beverage Products

Sulfur compounds possess extremely low odor and taste thresholds, making them very important in flavor and fragrance chemistry. The SCD is the detector of choice for analysis of sulfur compounds which affect product odor and flavor because of its sensitive and selective detection. As an example, sulfur chemiluminescence detection is useful for identification and quantitation of sulfur compounds in beer as shown in Chrom 3.

Gaseous Sulfur Emissions

The reliable measurement of atmospheric sulfur gases is important because these compounds are notorious for their foul odors when present at only part-per-billion levels and may be toxic at higher levels. There are numerous natural and industrial sources of sulfur gases. Natural sources of emission include: animals, vegetation, soils and volcanoes. Industrial sources include: refineries, smelters, kraft paper pulpers, food processors and power generators.

Technology



SIEVERS 355 SCD **Sulfur Chemiluminescence Detector** **Specifications**

Operating Specifications

| | |
|--|--|
| Sensitivity | Typical < 0.5 pg S/second (signal to noise of 3.3:1) |
| Typical Selectivity | g S/g C > 2 x 10 ⁷ |
| Linearity | > 10 ⁴ |
| Precision and Stability* | < 2% RSD 2 hours < 5% RSD 72 hours |
| Ozone flow through the Post-Ozone Restrictor | 20 to 30 mL/min at 3-6 psig |
| Reaction Cell Pressure | 3-10 Torr RV5 Oil-Sealed Pump 5-10 Torr Dry-Piston, Oil-Free Pump |
| Typical Burner Pressure | 250-400 Torr operating |
| Typical Burner Temperature | 800 °C |
| Typical Air Flow Rate | 65 SCCM recommended 3-10 SCCM recommended with FID adapter |
| Typical Hydrogen Flow Rate | Recommended 40 SCCM |
| Time to Reach 800 °C from Room Temperature | 10 minutes typical (120 VAC 60 Hz) |
| Safety Shroud Outside Temperature | < 65 °C at 800 °C burner temperature, typically |
| Signal Output | 0-1 V or 0-10 V |

* Based on thiophene in benzene at 1 ppm mass sulfur, 1 µL injection split 1:10; 30 m, 0.32 mm ID, 1 µm thick CP Wax (n=10 for 2.5 hours; n=42 for 72 hours)

Physical Specifications

| | |
|----------------------------|--|
| Power Requirements | |
| - Detector | 115 VAC, 50/60 Hz, 350 W 100 VAC, 50/60 Hz, 350 W 220-240 VAC, 50/60 Hz, 350 W |
| - Controller | 100-120 VAC, 50/60 Hz, 200 W 220-240 VAC, 50/60 Hz, 200 W |
| Dimensions | |
| - Detector | H: 40.6 cm (16.0 in); W: 23.4 cm (9.2 in); D: 55.3 cm (21.8 in) |
| - Burner | H: 31.2 cm (12.3 in); Diameter: 10.2 cm (4.0 in) |
| - Controller | H: 12.7 cm (5.0 in); W: 24.1 cm (9.5 in); D: 31.8 cm (12.5 in) |
| - Oil-Sealed Vacuum Pump | H: 26.1 cm (10.3 in); W: 15.2 cm (6.0 in); D: 43.0 cm (16.9 in) |
| - Dry-Piston Oil-Free Pump | H: 30.0 cm (12.0 in); W: 22.9 cm (9.0 in); D: 35.6 cm (14.0 in) |
| Weight | |
| - Detector | 17.0 kg (37.5 lb) |
| - Burner | 0.9 kg (1.9 lb) |
| - Controller | 4.5 kg (9.9 lb) |
| - Oil-Sealed Vacuum Pump | 21.5 kg (47.3 lb) |
| - Dry-Piston Oil-Free Pump | 13.6 kg (29.9 lb) |
| - Certifications | CE |

Gas Requirements

Air, Hydrogen, Oxygen (optional). All gases should be Ultrahigh-purity (UHP or 99.999%) grade. Traps for all gases are highly recommended.

Maintenance Schedule (based on continuous operation)

| | |
|------------------------|----------|
| Change Vacuum Pump Oil | 30 days |
| Replace Oil Filter | 90 days |
| Replace Chemical Trap | 30 days |
| Clean Reaction Cell | 6 months |



Ionic Instruments

Dual Plasma™ Technology



SIEVERS 255 NCD Nitrogen Chemiluminescence Detector

Overview

The NCD is a nitrogen specific detector that couples to major gas chromatograph (GC) brands. The NCD produces a linear and equimolar response to nitrogen compounds (except N₂ and nitrogen species with N-N₂ bonds), allowing analysis of complex samples without interference from other compounds in the sample matrix. In addition to detecting organic nitrogen compounds, the NCD responds to ammonia, hydrazine, hydrogen cyanide, and NO. Other NO_x species can be detected provided they can be chromatographed.

Key Features

- Nitrogen specific detection for Gas Chromatography (GC) or Supercritical Fluid Chromatography (SFC)
- Picogram detection limits
- No hydrocarbon quenching
- Linear, equimolar response to organic nitrogen compounds including nitrosamines
- Responds to ammonia, hydrazine, hydrogen cyanide, and NO_x

Major Applications

The 255 NCD may be used for the following applications: chemicals, environmental samples, food and beverages, fuels, gases, pesticides and herbicides, petrochemicals, polymers, and nitrosamines.


Principle of Operation

The effluent from a gas chromatography column enters the Dual Plasma Burner. A hydrogen and oxygen plasma

in the combustion tubes converts all nitrogen compounds to nitric oxide. A catalyst is used to prevent secondary nitrogen species from being formed and to ensure that all potential interferences are destroyed prior to detection.



Nitric oxide reacts with ozone to form electronically excited nitrogen dioxide. The excited nitrogen dioxide emits light in the red and infrared region of the spectrum (600–3200 nm) when it relaxes to its ground state. The light emitted is directly proportional to the amount of nitrogen in the sample.



355 SCD and 255 NCD Accessories and Options

FID Adapter

The Flame Ionization Detector (FID) Adapter eliminates the need to split the column effluent between two detectors operating at different pressures, and the associated problems with that approach. FID adapters are available for a number of major GCs.



Dry-Piston Oil-Free Pump

The optional Dry Piston Oil-Free Pump eliminates the need for vacuum pump oil, resulting in easier maintenance and lower operating costs.



SIEVERS 255 NCD

Nitrogen Chemiluminescence Detector

Specifications

Operating Specifications

| | |
|----------------------------------|---|
| Sensitivity | < 3 pg N/second (signal to noise of 3:1) in both N and Nitrosamine modes |
| Selectivity | g N/g C > 2 x 10 ⁷ in N mode (selectivity in Nitrosamine mode is matrix dependent) |
| Linearity | > 10 ⁴ |
| Repeatability* | < 1.5% RSD 8 hours (~ the same in both N and Nitrosamine modes) < 2% RSD 18 hours (~ 3% RSD in Nitrosamine mode over 21 hours) |
| Gas Flow through Ozone Generator | 20 to 30 mL/min at 3-6 psig (inlet pressure) |
| Reaction Cell Pressure | 5-10 Torr Oil-Sealed Pump 6-12 Torr Dry-Piston, Oil-Free Pump |
| Typical Burner Pressure | 130-150 Torr operating |
| Typical Burner Temperature | 950 °C |
| Typical Hydrogen Flow Rate | 6-10 SCCM |
| Typical Oxygen Flow Rate | 8-15 SCCM |
| Signal Output | 0-1 V or 0-10 V |

* Burner Temperature 950 °C, 11 SCCM oxygen and 6 SCCM hydrogen; 25 ppm N as nitrobenzene in toluene; 0.2 µL injection on column (HP 19095-121Z), n=7 for 3 hours; n=13 for 18 hours and n=10 n-dipropyl nitrosamine in toluene at 4 µL, 0.2 µL injection on column (HP 19095-121Z)

Physical Specifications

| | |
|----------------------------|--|
| Power Requirements | |
| - Detector | 115 VAC, 50/60 Hz, 350 W 100 VAC, 50/60 Hz, 350 W 220-240 VAC, 50/60 Hz, 350 W |
| - Controller | 100-120 VAC, 50/60 Hz, 200 W 220-240 VAC, 50/60 Hz, 200 W |
| Dimensions | |
| - Detector | H: 40.6 cm (16.0 in); W: 23.4 cm (9.2 in); D: 55.3 cm (21.8 in) |
| - Burner | H: 31.2 cm (12.3 in); Diameter: 10.2 cm (4.0 in) |
| - Controller | H: 12.7 cm (5.0 in); W: 24.1 cm (9.5 in); D: 31.8 cm (12.5 in) |
| - Oil-Sealed Vacuum Pump | H: 26.1 cm (10.3 in); W: 15.2 cm (6.0 in); D: 43.0 cm (16.9 in) |
| - Dry-Piston Oil-Free Pump | H: 30.0 cm (12.0 in); W: 22.9 cm (9.0 in); D: 35.6 cm (14.0 in) |
| Weight | |
| - Detector | 17.0 kg (37.5 lb) |
| - Burner | 0.9 kg (1.9 lb) |
| - Controller | 4.5 kg (9.9 lb) |
| - Oil-Sealed Vacuum Pump | 21.5 kg (47.3 lb) |
| - Dry-Piston Oil-Free Pump | 13.6 kg (29.9 lb) |



www.ionicsinstruments.com

Ionics Instruments

Headquarters: 6060 Spine Road • Boulder, Colorado USA 80301 • tel: 303.444.2009 / 800.255.6964 • fax: 303.444.9543 • Marketing@ionicsinstruments.com

Global Offices:

Ionics UK Ltd. — Unit 3 Mercury Way, Mercury Park Estate, Trafford Park, Manchester UK M41 7LY
tel: +44 (0) 161 866 9337 • fax: +44 (0) 161 866 9630 • e-mail: office@ionicsinstruments.co.uk

Ionics Instruments France — 83, avenue de Neuilly, 94120 Fontenay-sous-Bois, France
tel: +33 (0)1 45 14 44 84 • fax: +33 (0)1 45 14 44 89 • e-mail: IBGFrance@ionicsinstruments.com

Ionics Asia Pacific Pte. Ltd. — 59 Lorong 23 Geylang • Singapore 388383 • tel: (65) 6841 4388
fax: (65) 6841 4355 • e-mail: Julian@ionics.sg

For other sales or service locations, please visit the *Contact Us* section at www.ionicsinstruments.com.

Ionics Instruments is a subsidiary of Ionics, Incorporated, headquartered in Watertown, MA (www.ionics.com).



Local Sales Contact

