



DVLS³ Simply, Smart, Sensor



ensures safety when using Hydrogen as a carrier gas in GC systems

Reasons for using H₂ as a Carrier Gas



- **Fast Analysis:**Viscosity of H2
 - Low oven temp programming
- 2. High Efficiency:
 - H2 Flattest Van Deemter curve









Reasons for using H₂ as a Carrier Gas

3. Prolonged Column Life:

Due to lower oven temperature

4. Cost effective:

3x less expensive than its helium equivalent

Bottle or generator supply





Reluctant to use H₂ as a Carrier Gas?

 Hydrogen is an Explosive Gas: Undetected gas leaks can lead to an explosion in the GC oven



LEL of hydrogen in Air is at 4%





Safety measures

- Monitor hydrogen usage
- Existing measures in most GC models
 - Safety Shutdown in case of pressure loss
 - Flow Limiting Frit EPC
 - Explosion Test: GC designed to contain parts in case of explosion <u>6890exp.avi.MPG</u>
 - Valve box or column compartment need additional measures





Sensor principle of operation (1)

- Catalytic reaction on the surface of a pellet
- Exothermal oxidation of combustable gas $(2H_2 + O_2 \rightarrow 2H_2O(g) + heat)$
- A platinum coil inside the pellet
- The Platinum coil resistance is affected by temperature change







Sensor principle of operation (2)

 A second pellet, without catalyst, is used to compensate in a Wheatstone Bridge



 Current change between resistances is translated into an electronic signal









Gas Sensitivity Hydrogen Sensor







Drift specs vs sensor temperature







Typical long term zero offset drift



Elapsed time (days)





Hardware overview of a Hydrogen Sensor for GC applications







Four Sensors - One Controller







Multiple Sensor Setup







Automatic Carrier Gas Switch

After 1% Hydrogen detection, the carrier gas supply is switched to Nitrogen, ensuring continuous gas flow through the analytical columns







Automatic Stop Signal to GC

After 1% Hydrogen detection, a running sequence will stop, preventing waste of subsequent samples and analysis time







Software Control

New optional DVLS³ Sensor software controls:

- System configuration
- User management and creation of an audit trail
- Definition of alarm settings:
 - show a red status on dashboard
 - send an alarm message by:
 - 1. e-mail
 - 2. text (SMS)
 - 3. Growl
 - 4. Prowl

³ Sensor	Software - [Dash	nboard]							
Dashboa	ard Org	anization Configura	tion User N	Management	¦ ⊨ ∏ ≓ System Configu	iration	About	Exit	F
ard Deta	ils Audit								
_									
_	Sensor Status	Sensor Name	Sensor Type	Location	Min Value	Current Value	Max Value	Unit	
	-	Sensor6	H2 H2	Application	0	53230	5000	unit2	
•		Sensor5	H ₂ H2	GC	0	0	150	unit1	





Alarm Messages

After 1% H_2 detection the sensor alerts the staff by:

- Acoustic alarm
- **Optical** alarm (blinking display)
- Red Status of software dashboard
- Text (SMS) messages
- Growl
- Prowl







17

Calibration & Maintenance

- Zero Point Calibration Using Air
- Alarm level Calibrated using Hydrogen in Air Mixture
- Yearly frequency (or according to lab practices)
- Special Silicon-free calibration hose comes with the sensor









Features

- Safe technique to constantly monitor hydrogen leaks
- Automatic alarm transmission through:
 - 1. Optical and acoustic alarms
 - 2. Red status of Dashboard
 - 3. E-mail messages
 - 4. Text (SMS) messages
 - 5. Growl
 - 6. Prowl
- Dedicated optional software controls the alarm settings per sensor & user

1.0					11.4	1			
				.	1 - 1		(i)	2	1
Dashboard	Org	ganization Configur	ation User N	Management	System Configu	ration	About	Exit	H
oard Details	Audt								
Se	ensor Status	Sensor Name	Sensor Type	Location	Min Value	Current Value	Max Value	Unit	
		Sensor6	H, H2	Application	0	53230	5000	unit2	
•		Sensor5	H, H2	GC	0	0	150	unit1	





Features

- Automatic carrier gas switching from Hydrogen to Nitrogen
- Leak detection indicates analysis is compromised
- Automatic stop signal sent to the GC
- Easy calibration and installation
- Compatible with all GC brands
- Versatile design supports up to 4 sensors:
 - 1. Hydrogen
 - 2. Hydrogen as Leak Detector
 - 3. Temperature
 - 4. Barometric Pressure
 - 5. Level (liquid) weight









Questions?



www.davinci-ls.com