

CDSolutions


APPLICATIONS INFORMATION USING ADVANCED SAMPLE HANDLING TECHNOLOGY


Thermal Desorption and Pyrolysis to Study Soils

Natural organic matter (NOM) is from decaying plant and animal material, and it serves an important role in soils and sediments. It can help plants grow by improving uptake of nutrients and minerals, and improve the environment for microorganisms integral to plant growth. Its structure can determine the fate of toxic chemicals. A good understanding of the structure and components of organic matter may help us understand its properties.

Thermal desorption and pyrolysis is a good way to analyze both semivolatile and nonvolatile components of NOM. In pyrolysis, macromolecules normally unavailable for analysis by gas chromatography, are broken apart into smaller, more volatile components amenable to GC, producing information that can be used to characterize and help us understand NOM structure and origin.

Four sources of NOM were studied. Dry soil was added to a quartz sample tube so that it was about 3/4 full, anchored between quartz wool. Each sample was heated to 2 separate temperatures; 300 and 500°C. Two sources were commercial soils: spent mushroom compost and potting soil. The other two sources were household compost and soil from a wooded county park. Resulting chromatograms are shown in Figures 1 and 2. Phenols and methoxyphenols associated with lignin were seen at 300°C. Also, acetic acid, furans, and furaldehydes associated with cellulose and polysaccharides were seen at 300°C. More phenol and methoxyphenols were present, but the larger methoxyphenols were absent at 500°C, and nitrogen-containing compounds like pyrrole, and indolizine, usually associated with proteins, emerged at 500°C.





Model 5200 Conditions

Valve Oven: 325°C
 Transfer Line: 325°C
 Temperature: 300° or 500°C
 Time: 30 sec
 Interface Final: 300°C for 3 min

Trap: VoCarb
 Trap Rest: 50°C
 Trap Final: 300°C for 5 min

GC Conditions

Carrier: Helium
 Injector: 325°C
 Split: 50:1

Column: 35% phenyl (30m X 0.25mm)
 Detector: Quadrupole MSD
 Range: 35 - 550amu

GC Program:
 Initial: 40°C for 2 minutes
 Ramp: 10°C/min.
 Final: 300°C for 9.5 minutes

**FOR MORE INFORMATION
 CONCERNING THIS APPLICATION, WE RECOMMEND THE
 FOLLOWING READING:**

P. Buurman, F. Peterse & G. Almendros Martin, European Journal of Soil Science, 58(2007) 1330-1347.
 P. Buurman & R. Roscoe, European Journal of Soil Science, 62(2011) 253-266.
 P. Leinweber, H.R. Schulten, J. Anal. Appl. Pyrolysis, 49(1999) 359-383.
 Z. Parsi, N. Hartog, T. Górecki, J. Poerschmann, J. Anal. Appl. Pyrolysis 79(2007)9-15.
 D. M. White, D. S. Garland, L.R. Beyer, K. Yoshkawa, J. Anal. Appl. Pyrolysis 71(2004)107-118.
 D. M. White, D. S. Garland, C. Ping, B. Michaelson, Cold Regions Science and Technology 38(2004)63-

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.

CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto www.cdsanalytical.com.