

# Thermal Desorption:

## Solutions for 'seeing the stars'



# Seeing the stars.....


New Delhi, December 2018



Shanghai, Good day vs Bad day

[The WHO estimates](#) that 7 million people a year die prematurely from exposure to air pollution globally, with [the World Bank calculating](#) the cost to the world economy in lost labour as \$225bn.

## Why carryout air monitoring?



Ref: [WHO, 2018](#)

Where should  
we be  
monitoring?



Ref: [WHO](#),  
2018

# What do we need to monitor?



SULFUR OXIDES (SOX)

NITROGEN OXIDES (NOX)

CARBON MONOXIDE (CO)

CARBON DIOXIDE (CO<sub>2</sub>)

VOCS (VOLATILE ORGANIC  
COMPOUNDS)

PARTICULATE MATTER (PM)

HG IN GASEOUS FORM

RADIOACTIVE POLLUTANTS

AMMONIA (NH<sub>3</sub>)

# Markes focus



EUROPEAN COMMITTEE FOR STANDARDIZATION  
COMITÉ EUROPÉEN DE NORMALISATION  
EUROPÄISCHES KOMITEE FÜR NORMUNG

SULFUR OXIDES (SOX)

NITROGEN OXIDES (NOX)

CARBON MONOXIDE (CO)

CARBON DIOXIDE (CO<sub>2</sub>)

VOCS (VOLATILE ORGANIC  
COMPOUNDS)

PARTICULATE MATTER (PM)

HG IN GASEOUS FORM

RADIOACTIVE POLLUTANTS

AMMONIA (NH<sub>3</sub>)



# Air Monitoring of VOCs

## How?

- Environmental and public health  $\Rightarrow$  regulations  $\Rightarrow$  standard methods

- ASTM D6196 (pumped and passive tube)
- ASTM D5466 (canister)



- ISO 16017-1 (pumped tube)
- ISO 16017-2 (passive tube)

### Key Organisations

- CEN 14662-1 (Pumped tube benzene monitoring)
- CEN/TS 13649 (Stationary source emissions)



- US EPA Method TO-15 (canister)
- US EPA Method TO-17 (pumped tube)
- US EPA Method 325 (passive tube)
- US EPA Guidance on ozone precursor monitoring (on-line)

# Example legislative methods

## Key compound lists

- PAMS – (Photochemical Assessment Monitoring Stations)
  - Ozone precursors
  - On-line monitoring for VVOCs
- TO-15
  - Air toxics
  - Canister sampling
  - 65-component list
- TO-17
  - Air toxics
  - Tube sampling
  - Same as TO-15
- OVOC's
  - Oxygenated VOCs becoming more of an interest
  - Compounds in list fit into both categories





# Key monitoring activities and locations

- Environmental air monitoring
  - Ambient Air
  - Industrial Air
  - Fence-line monitoring
  - Odorous emissions
- Emissions from Products and Materials
  - Consumer goods
  - Air Fresheners
  - Toys
  - Construction products effecting indoor air
  - Medical device testing



## Simultaneous analysis of PAMS, TO-15 and OVOCs from humid air:

*117 compounds with no liquid nitrogen*

*The next step for air monitoring?*



# Simultaneous analysis of PAMS, TO-15 and OVOCs from humid air: *117 compounds with no liquid nitrogen*

- Automated canister analysis using CIA Advantage-xr with GC-MS



# Method Requirements

- **Target compounds**

- 109 compounds from PAMS & TO15 compound lists
- 8 OVOCs not in PAMS or TO15 compound lists
  1. Formaldehyde
  2. Acetaldehyde
  3. Crotonaldehyde
  4. Methacrylaldehyde
  5. Butyraldehyde
  6. Benzaldehyde
  7. Pentanal
  8. m -Tolualdehyde



- **Canister sampling**

- Manual sampling method requires the analysis of samples collected in canisters.
- Canisters must be inert coated stainless steel

- **Pre-concentration**

- To reach low detection limits large sample volumes are required
- The pre-concentration device must be capable of quantitatively retaining all compounds up to 600 mL, including very volatile C2 hydrocarbons.

- **Analysis**

- Separation, identification and quantitation of all target compounds by GC or GC - MS

# Successful PAMS, TO-15 & OVOC analysis....

This method requires pre-concentration and chromatographic separation of a high number of compounds ranging widely in volatility and polarity which presents significant challenges for many instruments:

- **Quantitative retention of very volatile to volatile organic compounds in a single analysis**
  - Trapping of the full compounds list up to 600 mL for maximum sensitivity
  - Fast desorption of all compounds for sharp peaks aiding GC separation
- **Automated unattended analysis**
  - With enough sample capacity to run all night without user intervention
- **Water removal with no loss of polar compounds**
  - Allows larger sample volumes for maximum sensitivity
  - Protects GC columns and detectors from wear due to water
- **Automated internal standard addition**
  - Provides an independent check of system performance with every sample
- **Ability to sample from pressurised or unpressurised sources**
  - No sample dilution saves time, reduces human error / risks of contamination and reduces reporting limits.
  - Allows the same instrumentation to be used for on-line or canister samples
- **Trapping and separation of 117 compounds with < 60 minute cycle times**
  - Ensures sampling is started at the same time every hour.

# Preconcentration: Introducing the CIA Advantage-xr

## UNITY-xr; Efficient cryogen-free trapping

Electrically-cooled focusing trap eliminates the cost of liquid nitrogen and ensures fast sample throughput.

- Trap all 117 compounds

## Outstanding productivity

Up to 27 samples at high and low concentrations can be analysed in one unattended sequence

- Automated, unattended analysis

## Platform-neutral

Compatible with all major makes of GC and GC-MS.

## Accurate MFC controlled sampling

- Sample from pressurised or unpressurised samples



## Unparalleled analyte range

Inert, optimised flow paths allow quantitative recovery of C<sub>2</sub> to C<sub>44</sub>, including reactive and thermally labile species .

## Dry-Focus3™: Trouble-free sampling of humid air

- Remove water with no loss of polar compounds
- Utilising Kori-xr

## Precise quantitative analysis

- Automated Internal standard capability is compliant with international methods


## Overlap mode

Prepare one sample while the previous analysis is running


- Fast sample – to sample cycle times

# Whole system configuration

**A** – Secondary column flow to FID



**B** – Primary column flow to MS





**Figure 4:** Dual-column GC-MS/FID instrument operation. ➔ = Analyte flow. ➔ = Gas flow. P = Carrier gas pressure supply.

# PAMS, TO-15 & OVOC in a single analysis with no liquid cryogen


- Dual column system with Deans switch
- Cryogen free pre-concentration AND cryogen free chromatography – oven start temp +35°C
- Excellent peak shape for VVOC at FID, as well as VOC & OVOCs *across the full compound range*
- <60 minute cycle time

Total Ion Chromatogram (TIC) of 400 mL of the 10 ppb, 100% RH standard





# PAMS, TO-15 & OVOC in a single analysis with no liquid cryogen




- **Highly reproducible data**
  - < 7.5% RSD on response across 10 replicates for all compounds
  - < 2.1% RSD for internal standard compounds
- **Very stable retention times**
  - < 0.17% RSD across 16 replicates for all compounds

*Example compounds covering the polarity and volatility range of the target list:  
10 replicate analysis of 10 ppb standard at 100% RH overlay perfectly for all compounds*

# PAMS, TO-15 & OVOC in a single analysis with no liquid cryogen

## Excellent linearity at 100% relative humidity

- 1.25 to 15 ppb equivalent
- All  $R^2$  values > 0.990
- Relative response factors highly reproducible
  - % RSD of RRF  $\leq 12$  % (method limit 30%)



## Low method detection limits

- All MDLs < 200 ppt
- Average MDL ~ 50 ppt

Compounds	MDL (ppt)
Toluene	8
Vinyl acetate	72
Naphthalene	26
<i>m</i> -Tolualdehyde	70
Formaldehyde	105
Dichlorodifluoromethane	22
Isopropanol	114
Chloroethene	47
Propane	22
N-Dodecane	73
Ethene	92
Acetylene	99
Ethanol	43
1,4-Dioxane	120

# Summary

## Markes cryogen-free CIA Advantage-xr solution provides:

- Analysis of VVOCs, OVOCs and VOCs simultaneously with no liquid nitrogen
- Excellent peak shape and retention of VVOCs, OVOCs and VOCs
- Low method detection limits
- Automated sequencing of up to 27 canister/on-line samples, standards and blanks
- Superior water removal: no loss of polar compounds
- Maximum productivity with sample overlap function
- Linearity from 1.25 to 15ppb and beyond
- High precision (low %RSD)
- Negligible carryover from high concentration samples
- Long-term retention time and response stability for high confidence data.



## From outdoor to indoor...


- Potential health risks that arise from poor indoor air quality have prompted regulation of a wide range of materials
- The average American spends 93% of their time indoors...
- 86% in houses, offices, etc...
- ... and 7% in automobiles.
- The majority of items which can effect the air quality in these locations are now regulated



# Standard methods

## The evolution of materials testing


- VIAQ – **Whole cabin testing**
  - (ISO12219-1, HJ 400, New UN specification - ECE/TRANS/WP.29/GRPE/2017/1, VDA)
- Material emission testing typically carried out using:
  - Assembly testing**
    - Small chambers (ISO12219-4 & 6) and Tedlar bags ~2000L for sampling assemblies (ISO 12219-9).
  - Screening**
    - Direct desorption (VDA278)
    - Micro-Scale Chambers (ISO12219-3)
    - ~20L Tedlar Bag Sampling (ISO12219-2)



# How to comply with standard methods

Emission testing methods based on the same principles/instrumentation


Whole Air testing (ISO 16000-6, EN16516, HJ400, ISO12219-1, UN (ECE/TRANS/WP.29/GRPE/2017/))



# How to comply with standard methods


Emission testing methods based on the same principles/instrumentation

Materials testing (ISO16000-9, ISO 12219- 2,3,4,6, EN16516, ASTM D7706-11, ASTM D7859)



# Product screening for manufacturers


**Bulk emissions testing** is valuable for profiling odours and emissions, and for testing of raw materials and foods.



Samples are placed straight into the chambers.

Vapours swept from the entire sample are collected.

**Surface emissions testing:** This approach is suitable for determining area-specific emission rates from flat samples.



Sprung spacers raise planar samples to the top of the chamber.

A seal forms when the lid is closed, so only vapours released from the sample's surface are collected.

## Why screen materials?

- Predict results of longer, more expensive tests
- To show uniformity across a range
- Trouble shoot customers complaints
- Compare products with competitors
- Determine new formulations for lower emitting products
- Odour profiling




# Breath Analysis




# The detrimental health effects of Air Pollution


## - early detection saves lives



# Analytical challenges in breath analysis




# What is breath analysis?



- VOC biomarkers in breath are relevant for a wide range of diseases
- Breath VOCs can also be used to monitor exposure to harmful VOCs
- Breath enables whole body blood sampling
- Non-invasive sampling
- No specialist medical training required
- Samples can easily be stored and transported

# Monitoring chemical exposure

- VOCs can be absorbed through the lungs and skin and are often metabolised by the liver.
- These VOCs and their metabolites can be measured on exhaled breath, providing an opportunity for non-invasive biomonitoring.



# Breath sampling

## How it works

- Alveolar ('end-tidal') breath collected and transferred onto sorbent tube
  - Provides information on VOC levels in blood (Pollutants and metabolites)

### ➤ Markes' Bio-VOC



### ➤ ReCIVA breath sampler (Owlstone Medical)



©

# Contact Markes



enquiries@markes.com



+44 (0)1443 230935



www.markes.com



@MarkesInt



www.linkedin.com/company/markes-international

