


Pyrolysis-GC/MS of Switchgrass

Switchgrass, like many other plant materials, is comprised largely of cellulose and lignin. Cellulose is a glucose polymer, and when pyrolyzed produces considerable char, water, carbon dioxide and many polar organic compounds. Significant among these are levoglucosan and furans such as furanone, furancarboxaldehyde and hydroxymethyl furancarboxaldehyde. Acetic acid is also produced.

Lignin is a complex, crosslinked aromatic biopolymer which can make up to 30% of plant material. It is essentially insoluble, and is responsible for much of the heat produced when plant materials are burned. When pyrolyzed, it produces mostly substituted phenolics, includ-

ing a series of methoxyphenols and dimethoxyphenols. These compounds are seen in the pyrograms of many wood and paper products including fiberboard, kraft paper and cardboard.

The figure below shows a pyrogram of switchgrass at 650°C. Since both cellulose and lignin are present, pyrolysis products from both biopolymers are seen in the pyrogram. From the cellulose, acetic acid and levoglucosan (eluting at about 17 minutes) are prominent, along with the characteristic furans. Lignin contributes the aromatics, including the methoxyphenols and dimethoxyphenols.



Equipment

This sample was analyzed using a CDS Model 5250 Pyroprobe Autosampler, interfaced to an Agilent6890/5975B gas chromatograph/mass spectrometer.

Model 5250 Conditions

Valve Oven: 300°C
Transfer Line: 325°C
Temperature: 650°
Time: 15 seconds
Sample Purge: 0.1 min
Equilibration: 0.1 min
Post Pyro Delay 0.1 min

GC Conditions

Carrier: Helium
Injector: 350°C
Split: 50:1
Column: HP-5MS (30m X 0.25mm)
Detector: 5975B MS
Range: 35 - 550

GC Program:

Initial: 40°C for 2 minutes
Ramp: 10°C/min.
Final: 300°C

FOR MORE INFORMATION
CONCERNING THIS APPLICATION,
WE RECOMMEND THE
FOLLOWING READING:

A. A. Boateng et al., Pyrolysis of switch-grass (*Panicum virgatum*) harvested at several stages of maturity, J. Anal. Appl. Pyrolysis 75 (2006) 55-64.

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.

CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto **www. cdsanalytical.com**.