



490 Micro GC Site Preparation Checklist

Thank you for purchasing an Agilent instrument. To get you started and to assure a successful and timely installation, please refer to this specification or set of requirements.

Correct site preparation is the key first step in ensuring that your instruments and software systems operate reliably over an extended lifetime. This document is an **information guide AND checklist** prepared for you that outlines the supplies, consumables, space and utility requirements for your equipment for your site.

For additional information about our solutions, please visit our web site at <http://www.chem.agilent.com/en-US/Pages/HomePage.aspx>

Customer Responsibilities

Make sure your site meets the following **prior to the installation date using the checklist below**. For details, see specific sections within this document, including:

- the necessary **laboratory or bench space is available**.
- the **environmental conditions for the lab** as well as laboratory gases, tubing,
- the **power requirements** related to the product (e.g. **number & location** of electrical outlets)
- the **required operating supplies** necessary for the product and installation
- please consult **Other/Special Requirements** section below for other product-specific information
- This procedure can be used for the CP4900, CP4900 PRO, and the 490 PRO Micro GC.
- Please refer to the Pre-installation Manual for additional information, part number CP501289.
- If Agilent is delivering installation and familiarization services, users of the instrument should be present throughout these services; otherwise, they will miss important operational, maintenance and safety information.**

Important Customer Information

1. If you have questions or problems in providing anything described as a *Customer Responsibilities* above, please contact your local Agilent or partner support/service organization for assistance prior to delivery. In addition, Agilent and/or its partners reserve the right to reschedule the installation dependent upon the readiness of your laboratory.
2. Should your site not be ready for whatever reasons, please contact Agilent as soon as possible to re-arrange any services that have been purchased.
3. Other optional services such as additional training, operational qualification (OQ) and consultation for user-specific applications may also be provided at the time of installation when ordered with the system, but should be contracted separately.



Dimensions and Weight

Identify the laboratory bench space before your system arrives based on the table below. Pay special attention to the total height and total weight requirements for all system components you have ordered and avoid bench space with overhanging shelves. Also pay special attention to the total weight of the modules you have ordered to ensure your laboratory bench can support this weight.

Instrument Description	Weight		Height		Depth		Width	
	Kg	lbs	cm	in	cm	in	cm	in
Micro GC 2-CH	6	14	28	11	30	12	16	6.5
Micro GC 4-CH	10	22	28	11	55	21.5	16	6.5
Power supply	1	2	6.4	2.5	17.8	7	9.5	4
Field case 2-CH	16	35	38	15	41	16	30	12
Field case 4-CH (with trolley)	30	68	47	18.5	73	28.5	38	15
Chromatography Workstation (computer with monitor, approximate values)	16	35	43	17	53	21	43	17




Environmental Conditions

Operating your instrument within the recommended temperature ranges insures optimum instrument performance and lifetime.

Special Notes:

1. Performance can be affected by sources of heat & cold e.g. direct sunlight, heating/cooling from air conditioning outlets, drafts and/or vibrations.
2. The site's ambient temperature conditions must be stable for optimum performance.

Instrument Description	Operating temp range °C (C)	Operating humidity range (%)
Micro GC 2/4CH	0 to +50	0 to 95





**490 Micro GC
Site Preparation Checklist**



Power Consumption

Special Notes:

1. If a computer system is supplied with your instrument, be sure to account for those electrical outlets.

Instrument Description	Line Voltage & Frequency (V, Hz)	Maximum Power Consumption (VA)	Maximum Power Consumption (W)
Micro GC 2/4CH	100 to 240V & 50-60Hz	130	130



Required Operating Supplies by Customer

Special Notes:

1. For information on Agilent consumables, accessories and laboratory operating supplies, please visit <http://www.chem.agilent.com/en-US/Products/consumables/Pages/default.aspx>
2. Installation kit to ensure successful installation.

Item Description	Part Number	Recommended Quantity
Installation kit for Micro GC	CP742355	1
Moisture filter	CP17971	1
Oxygen filter	CP17970	1



Other/Special Requirements

External gas supply

Gas cylinder provided with a proper working two-stage pressure assembly to adjust the carrier gas pressure to 550 Kpa \pm 10 % (80 PSI \pm 10%).

Gases in gas bottles must have a minimum purity of 99.999%

Internal gas supply of the field case

In order to fill the built-in carrier gas supply tank, a separate gas cylinder must be present having a pressure amply over 12000 kPa (120 bar, 1800 Psi).

Hydrogen must NEVER be used!


Recommended Gas clean filters

TCD

Carrier gas: He, N₂, Ar

Moisture filter: CP17971

Oxygen filter : CP17970



CP-Gas Clean filters are filled with nitrogen. If you are not using nitrogen as the carrier gas, flush filters and gas lines after installation of a new filter.



490 Micro GC Site Preparation Checklist

Gas samples

- Only gas samples should be supplied to “Sample-In”!
- Type of samples: non-condensing gas
- Samples other than non-condensing gases (aerosols, particles and polymers) must be filtered in advance (External filter kit part number: CP73672 Option: Genie filter)
- Sample condition: non-condensing gas ambient to 110°C
- Sample pressure systems: between 0 and 100 kPa (15PSI)
- Sample pressure Preconcentrator EDU-Varian: Atmospheric pressure
- Outlet of sample container must fit a stainless steel capillary of 1/16" outside diameter, provided with a Swagelok® female nut.

After storing Micro GC

Follow the procedure below if you're Micro GC has been stored for a long period.

- Put carrier gas on the Micro GC
- Switch the Micro GC ON
- Switch the TCD filament(s) OFF
- Set the column(s) temperature to its maximum (depending on the column module)
- Condition the column module, preferably overnight. This will ensure you that all the water has been removed from the column module and no damage will occur to the TCD filaments



**490 Micro GC
Site Preparation Checklist**

Document Control Logs

Revision Log:

Revision	Date	Reason For Update
7.1	2-22-2011	Christine Khainovsky
7.2	2-17-2012	Gabriele Tenconi

Approval Log:

Revision	Approver(s)	Title of Approver
7.1	Don Gage	Product support manager
7.2	Don Gage	Product support manager