

Agilent Solutions for Analyzing Polycyclic Aromatic Hydrocarbons in Seafood



Jim McCurry
October 2010



PAH analytical method summary

- **NOAA NMFS-NWFSC-59, 2004 GPC + GC/MS** is the current method specified for PAH analysis in seafood. It is a complex and time-consuming method so many customers are looking for an alternative.
- **QuEChERS + GC/MS** or **GC/MS/MS PAH analyzers** are the proposed solutions from Agilent to simplify sample preparation, reduce cycle time, and provide increased sensitivity. Backflush, Multi Mode Inlet, and MS/MS are the technologies contributing to the performance and sample through-put.
- **QuEChERS + HPLC/Fluorescence** is an excellent screening tool. This Agilent solution was chosen by the US Food Emergency Response Network.



Why QuEChERS for Sample Prep?



When Compared to traditional sample prep methods:

- 25-50%+ time savings
- Reduced solvent usage: 10-15 mL/sample
- No chlorinated solvents required
- Extract multiple families of compounds with one extraction method
- Does not require advanced sample preparation experience



QuEChERS Seafood Extraction Method:



Step1:Extraction

Finfish, shellfish

↓
Weigh 3 g into 50 mL tube,
add 2 ceramic homogenizers

↓
Add 12 mL of water, vortex 30 sec

↓
Add 15 mL of ACN (1% AA), vortex 1 min

↓
Vertically shake for
1 min, centrifuge at 4000 rpm for 5 min

↓
Aspirate and transfer 8 mL of extract
To dispersive SPE (fatty sample)

Step2: d-SPE (dispersive-SPE)

Transfer 8 mL of extract from Step 1
To d-SPE (fatty sample)

↓
Vortex 1 min

↓
Centrifuge 400 rpm, 5 min

↓
Aspirate extract , filter through 0.45 um
Nylon filter, transfer to GC or HPLC Vial

↓
Step 3: Analyze
GC/MS or GC/QQQ
LC/UV/FLD or LC/QQQ
(require dilution with water 1:4 or 1:5
prior to LC)



QuEChERS: PAH Determination in Fish

Amenable to both GC and LC, UV/FLD and MS detection

Note: Solvent exchange is not required for GC analysis





Figure 5. Overlay HPLC-FLD chromatograms of the spiked fish sample containing: 1. Nap. 2. Acy. 3. Ace. 4. Flu. 5. Phe. 6. Ant. 7. Fln. 8. Pyr. 9. BaA. 10. Chr. 11. BeP. 12. BeA. 13. BkF. 14. DahA. 15. BghiP. 16. InP. The spiking level for this sample was level 1 (see Table 3). The blue portion of the chromatogram used the following excitation/emission wavelengths: 260-nm/352-nm; the red portion 260-nm/420-nm; the light blue portion: 260-nm/440-nm. For acenaphthylene, UV detection at 230-nm was used. Chromatographic conditions are shown in Table 1.

50ppb EPA PAHs extracted from Swai fish using QuEChERS
DB-5ms 20m 0.18mm 0.18µmGC/MS SIM TIC



PAH Analysis: GC/MS with Column Backflush

-- Improved reliability and speed

Oven Program

50 °C for 0.8 min

then 70 °C/min to 180 °C for 0 min; then 7 °C/min to 230 °C for 1 min

then 40 °C/min to 280 °C for 1 min; then 25 °C/min to 335 °C for 3 min

Run Time

18.25 min


MM Inlet

Mode Pulsed Splitless Temperature: 320 °C

Column DB-EUPAH, 20 m x 180 µm x 0.14 µm

Column Flow constant flow at 1 mL/min (pressure = 25.885 psi)


MSD Transfer line : 320 C MS Source: 350 C MS Quad : 200 C



Internal Std	4	2-Methylnaphthalene	15	Phenanthrene	25	Benzo[k]fluoranthene	
1	Naphthalene-d8	5	Biphenyl	16	Anthracene	29	Benzo[a]pyrene
9	Acenaphthene-d10	6	2,6-dimethylnaphthalene	17	1-methylphenanthrene	30	Perylene
14	Phenanthrene-d10	7	HMB	18	Fluoranthene	32	Dibenz[a,h]anthracene
26	Benzo[a]pyrene-d12	8	Acenaphthylene	20	Pyrene	33	Indeno[1,2,3-cd]pyrene
		10	Acenaphthene	21	Benz[a]anthracene	34	Benzo[ghi]perylene
Target Compounds	11	2,3,5-trimethylnaphtha...	22	Triphenylene			
2	Naphthalene	12	Fluorene	23	Chrysene		
3	1-methylnaphthalene	13	Dibenzothiophene	24	Benzo[b]fluoranthene		



Low detection limit with GC/MS/MS



NIST Certified Reference Material 1974b Mussel sample

Solutions from Agilent for PAH Analysis in Seafood

-- for higher productivity, better performance on Day ONE

- **QuEChERS:** *simple and fast* extraction and sample clean-up
- **GC/MS with Backflush:** using innovative Capillary Flow Technology *shortens GC cycle time* while keeping column and source cleaner longer!
- **GC/MS/MS with Backflush:** highest *sensitivity* for PAH analysis .
- **PAH Optimized Column:** *optimized PAH separation* in shortest cycle time
- **Factory setup and tested analyzers:** Reduces method development time!
- **HPLC Fluorescence** - sensitive quick screening method

