


**sensitive***Most sensitive ELSD available***simple***Easiest to use and maintain***small***Integrates easily with any HPLC*

## Leading the Way in ELSD Technology for More than 20 Years

No one has more experience designing and manufacturing ELSDs than Alltech, now a part of Grace. Our customers' needs continue to drive the evolution of our ELSD design and development programs.



## ELSD Advantages

### Universal

- Have confidence you're seeing everything in your sample
- Improve quantitative measurements
- Estimate purity more accurately

### Sensitive

- Detection limits in the low nanogram range


### Gradient Compatible

- Maintain stable baselines during steep gradients

### Alltech® ELSDs simplify analysis of difficult compounds:


- pharmaceuticals
- impurities
- phospholipids
- fatty acids
- surfactants
- nutraceuticals
- excipients
- amino acids
- triglycerides
- polymers
- carbohydrates

## ELSD Detection Process:




### 1. Nebulization

Column effluent passes through a needle and mixes with nitrogen gas to form a dispersion of droplets.



### 2. Evaporation

Droplets pass through a heated "drift tube" where the mobile phase evaporates, leaving a fine mist of dried sample particles in solvent vapor.



### 3. Detection

The sample particles pass through a cell and scatter light from a laser beam. The scattered light is detected, generating a signal.

[www.elsd.com](http://www.elsd.com)

sensitive

simple


small

## ELSD Advantages Over Other LC Detectors


	ELSD	RI	UV	MS
<b>Sensitivity</b>	●	○	●	●
<b>Gradient Capability</b>	●	○	●	●
<b>Baseline Stability</b>	●	○	●	●
<b>Solvent Interference</b>	●	○	○	●
<b>Mass Balance</b>	●	●	○	○

**CHART KEY**  
**Excellent:** ●  
**Good:** ◐  
**Poor:** ○


**ELSD improves baseline stability and detection sensitivity compared to refractive index (RI)**




**ELSD shows what may be missing from your UV chromatogram**




**Use ELSD in parallel with UV/Vis and MS for maximum structural and concentration information**



**ELSD obtains a more accurate representation of sample mass than UV/Vis**






With reduced noise and greater baseline stability, the 3300 ELSD is more **SENSITIVE** than any other ELSD available




The 3300 ELSD's near ambient temperature evaporation maximizes **SENSITIVITY FOR SEMI-VOLATILES**

The 3300 ELSD is very **REPRODUCIBLE** as demonstrated by 10 identical injections of four ginkgolide components





A Windows®-based user interface makes navigating through menus intuitive and easy. All critical operating parameters and a scrolling chromatogram trace are shown on the large, full color LCD.



A simple, step-by-step method wizard guides you through detector setup

## Select from seven different languages for the user interface:

- English
- Español
- Français
- Deutsch
- Italiano
- Japanese
- Chinese



An informative help function is available in each screen

## Additional Features:

- Control and collect data digitally through Agilent ChemStation and EZChrom chromatography data software
- Ensure GLP compliance with password protected screen and error log
- Store up to 100 methods for fast set up
- Save gas and extend optics lifetime by programming unit to enter stand by mode at end of run

coming soon

www.elsd.com

sensitive

simple

small



## Maintenance Made Simple!



Parts requiring cleaning or service are accessible from the front



Remove the nebulizer with a simple twist and pull



Clean the drift tube easily with the included custom brush



Onboard diagnostics automate the troubleshooting process

## Full IQ/OQ/PQ Protocols Included with each ELSD

We make it easy for you to integrate the ELSD into a qualified HPLC system. The 3300 ELSD is thoroughly tested before shipment and calibration and performance documents are included with each instrument. Detailed IQ/OQ/PQ procedures are included with the detector to help you perform your own qualification. We also provide the option of onsite qualification performed by one of our service engineers.




[www.elsd.com](http://www.elsd.com)

**sensitive**


**simple**

**small**


## Stable Baselines with Fast, Steep Gradients




## Use ELSD to Trigger Collection of UV Transparent Compounds in Prep Applications




## Simplify Detection of Weak Absorbing Polymers




## Easily Detect Non-Chromophoric Drug Compounds



## Detect Difficult Samples without Derivatizing



## Improve Sensitivity for Carbohydrates



Alltech® 3300 ELSD*	
DESCRIPTION	PART NO.
3300 ELSD 110V – 240V	5135834
*Patent Pending	


Services for the Alltech® 3300 ELSD	
DESCRIPTION	PART NO.
Installation and Training	INSTALL
<b>Qualification</b>	
IQ/OQ at Customer Site	Q330010
PQ at Customer Site	Q330020
IQ/OQ/PQ at Customer Site	Q330030

Our dedicated applications chemists and technical support professionals are constantly adding new applications to our vast technical library. Visit [www.ELSD.com](http://www.ELSD.com) to search and download applications or one of our applications brochures.

#### Lipids



#### Carbohydrates



#### Pharmaceuticals



Alltech® 3300 ELSD Specifications	
<b>Light Source:</b>	Laser diode with collimating optics, 650nm, 30mW output, class IIIB
<b>Detector Element:</b>	Silicon photodiode
<b>Temperature Range:</b>	Ambient to 120°C in 1°C increments
<b>Nebulizer Gas:</b>	Nitrogen preferred; Calibrated to 4.0L/min
<b>Inlet Pressure:</b>	60 – 80psig
<b>Typical Operating Range:</b>	1.0 – 2.0L/min
<b>Flow Control:</b>	Digital Mass Flow Control
<b>Mobile Phase Flow Rate:</b>	To 3.0mL/min
<b>Analog Output:</b>	1V or 10mV full scale
<b>Communications:</b>	<i>Inputs:</i> TTL/contact closure for auto zero, gas shutoff, start, and standby <i>Outputs:</i> Contact Closure for pump shutdown in error conditions, general events <i>Data Ports:</i> USB, RS232, Ethernet
<b>User Interface:</b>	Full color Microsoft Windows® based Graphical LCD with alphanumeric keypad
<b>Power Requirements:</b>	120/240V, 50/60Hz
<b>Dimensions:</b>	11.6" H x 10.3" W x 19.5" D (29.5cm H x 26.2cm W x 49.5cm D)
<b>Weight:</b>	30lbs (16Kg)
<b>IEC/FDA Classification:</b>	Class 1 Laser Product; Complies with 21 CFR Subpart J
<b>Manufacturing Standards:</b>	Meets all applicable safety and EMC certifications; CE, UL, and CSA certified
<b>Warranty:</b>	1 year parts and labor

**See how the 3300 ELSD can benefit your analysis! Visit [www.elsd.com](http://www.elsd.com) and click on "Contacts". Ask about sample analysis or a product demonstration in your lab.**

ALLTECH®, ALTIMA™, DENALI®, FLEXIT™, GROM™, JONES CHROMATOGRAPHY™, MODCOL®, PLATINUM™, PREVAIL™, PROSPHERE™, and VYDAC® are trademarks, registered in the United States and/or other countries, of Alltech Associates, Inc. DAVISIL®, GRACE® and GRACE DAVISON DISCOVERY SCIENCES™ are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. MICROSOFT WINDOWS® and WINDOWS® are trademarks, registered in the United States and/or other countries, of Microsoft Corporation. ROBITUSSIN® is a trademark, registered in the United States and/or other countries, of American Home Products Corp. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership.

Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. © Copyright 2006 Alltech Associates, Inc. All rights reserved.

The information presented herein is derived from our testing and experience. It is offered, free of charge, for your consideration, investigation and verification. Since operating conditions vary significantly, and since they are not under our control, we disclaim any and all warranties on the results which might be obtained from the use of our products. You should make no assumption that all safety or environmental protection measures are indicated or that other measures may not be required. This product(s) may be covered by patents or patents pending. Grace reserves the right to change prices and/or specifications without prior notification. Printed in the USA.

# GRACE

#### Grace Davison Discovery Sciences Headquarters

2051 Waukegan Road • Deerfield, Illinois 60015-1899

Tel: 847.948.8600 • Fax: 847.948.1078 • Email: [contact.alltech@grace.com](mailto:contact.alltech@grace.com) • Web: [www.discoverysciences.com](http://www.discoverysciences.com)