



Analysis of Chocolate Products by Dynamic Headspace

Application Note

Food & Flavor

Author:

C. Zawodny

Thermal Desorption/Dynamic Headspace can be used to effectively characterize many food raw materials as well as finished products. The volatile components are purged from the sample chamber and collected onto a trap for transfer to the GC for analysis. Quantitation of profile components can also be done. This application note will focus on the analysis of milk chocolate (raw material) and a finished chocolate product (peppermint pattie).

Samples of ground milk chocolate and peppermint pattie (~500mg) were placed into separate dynamic headspace vessels, and desorbed. The chocolate sample was desorbed at 60°C and 90°C. The peppermint pattie was desorbed at 40°C.

Figure 1 shows the total ion chromatogram profile of milk chocolate desorbed at 60°C. Note the elution of Acetic(1), Butanoic (2), Hexanoic (3), and Benzoic acids (4) respectively. Figure 2 shows the same sample heated to 90°C for 10 minutes. Acid elution continues, but note the Vanillin peak (5) at approximately 14 minutes. Figure 3 shows the chromatogram of the ground peppermint pattie. Compound elutions include Acetic Acid (1), Eucalyptol (2) Menthone (3), t-Menthone (4), Menthol (5), and Menthyl Acetate (6).




Figure 1. Milk chocolate desorbed at 60°C.

Instrument Conditions

CDS Sample Concentrator

Valve Oven: 275°C
Transfer Line: 250°C
Tube Heat: 40°C, 60°C, 90°C 10 minutes
Trap Heat: 275°C 5.00 minutes

GC/MS

Column: 30m x 0.25 mm 5% phenyl
Carrier: Helium, 50:1 split
Injector: 350°C
Program: 40°C for 2 minutes,
10°/min to 295°C (2 min)




Figure 2. Milk chocolate desorbed at 90°C.




Figure 3. Peppermint patty desorbed at 40°C.