



Analysis of Dioxins and Dioxin-like PCBs in Feed and Food

Thierry Faye, Agilent Technologies

Jef Focant, University of Liege, Belgium





GC/MS/MS Analyzer for Analysis of Dioxins and Dioxin-like PCBs in Feed and Food

Thierry Faye

Market Development Manager

Agilent Technologies



Agilent Technologies

NEW Introducing the Agilent GC/MS/MS Dioxin Analyzer


⇒ Cutting edge technology enables fast & reliable Dioxin analysis at low levels:

- New standard for GC/MS/MS sensitivity with 7010 Triple Quadrupole EI source = An instrument with sensitivity of GC Sector MS.
- MultiMode Inlet (MMI) for effective cold split less injections and more.



⇒ Developed from successful collaborations with leading Dioxin Labs in Europe:

- The Agilent platform already validated according to new regulations in Europe for both food and feed (EC 589/2014, 709/2014).
- Refer to L'Homme/Focant's publication: *Journal of Chromatography A*, 1376 (2015) 149–158.
- Custom reporting with complete calculations have been developed and automated in MassHunter.



⇒ Ready for analysis:

- Pre-configured and pre-tested at our factory so installation in your lab is fast and efficient.
- **The RTL advantage:** guarantees the exact matching of our reference method on a new instrument.
- A service engineer runs a complete check out standard so validation can begin.
- Method never needs altering even when column maintenance is performed.


Validation of GC/MS/MS confirmatory method for the European official control of levels of dioxins, furans, and dioxin-like PCBs in foodstuffs

L'Homme B., Focant J.-F.



Organic and Biological Anal. Chem.
University of Liege, Belgium







Yu-Cheng



Yusho



95% Exposure by Food Consumption



Food and Feed control...



Dioxin EU Regulation for food-feed ‘started’ with the Belgian Dioxin Crisis in 1999...

EU Commission Documents


- Council Regulation 2001/102/EC
- Council Regulation (EC) No 2375/2001
- Commission Directive 2002/69/EC
- Commission Directive 2002/70/EC
- Commission Directive 2006/13/EC
- Commission Recommendation 2006/512/EC
- Commission Regulation (EC) No 1881/2006
- Commission Regulation (EC) No 1883/2006
- Commission Regulation (EC) No 152/2009
- Commission Regulation (EU) No 1259/2011
- Commission Recommendation 2011/516/EU
- Commission Regulation (EC) No 252/2012
- Commission Regulation (EC) No 277/2012
- Commission Regulation (EC) No 278/2012
- **Commission Regulation (EC) No 589/2014**
- **Commission Regulation (EC) No 709/2014**

<http://eur-lex.europa.eu>



EU Commission Strategy (food-feed)

- ✓ Continuous monitoring
- ✓ Maximum-Action-(Target) level strategy
- ✓ Screening-Confirmatory approach
- ✓ RASFF (high capacities)

 Based on state-of-the-art methods

Analytical Methods ?

Evolutionary Guidelines

HARMONISED QUALITY CRITERIA FOR CHEMICAL AND BIOASSAYS ANALYSES OF PCDDs/PCDFs IN FEED AND FOOD

PART 1: GENERAL CONSIDERATIONS, GC/MS METHODS

Rainer Malisch¹, Bert Baumann², Peter A. Behnisch³, Richard Canady⁴, Daniel Fraisse⁵,
Peter Fürst⁶, Douglas Hayward⁴, Ronald Hoogenboom⁷, Ronald Hoogerbrugge²,
Djien Liem², Olaf Pöpke⁸, Wim Traag⁷, Thomas Wiesmüller⁹

OHC 50 (2001) 53


HARMONISED QUALITY CRITERIA FOR CHEMICAL AND BIOASSAYS ANALYSES OF PCDDs/PCDFs IN FEED AND FOOD

PART 2: GENERAL CONSIDERATIONS, BIOASSAY METHODS


Peter A. Behnisch¹, Randy Allen², Jack Anderson³, Abraham Brouwer⁴,
David J. Brown⁵, T. Colin Campbell⁶, Leo Goeyens⁷, Robert O. Harrison⁸,
Ron Hoogenboom⁹, Ilse Van Overmeire⁷, Wim Traag⁷ and Rainer Malisch¹⁰

OHC 50 (2001) 59

www.dioxin20xx.org



PBMS



- GC-IDHRMS vs CALUX
- ISO17025
- Validation @ LOQs
(@1/5th level of interest)
- ≠Upper-Lowerbound < 20%
@1pg TEQ/g fat level
- Z-scores @ PTs
- Recovery rates...

Recent EN 16215

NEN-EN 16215:2012

EUROPEAN STANDARD

EN 16215

NORME EUROPÉENNE

EUROPÄISCHE NORM

April 2012

ICS 65.120

PBMS philosophy

English Version

**Animal feeding stuffs - Determination of dioxins and dioxin-like
PCBs by GC/HRMS and of indicator PCBs by GC/HRMS**

Aliments des animaux - Dosage des dioxines, des PCB de
type dioxine et des PCB indicateurs par GC/HRMS

Futtermittel - Bestimmung von Dioxinen und dioxin-
ähnlichen PCBs mittels GC/HRMS und von Indikator-PCBs
mittels GC/HRMS

This European Standard was approved by CEN on 9 March 2012.


9.1 General Xtraction

9.2.1 Feed
10 g

10.2
B I- PFE



Clean-up & Fractionation



Analysis



EU Com. Reg.
589 & 709

GC-QQQ

EU Com. Reg.
589 & 709

Confirmatory Tool

✓ HRMS sector or MS/MS instruments

ANALYTICAL CRITERIA FOR USE OF MS/MS FOR DETERMINATION OF DIOXINS AND DIOXIN-LIKE PCBS IN FEED AND FOOD

Kotz A¹, Malisch R^{1*}, Focant J², Eppe G², Cederberg TL³, Rantakokko P⁴, Fürst P⁵, Iadadis L⁶, Lovász C⁷, Scortichini G⁸, Diletti G⁸, di Domenico A⁹, Ingelido AM⁹, Kierkegaard A¹¹

¹European Union Reference Laboratory (EU-RL) for Dioxins and Furans, Institute of Food and Food, Freiburg, Germany; ²CART, University of Liege, Belgium; ³Swedish Institute for Food and Food, Kuopio, Finland; ⁵Chemisches und Veterinäruntersuchungsamt Münster-Lippe (CVUA-MEL), Münster, Germany; ⁶NCSR Demokritos, Agricultural Office, Budapest, Hungary; ⁸Istituto Zooprofilattico Sperimentale del Molise "G. Caporale, Teramo, Italy; ⁹Istituto Superiore di Sanità (ISS), Roma, Italy; ¹⁰Wageningen, The Netherlands; ¹¹The Food and Environment Research Agency (FERA), York, United Kingdom

EU RL Working group recommendations

Organohalogen compound 74 (2012), 156-159

EU Regulation 252/2012 amended by **589/2014** (2nd of June 2014)

EU Regulation 152/2009 amended by **709/2014** (20th of June 2014)

New EU Regulation

L 164/18

EN

Official Journal of the European Union

3.6.2014

COMMISSION REGULATION (EU) No 589/2014

of 2 June 2014

laying down methods of sampling and analysis for the control of levels of dioxins, dioxin-like PCBs and non-dioxin-like PCBs in certain foodstuffs and repealing Regulation (EU) No 252/2012

(Text with EEA relevance)

27.6.2014

EN

Official Journal of the European Union

L 188/1

COMMISSION REGULATION (EU) No 709/2014


of 20 June 2014

amending Regulation (EC) No 152/2009 as regards the determination of the levels of dioxins and polychlorinated biphenyls

(Text with EEA relevance)

- (9) In addition to the gas chromatography/high resolution mass spectrometry (GC-HRMS), technical progress and developments have shown that also gas chromatography/tandem mass spectrometry (GC-MS/MS) can be used as a confirmatory method for checking compliance with the maximum level (ML). Regulation (EU) No 252/2012 should therefore be replaced by a new Regulation providing for the use of gas chromatography/tandem mass spectrometry (GC-MS/MS) as an appropriate confirmatory method for checking compliance with the maximum level.

Criteria	PCDD/Fs and DL-PCBs GC-MS/MS (589/2014)	NDL-PCBs GC-MS/MS (589/2014)
Detectable quantity	-PCDD/F upper femtogram (10^{-15} g) -NO-PCB low picogram (10^{-12} g) -MO-PCB nanogram (10^{-9} g)	NDL-PCB nanogram (10^{-9} g)
Selectivity	-Chromatographic separation of 1,2,3,4,7,8-HxCDF and 1,2,3,6,7,8-HxCDF <25% valley peak to peak	Relative RT $\pm 0.25\%$ IS vs analyte
MRM tansitions	-Monitoring 2 specific precursors with each specific product ion transition for all labeled and unlabeled analytes -Relative ion intensities max $\pm 15\%$ -Resolution MS quadrupoles = unit	-Monitoring at least 1 precursor ion and 2 product ions -Tolerance ratio $\pm 20\%$ if rel. intens. >50% Tolerance ratio $\pm 25\%$ if rel. intens. 20-50% -Resolution MS quadrupoles = unit
Blank	-Used for LOQ calculation	-Used for LOQ calculation -Blank value <30% of maximul level ML
iLOQ	-iLOQ calculated from lowest cali. point -lowest concentration point on cali. must give acceptable and consistent deviation to the average RRF -Average RRF calculated for all points -Deviation to average RRF <30%	-ditto
LOQ	-LOQ caluclated from average blank level -LOQ < 1/5 of maximum level ML -Difference ub and lb levels <20%ML	-ditto -Diff. ub and lb for sum ind-PCB @ ML <20%
Accuracy Reproducibility	-Demonstrate performances at 0.5ML, ML, 2ML -Trueness (accuracy) $\pm 20\%$ -Within-lab reproducibility (RSD) <15%	-Demonstrate performances at 0.5ML, ML, 2ML -Trueness for sum ind-PCB @ ML $\pm 30\%$ -Within-lab reproducibility (RSD) <20%
Control	-QC chart for blanks -QC charts control sample	-QC chart for blanks -QC charts control sample
Recovery	-Individual internal std in range 60-120% -Out of range OK if contribu. to TEQ<10%	-Individual internal std in range 50-120% -Out of range OK if contribu. to sum ind-PCB<10%
Measurement uncertainty	-Expanded measurement uncertainty -Coverage factor = 2 (CL=95%) -If separate determination of congeners, make sum of separate uncertainty for sum of PCDD/F and DL-PCBs	-Expanded measurement uncertainty -Coverage factor = 2 (CL=95%)




From Regulation




Full validation

Validation for vegetable oil (feed)



**Agilent 7000C series
Triple Quad GC/MS**

PTV Injections




Dioxin fraction runtime = 42min

MO-PCBs fraction runtime = 28min

Chromatographic Profile


MO and NDL PCBs



PCDDs	RT
2378-TCDD	20,72
12378-PeCDD	24,31
123478-HxCDD	27,97
123678-HxCDD	28,10
123789-HxCDD	28,46
1234678-HpCDD	32,93
OCDD	39,35
PCDFs	
2378-TCDF	20,30
12378-PeCDF	23,26
23478-PeCDF	24,06
123478-HxCDF	27,02
123678-HxCDF	27,15
234678-HxCDF	27,79
123789-HxCDF	28,94
1234678-HpCDF	31,09
1234789-HpCDF	33,92
OCDF	39,76

MO-PCBs	RT
PCB-123	18,62
PCB-118	18,74
PCB-114	19,12
PCB-105	19,66
PCB-167	21,56
PCB-156	22,51
PCB-157	22,71
PCB-189	25,76
NDL-PCBs	
PCB-28	14,19
PCB-52	14,79
PCB-101	16,81
PCB-153	19,43
PCB-138	20,46
PCB-180	23,14



PCDDs, PCDFs, NO PCBs



NO-PCBs	RT
PCB 81	17,71
PCB 77	18,02
PCB 126	20,92
PCB 169	24,17

Chromatographic Separation


❖ Selectivity



✓ 25% valley separation HxCDF

Tandem *in-space* MS

7000 MS/MS System Optimized for gas chromatography



Compound name	ISTD?	Precursor ion	MS1 resolution	Product ion	MS2 resolution	Dwell	Collision energy
PePCB 126	☐	325.9	Unit ▼	255.9	Unit ▼	75	28
PePCB 126	☐	323.9	Unit ▼	253.9	Unit ▼	75	28
2378-TCDF	☐	305.9	Unit ▼	242.9	Unit ▼	75	33
2378-TCDF	☐	303.9	Unit ▼	240.9	Unit ▼	75	33
2378-TCDD	☐	321.9	Unit ▼	258.9	Unit ▼	75	24
2378-TCDD	☐	319.9	Unit ▼	256.9	Unit ▼	75	24
13C-TCDDs	☒	333.9	Unit ▼	269.9	Unit ▼	25	24
13C-TCDDs	☒	331.9	Unit ▼	267.9	Unit ▼	25	24
13C-PePCB 126	☒	337.9	Unit ▼	267.9	Unit ▼	25	28
13C-PePCB 126	☒	335.9	Unit ▼	265.9	Unit ▼	25	28
13C6-TCDD (RS)	☒	327.9	Unit ▼	264.9	Unit ▼	25	28
13C6-TCDD (RS)	☒	325.9	Unit ▼	262.9	Unit ▼	25	28
13C-2378-TCDF	☒	317.9	Unit ▼	253.9	Unit ▼	25	33
13C-2378-TCDF	☒	315.9	Unit ▼	251.9	Unit ▼	25	33

MRM transition Ratio

❖ Quant/Qual transition ratio

EU Reg 709/2014 says:

Tolerance

- ✓ PCDD/Fs, DL-PCBs: 2 specific precursors with each specific product ions

- 2378 TCDD: $319.9 > 256.9$
 $321.9 > 258.9$  Average ratio = 0.964

$\pm 15\%$
R(quad)=unit

- ✓ NDL-PCBs: at least 1 precursor and 2 product ions


- PCB 189: $393.8 > 323.8$
 $395.8 > 325.8$  Average ratio = 0.627

$\pm 20\%$
R(quad)=unit

In MassHunter

Agilent MassHunter Quantitative Analysis - Develop PCDD-F-NOPCB - PCDD-F_NO-PCB.batch.bin																							
File Edit View Analyze Method Update Report Tools Help																							
Analyze Batch Quantitate with Blank Subtraction Layout Restore Default Layout																							
Batch Table																							
Sample: QCG_9-95D Sample Type: Sample Compound: PCB 126 ISTD: 13C-PCB 126																							
Sample				PCB 126 Method						PCB 126 Results						Qualifier (325.9 -> 255.9) Results				Qualifier (337.9 -> 267.9) Results			
		Name	Type	Level	Vol.	Multiplier	Avg. RF	RSD	RT	ISTD Conc.	RT	MI	Calc. Conc.	Final Conc.	Ft/HM	Area	Ratio	MI	Area	Ratio	MI	Area	
		QCG_9-22D	Sample		5.00	0.14	0.9978	6.686823	20.949	360.00	20.927		13.8356	1.9999	0.042	28106	97.7		27451	93.9		701685	
		QCG_9-19D	Sample		5.00	0.14	0.9978	6.686823	20.949	360.00	20.916		15.5309	2.2271	0.050	3787	98.8		3742	95.1		85249	
		QCG_9-100D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		15.3375	2.2723	0.043	1190	90.9		1083	95.6		26136	
		QCG_9-76D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		13.8269	2.0349	0.042	1046	105.3		1102	90.1		26572	
		QCG_9-72D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		14.9119	2.1659	0.039	1340	98.9		1325	98.7		32019	
		QCG_9-70D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		19.7990	2.8893	0.040	1113	86.8		966	97.7		18716	
		QCG_9-44D	Sample		5.00	0.14	0.9978	6.686823	20.949	360.00	20.916		13.8320	1.9933	0.043	4048	91.1		3689	93.7		97650	
		QCG_9-5D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.927		14.2963	2.0991	0.044	3243	95.7		3103	94.3		77703	
		QCG_9-95D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	21.016		8.9093	1.3290	0.030	1273	22.3		285	95.5		30894	
		QCG_9-85D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		14.4958	2.1531	0.041	912	116.8		1066	98.8		24460	
		QCG_9-133D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		13.4025	1.9586	0.032	2945	90.5		2666	96.2		74059	
		QCG_9-11D	Sample		5.00	0.15	0.9978	6.686823	20.949	360.00	20.916		15.0040	2.1892	0.041	2692	88.6		2385	94.5		59299	
		QCG_9-82D	Sample		5.00	0.14	0.9978	6.686823	20.949	360.00	20.916		15.1048	2.1626	0.043	3461	86.8		3003	95.8		75531	

Tolerance:
 $95.1 \pm 15\%$



Proper Estimation of LOQs

❖ LOD/LOQ

ANALYTICAL CRITERIA FOR USE OF MS/MS FOR DETERMINATION OF DIOXINS AND DIOXIN-LIKE PCBS IN FEED AND FOOD

Kotz A¹, Malisch R^{1*}, Focant J², Eppe G², Cederberg TL³, Rantakokko P⁴, Fürst P⁵, Bernsmann T⁵, Leondiadis L⁶, Lovász C⁷, Scortichini G⁸, Diletti G⁸, di Domenico A⁹, Ingelido AM⁹, Traag W¹⁰, Smith F¹¹, Fernandes A¹¹

Approach 1: The LOQ can be calculated from the signal-to-noise ratio as already defined in the current regulations.

Approach 2: As an alternative approach, if the signal-to-noise ratio does not provide reliable results due to a very low, or no discernable noise level, the limit of quantification is based on the calibration curve. The limit of quantification is then defined as the lowest concentration point on a calibration curve that gives an acceptable ($\leq 30\%$) and consistent (measured at least at the start and at the end of an analytical series of samples) deviation to the average relative response factor calculated for all points on the calibration curve in each series of samples.

For conversion of the limit of quantification from the calibration curve to the sample, the recovery of the internal standards of the respective congener and the sample intake has to be taken into account.

Organohalogen compound 74 (2012), 156-159

- ✓ iLOQ = $10 \times \text{stdev}$ (8 replicate injections – cali point)
- ✓ LOQ = blank mean + $6 \times \text{stdev}$ (12 injections – blank)

Validation for vegetable oil (feed)

❖ iLOD/iLOQ

‘Acceptable and consistent deviation to the average RRF’

Quantifier						
Name	TS	RT	Transition	Scan	Type	
2378-TCDD	2	20.753	319.9 -> 256.9...	MRM	Target	
Calibration						
Level	Conc.	Response	Enable	RF	Level RSD	
CDC01	0.0160	424	✓	0.9341	2.534717	
CDC01	0.0160	459	✓	0.9707	2.534717	
CDC01	0.0160	450	✓	0.9257	2.534717	
CDC02	0.0400	968	✓	0.9107	7.259317	
CDC02	0.0400	958	✓	0.9098	7.259317	
CDC02	0.0400	944	✓	1.0297	7.259317	
CDC03	0.0800	2596	✓	0.9122	3.487891	
CDC03	0.0800	2710	✓	0.9728	3.487891	
CDC03	0.0800	2573	✓	0.9216	3.487891	
CDC04	0.4000	13128	✓	0.9350	2.407250	
CDC04	0.4000	13004	✓	0.9408	2.407250	
CDC04	0.4000	7104	✓	0.9772	2.407250	
CDC05	0.8000	28604	✓	0.9742	1.207290	
CDC05	0.8000	28255	✓	0.9670	1.207290	
CDC05	0.8000	28548	✓	0.9901	1.207290	
CDC06	2.8000	89229	✓	0.9997	1.368236	
CDC06	2.8000	89444	✓	1.0241	1.368236	
CDC06	2.8000	88721	✓	1.0235	1.368236	

1
2
3
4
5
6

< 15%

$$Av RRF_{1-6} = 0.9622$$

$$R^2(\text{linear fit}) = 0.9960$$

> 0.9900

$$Av RRF_1 = 0.9435$$




$$\text{Deviation} = -1.9\%$$


< 30%

Validation for vegetable oil (feed)

❖ iLOD/iLOQ


Calibration curve (lowest level)


 8 replicates

 $10 \times \text{stdev}$

❖ LOQ matrix

Procedure blanks

 10 independent injections

 Average + $6 \times \text{stdev.}$

Validation for vegetable oil (feed)

❖ (i)LOQ

Compound	Lowest cali point				Blanks		
	Avg Conc. pg/μL	Std. Dev.	iLOD pg/μL	iLOQ pg/μL	in blank?	LOQ MS/MS pg/g	LOQ HRMS pg/g
PCB 77	0,2917	0,0037	0,011	0,037	yes	49,66	64,59
PCB 81	0,2935	0,003	0,009	0,030	yes	3,67	9,53
2378-TCDD	0,0152	0,0018	0,005	0,018	no	0,02	0,06
2378-TCDF	0,0162	0,001	0,003	0,010	yes	0,10	0,12
PCB 126	0,2906	0,0077	0,023	0,077	yes	1,37	1,21
23478-PeCDF	0,0154	0,0021	0,006	0,021	yes	0,08	0,08
12378-PeCDD	0,0169	0,0029	0,009	0,029	no	0,01	0,06
PCB 169	0,3062	0,0071	0,021	0,071	no	0,02	0,27
12378-PeCDF	0,0142	0,0022	0,007	0,022	yes	0,56	0,09
123478-HxCDF	0,0164	0,0016	0,005	0,016	yes	0,07	0,10
123678-HxCDF	0,0169	0,0009	0,003	0,009	yes	0,06	0,06
234678-HxCDF	0,0166	0,0007	0,002	0,007	yes	0,08	0,09
123789-HxCDF	0,0156	0,002	0,006	0,020	yes	0,18	0,10
123789-HxCDD	0,0774	0,0061	0,018	0,062	yes	0,06	0,06
123478-HxCDD	0,017	0,0022	0,007	0,022	yes	0,02	0,06
123678-HxCDD	0,0426	0,0032	0,010	0,032	yes	0,14	0,08
OCDD	3,5201	0,0465	0,140	0,465	yes	1,99	2,02
OCDF	0,0158	0,0027	0,008	0,027	yes	0,60	0,77

Accuracy & Reproducibility

Spiked materials at 0.5 ML, ML, 2 ML (6 series, 3 days)


PCDD/Fs		Average ng WHO-TEQ/kg	Stdev	RSD %	Target ng WHO-TEQ/kg	Bias %
Spike level	ML/2	0,409	0,029	7,1	0,40	2,36
	ML	0,778	0,045	5,7	0,79	-1,54
	2ML	1,600	0,035	2,2	1,58	1,30
NO-PCBs		Average ng WHO-TEQ/kg	Stdev	RSD %	Target ng WHO-TEQ/kg	Bias %
Spike level	ML/2	0,307	0,028	9,0	0,33	-7,00
	ML	0,595	0,020	3,4	0,65	-8,53
	2ML	1,256	0,021	1,6	1,30	-3,42

- ✓ Bias < 20%
- ✓ Within lab reproducibility < 15%


Accuracy & Reproducibility

Proficiency test (PT)

	Reported pg/g TEQ	Target Value pg/g TEQ	Accuracy
Material 1			
PCDD/Fs	1.10±0.20	1.01	8,8%
DL-PCBs	0.80±0.19	0.89	-10,3%
Total TEQ	1.90±0.36	1.90	-0,1%
Material 2			
PCDD/Fs	0.55±0.11	0.48	15,6%
DL-PCBs	0.82±0.21	0.85	-3,0%
Total TEQ	1.38±0.26	1.33	3,7%



Z score = 0.80




Z score = 0.59

Wim Traag, Rikilt, NL

Within Lab Reproducibility

QC pork fat



MassHunter Report Generator

D:\MassHunter\GCMS\1\data\Develop PCDD-F-NOPCB\QuantResults\Laberca-PCDD-F-NO.batch.bin								
Batch Data Path								
Analysis Time	2014-06-12 10:33	Analyst	admin					
Report Time	2014-06-12 15:30	Reporter	admin					
	2014-05-12 16:56		Laberca-PCDD-F-NO.batch.bin					
Last Calib Update		Batch						
2014-06-	2014-06-	2014-05-12T16:56:08.3478554+02:00						
12T10:33:55.4621326+02:00	12T15:30:36.0656107							
0	+02:00							
Sample Name	14.4441_diox-1	Type	Sample	Vial	62			
Vol. [µl]	5	Comment	Diox					
Compound	RT [min]	Conc [ng/ml]	TEQ Conc [ng/ml]	LOQ	Upper Bound [ng/ml]	Medium Bound [ng/ml]	Lower Bound [ng/ml]	WHO-TEF 2005
2378-TCDD	20,94	<LOQ	<LOQ	0,0200	0,02000	0,01000	0,00000	1
12378-PeCDD	24,27	<LOQ	<LOQ	0,0300	0,03000	0,01500	0,00000	1
123478-HxCDD	28,02	0,0376	0,00376	0,0200	0,00376	0,00376	0,00376	0,1
123678-HxCDD	28,15	<LOQ	<LOQ	0,1400	0,01400	0,00700	0,00000	0,1
123789-HxCDD	28,51	<LOQ	<LOQ	0,0600	0,00600	0,00300	0,00000	0,1
1234678-HpCDD	33,08	<LOQ	<LOQ	0,4400	0,00440	0,00220	0,00000	0,01
OCDD	39,46	4,9895	0,00150	1,9900	0,00150	0,00150	0,00150	0,0003
2378-TCDF	19,90	<LOQ	<LOQ	0,1000	0,01000	0,00500	0,00000	0,1
12378-PeCDF	23,26	<LOQ	<LOQ	0,5600	0,01680	0,00840	0,00000	0,03
23478-PeCDF	24,08	<LOQ	<LOQ	0,0800	0,02400	0,01200	0,00000	0,3
123478-HxCDF	27,05	<LOQ	<LOQ	0,0700	0,00700	0,00350	0,00000	0,1
123678-HxCDF	27,19	<LOQ	<LOQ	0,0600	0,00600	0,00300	0,00000	0,1
123789-HxCDF	29,06	<LOQ	<LOQ	0,1800	0,01800	0,00900	0,00000	0,1
234678-HxCDF	27,84	<LOQ	<LOQ	0,0800	0,00800	0,00400	0,00000	0,1
1234678-HpCDF	31,29	<LOQ	<LOQ	0,5500	0,00550	0,00275	0,00000	0,01
1234789-HpCDF	34,01	<LOQ	<LOQ	0,0200	0,00020	0,00010	0,00000	0,01
OCDF	39,89	<LOQ	<LOQ	0,6000	0,00018	0,00009	0,00000	0,0003
PCB 77	18,01	<LOQ	<LOQ	49,6600	0,00497	0,00248	0,00000	0,0001
PCB 81	18,01	6,7691	0,00203	3,6700	0,00203	0,00203	0,00203	0,0003
PCB 126	21,04	<LOQ	<LOQ	1,3700	0,13700	0,06850	0,00000	0,1
PCB 169	24,22	<LOQ	<LOQ	0,0700	0,00210	0,00105	0,00000	0,03
Total TEQ					0,32144	0,16436	0,00729	
Sum TEQ PCDD/F					0,17534	0,09030	0,00526	
Sum TEQ PCBs					0,14610	0,07406	0,00203	
Sum TEQ PCDD/F					0,175 ±0,032			
Sum TEQ PCBs					0,146 ±0,035			

Take Home Message

- ✓ PTV-GC/MS/MS accepted as a confirmatory tool under EU Regs
- ✓ Full validation on challenging matrix
- ✓ MS/MS, but still dioxin analyses...
- ✓ MS/MS & sectors to be properly perceived