


Analysis of impurities in fine chemical octyl-dimethyl-4-aminobenzoate using the Agilent 1290 Infinity LC and ZORBAX RRHT and RRHD 1.8 μm columns

Application Note

Fine Chemicals

Author

Michael Woodman
Agilent Technologies, Inc.
Chemical Analysis Solutions
2850 Centerville Road
Wilmington, DE 19808
USA



Abstract

The Agilent 1290 Infinity LC has significant capabilities for a wide range of HPLC and UHPLC applications. It exhibits a broader power range (that is, the combination of pressure and flow capabilities) than any other commercially available system and has the flexibility to operate with a wide range of column dimensions and particle sizes. Additionally, advanced optical design in the diode array detector allows a wide dynamic range and high sensitivity, both of which are critical in the monitoring of small impurities in fine chemicals.

Introduction

The combined benefits are demonstrated by a separation of impurities found in a sample of octyl-dimethyl-4-aminobenzoate (Figure 1). The high pressure capability of the system allows the use of methanol, as well as acetonitrile, to explore the selectivity of the two solvents. At 1.5 mL/min, using a simple 2-min gradient and a 3.0 mm \times 50 mm 1.8 μm column, the analysis time is only 3 min. The separation of the main components is shown in Figure 2.



Agilent Technologies




Figure 1
Structure of the cited compound.

The speed, resolution and flexibility of the system are further demonstrated by a separation of the sample using methanol or acetonitrile with low solvent consumption 2.1 mm id, 1.8 μ m columns. The flow rate and gradient conditions are optimized for each solvent, to produce a gradient separation with maximum pressure of approximately 850 bar, a conservative setting for the 1200-bar capability of the Agilent 1290 Infinity LC. The separation of the main components, with the two organic solvents, is shown in Figure 3a (acetonitrile, top panel) and 3b (methanol, lower panel), where the chromatograms are zoomed to the region of peaks shown from approximately 1.2-2.5 min in Figure 2.

Configuration

- G4220A 1290 Infinity Binary Pump with Integrated Vacuum Degasser
- G4226A 1290 Infinity Autosampler
- G1316C 1290 Infinity Thermostatted Column Compartment
- G4212A 1290 Infinity Diode Array Detector

Conclusion

The combined high flow and high pressure capability of the system allows one to use high efficiency columns, producing rapid separations with remarkable resolution while conserving solvent over the use of 4.6 mm id columns. Impurity detection, due to high detector sensitivity and stability, is estimated to be < 0.01%.




Figure 2
Initial separation conditions showing a need for greater resolution and selectivity. Sample: Octyl dimethyl para-aminobenzoate, 1 mg/mL. Gradient: 1.5 mL/min, 40% to 90% ACN/water over 2 minutes. Up to 460 bar on ZORBAX StableBond RRHT C18, 3 mm \times 50 mm, 1.8 μ m, 40 $^{\circ}$ C. 0.75 minute retention: 4-amino-benzoic acid; 2.1 minute retention: Octyl dimethyl para-aminobenzoate.




Figure 3
Results using ACN vs. MeOH with the same gradient slope on the 1290 Infinity LC. Sample: ODPABA working standard, 1 mg/mL. Conditions: ACN gradient 0.6 mL/min, 40% to 90% ACN/water over 7.4 minutes. Up to 850 bar on ZORBAX StableBond RRHD C18 2.1 mm \times 150 mm, 1.8 μ m, 40 $^{\circ}$ C. Methanol gradient 0.52 mL/min, 50% to 100% MeOH/water over 5.7 minutes. Up to 850 bar on ZORBAX StableBond RRHD C18 2.1 mm \times 100 mm 1.8 μ m, 40 $^{\circ}$ C. The increased selectivity of methanol allowed a shorter column to be used, decreasing run time and avoiding the use of more expensive acetonitrile mobile phase.

www.agilent.com/chem/1290

© Agilent Technologies, Inc., 2009
October 1, 2009
Publication Number 5990-4880EN



Agilent Technologies