

NOTICE: Varian, Inc. was acquired by Agilent Technologies in May 2010. This document is provided as a courtesy but is no longer kept current and thus will contain historical references to Varian. For more information, go to www.agilent.com/chem.



Application Note SI-02091

Analysis of low ppb CO and CO₂ using the Low Level CO/CO₂ Analyzer

Coen Duvekot
Varian, Inc.

Introduction

The Varian 450-GC CO/CO₂ Analyzer is a standard turnkey solution for the analysis of carbon monoxide, methane and carbon dioxide. These components are separated from various gaseous hydrocarbon matrices. The instrument is specially suited for sub ppm levels of CO, CH₄ and CO₂ in hydrocarbon matrices such as ethylene or propylene.

To further enhance the detection limit a second configuration was developed using a 10 mL sample loop and a different column set able to accommodate large sample amounts. With this setup a detection limit of low ppb levels of CO can be achieved.

Instrumentation

The Varian Low Level CO/CO₂ Analyzer is 450-GC based and includes a gas sampling valve, two high performance packed columns, a methanizer to convert CO and CO₂ into methane and a Flame Ionization Detector (FID) (Figure 1).




Figure 1. Low level CO/CO₂ analyzer configuration.

The sample is injected onto the system via a gas sampling valve equipped with a sample loop. In the first column the volatile components through carbon dioxide are separated from the higher boiling components. These higher boiling components are back flushed to vent.

Conditions

Carrier Gas: Helium, settings are factory tuned
Column Oven: 50 °C
Valve Oven: 100 °C
Detector: Methanizer 400 °C, FID 200 °C

Table 1. Valve and event table.

Time (min)	Inj + Bfl to Vent	Sample
Initial	Fill + Bfl	Off
0.01	Fill + Bfl	On
1.80	Fill + Bfl	Off
2.00	Inject	Off
3.50	Fill + Bfl	Off

Results and Discussion

In order to show the improvement in detection limit 5 ppm CO, CH₄ and CO₂ in ethylene were analyzed with the standard analyzer and with the low level analyzer




Figure 2. 5 ppm CO and CO₂ in ethylene on two analyzer configurations.

In Table 2 the Limit of Detection (LOD) and Limit of Quantification (LOQ) is calculated. The LOD uses a signal to noise ratio of 3. The LOQ uses a signal to noise ratio of 5.

Table 2. LOD and LOQ calculations.

Standard configuration						
Component	Concentration (ppm)	Peak Height	Noise	S/N	S/N=3 LOD (ppb)	S/N=5 LOQ (ppb)
CO	5	1511	14.33	105.4	142	237
CH4	5	3251	14.33	226.9	66	110
CO2	5	893.9	14.30	62.5	240	400

Low level configuration						
Component	Concentration (ppm)	Peak Height	Noise	S/N	S/N=3 LOD (ppb)	S/N=5 LOQ (ppb)
CO	5	25767.6	15.16	1699.7	9	15
CH4	5	38366.2	15.16	2530.8	6	10
CO2	5	7720.8	15.16	509.3	29	49

Table 2 clearly shows the improved LOD and LOQ of the low level configuration. For example, the LOD for CO on a standard system is 142 ppb; on the low level configuration the LOD for CO is 9 ppb. This is an improvement by a factor of 15 compared to conventional approaches. Methane and CO₂ show somewhat less improvement as the peaks appear later in the chromatogram and thus have decreased peak heights. Nevertheless, the improvement here is also clear, by a factor of 8 to 11.

Another advantage of the low level configuration is that oxygen is separated from carbon monoxide. When oxygen passes through the methanizer, small amounts of an FID detectable gas are formed. With this low level configuration, separation is complete, making low level detection of carbon monoxide possible.

Figure 3 shows the analysis of CO, CH₄ and CO₂ in a propane sample.




Figure 3. Analysis of gases in a propane sample.

Figure 3 depicts the excellent separation of oxygen and carbon monoxide. Oxygen cannot be quantified as the conversion in the methanizer is not controlled and varies depending on the contamination of the methanizer. The analysis of another ethylene sample is shown in Figure 4. Again, excellent separation of all components and analysis at low ppm level are evident.




Figure 4. Analysis of gases in an ethylene sample.

The ethylene sample in Figure 4 was analyzed many times. Repeatability data are presented in Table 3 and Figure 5.

Table 3. Repeatability figures for an ethylene sample.

Run #	Area CO	Area CH4	Area CO2
1	37146.9	82992602.7	346293.8
2	36936.7	83168836	349159.6
3	36543.9	83355460.2	343262.5
4	36521.7	83498650.9	348612
5	36379.5	83662589.5	363218.5
6	36164.3	83667561.7	366011.7
7	35816.3	83724241.2	366008.4
8	35750	83647514.4	358569.8
9	35468.7	83446499.2	352556.4
10	35211.6	83428913.7	353663.9
11	35006.6	83162642.7	357820.8
12	34761.6	83021134.6	357907.7
13	34637.7	82922018.9	353037.4
14	34442.6	82751570.7	359377.6
15	34245.6	82549036.6	364962.9
Average	35668.9	83266618.2	356030.9
Std Dev	932	363234	7279
Rsd %	2.61	0.44	2.04

A relative standard deviation of 2.61 for CO, 0.44 for CH₄ and 2.04 for CO₂, clearly shows the suitability of the low level CO/CO₂ analyzer for these types of samples.




Figure 5. Repeatability data from Table 3 with 5% variation limits.

In Figure 5 the repeatability data of CO and CO₂ peak area is plotted with 5% variation limits. All values are nicely within these limits.

Figure 6 shows the analysis of a propylene sample. Again, separation is good. In this case low ppm range concentrations of CO, CH₄ and CO₂ are revealed.




Figure 6. Analysis of gases in a propylene sample.

Conclusion

The development of a specially configured 450-GC based CO/CO₂ analyzer achieves detection limits of low ppb in the analysis of carbon monoxide, methane and carbon dioxide in various gaseous hydrocarbon mixtures.

The improvement with respect to the standard CO/CO₂ analyzer varies by a factor of 8 - 15, depending on the component of interest.

Repeatability of such a system is very good, with a relative standard deviation of approximately 2%.

*These data represent typical results.
For further information, contact your local Varian Sales Office.*

*Varian and the Varian Logo are trademarks or registered trademarks of Varian, Inc. in the U.S. and other countries.
© Varian, Inc. 2009*

Application Note SI-02091

Varian, Inc.
www.varianinc.com
North America: 800.926.3000 – 925.939.2400
Europe: *The Netherlands*: 31.118.67.1000
Asia Pacific: *Australia*: 613.9560.7133
Latin America: *Brazil*: 55.11.3238.0400

