

Analysis of Synthetic Cannabinoids in Herbal Blends by GC Tandem Mass Spectrometry

Anthony Macherone*, Tom Gluodenis, Rafael Acosta, Ron Honnold
Agilent Technologies, Santa Clara, CA, USA

SOFT 2011
Poster # P077



Introduction

Synthetic cannabinoids represent a new and growing list of scheduled drugs. As these newly introduced compounds become regulated, there is a need for confirmatory analysis. However, obtaining reference standards can be problematic and libraries of mass spectra can be difficult to come by. Unlike other controlled substances, the matrix is not a "pure" drug or simple pill/powder. The botanical matrix that acts as a substrate creates significant data interpretation challenges and therefore a library of mass spectra for GC Gas Chromatography (GC) single quadrupole was developed and validated. To aid with data interpretation an automated Deconvolution Reporting Software (DRS) library was developed. This tool uses a mathematical algorithm to extract spectra, search potential components against a known spectral library and generate a quantitative report for the matched compounds in the sample. An alternative approach is to isolate the targets of interest via tandem mass spectrometry. The process of multiple reaction monitoring (MRM) offers high specificity for each analyte and lower detection limits in a complex matrix. The data from tandem mass spectrometry is significantly easier to interpret and offers higher confidence in the results.

Experimental

JWH-018, JWH-073, JWH-200, JWH-250, CP-47,497 (C7 analogue) and HU-210 were prepared as 5 µg/mL stock solutions. The matrix was prepared by grinding approximately 500 mg of material between two 5 in. x 5 in. sheets of sandpaper until a finely divided powder was developed. Both acidic and basic liquid-liquid extractions were employed. Derivatization was achieved via commercially available BSTFA in ethyl acetate. Specific MRM transitions and collision energies were determined through a series of experiments. GC-tandem mass spectrometry facilitated chromatographic resolution and molecular specificity required to identify and quantify each of the analytes.


Experimental

| GC Run Conditions | |
|---------------------------|--|
| Column 1 | HP-5MS UI (Agilent Santa Clara, CA) |
| Injection mode | Pulsed split-less |
| Inlet temperature | 300 °C |
| Injection volume | 1 microliter |
| Carrier gas | Helium, constant flow mode, 1.2 ml/min |
| Oven program | 80 °C (hold 0.17 min) then 30 °C/min to 300 °C (hold 0.5 min) then 5 °C/min to 340 °C (hold 5 min) |
| Transfer line temperature | 325 °C |
| MS conditions | |
| Tune | Autotune |
| Gain factor | 50 |
| Acquisition parameters | Electron Impact, multiple reaction monitoring |
| Collision gas | Nitrogen, 1.5 ml/min; Helium quench gas 2.25 ml/min |
| Solvent delay | 7.0 minutes |
| MS temperatures | Source 300 °C, Quadrupoles 150 °C |

MRM Transitions and collision energies


| Compound Name | Precursor Ion | Product Ion | Collision Energy |
|---------------|---------------|-------------|------------------|
| AM-694 | 435 | 232 | 27 |
| AM-694 | 435 | 220 | 13 |
| CP-47-497-C8 | 377 | 191 | 29 |
| CP-47-497-C8 | 377 | 167 | 33 |
| HU-211 | 530 | 446 | 13 |
| HU-211 | 446 | 299 | 21 |
| JWH-015 | 327 | 310 | 10 |
| JWH-015 | 310 | 268 | 23 |
| JWH-018 | 341 | 167 | 23 |
| JWH-018 | 324 | 254 | 23 |
| JWH-073 | 327 | 167 | 23 |
| JWH-073 | 310 | 254 | 23 |
| JWH-081 | 371 | 197 | 23 |
| JWH-081 | 354 | 269 | 31 |
| JWH-122 | 338 | 268 | 23 |
| JWH-122 | 298 | 181 | 12 |
| JWH-133 | 312 | 269 | 12 |
| JWH-133 | 269 | 93 | 23 |
| JWH-200 | 384 | 100 | 23 |
| JWH-200 | 100 | 56 | 17 |
| JWH-203 | 339 | 214 | 3 |
| JWH-203 | 214 | 144 | 17 |
| JWH-250 | 335 | 214 | 3 |
| JWH-250 | 214 | 144 | 17 |
| JWH-251 | 214 | 144 | 17 |
| JWH-251 | 144 | 116 | 12 |
| JWH-398 | 375 | 201 | 23 |
| JWH-398 | 318 | 189 | 23 |
| RCS-4 | 321 | 264 | 19 |
| RCS-4 | 264 | 135 | 17 |
| RCS-8 | 254 | 158 | 13 |
| RCS-8 | 254 | 144 | 19 |
| WIN55 212-3/2 | 100 | 70 | 13 |
| WIN55 212-3/2 | 100 | 56 | 15 |

Results and Discussion




Synthetic cannabinoids fall into the three structural categories as shown above. From right to left, one class possesses a structural scaffold similar to that of tetrahydrocannabinol (THC). The second type are synthetic naphthylindole analogues and finally, the third type posses a phenylcyclohexyl moiety. As is true for THC, a common motif inherent to most synthetic cannabinoids is a short aliphatic chain known to interact with the cannabinoid CB1 and CB2 receptors.

MRM TIC 100 ng/ml standard mix



MRM total ion chromatogram for 100 ng/mL of the standard mixture. All 17 synthetic cannabinoid standards were easily found.

Example Compound: JWH-018, 100 ng/ml



The ratio of the qualifying ion to that of the quantifier is used to further confirm positive identification of each analyte in matrix in addition to retention time information and precursor / product ion pairings. In this preliminary study, calibration curves were prepared over the range of 5 to 500 ng/mL and linearity was observed over this range.



The herbal blends analyzed in this study were: EX 565, K2 Blondie, K4 Purple Haze, K3 XXX, Lunar Diamond, Zombie, and K2 Diamond.


Lack of bulk sample homogeneity require some degree of mixing to ensure representative samples are tested.

- Matrix difficult to homogenize
- Mortar and pestle
- Hand-held herb grinder
- Electrical devices
- Sandpaper




Results and Discussion

Cal Curve, JWH-018 100ng/ml – 400 ng/ml



Black dots = Calibrators
Blue triangles = QC's
R² = 0.996
Externally standardized

Cal Curve, JWH-073 100ng/ml – 400 ng/ml



Black dots = Calibrators
Blue triangles = QC's
R² = 0.996
Externally standardized

Calibration curves were then constructed over the range of 100 – 400 ng/ml and replicate injections (n = 3) were made at each level. The average correlation coefficient of linearity for the all analytes in this study (r2) was 0.99 +/-0.012. The precision, measured as %RSD, was 13%, 7%, and 6% at 100 ppb, 200 ppb, and 400 ppb, respectively and levels of quantification as determined by a signal to noise ratio >= 10, were determined to range from 1 – 100 ng/ml in the heavy botanical matrix.

All of the blends analyzed contained two or more synthetic cannabinoids and several had as many as five different components. This was confirmed by the correct ratio of the qualifying ion to that of the quantifying ion and retention time accuracy JWH-073 and JWH-018 were detected in all of the blends analyzed at concentrations ranging from 50 to 150 ppb. Notably, K2 Blondie contained JWH-073 and JWH-018 at concentrations extrapolated to be as much as 1,000-fold higher than the other blends.

Herein, representative samples of synthetic cannabinoid herbal blends were analyzed to demonstrate the applicability of a GC/MS/MS analytical approach. Compared to selected ion monitoring (SIM), multiple reaction monitoring (MRM) made possible by GC/MS/MS systems offers significantly improved selectivity and sensitivity for the detection of trace-level synthetic cannabinoids complex matrices such as herbal incense blends. Furthermore, GC/MS/MS eliminates the need for post acquisition mass spectral deconvolution with subsequent library searching for identification and quantification. The utility of triple quadrupole MS cannot be overstated. Its ability to mitigate matrix effects and improve signal-to-noise significantly increases the analysts confidence.

Conclusions

Multiple reaction monitoring in GC tandem mass spectrometry offers a unique advantage over single quadrupole selected ion monitoring (SIM) in its ability to differentiate interfering matrix components from that of the analyte. The probability of an interfering ion having the same product ion spectrum is low, thus differentiation of matrix and analyte can be achieved. Through GC tandem mass spectrometry, highly selective and sensitive analytical methods in complex matrices can be developed, assuring positive identification and lowering detection limits far beyond those achieved in single quadrupole mass spectrometry. In this preliminary study, a selective and sensitive analytical method for the analysis of six common synthetic cannabinoids was developed.

Bibliography

Macherone A, Gluodenis TJ. Confirmation and Quantification of Synthetic Cannabinoids in Herbal Incense Blends by Gas Chromatography-Tandem Quadrupole Mass Spectrometry (GC/MS/MS). Agilent Application Note 5990-8987EN (2011)

Macherone A. Synthetic Cannabinoids: the Analytical Challenges. Agilent Live Webinar August 10, 2011