

Reveleris® X2 Flash System

RevealX™ Technology Can be Used in the Impurity Isolation of Valproic Acid from Cyclodextrin During Encapsulation

Introduction

Valproic acid is one of the major antiepileptic drugs used in the treatment of different kinds of epilepsy.¹ The drug has shown to be effective in treating cancer due to its antitumor properties. Although water is usually used as the reaction medium for enzymes, organic solvents or cyclodextrins have been added to improve the solubility of poorly soluble compounds. To increase its water solubility and bioavailability, the drug can be linked chemically by physical entrapment with polymers such as cyclodextrin or to polymeric carriers such as dextran.

Flash Chromatography Conditions:

Cartridge: Reveleris® C18 12g

Flow Rate: 36mL/min

Equilibration: 3.0 min

UV1 Wavelength: 220nm

UV2 Wavelength: 254nm

Solvent A: Water

Solvent B: Acetonitrile






Figure 1: In this application, the separation of valproic acid from its entrapment-carrier, γ -cyclodextrin, has been shown using the Reveleris® flash cartridge with the RevealX™ detection technology. Using a linear gradient with acetonitrile and water mobile phase, cyclodextrin is poorly retained on a reversed phase cartridge.

Gradient Table

Step	Time (min.)	%B
1	0.0	30
2	7.0	80



Valproic acid



γ -Cyclodextrin

Flash Chromatography Conditions:

Cartridge: Reveleris® amino 12g
Flow Rate: 36mL/min
Equilibration: 6.0 min
Run Length: 12.0 min
UV1 Wavelength: 212nm
UV2 Wavelength: 254nm
Solvent A: Water with 0.1% TFA
Solvent B: Acetonitrile




Figure 2: Compared to a reversed phase cartridge separation as shown where the retention is based on the hydrophobic interaction, the water-soluble cyclodextrin has been retained and separated from the free valproic acid using an amino phase media cartridge.

Gradient Table

Step	Time (min.)	%B
1	0.0	90
2	1.0	90
3	4.0	70
4	7.0	70

Conclusion

Polymer nanoparticles allow proper drug delivery and targeting with a wide range of properties under controlled conditions.² As polymer technology plays a key role in overcoming drug delivery challenges for those undeliverable molecules, nanospheres and nanocapsules designed as drug carriers can be synthesized with drugs linked through encapsulation or conjugation for better target delivery. There needs to be a better purification technique for isolation of these compounds when the drug is poorly soluble in water while its carrier is a polar molecule. Hence such a method helps medicinal chemists isolate water based encapsulated drug from its impurity either as free or degraded and unbound compounds from the reaction mixture.

References

1. Praveen, B., Shrivastava, P., Shrivastava, S. K.; In-Vitro release and pharmacological study of synthesized valproic acid-dextran conjugate; Acta Pharmaceutica Scientia, (2009), 51, pp. 169 – 176.
2. Vauthier, C. and Bouchemal, K.; Methods for the preparation and manufacture of polymeric nanoparticles; Pharmaceutical Research, (2008), 26, No. 5, pp. 1025 – 1058.

www.discoverysciences.com

ALLTECH® and REVELERIS® are trademarks, registered in the United States and/or other countries, of Alltech Associates, Inc. REVEALX™ is a trademark of Alltech Associates, Inc. GRACE® is a trademark, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. GRACE DAVISON DISCOVERY SCIENCES™ is a trademark of W. R. Grace & Co.-Conn. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. Grace Materials Technologies is a business segment of W. R. Grace & Co.-Conn., which includes products formerly sold under the Grace Davison brand. © Copyright 2013 Alltech Associates, Inc. All rights reserved. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification.

06/2013 M363a

Grace Headquarters:

W. R. Grace & Co.-Conn.
7500 Grace Drive, Columbia, MD 21044

IN THE AMERICAS:
2051 Waukegan Road
Deerfield, IL 60015-1899
Tel: 1 847 948 8600
Fax: 1 847 948 1078
Email: discoverysciences@grace.com

IN EUROPE:
Brandstraat 12
B-9160 Lokeren, Belgium
Tel: +32 (0)9-340-65 65
Fax: +32 (0)9 340 65 60
Email: discoverysciences.EU@grace.com

IN ASIA:
19th Floor, K.Wah Center
1010 Huai Hai Zhong Road
Shanghai, 200031 PRC
Tel: 86 21 5467 4678
Fax: 86 21 5405 1500
Email: dsbiz.asia@grace.com

IN INDIA:
17 Commerce Center
Opposite Krishna Hospital, Paud Road
Pune, 411038 India
Tel: +91 20 6644 9900
Fax: +91 20 2544 1740
Email: pune@grace.com

IN AUSTRALIA/NEW ZEALAND:
2 Kerr Court
Rowville, 3178
Victoria, Australia
Tel: +61 3 9237 6100
Fax: +61 3 9237 6101
Email: discoverysciences.AU@grace.com