

Mestranol Purification During Chemical Synthesis Using Reveleris® Navigator™ Software

Reveleris® X2 Flash Chromatography System

Introduction



Oestrogen compounds, along with progestogens, are the basis of combined oral contraceptives and hormone replacement therapy (HRT).¹ Among women, they are used to supplement natural oestrogen levels, suppress androgen formation in tumor growth of cancers, along with protection against osteoporosis, heart attacks, and possibly Alzheimer's disease.

Mestranol is one such compound which acts as a pro-drug that is formed from ethynylestradiol and methylsulfate. In this separation, mestranol can be purified using a simple method optimized by the Reveleris® Navigator™ software from the other two components that may be present in the chemical reaction as impurities.²

Experimental

Step 1

Initially two separate isocratic runs at different methanol concentrations were carried out on an HPLC system using a Reveleris® C18 analytical column. The retention times for the peaks to be resolved were obtained from the run. The chromatography conditions were as listed below.



Run Conditions

Cartridge: Reveleris® C18, 150 x 4.6 mm

Mobile phase: A: Water 20% B: Methanol 80%

Flow rate: 1.0 mL/min

Detector 1: UV @ 254 nm

Detector 2: ELSD (55° C, 1.5 LPM, Gain 1)

1. Methyl Sulfate (RT: 1.46 min, 1.26 min)
2. 17 α - Ethynylestradiol (RT: 7.04 min, 1.88 min)
3. Mestranol

Run Conditions

Cartridge: Reveleris® C18, 150 x 4.6 mm

Mobile phase: A: Water 0% B: Methanol 100%

Flow rate: 1.0 mL/min

Detector 1: UV @ 254 nm

Detector 2: ELSD (55° C, 1.5 LPM, Gain 1)

Experimental

Step 2

After selecting the appropriate HPLC conditions, the values for methanol solvent as percent B mobile phase concentration and the retention times of the peak of interest were entered. The preferred flash column and optimize for purity were selected. Only the retention times of the first two eluted peaks were used for flash chromatography optimization.

FlashNAV

TLC - Silica LC - C18 LC Transfer

Enter HPLC Data and select Reveleris parameters. Press "Calculate" to generate recommended gradient profile. [Learn More](#)

HPLC Conditions:

Column: **Reveleris LC-C18 (150 x 4.6mm)**

Flow Rate: mL/min

Solvents:

A: **Water**

B: **Methanol**

Lower %B

Higher %B

Flash Conditions:

Column: **Generic C18 12g**

Flow Rate: mL/min

Optimize for:

Speed


Purity

Continue gradient to 100%

Chromatograms showing peaks at retention times: t0, 1.25, t1, 1.46, t2, 7.04 (Lower %B) and t0, 1.25, t1, 1.26, t2, 1.88 (Higher %B).

Step 3

The Reveleris® Navigator™ software calculated the recommended gradient method to optimize for purity.



Results and Discussion

Run Conditions

Cartridge: Reveleris® C18 12g

Solvent A: Water

Solvent B: Methanol

Flow rate: 30 mL/min

UV1 wavelength: 220 nm

UV2 wavelength: 280 nm

Equilibration: 3.0 min

Injection type: Liquid

Gradient Method

Step	Time (min.)	%B
1	0	64
2	0.9	64
3	6	88
4	3.5	88

1. Methyl Sulfate 2.5mg
2. 17 α - Ethynylestradiol 2.5mg
3. Mestranol 2.5mg





Figure 1: Using a traditional flash chromatography system, chromatography method optimization is dependent on the user. Reveleris® Navigator™ software optimized the method purification by offering the appropriate gradient method for high purity isolation of all the three components that were present in the sample mixture.

Results and Discussion



GC-MS Conditions

AT-5ms (30 m, 0.25 mm, 0.25 μ m)

2 μ l injection at 250° C

30:1 split, 0.8 mL/min He

Figure 2: Analysis of mestranol and ethynylestradiol fractions show 100% purity using gas chromatography-mass spectrometry (GC-MS).

Conclusion

Using the Reveleris® Navigator™ software of the Reveleris® X2 flash chromatography system, mestranol in a sample mixture has been well resolved with high purity from the other two components without any flash method development. Such an optimization tool minimized any sample loss and eliminated any additional chromatographic runs using the flash system.

References

1. Medicinal natural products; a biosynthetic approach, 3rd edition; Dewick, P.; John Wiley & Sons, Inc., Hoboken, New Jersey, 2009.
2. Lesniowski, A.; McCreary, D.; Anderson, S.; Lawrence, K.; Automating gradient method development in flash chromatography provides productivity gains and minimizes solvent usage; ACS Fall 2010 poster, Boston, MA.

Americas
Deerfield, Illinois
Tel: +1 847.948.8600
discoverysciences@
grace.com

Europe
Lokeren, Belgium
Tel: +32 (0)9.340.65.65
discoverysciences.EU@
grace.com

India
Pune
Tel: +91 20.6644.9900
pune@grace.com

China
Shanghai
Tel: +86 21.5467.4678
dsbz.asia@grace.com

Australia
Rowville
Tel: +61 3.9237.6100
discoverysciences.AU@
grace.com

GRACE® is a trademark, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. REVELERIS® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. REVEALX™ and NAVIGATOR™ are trademarks of Alltech Associates, Inc. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Grace Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Materials Technologies is a business segment of W. R. Grace & Co.-Conn, which includes products formerly sold under the Grace Davison brand. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn, and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification. © Copyright 2014 W. R. Grace & Co.-Conn. All rights reserved.

1/2014 M348a

www.discoverysciences.com

GRACE