

# Fast and Efficient Purification of a Mixture of Glycerides, Mono-, Di-, and Tristearin

Reveleris® X2 Flash Chromatography System


## Introduction

Triglycerides are fatty acids present as esters in combination with glycerol.<sup>1</sup> Tristearin, a glyceryl ester of stearic acid, is a simple triglyceride having three identical acyl chains. It is found in both plants and animals and may be used as a drug delivery vehicle for target drug compounds.

This application demonstrates purification of a mixture of glycerides (mono-, di-, and tristearin) using C18 and amino phase chemistries. The C18 phase uses hydrophobic

interaction to determine separation, whereas the amino phase uses lipophilic interaction between the stationary phase and the fatty acyl chain of the analyte to determine separation. The amino phase cartridge can be used in a normal phase or a reversed phase mode, depending on the solvents employed. When used in normal phase mode, the less polar compounds elute first, followed by the elution of the polar ones selectively.

Tristearin



## Experimental

### Run Conditions

**Cartridge:** Reveleris® C18 12g (PN: 5152103)

**Load:** 0.4% mass load on column

**Flow rate:** 30 mL/min

**Equilibration:** 5.0 min

**Solvent A:** Acetonitrile

**Solvent B:** Methylene Chloride

**Detection:**

UV 1: 210 nm

UV 2: 254 nm

ELSD

**Run time:** 11 minutes

**Cartridge:** Reveleris® Amino 12g (PN: 5157331)

**Load:** 0.4% mass load on column

**Flow rate:** 30 mL/min

**Equilibration:** 5.0 min

**Solvent A:** Acetonitrile

**Solvent B:** Methylene Chloride

**Detection:**

UV 1: 210 nm

UV 2: 254 nm

ELSD

**Run time:** 13 minutes

### Gradient Method

Step	Time (min.)	%B
1	0	10
2	8	100
3	3	100

### Gradient Method

Step	Time (min.)	%B
1	0	100
2	2	100
3	9	10
4	2	10

### References


<sup>1)</sup> Dewick, Paul., *Medicinal Natural Products: A Biosynthetic Approach*, 3rd edition, (2009), John Wiley and Sons, UK, pp 40-44.

## Results and Discussion


Using the Reveleris® C18 cartridge and a methylene chloride/acetonitrile solvent gradient, all three components were separated in less than ten minutes (fig. 1). Due to the baseline drift from the methylene chloride and the non-chromophoric nature of the compounds, the UV detector failed to detect and collect the peak fractions. The evaporative light scattering detector (ELSD), an integral part of RevealX™ detection technology, detected all three peaks and triggered fraction collection during the separation. Purification of this mixture with UV-only detection would require collecting all the fractions based on time or volume, followed by subsequent thin layer chromatography (TLC)

confirmation of each fraction. This would take significantly longer and may result in sample loss.

Using the same solvent combination with the Reveleris® amino cartridge, the selectivity of these compounds was altered by using a stronger nonpolar solvent and then gradually increasing the polarity of the mobile phase with acetonitrile. This allowed the tristearin to elute earlier, while the polar components were retained longer (fig. 2). Such a change in selectivity may enable recovery of the most important compounds from a mixture in a shorter time with high purity.



**Figure 1: Separation of mono-, di-, and triglycerides using a Reveleris® C18 cartridge.**



**Figure 2: Separation of mono-, di-, and triglycerides using a Reveleris® amino cartridge. Note how altering polarity affects elution order.**

**Compound ID:**  
**1. Monostearin**  
**2. Distearin**  
**3. Tristearin**

## Conclusion

Using the Reveleris® X2 Flash Chromatography System with selective Reveleris® cartridge phases, method development can be optimized for efficiency and timesavings. Modifying gradients to further affect elution or retention of compounds provides additional

options to separate challenging compounds. The RevealX™ detection technology in the Reveleris® X2 Flash Chromatography System helps chemists to isolate and purify lipid-based compounds that are non-chromophoric with speed and high purity.

**Americas**  
Deerfield, Illinois  
Tel: +1 847.948.8600  
discoverysciences@grace.com

**Europe**  
Lokeren, Belgium  
Tel: +32 (0)9.340.65.65  
discoverysciences.EU@grace.com

**India**  
Pune  
Tel: +91 20.6644.9900  
pune@grace.com

**China**  
Shanghai  
Tel: +86 21.5467.4678  
dsbz.asia@grace.com

**Australia**  
Rowville  
Tel: +61 3.9237.6100  
discoverysciences.AU@grace.com

GRACE® is a trademark, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. REVELERIS® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. REVEALX™ is a trademark of Alltech Associates, Inc. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Grace Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification. © Copyright 2011 W. R. Grace & Co.-Conn. All rights reserved.

11/2013 M379a

[www.discoverysciences.com](http://www.discoverysciences.com)

**GRACE**