

Reveleris® X2 Flash System

RevealX™ Technology in Anti-malarial Drug Purification in Drug Discovery

Introduction

Cinchona and its alkaloids, such as quinine, have been used for the treatment of malaria, a disease caused by protozoa.¹ A range of synthetic anti-malarial drugs have been produced as alternatives to quinine. The emergence of Plasmodium falciparum strains resistant to the synthetic drugs has resulted in reintroduction of quinine, such as the Primaquine with 8-aminoquinoline, in drug discovery.

Using the RevealX™ technology of the Reveleris® X2 flash chromatography system and the Reveleris® reversed-phase cartridge, chromophoric and non-chromophoric target compounds can be separated and isolated in small amounts present in a chemical reaction by a single chromatography run.

Flash Chromatography Conditions:

Cartridge: Reveleris® C18 media 12g

Solvent A: Water

Solvent B: Methanol


Flow Rate: 36mL/min

UV1 Wavelength: 254nm



UV2 Wavelength: 220nm

Cartridge Equilibration: 5.0 min

Injection Type: Liquid



1. 2-Amino-5-Diethylamino Pentane 20mg
2. 8-Amino-6-Methoxyquinoline HBr 1mg
3. Primaquine 10mg



Gradient Table		
Step	Time (min.)	%B
1	0.0	30
2	5.0	100
3	7.0	100

Figure 1: Using traditional flash chromatography system, non-chromophoric compounds have poor sensitivity or are undetected during purification. RevealX™ technology of the Reveleris® X2 system can isolate such compounds using the ELSD in a reaction mixture during purification.

Conclusion

Using the RevealX™ technology, separation of both chromophoric and non-chromophoric compounds can improve sample recovery and reduce post-run analysis time in the drug discovery and development process.

References

1. Medicinal natural products; a biosynthetic approach, 3rd edition; Dewick, P.; John Wiley & Sons, Inc., Hoboken, New Jersey, 2009.

ALLTECH® and REVELERIS® are trademarks, registered in the United States and/or other countries, of Alltech Associates, Inc. REVEALX™ is a trademark of Alltech Associates, Inc. GRACE® is a trademark, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. GRACE DISCOVERY SCIENCES™ is a trademark of W. R. Grace & Co.-Conn. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. © Copyright 2010 Alltech Associates, Inc. All rights reserved. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification.

11/2013 M350a

W. R. Grace & Co.-Conn.
7500 Grace Drive, Columbia, MD 21044

Grace Discovery Sciences Regional Headquarters:

IN NORTH AMERICA
2051 Waukegan Road
Deerfield, IL 60015-1899
Tel: 1 847 948 8600
Fax: 1 847.948.1078
Email: discoverysciences@grace.com

IN EUROPE:
Brandstraat 12
B-9160 Lokeren, Belgium
Tel: +32 (0)9-340-65 65
Fax: +32 (0)9 340 65 60
Email: discoverysciences.EU@grace.com

IN ASIA:
19th Floor, K.Wah Center
1010 Huai Hai Zhong Road
Shanghai, 200031 PRC
Tel: 86 21 5467 4678
Fax: 86 21 5405 1500
Email: dsbiz.asia@grace.com

IN INDIA:
17 Commerce Center
Opposite Krishna Hospital, Paud Road
Pune, 411038 India
Tel: +91 20 6644 9900
Fax: +91 20 2544 1740
Email: pune@grace.com

IN AUSTRALIA/NEW ZEALAND:
2 Kerr Court
Rowville, 3178
Victoria, Australia
Tel: +61 3 9237 6100
Fax: +61 3 9237 6101
Email: discoverysciences.AU@grace.com