

# Fast and Efficient Purification of Sterols


## Reveleris® X2 Flash Chromatography System

### Introduction


Sterols are waxy insoluble substances or lipids synthesized from acetyl coenzyme A (CoA).<sup>1</sup> The most common sterol found among humans is cholesterol; though there are other non-cholesterol sterols derived from plants (phytosterol). Cholesteryl ester is the major transport and storage form of cholesterol in lipoprotein particles and most cells. Research work in sterols is of high interest to better understand human diseases such as atherosclerosis and in development of novel drugs to treat coronary heart disease.

Separation and identification of cholesterol and cholesterol-based compounds is challenging. They lack chromophores and have poor response due to higher background interference from certain organic solvents when detected at low wavelengths. Purification of such compounds may require additional method development with other solvents or a 'collect all' approach, adding time to the purification process. This application shows the benefit of RevealX™ detection technology for purifying a mixture of non-chromophoric sterols.


1) Ethyl Acetate



2) Cholesteryl Acetate



3) Cholesterol



### Experimental

#### Run Conditions

**Cartridge:** Reveleris® C18 12g (PN: 5152103)

**Load:** 2.2% mass load on column

**Flow rate:** 30 mL/min

**Equilibration:** 5.0 min

**Solvent A:** Acetonitrile

**Solvent B:** Methylene Chloride

**Detection:**

UV 1: 210 nm

UV 2: 254 nm

ELSD

**Run time:** 8 minutes

#### Gradient Method

Step	Time (min.)	%B
1	0	30
2	8	100


#### References

<sup>1</sup> Davidson, M. H.; Toth, P. P.; Maki, K. C.; *Therapeutic Lipidology: Phytosterolemia; Humana Press, Totowa, New Jersey, 2007, pp 291-319.*


## Results and Discussion

As sterols are hydrophobic and not easily soluble in polar solvents, the use of non-polar solvents such as methylene chloride is often required. When monitoring a purification by traditional UV detection at low wavelengths, the baseline signal drifts, giving poor sensitivity for fraction collection and completely masking peak number 3 (fig. 1). This results in extra time needed for post run work-up of all of the fractions.

Using RevealX™ detection technology, both UV and evaporative light scattering detection (ELSD) can be used to detect peaks and trigger fraction collection. In this run, ethyl acetate is detected by UV, while the cholesteryl acetate and cholesterol peaks are detected by ELSD (fig. 2). As there is no baseline drift in the ELSD signal, peaks can be collected based on the signal triggered from the ELSD while monitoring them with UV. This allows for separation and collection of all three components in a single run with high purity in less than six minutes.



**Figure 1: Fractions are collected continuously when using a UV detector only. Baseline drift masks peak 3.**



**Figure 2: With RevealX™ detection technology, the peaks can be monitored by UV detector while fractions are collected based on the ELSD signal. With ELSD, peak 3 is fully baseline resolved.**

**The integrated detection of RevealX™ detection technology gives confirmation of results by both UV and ELSD in a single run.**

**Compound ID:**  
**1. Ethyl Acetate**  
**2. Cholesteryl Acetate**  
**3. Cholesterol**

## Conclusion

The RevealX™ detection technology in the Reveleris® X2 Flash Chromatography System allows chemists to isolate and purify sterol-based compounds that are non-chromophoric with speed and high purity.

The separation is fast and efficient, resulting in three pure fractions within six minutes, meaning more time can be spent on discovery and less on purification.

**Americas**  
Deerfield, Illinois  
Tel: +1 847.948.8600  
discoverysciences@grace.com

**Europe**  
Lokeren, Belgium  
Tel: +32 (0)9.340.65.65  
discoverysciences.EU@grace.com

**India**  
Pune  
Tel: +91 20.6644.9900  
pune@grace.com

**China**  
Shanghai  
Tel: +86 21.5467.4678  
dsbiz.asia@grace.com

**Australia**  
Rowville  
Tel: +61 3.9237.6100  
discoverysciences.AU@grace.com

GRACE® is a trademark, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. REVELERIS® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. REVEALX™ is a trademark of Alltech Associates, Inc. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership or status. Grace Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. W. R. Grace & Co.-Conn. and its subsidiaries can not be held responsible for any damage or injury occurring as a result of improper installation or use of its products. Grace reserves the right to change prices and/or specifications without prior notification. © Copyright 2011 W. R. Grace & Co.-Conn. All rights reserved.

11/2013 M380a

[www.discoverysciences.com](http://www.discoverysciences.com)

**GRACE**