

**ASMS 2013**

**ThP-605**

**Fast Screening and  
Accurate Quantitation of  
Drugs of Abuse in  
Bioanalysis by CE-ESI-MS**


Isabelle Kohler<sup>1</sup>, Julie Schappler<sup>1</sup>,  
Serge Rudaz<sup>1</sup>, Tim Schlabach<sup>2</sup>, Martin  
Greiner<sup>3</sup>

<sup>1</sup>School of Pharmaceutical Sciences,  
University of Geneva, University of  
Lausanne, Switzerland; <sup>2</sup>Agilent  
Technologies Inc., Santa Clara, CA,  
USA ; <sup>3</sup>Agilent Technologies R&D  
and Mktg. GmbH & Co.KG,  
Waldbronn, Germany

## Introduction

The quantitation of xenobiotics in bodily fluids is of great importance in many fields, such as clinical and forensic toxicology, therapeutic drug monitoring, metabolism studies, workplace drug testing, or doping analysis. CE represents an alternative technique to GC and LC for a large range of clinical and toxicological applications with numerous advantages, such as high separation efficiency, short analysis time, and low solvent and sample consumption. UV/Vis is the most widely used detection technique with CE configuration but suffers from a lack of sensitivity due to the narrow optical path length afforded by the internal diameter of the capillary. This lack of sensitivity, combined with the relatively low selectivity of UV/Vis detection, is considered a challenging issue for the determination of potentially low concentration xenobiotics in bodily fluids. This issue can be circumvented by using highly selective and sensitive MS detection.

In this study, CE was hyphenated for quantitative purposes to a highly sensitive Agilent 6490 Triple Quadrupole LC/MS system equipped with Jet Stream and ion funnel technologies. Electrospray ionization (ESI) is the most widespread ionization source for coupling CE with MS and was used here with the sheath-flow configuration. The sheath-flow interface is characterized by an additional make-up liquid flowing through a so-called triple-tube ESI sprayer that mixes with the CE effluent at the capillary tip, providing electrical contact at the outlet end, and the appropriate flow rate ( $\mu\text{L}$  range) and solvent conditions for ionization of the analytes. Recently, a new sprayer, which presents an adapted design compared to the standard triple-tube sprayer, has been designed and was used here for the hyphenation of CE with QqQ.



**Fig. 1:** Workflow of CE-MS sample analysis, from real or spiked sample to final validated results. Screening done by CE-TOF and Quantitation using CE-QqQ

## Experimental

### Sample Preparation

Blank pooled urine was obtained from a pool of six healthy Caucasian non-drug consumers and stored after collection at  $-20^{\circ}\text{C}$ . Before analysis, the pooled urine was defrosted at ambient temperature, centrifuged at 10,000 rpm for 5 min and filtered through a  $0.45\text{-}\mu\text{m}$  nylon filter (BGB Analytik AG, Bökten, Switzerland). Stock standard solutions of the solid analytes were prepared by dissolving each compound in MeOH to obtain a concentration of  $1\text{ mg mL}^{-1}$  and stored at  $4^{\circ}\text{C}$  until use. Blank pooled urine was spiked daily at desired concentrations. Toxicological samples were received from the Laboratory of Clinical Chemistry (Geneva Hospitals, Geneva, Switzerland) and stored at  $-20^{\circ}\text{C}$  until use. Before analysis, samples were treated in the same manner in which the blank pooled urine was treated. For cocaine (COC) and methadone (MTD) quantitation, IS were spiked at  $50\text{ ng mL}^{-1}$  before dilution and injection. Two independent analyses were performed for each sample ( $n=2$ ).

### BGE

The BGE consisted of 1 M formic acid at pH 1.8. The pH value was measured with a Seven Multi pH meter (Mettler-Toledo, Schwerzenbach, Switzerland). The BGE was prepared every four days.

### Capillary Electrophoresis

CE experiments were performed with a G7100 CE system from Agilent Technologies (Waldbronn, Germany). Separation was performed using a fused-silica capillary (BGB Analytik AG, Bökten, Switzerland) with a total length of 80 cm and an internal diameter of  $50\text{ }\mu\text{m}$ .

For screening experiments, capillaries were coated with a commercial dynamic coating (CEofix™) compatible with MS. The capillary was conditioned daily with MeOH (5 min), water (5 min), 1 M NaOH (5 min), water (5 min), CEofix initiator (0.4 min with ESI source open), CEofix accelerator (0.4 min), and BGE (10 min) at 2 bar. Prior to each sample injection, the coated capillary was rinsed at 2 bar with BGE (3 min). The post-conditioning step was performed with water at 2 bar (2 min) and CEofix accelerator (1 min).

For quantitative experiments, an uncoated bare fused-silica capillary was used. Prior to each sample injection, the capillary was rinsed at 2 bar with BGE (3 min).

For screening and quantitation, a preplug of 7%  $\text{NH}_4\text{OH}$  ( $m/v$ ) was injected at 50 mbar for 10 s before hydrodynamic (HD) sample injection at 100 mbar for 150 s followed by a postplug injection of BGE at 50 mbar for 3 s. Experiments were carried out in positive polarity mode. A constant voltage of 30 kV with an initial ramping of  $1667\text{ V s}^{-1}$  (18 s) was applied, and the capillary temperature set at  $25^{\circ}\text{C}$ .

### Mass Spectrometry

For screening experiments, the CE instrument was coupled to a 6210 TOF LC/MS (Agilent Technologies, Santa Clara, CA, USA) via a coaxial sheath flow ESI interface with a standard triple-tube sprayer from Agilent Technologies.

For quantitative experiments, the CE instrument was coupled to a 6490 Triple Quadrupole LC/MS system (Agilent Technologies) via the coaxial sheath flow interface, and an Agilent Jet Stream (AJS) source.


**For details see:** I. Kohler, et al., Highly sensitive capillary electrophoresis-mass spectrometry for rapid screening and accurate quantitation of drugs of abuse in urine, *Anal. Chim. Acta* 780 (2013) 101-109.



## Results and Discussion

Validation criterion	COC	MTD	Validation criterion	COC	MTD
<b>Trueness</b>			<b>Accuracy</b>		
Relative bias [%]					
10 ng/mL	12.0	38.5	10 ng/mL	-1.1/25.1	-25.5/102.5
25 ng/mL	1.0	1.7	25 ng/mL	-10.5/12.4	-14.7/18.0
500 ng/mL	-3.1	-1.8	500 ng/mL	-12.5/6.3	-10.1/6.5
1,000 ng/mL	0.1	0.7	1,000 ng/mL	-7.4/7.6	-6.2/7.6
<b>Precision</b>			<b>LLOQ [ng/mL]</b>		
Repeatability/intermediate precision [RSD, in %]					
10 ng/mL	5.7/5.7	21.0/27.8	10	21	
25 ng/mL	5.0/5.0	7.1/7.1			
500 ng/mL	4.1/4.1	2.8/3.6			
1,000 ng/mL	3.0/3.3	3.0/3.0			

**Table 1:** Validation results for Cocaine and Methadone by CE-QqQ. LODs as low as 2 ng/mL were reached. LLOQ were determined as 10 and 21 ng/mL respectively.



**Fig. 4:** Evaluation of the method selectivity. (A) Electropherograms obtained for COC by injecting blank pool urine (CAL 00), pooled urine spiked with d<sub>3</sub>-COC at 50 ng mL<sup>-1</sup> (CAL 0) and pooled urine spiked with COC at 25 ng mL<sup>-1</sup> and d<sub>3</sub>-COC at 50 ng mL<sup>-1</sup> (VS). (B) Electropherograms obtained for MTD by injecting blank pooled urine (CAL 00), pooled urine spiked with d<sub>3</sub>-MTD at 50 ng mL<sup>-1</sup> (CAL 0) and pooled urine spiked with MTD at 25 ng mL<sup>-1</sup> and d<sub>3</sub>-MTD at 50 ng mL<sup>-1</sup> (VS).

## Conclusions

A fast and sensitive CE-ESI-MS two-step workflow was developed for the screening of drugs of abuse in urine samples prior to their quantitation. A CE-ESI-TOF/MS method was implemented for the screening step with a pH-mediated stacking procedure, which avoided a tedious off-line sample preparation with a simple urine dilution. The higher loading capacity (more than 20%) led to an increased sensitivity while maintaining strong efficiencies. The CE-ESI-TOF/MS method allowed for LODs as low as 2 ng mL<sup>-1</sup> for a varied set of common drugs of abuse and pharmaceutical compounds. The screening step was followed by compound quantitation by CE-ESI-MS/MS with a QqQ analyzer equipped with a new ESI source and a new triple-tube sprayer design, which did not show significant differences compared with the conventional ESI source and sprayer. The quantitative procedure was fully validated for COC and MTD according to reference guidelines based on selectivity, response function, trueness, precision, and accuracy. COC analysis was found to be accurate over the range of 10–1000 ng mL<sup>-1</sup>, with accuracy included within the ±30% tolerance limits, and MTD analysis was accurate in the concentration range of 21–1000 ng mL<sup>-1</sup>. The developed two-step strategy was eventually applied to the analysis of real cases and was found to be applicable for a fast and sensitive screening as well as for accurate quantitation in urine samples.