

# Automated Preparation of Simulated Distillation Samples for ASTM Methods D2887, D7213, D7398 and D6352 using a Dual Tower 7693A and Tray System

## Application Note

Hydrocarbon Processing

### Authors

Roger L. Firor  
Agilent Technologies, Inc.  
2850 Centerville Rd  
Wilmington, DE 19808  
USA

### Abstract

A dual tower 7693A and tray system installed on the 7890A Gas Chromatograph was used for preparation of hydrocarbon calibration standards, solvent blanks, and actual petroleum samples for the purpose of analysis by simulated distillation (SimDis). The front tower is equipped with a 5 or 10  $\mu\text{L}$  syringe while the back tower is equipped with a 250 or 500  $\mu\text{L}$  syringe. A 150 sample tray with heater and mixer/barcode reader is also used. Procedures are described for sample preparation for ASTM D2887, D7213, D7398 and D6352. The Multimode Inlet, G3510, operated in a temperature programmed split mode was used for all samples.



**Agilent Technologies**

## Introduction

Sample and calibration standard preparation for various simulated distillation methods is normally a manual process requiring dilution, mixing, and heating. Many procedures use volatile toxic solvents such as carbon disulfide. ASTM method D2887 commonly uses CS<sub>2</sub> for sample dilution while D6352 may use CS<sub>2</sub> or toluene for polywax calibration standard preparation. Sample heating is required for many of these procedures. Using the automation capabilities of the 7693A tower and tray system improves lab safety as well when working with CS<sub>2</sub> and other solvents by avoiding manual handling and uncontrolled heating of mixtures.

## Experimental

For all experiments, the 7890A GC was equipped with dual 7693A towers and tray. The front tower used a standard 5 or 10 µL syringe and the rear tower was equipped with the optional large syringe carriage with either a 250 or 500 µL syringe. Sample prep procedures were done on the rear tower and sample injection occurred on the front tower. The 7890A was configured with the multimode inlet operated in temperature programmed split mode. Detection was with FID. In addition, two 7890A oven systems were used. The first configuration used the conventional air bath oven and the second used the Low Thermal Mass (LTM) system. Instrumental parameters for various configurations are listed in Table 1.

Table 1. 7890A SimDis parameters

### LTM System for D2887

LTM module	5M × 0.32 mm × 0.50 µm DB1, 5 inch format
7890A oven	300 °C isothermal
Inlet	Multimode, 270 °C (0 min) to 355 °C at 200 °C/min
Liner	Single taper with glass wool, 5183-4647
Split ratio	20:1
Pressure program (Inlet)	8 psi (0 min) – 42 psi (0.9 min) at 14 psi/min
LTM program	40 °C (0 sec) to 350 °C (30 sec) at 100 °C/min

### Standard System for D2887

Column	10M × 0.53 mm × 3.0 µm D2887
Oven	40 °C (0 min) to 350 °C (5 min) at 15 °C/min
Inlet	Multimode, G3510, 50 °C (0 min) to 330 °C (4 min) at 200 °C/min
Liner	Single taper with glass wool, 5183-4647
Split	4 to 1
Flow	3.2 psig at 40 °C, constant flow mode

### 7890A system for D7213 and D7398 (Polywax 500 calibration)

#### LTM

Column	5M × 0.53 mm × 0.15 µm DB-HT SimDis 5-inch LTM format
Oven	LTM configuration, 7890A oven 325 °C isothermal, module 40 °C (0 min) to 400 °C (30 sec) at 50 °C/min
Inlet	Multimode, 270 °C (0 min) to 400 °C (3 min) at 300 °C/min
Split ratio	4 to 1 and 10 to 1
Pressure program	2.5 psi (0 min) to 9.5 psi (1.0 min) at 1 psi/min

#### Standard Air Bath Oven

Column	5M × 0.53 mm × 0.15 µm DB-HT SimDis
Oven program	40 °C (0 min) to 400 °C (5 min) at 15 °C/min
Inlet	Multimode, 210 °C (0 min) to 400 °C (10 min) at 200 °C/min
Split ratio	4 to 1
Flow	15 mL/min, constant flow mode

### 7890A system for D6352 (Polywax 655 calibration)

Column	5M × 0.53 mm × 0.15 µm DB-HT SimDis
Oven program	40 °C (0 min) to 430 °C (5 min) at 15 °C/min
Inlet	Multimode, 250 °C (0 min) to 430 °C (hold until end of run) at 200 °C/min
Split ratio	4 to 1
Flow	16 mL/min, constant flow mode

### 7693A System

Front tower	5 or 10 µL syringe, G4513A
Back tower	250 or 500 µL syringe, G4521A syringe carriage
Tray	150 sample capacity with heater and mixer/barcode reader, G4520A
Inlet	G3510 Multimode, CO <sub>2</sub> cooled
ChemStation	B.04.01
7890A firmware	A.01.10 or greater

## Discussion

A typical sample preparation program for D2887 setup is shown in Table 2. This illustrates just one way to program preparation of the calibration standard, reference gas oil (RGO), and blank that are necessary to set up a system for routine analyses. The commands can be assembled in other ways to produce the same end result. The following vials and tray locations are used with this program.

Tray position 1	Calibration mix, 0.5 µL of C5 to C40, Agilent part number 5080-8716
Tray position 2	1 mL RGO, Agilent part number 5060-9086
Tray position 3 to 5	Empty vials with 100 µL inserts, Agilent part number 5188-6592

When the procedure is complete, vial 3 will be the prepared RGO for injection, vial 4 will be the prepared calibration mix

for injection, and vial 5 will be a CS<sub>2</sub> blank. Next, a three-line sequence is set up that starts with vial 4 (calibration mix). Vial 4 is run with the ChemStation method set with this procedure active, then vial 3 (RGO) and vial 5 (CS<sub>2</sub> blank) are run using the same method but with the prep procedure inactive (unchecked in ChemStation's 7890A Injector Program pane under edit 7890A Parameters parameters menu because these samples are already prepared from the method in the first line of the sequence table). For all three samples, the core ChemStation method performs a sample preheat at 80 °C and a sample mix at 500 rpm for 20 seconds before injection. Lastly, the calibration, prepared RGO, and blank vials are fitted with 100 µL inserts so that the solvent amounts used for the procedure are minimized. Please note that when these inserts are used, mixing should be limited to speeds of approximately 500 rpm to avoid "spilling" liquid over the top of the insert into the bottom of the 2-mL vial.

Preparation of polywax standards for the higher temperature SimDis method is always challenging due to their low solubility. Solvents such as CS<sub>2</sub> and toluene are commonly used, and

Table 2. Sample prep procedure for D2887

Sampler program steps
Move vial from front sample vial offset by -3 vial(s) to back turret position #1
Dispense 750 µL from vial Wash A3 to vial Sample 1 on the Back tower
Move vial from back turret position #1 to front sample vial offset by -3 vial(s)
Move vial from front sample vial offset by -1 vial(s) to back turret position #3
Move vial from front sample vial offset by 0 vial(s) to back turret position #2
Load 150 µL from vial Wash A1 with 0 µL airgap
Load 50 µL from vial Sample 3 with 0 µL airgap
Load 0 µL from vial Waste A1 with 0 µL airgap
Load 150 µL from vial Wash A1 with 0 µL airgap
Load 0 µL from vial Sample 2 with 0 µL airgap
Move vial from front sample vial offset by -3 vial(s) to heater
Heat vial at 80 degrees C for 300 seconds
Move vial from heater to back turret position #1
Load 5 µL from vial Sample 1 with 0 µL airgap
Load 0 µL from vial Sample 2 with 0 µL airgap
Load 150 µL from vial Wash A2 with 0 µL airgap
Wait for 1 minutes
Load 0 µL from vial Waste A3 with 0 µL airgap
Dispense 150 µL from vial Wash A3 to vial Waste A1 on the Back tower
Wash syringe in Back tower, drawing from Wash A2 dispensing into Waste B1 3 times
Move vial from back turret position #1 to front sample vial offset by -3 vial(s)
Move vial from back turret position #2 to front sample vial offset by 0 vial(s)
Move vial from front sample vial offset by -2 vial(s) to back turret position #1
Dispense 20 µL from vial Sample 1 to vial Sample 3 on the Back tower
Move vial from back turret position #3 to front sample vial offset by -1 vial(s)
Move vial from back turret position #1 to front sample vial offset by -2 vial(s)
Move vial from front sample vial offset by 1 vial(s) to back turret position #1
Wash syringe in Back tower, drawing from Wash B3 dispensing into Waste B2 3 times
Dispense 150 µL from vial Wash A3 to vial Sample 1 on the Back tower
Move vial from back turret position #1 to front sample vial offset by 1 vial(s)
Wash syringe in Back tower, drawing from Wash A1 dispensing into Waste A1 2 times
Wash syringe in Front tower, drawing from Wash A1 dispensing into Waste A1 2 times

heating of the solvent/polywax vial is required just prior to injection. This entire procedure can be automated with the 7693A tower and tray system. The basic procedure for Polywax 500 is as follows:

- Place approximately 80–100 mg of Polywax 500 in a 2-mL vial and seal
- Add 125  $\mu\text{L}$  of a C20/toluene solution to the polywax vial
- Add 1.25 mL of toluene to the polywax-C20 vial
- Mix the vial
- Heat the vial at 80  $^{\circ}\text{C}$  for 4 min
- Return to tray

- Heat one final time (3 min. typical) just prior to injection

Table 3 shows the basic prep procedure using a dual tower/tray system automating the steps shown above. The only manual step is adding the solid polywax to Vial 1. Vial 2 contains a C20/toluene mixture. Preparation of this sample could be automated as well. This procedure is applicable to D7213 SimDis and D7398 (Boiling Range Distribution of Fatty Acid Methyl Esters).

A resulting chromatogram from injection of the prepared Polywax 500 vial (vial 1) is shown in Figure 1. A symmetric distribution of the polywax fragments with good resolution to C80 can be seen.

Table 3. Preparation of Polywax 500

Sampler program steps
Wash syringe in Back tower, drawing from Wash A1 dispensing into Waste A1 2 times
Move vial from tray vial #1 to back turret position #1
Move vial from tray vial #2 to back turret position #2
Dispense 1000 $\mu\text{L}$ from vial Wash A2 to vial Sample 1 on the Back tower
Dispense 5 $\mu\text{L}$ from vial Sample 2 to vial Sample 1 on the Back tower
Move vial from back turret position #1 to mixer
Move vial from back turret position #2 to tray vial #2
Mix at 2000 rpm 2 times for 10 seconds
Move vial from mixer to heater
Heat vial at 80 degrees C for 240 seconds
Move vial from heater to tray vial #1
Wash syringe in Back tower, drawing from Wash A2 dispensing into Waste B1 3 times




Figure 1. Polywax 500 with C20 marker. Multimode inlet with 7890A oven.

The preparation program for Polywax 655 is essentially the same as shown above for Polywax 500 except that heating is extended to 6 minutes, for better dissolution. Then just prior to injection, the prepared vial is heated for another 3 minutes. In the chromatogram shown below in Figure 2, a small amount (5  $\mu$ L) of C5-C18 mix was added to the Polywax 655/toluene solution as part of the automated procedure.




Figure 2. Chromatogram of Polywax 655.

The chromatogram was produced with the multimode inlet used in temperature-programmed split mode. Good definition of polyethylene fragmented to C110 is shown in Figure 3 where the last 5 minutes of the chromatogram are enlarged to show detail. Producing this detail out to C110 is extremely difficult for most chromatographic systems. The 7890A/7693A system produces excellent results with this sample.




Figure 3. Polywax 655 to C110. Multimode inlet program: 150 °C (0 min) to 430 °C (hold until end of run) at 200 °C/min. 7890A oven: 40 °C (0 min) to 430 °C (5 min) at 15 °C/min. 3 µL injection. Solvent is toluene.

Reproducibility of the sample preparation steps is excellent as seen in Figure 4, for the dilution of a heavy vacuum gas oil sample (HVGO). The program steps that were followed to produce these chromatograms are given in Table 4. The back tower equipped with a 500-µL syringe, was used for sample preparation and the front tower with a 5-µL syringe was used

for sample injection. Carbon disulfide was used for sample dilution. This program assumes a sequence is run using vial 2. Vial 1 is the stock HVGO sample that is first prepared by adding 0.5 g of the oil to a 2-mL vial. This material is extremely viscous and cannot be drawn into a syringe. Therefore the program performs a fully automated two-stage dilution prior to injection.




Figure 4. Overlay of 11 runs of HVGO, each prepared using 7693A towers and tray.

Table 4. Preparation of HVGO for injection. CS<sub>2</sub> is used as the solvent.

Sampler program steps
Move vial from front sample vial offset by -1 vial(s) to back turret position #1
Dispense 600 µL from vial Wash A3 to vial Sample 1 on the Back tower
Move vial from back turret position #1 to heater
Move vial from front sample vial offset by 0 vial(s) to back turret position #1
Heat vial at 60 degrees C for 180 seconds
Mix at 3000 rpm 3 times for 20 seconds
Move vial from heater to back turret position #2
Dispense 250 µL from vial Sample 2 to vial Sample 1 on the Back tower
Wash syringe in Back tower, drawing from Wash A3 dispensing into Waste A1 3 times
Dispense 1000 µL from vial Wash A1 to vial Sample 1 on the Back tower
Move vial from back turret position #2 to front sample vial offset by -1 vial(s)
Move vial from back turret position #1 to mixer
Mix at 3000 rpm 4 times for 20 seconds
Move vial from mixer to front sample vial offset by 0 vial(s)

## Conclusions

Difficult sample preparation procedures that are commonly used for petroleum and fuel samples can be easily automated with the 7693A tower and tray system for the 7890A and the 6890A. The system is particularly well suited for preparation of polywax calibration samples that are used for higher temperature methods. Tasks such as mixing, solid dissolution, dilution, heating, and internal standard addition are easily accomplished.

Chromatographic performance is enhanced through use of the multimode inlet. Using standard split injection liners, good sample capacity without carryover and with minimal discrimination of wide boiling samples is seen. The inlet was used in the temperature-programmed split mode for this work. Cryo cooling was not used, however, cryo can be used optionally to shorten inlet cool down between runs if desired.

For samples that fall within the boiling point range of D2887, D7213, and D7398, the Low Thermal Mass (LTM) system can be used to shorten typical analysis cycle times by 30 to 50% [1]. The high temperature method D6352 requires the standard 7890A oven.

The sample prep procedures listed here represent just one way of accomplishing a given task. Given the commands available with the system, there are many variants that will lead to the same end result.

## Reference

1. C. Wang, R. Firor, and P. Tripp, "Fast Hydrocarbon and Sulfur Simulated Distillation Using the Agilent Low Thermal Mass (LTM) System on the 7890A GC and 355 Sulfur Chemiluminescence Detector," November 2008, Agilent Technologies publication 5990-3174EN.

## For More Information

For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2009  
Printed in the USA  
March 26, 2009  
5990-3778EN



**Agilent Technologies**