# Azura

## ▶ Degasser DG 2.1S User Manual

V6880

# HPLC

# Table of Contents

**Note:** For your own safety, read the manual and always observe the warnings and safety information on the device and in the manual!

<b>Intended Use</b> .....	<b>5</b>
Device Overview .....	5
Features .....	7
Degassing Principle .....	7
Eluents .....	8
<b>Scope of Delivery</b> .....	<b>9</b>
<b>Safety for the User</b> .....	<b>9</b>
Definition of Personal and Material Damage .....	12
Decontamination .....	13
<b>Symbols and Signs</b> .....	<b>14</b>
<b>Unpacking and Setup</b> .....	<b>14</b>
Contacting the Technical Support .....	14
Location Requirements .....	14
Unpacking .....	15
<b>Initial Startup</b> .....	<b>16</b>
Screw Fitting of the Solvent Tube .....	16
Solvent Tube Connection .....	16
Connecting Degasser and Pump .....	16
Switching On the Degasser .....	17
Purging .....	18
Leak Test .....	18
Switching Off the Degasser .....	18
Short-Term Switch-Off .....	19
Long-Term Switch-Off .....	19
<b>Troubleshooting</b> .....	<b>19</b>
Possible Problems and Rectifications .....	20


<b>Maintenance and Care</b> . . . . .	<b>21</b>
Maintenance Contract . . . . .	21
Cleaning and Caring for the Device . . . . .	21
What is to be done if the degasser leaks? . . . . .	21
Tightening Screw Fittings . . . . .	21
Preventative Measures . . . . .	22
<b>Technical Data.</b> . . . . .	<b>22</b>
<b>Accessories and Spare Parts</b> . . . . .	<b>23</b>
<b>Legal Information</b> . . . . .	<b>24</b>
Warranty Conditions . . . . .	24
Transport Damage . . . . .	24
<b>Disposal</b> . . . . .	<b>25</b>
<b>Index</b> . . . . .	<b>26</b>
<b>Declaration of Conformity</b> . . . . .	<b>27</b>

# Intended Use

**Note:** Only use the device for applications that fall within the range of the intended use. Otherwise, the protective and safety equipment of the device could fail.

## Device Overview


The degasser is a device for degassing eluents. On the front of the device, there are 2 analytical degasser chambers with inlet and outlet side, one bicolor LED and the exhaust port. Air and possibly solvent vapors will be pumped out through the exhaust port. If necessary, it can be combined with a fume hood or similar extraction devices.

The ports for ground and mains connection are located on the rear side of the degasser.

**Legend**

- |                 |   |
|-----------------|---|
| ① CE mark       | ⑤ Ground  |
| ② Serial number | ⑥ Mains connection for power supply with power switch |
| ③ Warning 1     |   |
| ④ Warning 2     |   |

**Fig. 2 Analytical 2-channel degasser, rear view***Operating range*

Gases that have dissolved into the solvent can cause bubbles in pumps and detectors. Good chromatographic separation therefore requires degassing the solvent. The analytical 2-channel degasser DG 2.1S is equipped with two degassing chambers and can thus simultaneously degas two solvents.

The device can be used in the following areas:

- Biochemistry analyses
- Chemical analyses
- Food analyses
- Pharmaceutical analyses
- Environmental analyses
- Chiral analyses

The degasser is, e. g., used at universities, research institutions, and routine laboratories.

## Features

- Reliable and convenient degassing of up to two eluents at the same time.
- Teflon AF<sup>®</sup> allows high gas diffusion rates
- Constant performance vacuum pump for high baseline stability
- Chamber volumes lower than 0.5 ml, thus fast eluent changes are possible
- Very good physical and chemical stability

## Degassing Principle

While the solvent flows through the short Teflon AF<sup>®</sup> tubes, the dissolved gases are removed from the solvent by means of a vacuum on the other side of the membrane. An integrated vacuum pump maintains a constant vacuum.

- Structural properties of the Teflon AF<sup>®</sup> membrane accelerate the degassing process.
- Thanks to the molecular structure of the Teflon AF<sup>®</sup> membrane, only a very small surface is required for degassing the solvent. Thus the risk of carryover during a solvent change is reduced.
- The solvophobic and hydrophobic properties of Teflon AF<sup>®</sup> reduces the risk of contamination when exchanging the solvent.

The special pump design makes it possible to flush the pump head with a small blast of air, to remove solvent fumes, which could escape from the vacuum chamber.


Fig. 3 Degassing principle

## Eluents

Even small quantities of other substances, such as additives, modifiers, or salts can influence the durability of the materials.

**Note:** The list of selected solvents was compiled based on research in the pertinent literature and is only a recommendation by the manufacturer. If there is any doubt, contact the technical support of the manufacturer.

### Not suitable eluents

- Mineral and organic acids (except in buffer solutions)
- Bases (except in buffer solutions)
- Eluents containing particles
- Perfluorinated eluents, e.g. Fluorinert® FC-75, FC-40
- Perfluorinated polyether, e.g. Fomblin®
- Halogenated hydrocarbons, e.g. Freon®
- Hexane/heptane

### Less suitable eluents

- Methylene chloride
- Tetrahydrofuran (THF)
- Dimethyl sulfoxide (DMSO)
- Slightly volatile eluents

**Suitable eluents**

- Acetone
- Acetonitrile
- Ammonia
- Benzene
- Chloroform
- Acetic acid (10 – 50 %), at 25 °C
- Ethyl acetate
- Ethanol
- Isopropanol
- Carbon dioxide (liquid 99.999% CO<sub>2</sub>)
- Methanol
- Sodium hydroxide (1 M)
- Phosphate buffer solutions (0.5 M)
- Phosphoric acid
- Toluol
- Water

## Scope of Delivery

**Note:** Only use original parts and accessories made by the manufacturer or a company authorized by the manufacturer.

Degasser DG 2.1S	<input type="checkbox"/>
User manual	<input type="checkbox"/>
Power unit, 24 V , 60 W	<input type="checkbox"/>
Accessory kit AZURA small devices <sup>1</sup>	<input type="checkbox"/>
Accessory kit DG 2.1S	<input type="checkbox"/>

## Safety for the User

### *Professional group*

The user manual addresses persons who are qualified as chemical laboratory technicians or have completed comparable vocational training.

The following knowledge is required:

- Fundamental knowledge of liquid chromatography
- Knowledge regarding substances that are suitable only to a limited extent for use in liquid chromatography
- Knowledge regarding the health risks of chemicals
- Participation during an installation of a device or a training by the company KNAUER or an authorized company.

If you do not belong to this or a comparable professional group, you may not perform the work described in this user manual under any circumstances. In this case, please contact your superior.

### *Safety equipment*

When working with the device, take measures according to lab regulations and wear protective clothing:

- Safety glasses with side protection
- Protective gloves
- Lab coat

*What must be taken into account?*

- All safety instructions in the user manual
- The environmental, installation, and connection specifications in the user manual
- National and international regulations pertaining to laboratory work
- Original spare parts, tools, and solvents made or recommended by KNAUER
- Good Laboratory Practice (GLP)
- Accident prevention regulations published by the accident insurance companies for laboratory work
- Filtration of substances under analysis
- Use of inline filters
- Once they have been used, never re-use capillaries in other areas of the HPLC system.
- Only use a given PEEK fitting for one specific port and never re-use it for other ports. Always install new PEEK fittings on each separate port.
- Follow KNAUER or manufacturer's instructions on caring for the columns

More safety-relevant information is listed in alphabetical order in the following table:

Topic	Explanations
Flammability	Organic solvents are highly flammable. Since capillaries can detach from their screw fittings and allow solvent to escape, it is prohibited to have any open flames near the analytical system.
Solvent tray	Risk of electrical shock or short circuit if liquids get into the device's interior. For this reason, place all bottles in a solvent tray.
Solvent lines	Install capillaries and tubing in such a way that liquids cannot get into the interior in case of a leak.
Leaks	Regularly check if any system components are leaking.

Topic	Explanations
Power cable	Defective power cables are not to be used to connect the device and the power supply system.
Self-ignition point	Only use eluents that have a self-ignition point higher than 150 °C under normal ambient conditions.
Power strip	If several devices are connected to one power strip, always consider the maximum power consumption of each device.
Power supply	Only connect devices to voltage sources, whose voltage equals the device's voltage.
Toxicity	Organic eluents are toxic above a certain concentration. Ensure that work areas are always well-ventilated! Wear protective gloves and safety glasses when working on the device!

*Where is use of the device prohibited?*

Never use the system in potentially explosive atmospheres without appropriate protective equipment. For further information, contact the Technical Support of KNAUER.

*Taking the device out of operation*

At any time, take the device completely out of operation by either switching off the power switch or by pulling the power plug.

*Opening the device*

The device may be opened by the KNAUER Technical Support or any company authorized by KNAUER only.

## Definition of Personal and Material Damage

Possible dangers related to the device are divided into personal and material damage in this user manual.

Category	Explanations
DANGER!	Lethal or very serious injuries can occur.
WARNING!	Serious or moderate injuries can occur.
CAUTION!	Moderate injuries can occur. Device defects can occur.

## Decontamination

Contamination of devices with toxic, infectious or radioactive substances poses a hazard for all persons during operation, repair, sale, and disposal of a device.

### **DANGER!**

**Danger if getting in contact with toxic, infectious or radio-active substances.**

**Before disposing off or sending away contaminated devices, commission an expert with the decontamination.**

All contaminated devices must be properly decontaminated by a specialist company or the operating company before they can be recommissioned, repaired, sold, or disposed of.

All materials or fluids used for decontamination must be collected separately and disposed of properly.

# Symbols and Signs

The following table explains those symbols and labels which are used on the device, in the software or in the manual:

Symbol	Explanation
	Electrostatic-discharge hazard.
	General warning sign that moderate injuries or device defects can occur
	Device fulfills the requirements of the Confor- mité Européenne, which is confirmed by the Declaration of Conformity.

## Unpacking and Setup

### Contacting the Technical Support

You have various options to contact the Technical Support:

Phone: +49 30 809727-0

Fax: +49 30 8015010

E-mail: [support@knauer.net](mailto:support@knauer.net)

You can make your requests in English and German.

### Location Requirements

#### *Requirements*

The location for the device must meet the following requirements:

- Weight 2.3 kg

- (Width × Height × Depth) 121 × 138 × 190 mm
- Power supply 85-265 V
- Air humidity < 90 %, non-condensing
- Temperature 4-40 °C  
39.2-104 °F

**Caution!**

**Defect of the device due to overheating possible!  
Protect the device against exposure to direct sunlight.  
Keep clear at least 15 cm at the rear and 5–10 cm at each  
side for air circulation.**

## Unpacking

Store all packing materials. Packing list should be kept for repeat orders.

### *Tools*

Utility knife


### *Procedure*

1. Check for damages caused during transportation. In case you notice any damage, contact the technical support and the forwarder company.
2. Setup the delivery so the label is in the correct position. Using the utility knife, cut the adhesive tape. Open the delivery.
3. Remove the foam insert. Take out the accessories kit and the manual.
4. Open the accessories kit and take out all accessories. Check the scope of delivery. In case any parts are missing, contact the technical support.
5. Grip the device at its side panels and lift it out of the packaging.
6. Remove the foam inserts from the device.
7. Check for damages caused during transportation. In case you notice any damage, contact the technical support.
8. Set-up the device in its location.

## Initial Startup

### Screw Fitting of the Solvent Tube

Make sure that the tapered side of the seal ring ① faces the screw fitting ② of the solvent tube.

### Solvent Tube Connection

Appropriate solvent tubes are included in the device accessories.

- Make sure that the tube has a straight ending.
- Only tighten screw fittings by hand.
- Use the supplied blind fittings to seal off all unused inlets and outlets.

### Connecting Degasser and Pump

Each degassing chamber features an inlet and an outlet.

#### Caution!

**Very high pressures damage the degasser membrane. The membrane can withstand a maximum pressure of 7 bar.**

**Do not connect the degasser to the pump outlet.**

*Procedure*

1. Connection of the chambers in the direction of the arrow.
2. Remove the blind fittings from the degassing chamber inlets and outlets that are to be used. Keep the filler caps.
3. Unused in and outlets have to be closed with blind fittings to prevent contamination of the disconnected chambers.
4. The inlets of the degasser chambers are connected to the solvent containers by means of Teflon tubes.
5. Teflon tubes are also connected to the outlets of the degassing chambers, with screw fittings and seal rings. These tubes are then connected to the respective pump.

**Switching On the Degasser***Requirements*

- All lines have been connected correctly.
- Unused degasser chambers have been sealed off with blind fittings.

**Caution!**

**Possible damage to the device caused by condensed water!**

**Allow device to acclimate for 3 h, before it is connected to power supply and taken into operation.**

**Note:** Purge the degasser when changing solvents

**Note:** Regularly check the capillaries for leaks.

*Procedure*

1. Switch off the power supply.
2. Connect device with power supply
3. Connect the power supply of the device to the mains.
4. Switch on the power supply. During the start procedure, the LED flashes yellow for approx. 30 seconds. After the specified vacuum value has been reached, the LED lights up green permanently. The device is now ready for operation.
5. Switch on the HPLC pump.

When you switch off the HPLC pump, also switch off the degasser.

## Purging

Flush the degasser and the connecting tubes with the appropriate solvent prior to use. Thanks to the low chamber volume, only small amounts of the solvents are needed for flushing.

### *Prerequisite*

The degasser is connected to the pump.

**Note:**        **Prior to a solvent change in the degasser, purge with a liquid in which both solvents are soluble.**

### *Procedure*

1. Configure the pump: Flow rate 2 ml/min.
2. Start the pump and flush for 1 to 2 minutes.

## Leak Test

The vacuum is continuously checked by means of a built-in microprocessor and the pump output adjusted. Thus it is possible to detect leaks.

### *Procedure in the event of a detected leak*

- If a leak is detected, the microprocessor of the degasser automatically increases the pump speed to maintain the vacuum. If the vacuum still falls below the target value, the yellow LED flashes.
- If it is not possible to achieve the required vacuum within 30 minutes by means of the increased pump output, the LED lights up yellow at the degasser to indicate that there is possibly a leak in the system.
- The degasser will be switched off automatically.

## Switching Off the Degasser

Here a differentiation is made between short-term and long-term switch-off.

**Note:**        **To separate the degasser from the mains, remove the switched-mode power supply from the rear of the device.**

## Short-Term Switch-Off

**Note:**        **Damage to the degasser that results from using buffer-containing solvents are not covered by the warranty.**

### *Procedure*

1. Remove critical solvent from the degassing chambers and other system components by purging with for example isopropanol (2-propanol) .
2. Subsequently purge with water.
3. If the degasser is switched off for several days (e.g. over the weekend), purge the device with isopropanol prior to ending operation.
4. Switch off the degasser at the power supply.

## Long-Term Switch-Off

### *Procedure*

1. Remove critical solvent from the degassing chambers and other system components by purging with a suitable solvent.
2. First purge the degasser with water and subsequently with isopropanol.
3. Switch off the degasser and remove the connection tubes to the degasser.
4. The degasser chambers can be dried by means of laboratory gas.
5. Seal all inlets and outlets of the degasser with blind fittings and subsequently store the device at a safe and dry location.
6. Prior to the reuse of the degasser, purge the device briefly with the solvent which will be used afterwards.

# Troubleshooting

### *First measures*

1. Check all cabling.
2. Check all screw fittings.
3. Check whether air has gotten into the supply lines.
4. Check device for leaks.

### Further measures:

- Check errors against error list
- Contact the technical support hotline of the manufacturer

## Possible Problems and Rectifications

Problem	Cause	Solution
The device is on, but the LED is not light and the device is not energized.	<ul style="list-style-type: none"> <li>▪ Power supply or plug is not plugged</li> <li>▪ Fuse blown</li> <li>▪ Power supply is defective</li> </ul>	Contact the service department.
The device is running but you can not/hardly hear the pump.	At low rotational speeds, the pump can hardly be heard, although the vacuum is good and the degassing process is taking place.	Compare the baseline of non-degassed ethanol at 251 nm with the baseline of degassed ethanol. If the degasser is working properly, the noise of the baseline of the non-degassed ethanol should be significantly higher.
Air bubbles occur in the solvent tubes of the degasser outlet.	Loose or damaged screw fittings	Check the inlet and outlet fittings. Replace old fittings, if necessary.
No solvent flow	Air in the degasser chambers	Purge degasser chambers
	If a buffer solution is left in the degasser for an extended period, it can cause clogging in the degasser.	Use another degasser chamber. Carefully let water flow into the clogged degasser chamber to dissolve the buffer. If this does not help, contact the service department.

### *Further measures*

Inform the Technical Support of the manufacturer.

# Maintenance and Care

Proper maintenance of your HPLC device will ensure successful analyses and reproducible results.

## Maintenance Contract

The following maintenance work on the device may only be performed by the manufacturer or a company authorized by the manufacturer and is covered by a separate maintenance contract:

- Opening the device
- Removing the hood or the side panels.

## Cleaning and Caring for the Device

### Caution!

**Intruding liquids can cause damage to the device!**  
**Place solvent bottles next to the device or in a solvent tray.**

**Moisten the cleaning cloth only slightly.**

All smooth surfaces of the device can be cleaned with a mild, commercially available cleaning solution, or with isopropanol.

## What is to be done if the degasser leaks?

If the degasser leaks, solvent can reach the motor and from there reach the degasser exhaust and flow out. In this case, the automatic pump control will switch the device off.

- Inform the Technical Support of KNAUER.

## Tightening Screw Fittings

1. If the screw fitting of the solvent tube is untight, carefully tighten the screw fitting.
2. If it still leaks, loosen and check the screw fitting.
3. If no damage or deformation is visible at the seal ring and at the screw fitting, repeat step 1.
4. If it still leaks, replace the screw fitting and seal ring.

## Preventative Measures

To ensure that the degasser can constantly run at maximum capacity, pay attention to the following:

- Only use ultrapure solvent as mobile phase. If water is used as solvent, it also has to fulfill this standard or be filtered and deionized.
- Filter all solvents to prevent clogging.
- To dry the degasser, only use ultrapure gas.
- Observe the information on switching off the degasser.

## Technical Data

Channels	2
Maximum flow rate	10 ml/min
Degassing method	gas permeation through a fluoropolymer membrane
Dead volume	approx. 285 µl per channel
Solvent applicability	universal, except hydrochloric acid and halogenated hydrocarbons especially hexafluoro isopropanol (HFIP), hexane, heptane
Pressure decline	1.37 mm (Hg/ml/min)
Maximum pressure stability	7 bar
Wetted materials	PEEK, glass-filled PTFE, Teflon AF®
Display	1 LED
Dimensions	121 x 138 x 190 mm
Weight	2.3 kg
Power supply	85–265 V, 50–60 Hz, 20 W
Ambient conditions	temperature range: 4–40°C; 39.2–104 F air humidity: Below 90 %, non-condensing

## Accessories and Spare Parts

For repeat orders of spare parts please contact the Technical Support of the manufacturer or conform to the packing list which is enclosed when the device is delivered.

**Note:** For repeat orders of spare parts use the enclosed packing list. Contact the Technical Support in case there are any questions on spare parts or accessories.

Name	Order number
Degasser DG 2.1S	AZE02
Accessory kit AZURA small devices	FZA01
Accessory kit DG 2.1S	FZE02
User manual	V6880

# Legal Information

## Warranty Conditions

The factory warranty for the device is valid for 12 months after the date of dispatch. All warranty claims shall expire in the event that any unauthorized changes are made to the device.

During the warranty period, any components with material or design-related defects will be replaced or repaired by the manufacturer free of charge.

This warranty excludes the following:

- Accidental or willful damage
- Damage or errors caused by third parties that are not contractually related to the manufacturer at the time the damage occurs
- Wear parts, fuses, glass parts, columns, light sources, cuvettes and other optical components
- Damage caused by negligence or improper operation of the device and damage caused by clogged capillary
- Packaging and transport damage

In the event of device malfunctions, directly contact the manufacturer.

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Internet:www.knauer.net

## Transport Damage

The packaging of our devices provides the best possible protection against transport damage. Check the devices for signs of transport damage. In case you notice any damage, contact the Technical Support and the forwarder company within three workdays.

# Disposal

Hand in old devices or disassembled old components at a certified waste facility, where they will be disposed of properly.

## *AVV marking in Germany*

According to the German "Abfallverzeichnisverordnung" (AVV) (January, 2001), old devices manufactured by KNAUER are marked as waste electrical and electronic equipment: 160214

## *WEEE registration*

KNAUER as a company is registered by the WEEE number DE 34642789 in the German "Elektroaltgeräteregister" (EAR). The number belongs to category 8.

All distributors and importers are responsible for the disposal of old devices, as defined by the WEEE directive . End-users can send their old devices manufactured by KNAUER back to the distributor, the importer, or the company free of charge, but would be charged for the disposal.

## *Solvents and other operating materials*

All solvents and other operating materials must be collected separately and disposed of properly.

All wetted components of a device, e. g. flow cells of detectors or pump heads and pressure sensors for pumps, have to be flushed with isopropanol first and water afterwards before being maintained, disassembled or disposed.

# Index

## A

Accessories 9  
Acids 9  
Additives 8  
Ambient conditions 22  
AVV marking 25

## B

Buffer solutions 8

## C

CE marking, see Declaration of conformity 27

## D

Declaration of Conformity 27  
Decontamination 13  
Degasser  
  Front view 5  
Degassing principle 7

## E

Eluent 7  
Environmental protection 25

## F

Flushing the device  
  Flushing the degasser 18

## I

Installation location, see location 14

## L

Leak 11  
Leak test 18  
Leakage 21  
Location 14  
Long-term switch-off 19

## M

Modifiers 8

## O

Order numbers 23

## P

Power supply 12  
  cable 12  
  strip 12  
Professional group 10

## S

Safety equipment 10  
Salts 8  
Screw fittings  
  tightening 21  
Short-term switch-off 19  
Solvent  
  flammability 11  
  line 11  
  self-ignition point 12  
  toxicity 12  
  tray 11  
Solvent exchange 18  
Solvent tube  
  Screw fitting 16  
Spare Parts 9  
Spare parts 23  
Switch-off 18  
Switch-on 17

## T

Technical Support 14  
Transport damage 24

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## *Product:*

Degasser DG 2.1S: EZE02

The product complies with the following requirements and specifications:

- DIRECTIVE 2006/42/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 17 May 2006 on machinery, and amending Directive 95/16/EC (recast)
- IEC 60799 (1998) Electrical accessories – Cord sets and interconnection cord sets
- EN 61000-3-2 (2005 + A1:2008 + A2:2009) Electromagnetic compatibility (EMC) Part 3-2
  - EMC standards (2004/108/EC)
- IEC 61010-1 (2010 + Corrigendum: 2011) Safety requirements for electrical equipment for measurement, control and laboratory use
  - Low voltage directive (2006/95/EC)
- EN 61326-1 (2006) Electrical equipment for measurement, control and laboratory use – EMC requirements
  - EN 61326-1 Corrigendum 2 (2010)
- Environmentally sound use of electrical and electronic equipment:
  - RoHS directives 2002/95/EC (2003) and 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment
  - WEEE directive 2002/96/EC (2003) on waste electrical and electronic equipment

The device was tested with a typical configuration. The mark of conformity has been applied to the rear panel of the device:

Berlin, 2012-03-23

Alexandra Knauer (CEO)

© KNAUER Wissenschaftliche Geräte GmbH  
All rights reserved.  
Technical modifications are subject to change.  
Translation of the original German editions of  
this manual, version 2.0 from 2015-01-28.  
Printed in Germany on environmentally  
friendly paper from sustainable forests.

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® AZURA is a registered trademark of