

CDSolutions

APPLICATIONS INFORMATION USING ADVANCED SAMPLE HANDLING TECHNOLOGY

CAPILLARY GC OF POLYOLEFIN PYROLYSATES





Pyrolysis gas chromatography is a simple technique which has been utilized by scientists for almost three decades. The process of pyrolysis involves the liberation of organic material from a relatively nonvolatile matrix. Natural and synthetic polymers can be fragmented to form smaller, more volatile compounds which can be analyzed by capillary gas chromatog-

raphy. Polymer samples can be pyrolyzed as solids in a quartz tube, or as polymer solutions dissolved in a solvent and deposited on a ribbon as a thin film.

Polyolefins are essentially very high molecular weight hydrocarbons. Upon pyrolysis, smaller hydrocarbons are produced which are characteristic of the polyolefin composition.

Polyolefins undergo a process of random scission producing olefins, alkanes and diolefinic compounds. The presence of isoalkanes indicates branching in the polymer chain. The product distribution in pyrolysis is dependent on the pyrolysis temperature. Pyrolyzing polyolefins at higher temperatures will cause an increase in the abundance of

Continued on back



Continued from front

pyrolysates at the earlier retention regions due to increased fragmentation.

Pyrolysis of polyethylene produces primarily straight chain alkanes and alkenes. The nature of the polyethylene chain produces oligomeric groups of pyrolysis fragments differing from each other by one carbon unit. This pattern is predictable and reproducible for polyethylene. Polypropylene (Figure 2), when pyrolyzed, produces fragments which are more highly substituted. The oligomeric groups in polypropylene pyrolysis vary by three carbon units. Structural defects may be determined in polyolefins by the degree of substituted alkanes produced by pyrolysis. Polyisobutylene and polytransisoprene (Figures 3 and 4), give more complex chromatograms due to the increased substitution of the polymer backbone. Polyisobutylene can be identified by

alternate carbons in the polymer chain being disubstituted with methyl groups, causing a greatly different pattern than that of polyethylene and polypropylene. The primary pyrolysis products of polytransisoprene are isoprene monomer and dipentene. Pyrolysis can thus be used in a qualitative manner to identify polymers, but also in a quantitative sense in the analysis of copolymer blends.

EQUIPMENT PYROLYSIS

CDS Model 120 Pyroprobe

Pyrolysis temperature:

750°C

Interface temperature: 285°C

GAS CHROMATOGRAPHY

Column: 50m x 0.25mm, SE-

54 capillary

Initial temperature: 50°C for 3 minutes

Program rate: 8°C/min to

285°C

For more information on this and related applications, we recommend the following readings:

Levy, E. J. and T. P. Wampler. "Effects of Slow Heating Rates on Products of Polyethylene Pyrolysis." *Analyst*, Vol. III, (1986), pp. 1065-1067.

Nagaya, T. et al. "Microstructural Characterization of Polypropylenes by High Resolution Pyrolysis-Hydrogenation glass capillary gas chromatography." *Macromolecules*, Vol. 13, (1980).

Additional literature may be obtained from your Chemical Data Systems representative, or by writing to the CDS Applications Lab.

ABOUT CDS

CDS Analytical, Inc. is a leader in the design and manufacture of laboratory instruments for sample preparation and analysis. With 20 years experience in the field, CDS is dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of analytical pyrolysis, CDS manufactures the Pyroprobe 1000 and 2000 for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of purge and trap instruments for the analysis of volatile organic compounds in the environmental, food and pharmaceutical areas, as well as custom systems for complex, multicomponent materials investigation. Our customers, their requirements and applications are important to us. To help meet their needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, or call us at 1 800 541 6593.