



# Agilent's GC Analyzer Solutions: Introduction to the Controlled Substances Analyzer

---

Get your laboratory on the  
Analytical ***“FAST TRACK”***

December 2013

# Today's Challenges

## Expecting More from Less

### Business Challenges


- Increasing cost of ownership
- Trace analysis in complex matrices

### Resource Allocation

- Less time for method development
- Limited technical experience

### Optimizing Efficiency

- Demands for higher throughput
- Quicker return on capital investment



# Analyzer Portfolio

Resolving the Your Search for Value

**GO.**


**MOVE BEYOND  
INSTRUMENTATION  
AND GET YOUR LAB ON THE  
"FAST TRACK" TO RESULTS**

**AGILENT ANALYZERS  
AND APPLICATION KITS**

The Measure of Confidence


# What is an "Analyzer"?

*An Analyzer is a factory-configured GC or GC/MSD system constructed, tuned and tested for the application prior to delivery to the customer's laboratory.*




# Pre-Configured Analyzers

## Increasing the Value Proposition



**Focus your team on analyses; not method development!**


# Analytical Requirements





# Intrinsic Value of Analyzers and Application Kits

Reduce the time required for system deployment



...Faster Application Startup with a Guaranteed Method

# GO.



## **Controlled Substances Analyzer** ... IDENTIFY Drugs of Abuse Targets with Confidence!

*Developed in Collaboration with:*

***Sarah Keeling and Fran Diamond***  
***NMS Labs, Willow Grove, PA***





# Analysis of Controlled Substances: Challenges

## Routine Seized Drug Analysis

### Multiple sample types:

- Pharmaceuticals. examples: oxycodone, hydrocodone
- Illicit drugs. examples: heroin, cocaine
- Botanical materials. examples: marijuana, synthetic marijuana (Spice)
- Syringe residues. examples: heroin, methamphetamine
- Charred pipe residues. examples: crack cocaine, hashish

### Analyte list constantly growing

- Synthetic cannabinoids. examples: JWH-018, PB-22
- Designer stimulants. examples: MDPV, BZP
- Novel psychoactive substances



# Analysis of Controlled Substances: Challenges

## Routine Seized Drug Analysis

### Laboratory Needs:

- Analyze complex mixtures by GC/MS
- Efficiency of testing
- Speed of analysis
- High quality results
- Time efficient data processing
- Simple sample prep



### Analytical Challenges

- Samples of varying nature
- Identify closely related analytes: isomers, analogs, etc.
- Variable concentrations: column overload to trace concentrations
- Complex matrix: cutting agents, adulterants, non-volatile compounds

# Current Situation

## Routine Seized Drug Analysis

Powders &



Expanding list of new drug targets

Long Cycle Times (23.0 Min)

Carry-over interferences between samples

PBM Based Searching

Manual Data Review

- Manual Spectrum subtraction
- Manual Baseline Review

Methods vary from lab-to-lab

- In place validation per lab SOP
- Limited standardization of method

Inconsistent spectral libraries and retention times

Cost of Analysis with Helium Carrier

# Wanted to Create a New Configuration



## Rapid Identification of Seized Drug Samples

- **Reduce operating expense:** Use hydrogen carrier gas
- **High speed analysis:** from 23.0 to 13.5 min run-to-run for analysis to Vardenafil
- **Backflushing:** Minimize column trimming and eliminate sample carry-over
- **Retention Time Locking (RTL):** Precise RT matching column-to-column, instrument-to-instrument, and lab-to-lab
- **Deconvolution Reporting Software:** Rapid reliable compound identification
- **New 460 compound database:** Comprehensive list of controlled substances including new psycho active substances (synthetic cannabinoid, and “Bath Salt” compounds)
- **Ultra Inert Flow Path: Inlet to Detector** Minimizes sample decomposition

# Configuration of Controlled Substance Analyzer

Post column Capillary Flow Technology (PUU) device provides:

- A) Backflushing
- B) No vent column changing



Purged Ultimate Union (PUU)

# Method Parameters for Hydrogen Carrier Controlled Substance Analyzer

Ramp	'C/min	'C	Hold min
Initial		50	1
Ramp 1	60	325	4
Runtime	9.58 min		
Backflush	325 'C/min for 1.0 min		

Inlet	Split/Splitless
Temp	265 'C
Mode	Pulsed Splitless, Constant Flow
Flow	1.4 (adj to lock)
Pulse Press	11 psi
Pulse time	0.3 min
Purge time	1 min
Purge Flow	50 mL/min

Column	DB-5MSUI part # (100-2000) 10m x 0.25 mm id x 0.25 um film
Outlet Pressure	Programmed for constant MSD restrictor flow (2.81 psig initial)

Injection volume 1 uL  
RT locked to Proadifen at 5.039 min

MSD	Agilent 5977
Solvent Delay	2.3 min
Acquisition Mode	Scan
Scan Range	40 to 570
Threshold	0
Sampling	2
TID	ON
Quad Temp	200 'C
Source Temp	300 'C
Transfer Line	300 'C
Tune	Atune, Gain Normalized 1X

Backflush Device Purged Ultimate Union  
with overpressure vent

MSD Restrictor	0.81 m x 0.12 mm id Inert fused silica tubing
Restrictor Flow	2.0 mL/min Constant Flow
Backflush Time	1.0 min
Backflush Temp	325
Backflush Press	40

# GO.



## Controlled Substances Analyzer

... Development Challenges using  
Hydrogen Carrier Gas



# Incompatibility of O-Rings with H2 Carrier

## Restek Publication

**SeparationScience**  
PREMIER LEARNING FOR ANALYTICAL CHEMISTS

North America  
Volume 5, Issue 5  
2013



**Identification of Allergens in Common  
Snack foods using Transportable  
Single-quadrupole GC-MS**

**Ghost Peaks in Gas Chromatography Part 2:  
The Injection Port**

*Jaap de Zeeuw*

**Identification of Flavour and Fragrance Allergens in  
Common Snack foods using Transportable  
Single-quadrupole GC-MS**

*Suli Zhao*

**Comparison of Regulated Methods Using Ion  
Chromatography for Bromate Analysis in Drinking  
Water**

*Richard F. Jack*

**Selective Therapeutic Monitoring of Drugs and Their  
Pharmacologically Active Metabolites by Liquid  
Chromatography**

*Gellért Karvaly, Judit Szaniszló and József Fűrész*

---

**Analysis of Beer by Comprehensive 2D-LC**

**Transition metals in feedwater using post-column reaction and  
subsequent UV/VIS detection**

**High-performance screening of pollutants in water using GCxGC/TOF  
MS**

**Automated Extraction and GC/MS Determination of Phthalates in  
Consumer Products**

**Quantitation and Identification of Pharmaceuticals and Personal  
Care Products (PPCP) in Environmental Samples using Advanced  
TripleTOF® MS/MS Technology**

[Click here to view more>>](#)



# Incompatibility of O-Rings with H<sub>2</sub> Carrier

Restek Publication



Figure 7: Deformation of "O"-ring when used with hydrogen as the carrier gas. Helium did not show this effect.

## "O" – Ring Contamination

Using an Agilent gas chromatograph, the injection port also contains a "O"-ring. This ring prevents the sample from re-entering the liner after injection. The "O" ring is a challenge as it is continuously at high temperature and it is a polymeric material that has to make a good seal. We received customer input that shared issues with deformed rings in split/splitless injection systems (see Figure 7) which caused injection issues. In this case the "O"-ring deformed and hardened using hydrogen as the carrier gas. It did not happen when helium was used.

[http://www.sepscience.com/docs/Bespoke/Editions/SepSci/Sepsci0513na.pdf?utm\\_source=iContact&utm\\_medium=email&utm\\_campaign=eLearning%20Newsletter%20-%20North%20America&utm\\_content=](http://www.sepscience.com/docs/Bespoke/Editions/SepSci/Sepsci0513na.pdf?utm_source=iContact&utm_medium=email&utm_campaign=eLearning%20Newsletter%20-%20North%20America&utm_content=)

# Incompatibility of O-Rings with H2 Carrier

Restek Publication




Figure 8: Ghost peaks produced by a dirty injection port. Showing peaks for triphenyl phosphine oxide ( $m/z$  279), phthalates ( $m/z$  149, 167 and 277) and PAH ( $m/z$  278). Full story, see ref. [4].

[http://www.sepscience.com/docs/Bespoke/Editions/SepSci/Sepsci0513na.pdf?utm\\_source=iContact&utm\\_medium=email&utm\\_campaign=eLearning%20Newsletter%20-%20North%20America&utm\\_content=](http://www.sepscience.com/docs/Bespoke/Editions/SepSci/Sepsci0513na.pdf?utm_source=iContact&utm_medium=email&utm_campaign=eLearning%20Newsletter%20-%20North%20America&utm_content=)

The “O”-ring itself can also produce ghost peaks. Figure 8 shows a peak that typically shows up as a “ghost-peak” generated by the “O”-ring. The component, triphenyl phosphine oxide, will accumulate on the column when it’s set at low temperature. It has a  $m/z$  of 266. The same chromatogram also shows phthalate peaks that are used as septum-softeners.

It is important to make sure that “O”-rings are used with sufficient thermal stability. There are many brands of “O”-rings and problems can result — not only from the formation of ghost peaks but also from the challenge sometimes of getting liners out of injection systems, when maintenance is due.

# System Development: O-Ring Issues


H<sub>2</sub> Carrier for Controlled Substances



- Cotton swabs used to clean injection port metal.
- Swabs with dirt used to clean top of the weldment assembly.

# System Development: O-Ring Issues

H<sub>2</sub> Carrier for Controlled Substances




O-ring suspected of leaving a deposit (breakdown product) around the top of the weldment



Figure 7: Deformation of "O"-ring when used with hydrogen as the carrier gas. Helium did not show this effect.

# System Development: O-Ring Issues

## H<sub>2</sub> Carrier for Controlled Substances



Residue on top of  
inlet weldment

### Problem Resolution:

- Convert to 5190-2293 Liner, UI, splitless, single taper, glasswool with o-ring
- 5188-5365 Liner O-Ring, Non-Stick 10PK
- 5190-6145 Gold Seal, Ultra Inert, w/Washer, 10/PK
- Provided UI Inlet Weldment

# New Inert Flow Path for the 7890B GC / 5977A MSD

Reliability, Durability, Speed and Ease of Use



From inlet-to-detector ... **COMPLETE DEACTIVATION** of sample path!

# GO.



## Controlled Substances Analyzer

### Analytical Performance





# Test Mix Chromatography

## Original Helium Method and New Analyzer Method

- 1 Amphetamine
- 2 Phentermine
- 3 Methamphetamine
- 4 N-Propylamphetamine
- 5 Pentobarbital
- 6 10,11-Dihydrodibenz(b,f)(1,4)oxazepin-11-one

- 7 Cocaine
- 8 Oxycodone
- 9 Alprazolam
- 10 Strychnine






# Example: Street Heroin

## Controlled Substances Analyzer

- |   |   |    |                      |
|---|---|----|----------------------|
| 1 | N-Propylamphetamine (ISTD)                          | 7  | 6-Monoacetylmorphine |
| 2 | Benzocaine  | 8  | Heroin               |
| 3 | Caffeine  | 9  | Papaverine           |
| 4 | Lidocaine   | 10 | Noscapine            |
| 5 | 10,11-Dihydrodibenz(b,f)(1,4)oxazepin-11-one (ISTD) |    |                      |
| 6 | Acetylcodeine                                       |    |                      |



# DRS Report for Heroin Sample

## Controlled Substances Analyzer

MSD Deconvolution Report

Sample Name: ALA33

Data File: D:\MassHunter\GCMS\1\data\Heroin\_Sample.D

Date/Time: 15:14 Wednesday, Oct 2 2013

Adjacent Peak Subtraction = 1

Resolution = High

Sensitivity = Medium

Shape Requirements = Medium


The NIST library was searched for the components that were found in the AMDIS target library.

Retention Time (R.T.)	Cas #	Compound Name	Amount (ng)		AMDIS		NIST	
			Chemstation	AMDIS	Match	R.T. Diff	Reverse Match	Hit Number
3.285	51799327	N-Propylamphetamine			99	-0.4	90	1
3.8059	94097	Benzocaine			100	-0.5	95	2
4.2958	58082	Caffeine			100	-0.4	91	1
4.3613	137586	Lidocaine			99	-0.4	92	1
4.7037	3158858	10,11-Dihydrodibenz(b,f)(1,4)oxazepin-11-one			97	0.2	93	1
5.3308	6703271	Acetylcodeine			100	-0.2	95	1
5.3687	2784738	6-Monoacetylmorphine			100	0.2	95	1
5.4999	561273	Heroin			100	0.9	97	1
5.7215	58742	Papaverine			98	-0.2	94	1
6.2532	128621	Noscapine			92	-0.5	90	1

# Example: Marijuana Sample

## Controlled Substances Analyzer

- 1 N-Propylamphetamine (ISTD)
- 2 10,11-Dihydrodibenz(b,f)(1,4)oxazepin-11-one (ISTD)
- 3 Delta-9-tetrahydrocannabinol (THC)
- 4 Cannabinol
- 5 JWH-250
- 6 JWH-018



# DRS Report for Marijuana Sample

## Controlled Substances Analyzer

MSD Deconvolution Report  
 Sample Name: KCJ06k  
 Data File: D:\MassHunter\GCMS\1\data\Agilent  
 helium to hydrogen\Case 07.D  
 Date/Time: 15:17 Friday, Nov 1 2013

Adjacent Peak Subtraction = 1  
 Resolution = High  
 Sensitivity = Medium  
 Shape Requirements = Medium

The NIST library was searched for the components that were found in the AMDIS target library.



Retention Time (R.T.)	Cas #	Compound Name	Amount (ng)		AMDIS Match	R.T. Diff (sec)	NIST	
			Chemstation	AMDIS			Reverse Match	Hit Num.
3.2713	51799327	N-Propylamphetamine			100	-1.2	89	1
4.6926	3158858	10,11-Dihydrodibenz(b,f)(1,4)oxazepin-11-one			98	-0.4	93	1
5.2699	1972083	Delta-9-tetrahydrocannabinol (THC)			100	-0.1	95	1
5.3508	521357	Cannabinol			100	0	93	1
5.8397	864445432	JWH-250			99	-0.8	89	1
6.2516	209414073	JWH-018			96	-1	89	1



Synthetic Cannabinoids

# Analysis of Mushroom Drug Sample

## Original Helium Method




**Heavy matrix compounds can carry over to following analysis, necessitating very long post run bake out**


# Analysis of Blank Following Mushroom Sample

H<sub>2</sub> Carrier Instrument

Blank following mushroom sample,  
0.5 min backflush  
**(too short)**



Blank following mushroom sample,  
1.0 min backflush  
**(correct length)**




**Using a 1.0 minute backflush prevents heavy matrix compounds from carrying over to next run**

# Deconvolution Training Standard

## Deconvolution Reporting Software (DRS) Familiarization


**What is this peak?**  
**First, use conventional approach of average spectrum and PBM search to identify it**




# PBM Identifies THC With 96 Match

## Conventional Approach

### Average Spectrum



### THC, 96 Match



- First hit is THC with 96 PBM match.
- Second hit is hydrocodone with only 35 match

PBM Search Results: C:\DATABASE\FORCHEM.C.L


Rank	Name	Ref No.	MW	Qual
1	Delta-9-tetrahydrocannabinol (THC)	149	314	96
2	Hydrocodone	80	299	35
3	Codeine	32	299	16

Difference




# Peak Is THC, Right?


## Average Spectrum



## Hit 1: THC, 96 Match PBM




## Hit 2: Hydrocodone, 35 Match PBM



# Overlap of THC and Hydrocodone

As viewed in AMDIS




- Deconvolution finds 2 compounds closely merged together.
- Peaks are separated by only **0.002 min.** (Less than 1 scan)


# First Peak Found By AMDIS

## THC

### Deconvoluted Spectrum at 5.281 min



### THC, 95 Net Match




**AMDIS deconvolution confirms that the first peak is THC**


# Second Peak Found By AMDIS

Hydrocodone

## Deconvoluted Spectrum at 5.283 min



## Hydrocodone, 95 Net Match



This peak would probably have been missed without AMDIS deconvolution

# DRS Report For Training Sample

## Automated Deconvolution and Reporting

### MSD Deconvolution Report

Sample Name: THC/HC

Data File: D:\MassHunter\GCMS\1\data\HCOD\_THC.D

Date/Time: 17:18 Wednesday, Oct 2 2013

Adjacent Peak Subtraction = 1

Resolution = High

Sensitivity = Medium

Shape Requirements = Medium

The NIST library was searched for the components that were found in the AMDIS target library.

Retention Time (R.T.)	Cas #	Compound Name	Amount (ng)		AMDIS		NIST	
			Chemstation	AMDIS	Match	R.T. Diff	Reverse Match	Hit Number
3.2702	51799327	N-Propylamphetamine			99	-1.3	91	1
4.6881	3158858	10,11-Dihydrodibenz(b,f)(1,4)oxazepin-11-one			97	-0.7	90	1
5.2479	521379	Cannabidiol			68	5.6	64	3
5.281	1972083	Delta-9-tetrahydrocannabinol (THC)			95	0.5	85	1
5.283	125291	Hydrocodone			95	0.2	88	1
5.3145	466999	Hydromorphone			78	-0.1	67	2
5.3637	521357	Cannabinol			67	0.8	55	4

Automated report generation takes about ~ 60 sec!



# New Controlled Substances Analyzer

## Rapid Identifying Seized Drug Samples

- ✓ • **Reduce operating expense:** Use hydrogen carrier gas
- ✓ • **High speed analysis:** from 23.0 to 13.5 min run-to-run
- ✓ • **Backflushing:** Minimize column trimming and eliminate sample carry-over
- ✓ • **Retention Time Locking (RTL):** Precise RT matching column-to-column, instrument-to-instrument, and lab-to-lab
- ✓ • **Deconvolution Reporting Software:** Automated, rapid and reliable compound identification
- ✓ • **New 460 compound database:** Comprehensive list of controlled substances including new psycho active substances (synthetic cannabinoid, and “Bath Salt” compounds)
- ✓ • **Ultra Inert Flow Path:** “Designed-in” from **Inlet-to-Detector** to minimize sample decomposition

# GO.



## Controlled Substances Analyzer

How do I buy it?





# Controlled Substances Analyzer

Helium Carrier Gas

**G3445B**  
**Option #472**

- 7890 GC System, Column, c/o Sample
- Opt 002/003 Fast Oven Option
- Opt 114 Inert Capillary S/SL EPC
- G7038A 5977A Inert EI MS Turbo
- G1679AA Controlled Subs. RTL DB
- G1716AA Deconv. & Reporting SW
- G1033A MS Library Bundle

**Contact your Agilent Account Manager for full ordering details**



# Controlled Substances Analyzer

Hydrogen Carrier Gas

**G3445B**  
**Option #476**

7890 GC System, Column, c/o Sample

Opt 002/003 Fast Oven Option

Opt 114 Inert Capillary S/SL EPC

G7038A 5977A Inert EI MS Turbo

G1679AA Controlled Subs. RTL DB


G1716AA Deconv. & Reporting SW

G1033A MS Library Bundle

**Contact your Agilent Account Manager for full ordering details**

# Analyzer Value to Customers

Helping your customer enhance their competitive advantage





# Questions or Comments






# APPENDIX

# Codeine, H<sub>2</sub> Instrument

**Some hydrocodone is formed when codeine is injected.  
About 6% is converted.**



# Caffeine, H<sub>2</sub> Instrument

**Severely overloaded caffeine peak has RT shifted outside recognition window. Use alternate window setting in AMDIS when these are encountered.**

