

CDS Solutions


APPLICATIONS INFORMATION USING ADVANCED GC SAMPLE HANDLING TECHNOLOGY

DEGRADATION MECHANISMS - RANDOM SCISSION


During pyrolysis, polymeric materials may degrade via a number of mechanisms which are generally grouped into three classes - random scission, depolymerization, and side group elimination. Random scission results from the production of free radicals along the backbone of the polymer which causes the macromolecule to be fragmented into smaller molecules of varying chain lengths. These fragments chromatograph to reveal a repeating series of oligomers frequently differing in chain length by the number of carbons in the original monomer.

The polyolefins generally degrade through a random scission mechanism, and polyethylene is a good example of this behavior. When a free radical is formed along the chain of polyethylene, chain scission occurs, producing a molecule with an unsaturated end and another with a terminal free radical. This free radical may take a hydrogen from a neighboring carbon, producing a saturated end and a new radical, or combine with another free radical to form an alkane. Multiple cleavages produce molecules small enough to be volatile, with double bonds at both ends, one end, or neither. Since the scission was random, molecules are made with a wide variety of


RANDOM SCISSION



PRODUCING



PYROLYSIS OF POLYETHYLENE 750°C FOR TEN SECONDS



1. Decadiene
2. Decene
3. Decane

chain lengths. These appear in the pyrogram as a series of triplet peaks. Each triplet consists of an alkane, an alkene and a diene of a specific chain length. The hydrocarbons in each triplet have one more carbon than the molecules in the triplet which eluted just prior to it.

The accompanying chromatogram resulted from the pyrolysis of polyethylene at 750°C showing oligomers containing up to 30 carbons, with the C₁₀ fragments marked.

EQUIPMENT PYROLYSIS

CDS Model 120 Pyroprobe,
coil probe with quartz tube
Temperature: 750°C for ten
seconds

Interface temperature: 280°C

GAS CHROMATOGRAPHY

Column: 25m x 0.25mm fused
silica capillary, SE-54

Detector: Flame ionization

Initial temperature: 50°C for 2
minutes

Rate: 8°C/min

Final temperature: 300°C for
10 minutes

Chart speed: 1 cm/min

Split ratio: 75:1

Carrier gas: Helium

For more information on this
and related applications, we

recommend the following
readings:

Irwin, William J. *Analytical
Pyrolysis: A Comprehensive
Guide*. Marcel Dekker,
publisher.

Levy, E. J. and S. A. Liebman.
*Pyrolysis and GC in Polymer
Analysis*. Marcel Dekker,
publisher.

Levy, E. J. and T. P. Wampler.
"Effects of Slow Heating Rates
on Products of Polyethylene
Pyrolysis." *Analyst*, Vol. III,
(1986), pp. 1065-1067.

Additional literature may be
obtained from your Chemical
Data Systems representative, or
by writing to the CDS Applica-
tions Lab.

ABOUT CDS

CDS Instruments, a unit of Autoclave Engineers Inc., is a worldwide leader in the manufacture of instruments and equipment for research applications. With over 15 years of service to the research community, CDS is dedicated to the development of instrumentation for use in the preparation, separation and analysis of complex organics. We are also dedicated to providing researchers with the technical and applications information they need to perform their work efficiently. The purpose of CDSolutions is to supply the researcher with a constant stream of applications solutions for advanced GC sample handling. The article you just read is one of many available at no charge to our customers. For a complete list of topics, please write to our Technical Information Dept. in Oxford, Pennsylvania.



INSTRUMENTS
A unit of Autoclave Engineers

• 7000 Limestone Rd., Oxford, PA 19363, U.S.A. • Phone: (215) 932-3636 • Telex: 83-5308