


Environmental and Food Safety Analyses by UHPLC and LC/MS Using the New Agilent 1290 Infinity LC System and Agilent Mass Spectrometers

January, 2010

Mike Woodman, Jerry Zweigenbaum
Applications Development
Agilent Technologies
Wilmington, Delaware



Topics of Discussion

UHPLC theory and advantages

Challenges in Instrument design

- Pump: flow and pressure range, mixing requirements
- Injector: injection range, principle, carryover
- Column availability, stability, variety
- Column management: temperature control, column dimensions and valve configurations
- UV/VIS detection: balancing dispersion and sensitivity
- MS Considerations: acquisition rate, optimal flow rate

Resolution Equation

Although we may be aiming for faster separations, we normally also need to maintain or increase resolution

$$R_s = \frac{\sqrt{N}}{4} \frac{(\alpha-1)}{\alpha} \frac{k}{(k+1)}$$


α = **Selectivity** – influenced by mobile and stationary phase

N = **Column Efficiency** – influenced by length and particle size

k = **Capacity Factor** (retention) – influenced by stationary and mobile phase, gradient slope and dwell volume (gradients)

Method Transitions to UHPLC sub 2 micron materials

Maintain Acceptable Flow Rates at High Linear Velocity by Progressively Decreasing Column Diameter



Lin. Vel.	mm/sec	1	2	3	4	5
4.6 mm	ml/min	0.7	1.4	2.1	2.8	3.5
3 mm	ml/min	0.3	0.6	0.9	1.2	1.5
2.1 mm	ml/min	0.14	0.29	0.44	0.58	0.73
1 mm	ml/min	0.033	0.066	0.1	0.133	0.166

Increased Pressure Requirements and Enhanced Speed and Performance Associated with Reduced Particle Size


d_p (μm)	ΔP (bar)	t_R (min)	H_{MIN} (μm)	N
5.0	14.5	34.7	10.0	25,000
3.0	66.9	20.8	6.0	41,000
1.5	531	10.4	3.0	83,000
1.0	1800	6.9	2.0	125,000
0.75	4270	5.2	1.5	166,000
0.50	14400	3.5	1.0	250,000

With a 5.0 to 1.5 μm particle size reduction

- t_R decreases >3-fold
- N increases >3-fold


But . . . ΔP increases >35-fold

Maximum Power Range




Advanced UHPLC Binary Pump

- Flow path and integration of new materials




What to Look for in Gradient Performance

- Steep 0.5ml (25%/cv) Gradient, alkylphenones, system 38




System to System Performance


- systems 38 and 0x67, same column installed, 2ml gradient



Resolution – effect of increasing capacity factor k' in isocratic Resolution Equation



Resolution – effect of increasing capacity factor k^* in gradient peak capacity equations



for a 1ml/min gradient, 30-90% ACN, 3x50mm 1.8 μ m, 11ul disp., avg. analyte mc.wt. 162g/m

My Lab Environment?



1290 Infinity LC

Fluorescence

Diode Array

Outboard 12 position solvent
select valves

Column Compartment with 2
position 6 port valve

Autosampler

Binary pump with integral
vacuum degasser and solvent
select valves

What Kind of Lab Environment is Needed?





Design Considerations for UHPLC Sample Introduction

- Injection principle (variable loop, fixed loop designs)
- Usable injection volume range
- Speed of injection cycle
- Minimal dispersion
- Minimal carryover
- Minimal sample waste

Design Considerations for UHPLC Samplers

- Injection principle (variable loop, fixed loop designs)



Fixed loop designs leave flush solution in loop, which can distort peak shape or alter selectivity in scouting separations

Tips for minimizing dispersion in LC systems

- Keep injection volumes small, ensure weak diluent when possible
- Ensure autosampler flush solvent is matched with starting conditions





Figure 2. Effect of Strong Diluent Injections – Bandbroadening

Methanol diluted sample, gradient from 30% ACN initial condition, 3x50mm 1.8µm column

Flexibility in injection options – more freedom of choice - *extending the capabilities*



Applications:



- a) Automated needle seat back-flushing for lower carry over*
- b) Fixed loop injection mode for fastest cycle times*
- c) By-pass metering device for lower delay volume*
- d) ...and more*

Carry-over comparison

- *two system designs*

0.021% on popular UHPLC system (PLUNO mode)

0.002% to 0.004% on another UHPLC without needle seat backflushing



Column: 2.1x50mm BEH C 18

Design Considerations for Column Compartments

- Column management: use any analytical column size or hardware
- Use fully porous, superficially porous or monolithic materials
- Temperature control: sub-ambient to 90C for all application types with independent temperature zones for each column compartment
- Valve configurations: 6-port, 10-port, two position valves for 2 column selection, including alternating column regeneration (10 port only)
- With 2 column compartment modules, 9-port eight position column selection valves, you may expand column selection to 4 columns with independent thermal zones, or up to 8 columns total
- Change valve type in minutes using universal drive valves and RFID identified liquid ends

UHPLC Therm. Column Compartment

- look for flexibility:



RFID-tag, stores
type, pressure
range, valve
switches

New 2PS/6PT, 2PT/10PT and 9PT/8PS high pressure valves





Design Considerations for UHPLC Optical Detection

- Maximizing sensitivity without excessive cell volume
- Flow rate and gradient compositional stability
- Refractive index correction approaches
 - - direct correction through optical design,
 - - acquisition with reference wavelength,
 - - baseline subtraction from reference chromatogram

Long Path Cartridge Cell

- *Optofluidic waveguides*




High Light Transmission due to Total-Internal Reflection (TIR) principle (~ 100 % Light efficiency) supports longer path length cell designs

Benefits:

- *High sensitivity (S/N) with small cell volumes (dispersion effects)*
- *Coating-free fused silica (no special care instructions or smiling baseline effects)*
- *Easy cell selection (one cell for all major applications)*
- *Cartridge design for ease of use*

Effect of Data Acquisition Rate (time constant) Peak Width, Resolution and Peak Capacity in Ultra-Fast LC



80Hz **20Hz**

- 30% Peak Width
- + 30% Resolution
- + 40% Peak Capacity
- + 70% Apparent Column Efficiency

80Hz **10Hz**

- 55% Peak Width
- + 90% Resolution
- + 120% Peak Capacity
- + 260% Apparent Column Efficiency

Conclusion –UHPLC Design Requirements

- Look for best overall performance and feature set according to your application requirements (routine vs. method development, etc.)
- Choose flexible injection designs that allow variable or fixed loop capability, should usage priorities change
- Versatile column compartments support a wide range of applications on one common platform
- Understand available column technologies to ensure a good mix of performance and chemical selectivity

Pharmaceuticals in Water

- Rising concern about pharmaceuticals in surface water
- Wide range of analytes at low levels
- High speed MS sampling and mass resolution needed



Choose the Appropriate MS Technology

MS/MS

- Target compound analysis
- MRM
- Quantitation in complex matrices

TOF

- Identification of unknowns
- High mass accuracy
- Fast acquisition
- Full scan

Note: MRM with two transitions is normally required

UHPLC is Important for Environmental Chemistry

Many analytes
in complex
matrices




Maximum peak
capacity and
resolution




1.8 micron
columns and
high pressure

Agilent ZORBAX
RR-HT columns – 600 bar
RR-HD columns – 1200 bar

Improve Group 1 Analysis with 1.8 Micron




EPA method
2.1 x 100 mm
3.5 micron
22 minutes



Modified method
2.1 x 100 mm
1.8 micron
Eclipse Plus C-18
750 bar
11 minutes


Higher pressure with reduced particle size – more power

Improve Further with Gradient Change



2.1 x **100 mm**
 1.8 micron
 Eclipse Plus C-18
750 bar
 11 minutes

Time	Composition
0.0	10% ACN
3.4	10% ACN
20.0	100% ACN
20.3	100% ACN




2.1x **50 mm**
 1.8 micron
 Eclipse Plus C-18
375 bar
 6 minutes

Time	Composition
0.0	10% ACN
1.7	10% ACN
10.0	100% ACN
10.3	100% ACN

shorter column further increases throughput

Fast Resolution of 12 Tetracyclines




2.1x 50 mm
1.8 micron
Eclipse Plus C-18
375 bar
6 min

Time	Composition
0.0	10% ACN
1.7	10% ACN
10.0	100% ACN
10.3	100% ACN

Good resolution with shallow gradient on short column

Peak Capacity for 90 Pharmaceuticals in 20 min




2.1 x 50 mm, 1.8 micron Eclipse Plus C-18,
750 bar observed operating pressure at 0.6 ml/min
1290 Infinity LC

Time	Composition
0.0	10% ACN
3.4	10% ACN
20.0	100% ACN
20.3	100% ACN

“Infinity” resolution for complex mixtures

Ultrafast Chromatography for 15 Analyte Subset




2.1x 50 mm
1.8 micron
Eclipse Plus C-18
750 bar
1 minute

Time	Composition
0.0	10% ACN
1.5	100% ACN

Ultimate speed on a short column with ballistic gradient

Rapid Analysis of Boulder Wastewater Effluent



2.1x 50 mm, 1.8 micron Eclipse Plus C-18,
running at 750 bar

Time	Composition
0.0	10% ACN
1.5	100% ACN

It works on real samples, too!!

Pesticides

- Number of Compounds limited with Fast Target Analyses
- Rapid screening of many pesticides needed
- Complex matrices make identification difficult



Find These Compounds in a Sample

Methamidophos	Carbofuran	Forchlorfenuron	-Mefenacet	Phenthoate	Triflumuron	Fenpyroximate(Z)
1-Naphthol,	Formetanate	-Linuron	-Ditalimfos	Iprovalicarb	Isoxaflutole	Dioxathion
Methomyl	Acetamidrid	-Fludioxonil	Desmedipham	Pyriproxyfen	Clethodim	Fluoxastrobin
Fenuron	Aldicarb-sulfone	-Clothianidin	Phenmedipham	Sulfotep	toxazole	Fluazinam
Acephate	Butoxycarboxim	-Thiacloprid	Bifenazate	Sulprofos	Flufenacet	Temephos
Fuberidazole	Mexacarbate	-Imidacloprid	-Fenhexamid	EPN	Benzoximate	Butafenacil
Propamocarb	Monocrotophos	Thiobencarb	-Flutriafol	Flutolanil	Topramezone	Hydramethylnon
Tricyclazole	Bendiocarb	Metobromuron	-Fenoxycarb	Cyazofamid	Pyridaben	Lufenuron
Butocarboxim	Dioxacarb	-Demeton-S Phorate	-Furalaxyl	Famphur	Picoxystrobin	Indoxacarb
Carbendazim	Mepanipyrim	Diethofencarb	-Methidathion	Diniconazole (E)	Methoxyfenozide	Moxidectin
Isoprocarb	Mevinphos	-Mepronil	-Clofentezin	Diniconazole (Z)	Spiromesifen	Spinosyn A
Cymoxanil	Ethiofencarb	-Nitenpyram	-Fenamiphos	Benalaxyl	Tetraconazol	Spinosyn D
Cycluron	Methiocarb	Methoprotryne	Fenpropimorph	Dimoxystrobin	Pyrazophos	Avermectin B1b
Pyrimethanil	Cyprodinil	-Neburon	-Diazinon	Etaconazole	Famoxadon	Emamectin B1b
Diamidafos	Prometon	-Fenthion	-Pirimifos-methyl	Diclobutrazol	Bromuconazole	Avermectin B1a
Thiabendazole	Secbumeton	-Oxadixyl	-Buprofezin	Fenarimol	Fluquinconazole	Ivermectin B1a
Carbaryl	Terbumeton	-Metalaxyl	-Fenazaquin	Malathion	Prochloraz	Emamectin B1a
Dinotefuran	Ametryn	-Propetamphos	-Quinoxifen	Pirimifos-ethyl	Furathiocarb	Doramectin
Aldicarb-sulfoxide	Tebuthiuron	Penconazole	-Tebuconazole	Ipconazole	Ethion	Eprinomectin B1a
Isoproturon	Dimethoate	Ethofumesate	Fensulfothion	Tebufenpyrad	Sulfentrazone	
Promecarb	Fonicamid	Vamidothion	-Ediphenphos	Zoxamide	Pyraclostrobin	
Aminocarb	Diuron	-Terbufos	-Diflubenzuron	Fenbuconazole	Dimethomorph(E)	
Propoxur	Fluometuron	- Iprobenfos	-Fenamidone	Bitertanol	Dimethomorph(Z)	
Chlortoluron	Siduron	-Myclobutanil	Triazophos	Piperonyl butoxide	Rotenone Ethiprole	
Omethoate	Carbetamide	Chloroxuron	Hexaconazole	Pyridafenthion	Flucarbazone	
Simetryn	Dicrotophos	Thiamethoxam	Kresoxim methyl	Propiconazole	Alanycarb	
Monolinuron	Carbofuran,	-Parathion	Nuarimol	Thiophanate-methyl	Pinoxaden	
Metribuzin	-3 hydroxy	Cyproconazole	Tribufos	Prothioconazole	Azoxystrobin	
Pymetrozine	Pirimicarb	Uniconazole	Flusilazol	Azinphos-ethyl	Difenconazole	
Pyracarbolid	Prometryn	-Amitraz	Bupirimate	Triflumizol	Trifloxystrobin	
Thiofanox	-Terbutryn	-Triadimenol	Azinphos-methyl	Isofenphos	Benfuracarb	
Oxamyl	-Ethoprop	- Imazalil	Triticonazole	Propargite	Carfentrazone-ethyl	
Thidiazuron	-Cyanophos	-Spiroxamine	Tebupirimfos	Hexythiazax	Mandipropamid	
Methabenzthiazuron	-Fonophos	-Quinalphos	Metconazole	Tebufenoxide	Fenpyroximate(E)	

How Did We Do It?



1290 Infinity LC

High pressure



**Zorbax Eclipse Plus C18
column**

High resolution



6230 TOF with JetStream


High sensitivity



Pesticide database


Screening

Is 3 Minutes Fast Enough?



Remarkable specificity, speed and sensitivity

TOF Can Identify More Compounds in Less Time



2.1 x 50 mm x 1.8 micron
Eclipse Plus C-18
900 bar
1.5 mL/min
1290 Infinity

Time	Composition
0.0	10% ACN
1.5	100% ACN

TOF fast acquisition rates ensure maximum throughput

High Accuracy Confirmed by Database

3-minute 1290/6230 TOF Screening (224 +ve Pesticide suite)

Analyte amount (on-column)	# of compounds identified (%)	Average Mass Accuracy (ppm)
500 fg	45 (20%)	0.24
1.25 pg	96 (43%)	0.51
2.5 pg	124 (55%)	0.49
5 pg	163 (73%)	0.12
25 pg	202 (90%)	0.70
50 pg	217 (97%)	0.06
125 pg	224 (100%)	0.60

More than 70% confirmed with only 5pg on column

TOF/Q-TOF


Capabilities

- Full-Spectral Acquisition
- Accurate mass
- Fast acquisitions
- High Spectral Resolution

New 1290 Infinity LC with 6540 QTOF LC/MS




ECC of 100 Pesticides in 3 min using new *Agilent 1290 Infinity LC* with the new *Agilent 6540 QTOF*



With FULL SPECTRA Benefit from FAST RUN TIMES

Example mass spectrum from data on 3 min run with 1290 Infinity LC and 6540 QTOF.




Note the mass resolution at 10 spectra per second

How Do We Do It with a Triple Quad

Dynamic MRM Required for Overlapping Peaks

Necessary to have at least 20 points across the peak



once again, true UHPLC demands fast data acquisition

Dynamic MRM Acquisition Parameters – Customized Methods from Database


The screenshot displays the Agilent software interface for Dynamic MRM acquisition. The 'Acquisition' tab is active, showing a table of scan segments. The table includes columns for Compound Name, ISTD?, Precursor Ion, MS1 Res, Product Ion, MS2 Res, Fragmentor, Collision Energy, Ret Time (min), and Delta Ret Time. The table lists six compounds: Bentazon, 2,4,5-T, Silvex, Acifluorfen, Dinoseb, and Hexaflumuron. A context menu is open over the table, with 'Import from optimizer...' selected.

Compound Name	ISTD?	Precursor Ion	MS1 Res	Product Ion	MS2 Res	Fragmentor	Collision Energy	Ret Time (min)	Delta Ret Time
Bentazon	<input type="checkbox"/>	239.1	Unit	132	Unit	80	32	7.117	1
2,4,5-T	<input type="checkbox"/>	252.9	Unit	194.8	Unit	76	9	8.565	1
Silvex	<input type="checkbox"/>	266.9	Unit	194.9	Unit	90	5	9.325	1
Acifluorfen	<input type="checkbox"/>	360	Unit	315.9	Unit	78	5	10.119	1
Dinoseb	<input type="checkbox"/>	239.1	Unit	207	Unit	154	21	11.226	1
Hexaflumuron	<input type="checkbox"/>	459	Unit	438.9	Unit	102	5	11.423	1

*** Note that, with Dynamic MRM, the user sets the cycle time and dwell times are automatically managed.**

Dynamic MRM 11 Compounds (22 MRMs) eluting in about 1 min window


Compound Name	Precursor Ion	Product Ion	Retention Time
Cinosulfuron	414.1	183	5.579
Cinosulfuron (Q)	414.1	157	5.579
Chlorotoluron	213.1	72	5.642
Chlorotoluron (Q)	213.1	140	5.642
Atrazine	216.1	174	5.682
Atrazine (Q)	216.1	132	5.682
Carbaryl	202.1	145	5.736
Carbaryl (Q)	202.1	117	5.734
Carboxin	236.1	143	5.836
Carboxin (Q)	236.1	87	5.836
Chlorsulfuron	358.0	167	5.896
Chlorsulfuron (Q)	358.0	141	5.896
Ethiofencarb	226.1	107	5.937
Ethiofencarb (Q)	226.1	164	5.936
Dodemorph	282.3	116	6.073
Dodemorph (Q)	282.3	98	6.074
Diuron	233.0	72	6.101
Diuron (Q)	233.0	160	6.101
Cyprodinil	226.1	108	6.245
Cyprodinil (Q)	226.1	93	6.246
Difenoxurone	287.1	123	6.509
Difenoxurone (Q)	287.1	72	6.509




Typical results: 10 pg Atrazine on-column using Dynamic MRM 6460A QQQ with Jet Stream Technology

- *Avg. Signal Height: 15,650*
- *Avg. Signal Area: 50,966*
- *RSD: 3.2 %*
- *Estimated LOQ: 100 fg or less*

- *6-7 data points above FWHM*
- *3 sec FWHM, 6 sec @10% valley*
- *20 data points baseline-to-baseline*




Dynamic MRM of 300 Pesticides with Two Transitions Each



Run with 1200-SL LC 2.1 x 100 C18 Eclipse PLUS 1.8 u

Pesticide Mix run with Agilent 1290 Infinity LC : Benefit of *HIGHER PEAK CAPACITY*



LC 2.1 x 150 C18 Eclipse PLUS HD 1.8 μ , 0.6 mL/min, 800 bar

Conclusions

- **For maximum peak capacity, choose high pressure rated columns and optimize effective gradient slope (k^*) vs. available system pressure and preferred analysis time. Use the Agilent Method Translator for modeling**
- **Use MS conditions that allow 20 cycles/second for sharp peaks of 1-3 seconds, either
LC/MS/MS by triple quadrupole with fast dwell times of 5 msec*, or
LC/TOF-MS or Q/TOF-MS and 20 cycles per second (i.e.~500 transients per second)**

LC Systems and Method Translator: <http://www.chem.agilent.com/en-US/Products/Instruments/lc/systems/1290infinitylcsystem/Pages/default.aspx>

MS systems:

<http://www.chem.agilent.com/en-US/products/instruments/ms/Pages/default.aspx>