

Reveleris™ High Resolution Flash Cartridges

Recommended Loading and Operating Conditions

- Up to 3x Increase in Recovery of Pure Compound Compared to Leading Cartridges
- Up to 2x Increase in Loading Capacity Compared to Leading Cartridges



Reveleris™ Silica Cartridges - Luer Endfittings

Column Size (g)	Suggested Loading Amount (mg)	Suggested Flow Rate* (mL/min)	Expected Back Pressure* (psi)	Single Column Volume (mL)	Part No.
4	4-800	18	30	5	5146130
12	12-2,400	36	30	17	5146131
40	40-8,000	40	30	48	5146132
80	80-16,000	60	55	125	5146133
120	120-24,000	80	60	195	5146134
330	330-66,000	120	40	445	5146135

*Based on a solvent system of 80% Hexane, 20% Ethyl Acetate



CAUTION

Do not operate Reveleris™ Flash Cartridges above the rated pressure listed on the Certificate of Analysis. Operation above the rated pressure may result in column failure and solvent leakage.

www.discoverysciences.com

GRACE® and GRACE DAVISON® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. GRACE DAVISON DISCOVERY SCIENCES™ is a trademark of W. R. Grace & Co.-Conn. REVELERIS™ is a trademark of Alltech Associates, Inc. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership. Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn. © Copyright 2008. All rights reserved.

The information presented herein is derived from our testing and experience. It is offered for your consideration and verification. Since operating conditions vary significantly, and are not under our control, we disclaim all warranties on the results that may be obtained from the use of our products. Grace reserves the right to change prices and/or specifications without prior notification.

0908 F401

GRACE