



CATALYST FOR SUCCESS

➔ EXPANDED HPLC METHOD FOR N-METHYL CARBAMATES

POST-COLUMN ANALYSIS OF PESTICIDES IN FOOD AND POTABLE WATER SAMPLES

Carbamate pesticides are widely used around the world to protect crops. In addition, they are used as biocides for industrial or other applications and in household products. Though carbamates are biodegradable compounds and bioaccumulation usually happens only to a low extent it is important to monitor produce to make sure enough time has elapsed between harvest and applying pesticides. Also, because of their high solubility, carbamates can leach into ground waters in porous soils and consequently find their way into drinking water supplies.

As part of FDA's pesticide monitoring program individual lots of domestic and imported foods and feeds are sampled and tested for pesticide residues to enforce the tolerances set by EPA. There are 11 compounds mandated by USEPA Method 531.2 for drinking water but they represent only a fraction of the carbamates that require monitoring in domestic and imported products. Methyl carbamates are separated using a reversed-phase column and then readily react with *o*-Phthalaldehyde and a mercaptan after hydrolysis to form highly fluorescence compounds. This post-column reaction is the basis for official EPA Method 531.2 and AOAC Method 985.23.

This new expanded method is suitable for detecting a wide range of carbamates; post-column derivatization with fluorescence detection is a sensitive and selective method for residue analysis in water, food and feed samples. This method employs the same HPLC and post-column equipment and chemicals as USEPA Method 531.2 and will allow laboratories to increase the range of tested compounds.

The separation is achieved on a C₈ stationary phase with a water/Methanol gradient. Differences in selectivity of a water/Acetonitrile gradient may be used for confirmation.

METHOD

Analytical Conditions

Column: Expanded resolution C₈ analytical column 4x250 mm,
P/N 0840250 Guard column P/N 18ECG001

Flow Rate: 0.8 mL/min

Column Temperature: 50 °C

Mobile Phase: Water/Methanol

Post-Column Conditions

Post-Column System: Pinnacle PCX or Vector PCX

Reactor Volume: 0.5 mL

Reactor Temperature: 100 °C

Reagent 1: CB130

Reagent 2: 100 mg of OPA and 2 g of Thiofluor in 950 mL
of CB910 Diluent

Reagents Flow Rate: 0.3 mL/min

Detection:

Fluorescence Detector

λ_{ex} = 330 nm, λ_{em} = 465 nm




HPLC Gradient:

TIME	WATER %	METHANOL %
0	85	15
2	85	15
42	30	70
46	30	70
46.1	0	100
50	0	100

Equilibration: 10 min

ACKNOWLEDGEMENTS:

Pickering Laboratories would like to thank John Casanova and Frank Schenck (FDA) for their help with this project.



List of Carbamates

1 Aldicarb sulfoxide	13 3-Hydroxycarbofuran	25 Banol
2 Butoxycarboxim	14 Butocarboxim	26 Isoprocarb (MIPC)
3 Aldicarb sulfone	15 Aldicarb	27 2,3,5-Trimethacarb
4 Oxamyl	16 Metolcarb	28 3,4,5-Trimethacarb
5 Methomyl	17 Cloethacarb	29 Fenobucarb (BPMC)
6 Ethiofencarb sulfone	18 Bendiocarb	30 Methiocarb
7 Ethidimuron	19 Carbetamide	31 BDMC
8 Thiofanox sulfoxide	20 Propoxur	32 Bufencarb
9 Thiofanox sulfone	21 Carbofuran	33 Promecarb
10 Formetanate HCl	22 Carbaryl	34 Zectran
11 Ethiofencarb sulfoxide	23 Ethiofencarb	35 Aminocarb
12 Dioxacarb	24 Thiofanox	