

Fast determination of residual glycerol and glycerides in biodiesel by SFC/MS using the Agilent 1260 Infinity Analytical SFC System

Application Note

Petrochemical

Authors

Maria Rambla Alegre,
Melissa N. Dunkle, Gerd Vanhoenacker,
Frank David, and Pat Sandra
Research Institute for Chromatography
Kennedypark 26
B-8500 Kortrijk, Belgium

Martin Vollmer
Agilent Technologies Inc.
Waldbronn, Germany



Abstract

Supercritical fluid chromatography (SFC) in combination with mass spectrometric detection (MS) is a valuable technique for the determination of glycerol and mono-, di-, and triglyceride impurities at trace levels in biodiesel. The diluted biodiesel sample is analyzed without derivatization on a cyanopropyl silica column. Using MS detection in fast polarity switching positive and negative ion detection mode, detection limits of 0.02% (w/w) were reached for glycerol and all glycerides (mono-, di-, and tri-). Good repeatability was obtained, which allowed the system to be used for qualitative as well as for quantitative analysis, with an analysis time of less than 6 minutes.



Agilent Technologies

Introduction

The European Norm EN 14214 specifies requirements and test methods for fatty acid methyl esters (FAMES) to be used as automotive fuel for diesel engines¹. The FAME mixtures, referred to as B100 biodiesel, are typically produced from vegetable oils, mainly consisting of triglycerides. After conversion, the B100 biodiesel should be free from lipid contaminants. The maximum residue levels are respectively 0.8% (w/w) for monoglycerides (MG), 0.2% for diglycerides (DG) and triglycerides (TG), and 0.02% for free glycerol (Gly). The reference method for quality control EN 14105² is based on derivatization of glycerol and mono- and diglycerides into silyl derivatives in the presence of pyridine and N-trimethylsilyltrifluoroacetamide (MSTFA), followed by high temperature gas chromatography on a short capillary column with thin film thickness, using on-column injection (or equivalent), and flame ionization detection².

SFC is a valuable alternative to the GC-FID method. Previously, it was demonstrated that the analysis of mono-, di-, triglycerides and glycerol in biodiesel can be performed on a cyano column in 6 minutes using the Agilent 1260 Infinity Analytical SFC System in combination with evaporative light scattering (ELSD) detection³. The SFC-ELSD method performs well for fast quality screening and glycerides can be detected at levels down to 0.1% (w/w). Using ELSD, glycerol cannot be detected below 0.1% level and differentiation of different glycerides (with, for instance, palmitin, olein, or linolein fatty acids) is not possible.

This Application Note used SFC in combination with mass spectrometric detection (MS) with a single quadrupole system operated in fast polarity switching mode. The nonderivatized compounds can be detected at trace concentration and the mass spectra can be used for additional confirmation and identification of the impurities.

Experimental

Solutions

Stock solutions of monopalmitin (P), dipalmitin (PP), and tripalmitin (PPP) were prepared from fine chemicals (Sigma-Aldrich) dissolved in chloroform (20 mg/mL). A glycerol stock solution was prepared in ethanol (20 mg/mL). Stock solutions of mono-olein (O), diolein (OO) and triolein (OOO) in pyridine (5 mg/mL) were obtained from Sigma-Aldrich.

Aliquots of these stock solutions were spiked in a B100 biodiesel sample diluted at 10 mg/mL in ethanol to obtain two series of calibration samples: Series 1 contained P, PP, PPP, and glycerol at different levels from 0.05–5% (mass versus biodiesel mass) in the biodiesel solution, and Series 2 contained O, OO, OOO, and glycerol (0.05–5% mass versus biodiesel mass) in biodiesel solution.

Instrumentation

Analyses were performed on an Agilent 1260 Infinity Analytical SFC System coupled to an Agilent Agilent 6130 Single Quadrupole LC/MS (G6130A) equipped with an APCI source. The MSD detector was coupled to the SFC module by means of a heating device (Caloratherm) prior to the APCI ionization source. Additionally, a make-up flow was added between the UV detector and the BPR, through an Agilent Zero dead volume T-piece, as previously described⁴. Table 1 summarizes the system components, and Figure 1 shows the instrumental configuration.

Part number	Description
G4309A	Agilent 1260 Infinity Analytical SFC System
G1310B	Agilent 1260 Infinity Isocratic Pump (make-up flow)
G6130A	Agilent 6130 Single Quadrupole LC/MS
AG1	Caloratherm Available through RIC ^[5]
AG004	Pre-heater Available through RIC ^[5]

^[5] Contact info@richrom.com for more information

Table 1
Agilent 1260 Infinity Analytical SFC System – System modules.

Experimental

Separation is done on an Agilent ZORBAX SB-CN, 4.6 × 250 mm, 5- μ m column (p/n 880975-905) using methanol as modifier. Detection was performed by MS in scan/SIM mode using fast polarity switching. Three signals were acquired: scan positive mode is used for detection and identification of triglycerides, scan negative mode is used for detection and identification of mono- and diglycerides, and the SIM negative ion detection mode is used for the selective and sensitive detection of free glycerol. This compound is measured as the $[M-CH_3OH]^-$ adduct (m/e 124). The experimental conditions are summarized in Table 2.




Figure 1
1260 Infinity Analytical SFC/6130 MSD configuration.

Conditions

Column:	Agilent ZORBAX SB-CN, 4.6 × 250 mm, 5 μ m
Supercritical fluid:	CO ₂
Modifier:	MeOH
Outlet pressure:	150 bar
Flow rate:	3 mL/min
Gradient:	0–1.3 minutes: 2% and 1.3 to 3.2 minutes: 2%–10%
Temperature:	60 °C
Injection volume:	5 μ L
Make-up:	IPA at 0.6 mL/min (before BPR)
Caloratherm:	60 °C
APCI:	Capillary V \pm 4000 V Corona I: 4.0 μ A (+), 20 μ A (–) Drying gas: 6.0 L/min at 325 °C Nebulizer: 55 psig Vaporizer: 350 °C
APCI:	Scan pos (200–1,000) Frag 70 Scan neg (50–800) Frag 70 SIM neg (ion 124.2) Frag 70

Table 2
Experimental conditions.

Results and discussion

It has been demonstrated that, by using SFC on a cyano column, a very interesting separation is obtained with the FAMES (matrix) eluting first, following by triglycerides, diglycerides, monoglycerides, and finally, free glycerol³. Column choice and temperature are important parameters to obtain this separation. At 60 °C, optimum separation was obtained between the FAME fraction and the first eluting triglycerides (tripalmitin, PPP). Figure 2 illustrates this, showing the separation of a biodiesel sample spiked at 1% (w/w) with tripalmitin (PPP), dipalmitin (PP), monopalmitin (P), and glycerol. The highest response for triglycerides is obtained in positive ion detection mode, while di- and mono-glycerides give higher response in negative ion detection mode. Glycerol is clearly detected at 4.9 minutes in the SIM trace. Interestingly, the FAME matrix is hardly detected in any of these traces. Similar profiles are obtained for biodiesel spiked with triolein (OOO), diolein (OO), mono-olein (O), and glycerol, as illustrated in Figure 3.




Figure 2
Biodiesel (10 mg/mL) spiked with P, PP, PPP, and glycerol at 1% (w/w). (a) Positive ion TIC, (b) Negative ion TIC and (c) EIC from SIM NEG (m/e 124.2 for Glycerol).




Figure 3
Biodiesel (10 mg/mL) spiked with O, OO, OOO, and glycerol at 1% (w/w). (a) Positive ion TIC, (b) Negative ion TIC and (c) EIC from SIM NEG (m/e 124.2 for Glycerol).

The performance of the SFC-MS method was evaluated by analyzing dilution series and performing a repeatability test. Table 3 summarizes the results. Good linearity was obtained in a concentration range from 0.05 to 5% spiked level. Repeatability is in the order of 5% RSD (obtained on extracted ion signals from scan TICs) for the glycerides and is better than 2% for glycerol (SIM acquisition). The detection limit was below 0.05% for all mono-, di-, and triglycerides, and 0.02% for glycerol. Figure 4 illustrates the sensitivity obtained for trace detection of free glycerol showing the analysis of a blank (glycerol free) biodiesel and a sample spiked at 0.05% level (normative level). The glycerol peak can easily be detected and quantified.

	Linearity (R ²) ⁽¹⁾	Range % (w/w)	Repr. (% RSD)
P	0.9986	0.05-2	5.2 ⁽²⁾
PP	0.9952	0.05-2	4.5 ⁽²⁾
PPP	0.9919	0.05-2	6.2 ⁽²⁾
O	0.9924	0.05-5	5.6 ⁽²⁾
OO	0.9972	0.05-2	7.2 ⁽²⁾
OOO	0.9930	0.05-2	5.1 ⁽²⁾
Glycerol	0.9986	0.05-5	1.4 ⁽³⁾

⁽¹⁾ 0.05, 0.1, 0.2, 0.5, 1, 2, and 5% spiked level on biodiesel sample, one injection/level

⁽²⁾ Six consecutive injections of 1% (Scan)

⁽³⁾ Six consecutive injections of 0.2% (SIM)

Table 3
Method performance data.




Figure 4
SIM NEG trace (m/e 124) of (a) biodiesel blank and (b) biodiesel spiked glycerol at 0.05% (w/w) in biodiesel (10 mg/mL).

As seen in Figures 2 and 3, complete separation of different triglycerides, diglycerides and monoglycerides is not achieved under these conditions. Triolein and diolein, for instance, co-elute (Figure 3). For quality control of biodiesel, this separation is, however, less important compared to the separation between FAMES and glycerides. Moreover, since detection is done by MS in scan mode, further elucidation and characterization of glyceride contaminants is possible. Figure 5 illustrates this, showing the extracted ion chromatograms from the analysis of the biodiesel sample spiked with olein-glycerides. Using selective ions, triolein, diolein, and monoolein can easily be differentiated.

Conclusion

The Agilent 1260 Infinity Analytical SFC System, in combination with the Agilent 6130 LC/MS Single Quadrupole System was successfully used for the analysis of impurities in biodiesel. Trace levels of free glycerol, mono-, di-, and triglycerides can be determined in biodiesel samples without the need for derivatization. Good sensitivities, complying with EU regulation, are obtained for all solutes using scan/SIM mode and fast polarity switching. Both quantification and identification of the compounds is possible within an analysis time of 6 minutes.




Figure 5
Biodiesel (10 mg/mL) spiked with 0, 00, 000, and glycerol at 1% (w/w). (a) EIC from TIC POS for 000 detection, (b) EIC from TIC NEG for 00 detection, (c) EIC from TIC NEG for 0 detection.

References

1.
Automotive fuels-fatty acid methyl esters (FAME) for diesel engines - Requirements and test methods (per EN 14214:2002 E) (Brussels, Belgium).
2.
Fat and oil derivatives – Fatty Acid Methyl Esters (FAME) - Determination of free and total glycerol and mono-, di-, triglyceride content (Reference method). Norm EN 14105:**2011**.
3.
M. Rambla Alegre, M.N. Dunkle, F. David, P. Sandra, Fast screening of impurities in biodiesel using the Agilent 1260 Infinity SFC in combination with ELSD, Agilent Technologies application note XXX, **2012**.
4.
M. Dunkle, G. Vanhoenacker, F. David, P. Sandra, M. Vollmer, The Agilent 1260 Infinity SFC/MS Solution, Agilent Technologies Publication Number 5990-7972EN, **2011**.
5.
The Caloratherm (p/n AG1) and pre-heater (p/n) are products from SandraSelerity Technologies and are available through RIC. For more information, contact info@richrom.com.

www.agilent.com/chem/sfc

© Agilent Technologies, Inc., 2012
Published in the USA, October 1, 2012
5991-1143EN



Agilent Technologies