

High-Throughput Analysis of Epigenetic Targets with Agilent RapidFire/MS Systems: p300 Histone Acetyltransferase (HAT)

Application Note

Author

Peter Rye, Lauren Frick,
and William LaMarr
Agilent Technologies, Inc.
Wakefield, MA USA

Introduction

The field of epigenetics focuses on the investigation of enzymes that alter gene expression through modification of their target substrates, frequently by adding or removing acetyl groups. Historically, studying the activity of such enzymes in a high-throughput manner has proved challenging due to the relatively small molecular alterations as well as the possibility of sequential modifications leading to multiple analytes of interest. The p300 protein, for example, can serially acetylate multiple lysine residues on the C-terminus of the p53 protein. As such, high-throughput bioassays that allow the direct, concurrent quantification of multiple modification states are advantageous.^{1,2}



Agilent Technologies

Rapid, Direct Measurement of Multiple Reaction Products

In Figure 2, we demonstrate the p300-mediated sequential acetylation of a p53 peptide (Amino acids 368-386) by RapidFire/MS/MS. Analysis by RapidFire/MS/MS enabled the fast and direct measurement of multiple acetylated species concurrently.

In addition, we analyzed whole histone proteins by RapidFire/TOF to illustrate a range of options for measuring HAT activity by high-throughput mass spectrometry. Figure 3 shows the inhibition of p300 by garcinol using whole histones as the substrate. Deconvoluted RapidFire/TOF data were acquired for 18 reactions (endpoint IC_{50} for 6 concentrations in triplicate) and indicated increasing acetylation with decreasing inhibitor concentration.




Figure 2. Sequential acetylation of a p53 (368-386) peptide by p300 over time.




Figure 3. Inhibition of p300 by garcinol when using whole histone proteins as a substrate. Deconvoluted RapidFire/TOF data show increasing acetylation (going towards back) with decreasing inhibitor concentration.

Further evidence of the benefit of the RapidFire/MS technique is demonstrated in Figure 4. Here the data shows the inhibition of p300 acetylation by curcumin, garcinol, and anacardic acid as determined by RapidFire/MS. The calculated IC_{50} values are in good agreement with those previously reported.^{4,5}

Conclusions

The RapidFire/MS/MS and RapidFire/TOF systems were used to analyze p300 reactions at a sustained rate of ~7 seconds per sample. The label-free methods circumvented the need for non-native surrogate substrates and radioactive methodologies.

RapidFire/MS/MS enabled direct detection of peptide reactants while RapidFire/TOF enabled direct detection of whole protein species. Analysis of the reactions established kinetic parameters suitable for determining IC_{50} values that agreed well with those previously reported.^{2,3}

Epigenetic changes which involve sequential modifications are an ideal application for the RapidFire System because RapidFire/MS can directly detect multiple states of peptides and proteins discretely.




Figure 4. Inhibition of p300 by curcumin, garcinol, and anacardic acid as determined by RapidFire/MS techniques.

References

1. Rye, PT, *et al.* Advances in Label-Free Screening Approaches for Studying Sirtuin-Mediated Deacetylation. *J Biomol Screen.*, **2011**, Sep 12. [Epub ahead of print].
2. Rye, PT, *et al.* Advances in Label-Free Screening Approaches for Studying Histone Acetyltransferases. *J Biomol Screen.*, **2011**, Sep 9. [Epub ahead of print].
3. Prives, C, and Manley, JL. Why is p53 acetylated? *Cell*, **2001**, 107(7):815-8.
4. Balasubramanyam, K, *et al.* Curcumin, a novel p300/CREB-binding protein-specific inhibitor of acetyltransferase, represses the acetylation of histone/nonhistone proteins and histone acetyltransferase-dependent chromatin transcription. *J Biol Chem.*, **2004**, 279(49):51163-71.
5. Reaction Biology Corp. (2010). p300 Datasheet. Retrieved April 26, 2011, from <http://www.reactionbiology.com/datasheets2010/p300.pdf>.

www.agilent.com/lifesciences/rapidfire

For research purposes only and not for use in diagnostic procedures. The information described here is intended for reference and research purposes only. Agilent Technologies offers no guarantee as to the quality or suitability of this data for your specific application.

Information, descriptions, and specifications in this publication are subject to change without notice.

©Agilent Technologies Inc., 2011
Published in the USA, November 21, 2011
5990-9343EN



Agilent Technologies