

## Headspace Analysis of Tea

### Introduction

Tea is one of the most common drinks worldwide.

The analysis of the composition of the volatile organic compounds in the headspace of tea allows to classify the tea samples with respect to different criteria like storage age, quality grade or heritage.

Relevant olfactory volatiles in tea are e.g. 3-hexenal, 2,3,-butandione, 2-nonenal and monoterpenes. As the IMS technology (IMS = ion mobility spectrometry) is very sensitive for these kind of compounds, the FlavourSpec® (figure 1) was used for the analysis with respect to the differentiation of tea samples. Data analysis was performed using G.A.S. the software suite.



Figure 1: FlavourSpec® made by G.A.S. mbH

### Experimental

All measurements were carried out with the FlavourSpec®, a GC-IMS equipped with an autosampler with headspace option. For an efficient separation an isothermal heated multi capillary column is used (OV-5, 5% diphenyl - 95% dimethyl polysiloxane, 20 cm, 0.2 µm). Ionisation source is Tritium (<sup>3</sup>H). The activity (300 MBq) is below the threshold of 1 GBq so that no licence is required in all EURATOM countries. The measurement parameters are listed in table 1.



Figure 2: Different Samples of Tea

For analysis purpose 1g of each tea was transferred into a 20ml headspace vials. Each sample (Table 2) was put into two vials for reproducibility testing (index a and b). The vials were then placed onto the FlavourSpec® for automated processing and analysis. The instrument's analysis parameters were set to default (Table 1).




Figure 3: Sample preparation


**Table 1: Experimental Parameters**

|                     |                                    |   |
|---------------------|------------------------------------|---|
| <b>FlavourSpec®</b> |                                    |   |
|                     | Polarity IMS                       | <b>Positive</b>   |
|                     | T1 (IMS)                           | <b>45 °C</b>  |
|                     | T2 (MCC-OV5, 20cm)                 | <b>40 °C</b>  |
|                     | T3 (Injektor)                      | <b>80 °C</b>  |
|                     | E1 (Driftgas flow rate)            | <b>150 ml/min<br/>Nitrogen 5.0</b>  |
|                     | E2 (Carriergas flow rate, dynamic) | <b>25 ml/min (step<br/>to 150ml/min after<br/>10min)<br/>Nitrogen 5.0</b> |
|                     | Run time                           | <b>20 min</b>   |
|                     | Average                            | <b>6</b>  |
| <b>Agitator</b>     |                                    |   |
|                     | Incubation temperature             | <b>90 °C</b>  |
|                     | Incubation time                    | <b>10 min</b>   |

**Results**

Figure 4 exemplary shows the IMS chromatograms of two samples (I and II) of the analysed tea.

Every tea sample gives a characteristic peak pattern.


**Table 2: List of analysed Tea Samples**

| Tea samples |                           |               |                 |
|-------------|---------------------------|---------------|-----------------|
| NO.         | Sample name               | Grade         | Category number |
| 1           | 2012 Sichuan® Green Tea   | Special Grade | G-1             |
| 2           | 2013 Ningde Green Tea     | Special Grade | G-4             |
| 3           | 2013 Anhui Green Tea      | Special Grade | G-5             |
| 4           | 2013 Fuding Jasmine Tea   | Special Grade | S-1             |
| 5           | 2013 Wuyi Black Tea       | Special Grade | B-1             |
| 6           | 2013 Jianou Narcissus Tea | Level A       | O-2             |
| 7           | 2013 Anxi Tie Guanyin Tea | Special Grade | O-12            |
| 8           | 2013 Anxi Tie Guanyin Tea | Level A       | O-13            |
| 9           | 2013 Anxi Tie Guanyin Tea | Level B       | O-14            |
| 10          | Pu'er Tea-20 years        | Level A       | D-4             |
| 11          | Pu'er Tea-15 years        | Level B       | D-5             |
| 12          | Pu'er Tea-8 years         | Level C       | D-6             |
| 13          | Pu'er Tea-5 years         | Level D       | D-7             |
| 14          | White Tea                 | Level A       | W-1             |

In order analyse/classify the samples the intensities of 101 signals were selected by marking evaluation areas on characteristic signals. The criterion for selection is the variation of individual signals in the chromatograms of the different samples. Thus only those compound signals are selected that vary in between the samples (Figure 5).




Figure 5: Exemplary clipping of chromatogram of sample D7 with indicated evaluation areas

The following figure gives an overview on the variation of the selected signals by plotting the selected evaluation areas for all samples using the LAV plug-in 'Gallery-Plot'. Note that for each tea two equal samples are analysed. The tea classes are obvious. Each of the reproducibility measurements validates the high reproducibility of the system (Figure 7).




Figure 6: 'Gallery-Plot' of the evaluation area selection: Samples in rows, areas in columns

The intensities of all signal peak areas automatically can be determined for all samples. Due to the huge number of areas and samples the presentation of the result has to be omitted here. Exemplarily, the intensities of the signal areas for 5 selected peaks (G28, S2, B8, D9 and W3) are presented in table 3.

Table 3: Automatically retrieved intensities of five selected evaluation areas (values given are in mV)

|        | "G28" | "S2" | "B8" | "D9" | W3   |  |
|--------|-------|------|------|------|------|--|
| G-1-a  | 171   | 69   | 102  | 184  | 363  |  |
| G-1-b  | 165   | 69   | 91   | 194  | 409  |  |
| G-4-a  | 251   | 63   | 107  | 181  | 363  |  |
| G-4-b  | 224   | 79   | 104  | 194  | 348  |  |
| G-5-a  | 340   | 63   | 113  | 184  | 348  |  |
| G-5-b  | 340   | 62   | 102  | 184  | 348  |  |
| D-4-a  | 78    | 71   | 93   | 828  | 297  |  |
| D-4-b  | 78    | 68   | 98   | 786  | 323  |  |
| D-5-a  | 76    | 58   | 98   | 1000 | 297  |  |
| D-5-b  | 79    | 61   | 100  | 847  | 313  |  |
| D-6-a  | 90    | 64   | 100  | 652  | 318  |  |
| D-6-b  | 94    | 64   | 93   | 630  | 323  |  |
| D-7-a  | 106   | 61   | 93   | 557  | 297  |  |
| D-7-b  | 99    | 73   | 93   | 595  | 292  |  |
| O-12-a | 914   | 55   | 139  | 238  | 681  |  |
| O-12-b | 1000  | 58   | 123  | 213  | 661  |  |
| O-13-a | 822   | 58   | 111  | 216  | 601  |  |
| O-13-b | 865   | 57   | 114  | 222  | 601  |  |
| O-14-a | 818   | 54   | 114  | 216  | 469  |  |
| O-14-b | 755   | 59   | 113  | 219  | 489  |  |
| O-2-a  | 992   | 52   | 176  | 216  | 560  |  |
| O-2-b  | 998   | 55   | 173  | 216  | 575  |  |
| W-1-a  | 39    | 58   | 160  | 267  | 1000 |  |
| W-1-b  | 40    | 62   | 141  | 248  | 979  |  |
| B-1-a  | 208   | 59   | 1000 | 229  | 575  |  |
| B-1-b  | 190   | 61   | 964  | 222  | 595  |  |
| S-1-a  | 116   | 750  | 134  | 194  | 590  |  |
| S-1-b  | 130   | 1000 | 132  | 197  | 606  |  |




Figure 7: Principal Component Analysis (PCA) of selected samples based on the signal intensities allows a classification.




Figure 8: Gallery plot and PCA of sample set

### Summary

The headspace composition of all tea samples was successfully analysed using a FlavourSpec®. The individual compounds are 2-dimensionally separated by gas-chromatography plus ion-mobility-spectrometry. A set of 101 individual signals representing headspace compounds was analysed in order to successfully determine variations between the tea samples.

A more detailed analysis was carried out within the different samples groups:

- Different habitants (green teas)
- Different qualities (Anxite guanyin teas)
- Different storage ages (Pu'er teas)

Therefore, by visual data mining (using the G.A.S. 'Dynamic PCA' software) those areas are selected, that exhibit a significant difference in intensity.

The results are displayed in figure 8. The tea classes can easily be separated by occurrence/absence and concentration, resp., of specific volatile compounds.

- The green teas (2012 Sichuan Green Tea, 2013 Ningde Green Tea and 2013 Anhui Green Tea) from different habitants can easily be differentiated based on the FlavourSpec® analysis (figure 8A).

- The quality grade of the 2013 Anxi Tie Guanyin Tea can precisely be determined by FlavourSpec® analysis (figure 8B).

- The storage age of the Pu'er Teas can clearly be determined by FlavourSpec® headspace analysis (figure 8C).