

# Agilent MassHunter Workstation – Data Acquisition for 6400 Series Triple Quadrupole LC/MS

## Familiarization Guide

Before you begin 3

Prepare your system 3

Prepare to acquire data 4

Exercise 1 – Develop an acquisition method 6

Task 1. Enter acquisition parameters and acquire data 6

Task 2. Determine precursor ion masses 11

Task 3. Find optimum fragmentor voltage for maximum response 14

Task 4. Determine product ion masses 24

Task 5. Find optimum collision energy for MRM acquisition 30

Exercise 2 – Develop a Dynamic MRM acquisition method from an MRM acquisition data file or an MRM method 33

Task 1. Create a batch file from an existing MRM data file 33

Task 2. Print a report in the Quantitative Analysis program 36

Task 3. Create a Dynamic MRM method using Update dMRM 37

Task 4. Create a Dynamic MRM method from an MRM method 39

Exercise 3 – Create a Triggered Dynamic MRM acquisition method 40

Task 1. Create a Triggered Dynamic MRM method from a Dynamic MRM method manually 40

Task 2. Add/Modify compounds in an existing database 42

Task 3. Create a Triggered Dynamic MRM method from an existing database 52

Exercise 4 – Optimize Acquisition parameters 56

Task 1. Use the Optimizer Software to optimize acquisition parameters 56

Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters 63



Use the exercises in this guide to learn how to use the Agilent 6400 Series Triple Quad LC/MS. You can do these exercises with the demo data files, SulfaDrugs, shipped with the system (in the **Data** folder of your Qualitative Analysis installation disk), or with data you acquire.

In Exercise 1, you learn how to determine the best acquisition settings for analyzing your compounds of interest. These instructions help you understand not only how to set up a worklist to optimize instrument parameters for best sensitivity in acquisition, but also how to use the Qualitative Analysis program to identify parameter values producing optimum signal response. You can also learn about the Qualitative Analysis program by using the *Qualitative Analysis Familiarization Guide* or the Qualitative Analysis online Help.

In Exercise 2, you learn how to use either an acquired data file or the Quantitative Analysis report results to update a dynamic MRM method. This method allows you to easily set up a dynamic MRM method.

In Exercise 3, you learn how to create a triggered dynamic MRM method.

In Exercise 4, you learn how to use two programs to optimize parameters. The Optimizer Software helps you optimize acquisition parameters. Specifically, it automates the selection of the best precursor ion and the fragmentor voltage for the most abundant precursor ion, selection of the best product ions, and optimization of collision energy values for each transition for a list of compounds you specify. The “Source and iFunnel Optimizer” program helps you to find the optimal source and iFunnel parameters.

## NOTE

See the *Concepts Guide* to learn more about how the triple quadrupole mass spectrometer works and why the fragmentor and collision energy voltages are important. For background information, see Chapter 3, “Agilent Triple Quad MS and Sensitivity”, in the *Concepts Guide*. See the online Help for detailed information on how the program works.

---

Each task is presented in a table with three columns:

- Steps – Use these general instructions to proceed on your own to explore the program.
- Detailed Instructions – Use these if you need help or prefer to use a step-by-step learning process.
- Comments – Read these to learn tips and additional information about each step in the exercise.

## Before you begin

Before you begin, you need to check that your system is ready. If you plan to acquire data, you also need to set up the instrument.

### Prepare your system

- 1 Check that:
  - The Data Acquisition program has been installed.
  - The LC modules and the 6400 Series Triple Quad LC/MS have been configured.
  - The performance has been verified.
  - The system has been turned on.

If these actions have not yet been done, see the *Installation Guide* for your instrument.

- 2 Copy the data files to your PC.

Copy the folder named **SulfaDrugs** in the **Data** folder on your Qualitative Analysis installation disk to any location on your hard disk. This folder contains all the data files needed for this exercise.

#### NOTE

Do not re-use the sulfa drug data files already on your system unless you know that you copied them from the originals on the disk and you are the only one using them. Data files that are already on the system may contain processed results, leading to different behavior during the exercises in this guide.

---

## Prepare to acquire data

If you do not intend to acquire data but want to learn how to use the Qualitative Analysis program for method development, you can skip this step, which tells you how to prepare the demo sample. You then do those tasks that show you how to use the Qualitative Analysis program with the sulfa drug data files shipped with the system.

**Parts List** The exercise in this guide uses this equipment and materials:

- Agilent 1200, 1260 Infinity or 1290 Infinity LC modules: well-plate sampler, binary pump, thermostatted column compartment, DAD
- Zorbax column (see [Table 1](#) on page 4)
- A 1 ng/μL concentration of the sulfa mix sample (prepared in this step)

**Table 1** Zorbax columns

Triple Quadrupole	Column Description	Film Thickness	Pore Size	Part Number
6410B, 6420, 6430, 6460 and 6490	SB-C18 2.1mm x 50mm	1.8 μm	80Å	822700-902

**1** Prepare the LC solvent.

In 1-liter reservoirs of HPLC-grade water and acetonitrile (ACN), add 1 mL of 5M ammonium formate each to make 5mM ammonium formate in water and ACN and use for the A and B channels, respectively.

**2** Prepare the sample.

**a** Add 10 μL of the sulfa mix from one of the ampoules (500 μL) to 990 μL of solvent A in a 2 mL glass sample vial so that the final concentration is 1 ng/μL.

**b** Cap the vial and place in a sample location in the autosampler.


**3** Set up the LC column.

Use the Agilent column from [Table 1](#). Other columns and instrument parameters may be used in these exercises, but some parameters may need adjustment, and the results will differ.


**4** Set the column temperature to 60°C. Lower temperatures may be used; however, the retention times will be longer, and the pump pressure may exceed the limit of some LC systems.

The Electrospray LC Demo Sample (P/N 59987-20033) contains five ampoules with 100 ng/ $\mu\text{L}$  each of sulfamethizole ( $\text{M}+\text{H}^+$  = 271), sulfamethazine ( $\text{M}+\text{H}^+$  = 279), sulfachloropyridazine ( $\text{M}+\text{H}^+$  = 285), and sulfadimethoxine ( $\text{M}+\text{H}^+$  = 311).


**Sulfamethizole**




**Sulfamethazine**



**Sulfachloropyridazine**



**Sulfadimethoxine**



**NOTE**

Determining optimal parameter values for acquiring sample compound data requires that the Agilent Triple Quad instrument already be tuned on the Tuning Mix calibrant ions. Before proceeding with this exercise, make sure you have used Checktune or Autotune to verify that calibrant ions each have the proper mass assignment, peak width, and signal intensity.

See the *Quick Start Guide*, *Installation Guide* or online Help for instructions on tuning the instrument.

## Exercise 1 – Develop an acquisition method

### Task 1. Enter acquisition parameters and acquire data

## Exercise 1 – Develop an acquisition method

For this exercise you analyze a mixture of four sulfonamide compounds.

### Task 1. Enter acquisition parameters and acquire data

In this exercise, you enter the conditions for the analysis of the sulfa drug mix.

Steps	Detailed Instructions	Comments
1 Enter LC parameters appropriate for sulfa drug mix.  See <a href="#">Table 2</a> .	<p><b>a</b> Double-click the <b>Data Acquisition</b> icon.</p> <p><b>b</b> Make sure that Acquisition appears as the selection in the <b>Context</b> text box. If Tune is the selection, click <b>Acquisition</b> from the <b>Context</b> dropdown menu in the Combo bar.</p> <p><b>c</b> Enter the LC parameters listed in the <a href="#">Table 2</a>.</p>	<ul style="list-style-type: none"><li>• The Data Acquisition window appears. See <a href="#">Figure 1</a>.</li></ul>

**Table 2** LC parameters for sulfa drug mix

Parameter	LC Parameter
<b>PUMP</b>	
• Flowrate	800 µL/min
• Solvent A	5 mM ammonium formate in water
• Solvent B	5 mM ammonium formate in 90:10 acetonitrile:water
• Gradient (min - %B)	0 min - 13% 1.80 min - 60% 2 min - 60%
• Stop Time	2.5 min
• Post Time	3.0 min
<b>INJECTOR</b>	
• Inj. Vol.	2.0 µL

**Table 2** LC parameters for sulfa drug mix (continued)

<b>Parameter</b>	<b>LC Parameter</b>
• Injection	Standard
• Draw Position	0.0 mm
<b>UV DETECTOR</b>	
• Ch A	254 nm (4 nm BW on DAD)
• REF A (DAD only)	400 nm (80 nm BW)
<b>COL THERM</b>	
• Temp	60 °C for the 6460 and 6490 with Agilent Jet Stream Technology 40 °C for other instruments

## Exercise 1 – Develop an acquisition method

### Task 1. Enter acquisition parameters and acquire data





Figure 1 Agilent MassHunter Workstation Software – Data Acquisition window

Steps	Detailed Instructions	Comments
2	<p>Enter MS parameters appropriate for sulfa drug mix and save the method as <i>iiMS2Scantest.m</i>, where <i>iii</i> are your initials.</p> <p>See Table 3.</p>	<p>a Click the <b>QQQ</b> tab in the <b>Method Editor</b> window.</p> <p>b Select <b>MS2Scan</b> from the <b>Scan Type</b> list in the Time Segments table.</p> <p>c Enter the other MS parameters as listed in Table 3. These parameters are in either the Acquisition or the Source tabs.</p> <p>d Save the method as <i>iiMS2Scantest.m</i>, where <i>iii</i> are your initials.</p>

**Exercise 1 – Develop an acquisition method**  
**Task 1. Enter acquisition parameters and acquire data**

**Table 3** MS parameters for sulfa drug mix

Parameter	Value (ESI)	Value (AJS ESI)
• Inlet	ESI (positive polarity)	AJS ESI (positive polarity)
• Scan Type	MS2Scan	MS2Scan
• Delta EMV pos	400 V	200 V
• Mass Range	100 to 400	100 to 400
• Cell Acceleration Voltage	7 V	7 V
• Gas Temp	350 °C 250 °C for Agilent 6490	350 °C 250 °C for Agilent 6490
• Gas Flow	12 L/min 14 L/min for Agilent 6490	10 L/min 14 L/min for Agilent 6490
• Nebulizer	50 psi	35 psi
• Sheath Gas Temperature	not applicable	400 °C
• Sheath Gas Flow	not applicable	12 L/min
• Nozzle Voltage	not applicable	0 V
• Capillary Voltage positive	4000 V	4000 V
• Fragmentor	100 V (not adjustable on 6490, comes from the Tune file)	100 V (not adjustable on 6490, comes from the Tune file)



**Figure 2** Select **Scan Type** of MS2 Scan in the QQQ tab

## Exercise 1 – Develop an acquisition method

### Task 1. Enter acquisition parameters and acquire data

Steps	Detailed Instructions	Comments
<p><b>3</b> Acquire data (optional).</p> <ul style="list-style-type: none"><li>• Set up a one-line worklist with the method you just created.</li><li>• Name the data file <b>iiisulfamix01.d</b>, where <b>iii</b> are your initials.</li><li>• Designate a directory path to hold your data files and method.</li></ul>	<p><b>a</b> If necessary, click <b>View &gt; Worklist</b> to display the Worklist window.</p> <p><b>b</b> Click <b>Worklist &gt; Worklist Run Parameters</b>. Verify that the parameters are set properly. Click <b>OK</b>.</p> <p><b>c</b> Click <b>Worklist &gt; Add Multiple Samples</b>.</p> <p><b>d</b> Type <b>iiisulfamix01.d</b> as the data file name</p> <p><b>e</b> Select <b>iiiMS2Scantest.m</b> as the method name.</p> <p><b>f</b> Click the <b>Sample Position</b> tab.</p> <p><b>g</b> Select the Autosampler, Well-plate or Vial Tray.</p> <p><b>h</b> In the graphic, select a single position. Click <b>OK</b>.</p> <p><b>i</b> In the Worklist window, mark the check box to the left of the sample.</p>	<ul style="list-style-type: none"><li>• The Worklist window is tabbed with the Method Editor window by default. Click the <b>Worklist</b> tab to show the Worklist window.</li><li>• The <b>Number of samples</b> is set to 1.</li><li>• You have just acquired a full scan MS data file to see what ions are being formed from the sample.</li><li>• This step is optional because you can perform the next step with an example data file that comes with the program. If you prefer, you can create your own data file as described in this step.</li></ul>
	<p><b>j</b> Click the <b>Start Worklist Run</b> icon in the main toolbar, the <b>Run Worklist</b> icon in the Worklist toolbar or click the <b>Worklist &gt; Run</b> command.</p>	




	Sample Name	Sample Position	Method	Data File
<input checked="" type="checkbox"/>	Sample 1	Vial 1	pthMS2Scantest.m	pthulfamix01.d

## Task 2. Determine precursor ion masses

In this exercise, you determine the precursor ions for each of the sulfa drugs in the acquired data file.

Steps	Detailed Instructions	Comments
<p>1 Open the acquired data file.</p> <ul style="list-style-type: none"> <li>In the Qualitative Analysis program, open either the example file, <b>sulfamix01.d</b>, or the data file you created in “<a href="#">Task 1. Enter acquisition parameters and acquire data</a>” on page 6.</li> </ul>	<p>a Double-click the <b>Qualitative Analysis</b> icon. </p> <p>The program displays the “Open Data File” dialog box.</p>	<ul style="list-style-type: none"> <li>When you open the sulfa drug directory after installation, the <b>Load result data</b> (lower left corner) check box is grayed out.</li> <li>If you see the check box marked, this means that the data file(s) already contains results. <b>Clear this check box before opening the file.</b></li> </ul>



## Exercise 1 – Develop an acquisition method

### Task 2. Determine precursor ion masses


Steps	Detailed Instructions	Comments
	<p><b>b</b> Do one of the following:</p> <ul style="list-style-type: none"><li>• Select the example data file <b>sulfamix01.d</b>, and click <b>Open</b>.</li><li>• Select the data file you created in “Task 1. Enter acquisition parameters and acquire data” on page 6, and click <b>Open</b>.</li></ul> <p>By default, the system displays the Total Ion Chromatogram (TIC).</p>	<ul style="list-style-type: none"><li>• The figure below shows the default layout.</li><li>• The Qualitative Analysis program displays a newly opened data file with the same layout and display settings used for the previous data file. Therefore, you <b>MUST</b> make sure to return to the default settings for this exercise.</li></ul>

Before you begin, make sure that all previous settings are returned to their default values:

- Restore default layouts
  - Click Configuration > Window Layouts > Restore Default Layout.
- Make sure the method is default.m. (see title bar)
  - Click Method > Open.
  - Select default.m, and click Open.
- Return display options to default settings.
  - In the Configuration menu, click each of the Display Options commands.
  - Click Default, and then OK.

Or...

- Restore the General layout.
  - Click Configuration > Configure for Workflow > General.
  - Click OK.
  - (optional) You may be asked to save method changes.
- Return display options to default settings.
  - In the Configuration menu, click each of the Display Options commands.



## Exercise 1 – Develop an acquisition method

### Task 2. Determine precursor ion masses

Steps	Detailed Instructions	Comments
-------	-----------------------	----------

2 Determine precursor ion masses for all four peaks.


- You have determined them correctly if you find the values are similar to those shown in this table:

Compound	RT	m/z
Sulfamethizole	0.47	270.9
Sulfachloropyridazine	0.88	284.9
Sulfamethazine	1.20	279.0
Sulfadimethoxine	2.23	311.0

- If you acquired the data file using the Agilent Jet Stream Technology, the retention times may be different.
- The sulfamix01.d data file was acquired with a different column so your retention times are different.
- Close the data file after finding the precursor ion masses.

- In the Chromatogram Results window, make sure that the Range Select icon in the toolbar  is on.
- Click the left mouse button and drag the cursor across the first peak to produce a shaded region, as in the figure below.
- Right-click the shaded area, and click **Extract MS Spectrum** from the shortcut menu.

- The system displays an averaged spectrum across the peak in the MS Spectrum Results window.
- The precursor mass of the first compound, sulfamethizole, is determined to be m/z 270.9.
- To obtain a single scan, double-click the apex of the peak.



- Repeat **step a** through **step c** for the other compounds. The precursor ion masses should match those in the table in step 2.
- Click **File > Close Data File**.
- When asked if you want to save the results, click **No**.

- Some compounds form sodium (Na) and/or potassium (K) adducts as well, corresponding to M + 23 and M + 39 masses respectively. Seeing these masses along with the M + H can make for an easy confirmation of which ion is the pseudo-molecular ion (M + H)+.


## Exercise 1 – Develop an acquisition method

### Task 3. Find optimum fragmentor voltage for maximum response

## Task 3. Find optimum fragmentor voltage for maximum response

Task 3 shows you how to carry out the optimization for fragmentor voltage by creating selected ion-monitoring experiments for each compound within a method and setting up multiple methods with varying fragmentor voltages.

The Fragmentor Voltage for the 6490 is set automatically during Autotune, and it cannot be set in the Data Acquisition program. If your instrument is a 6490, skip to “Task 4. Determine product ion masses”. You can do the Qualitative Analysis part of this task by using the data files that were shipped with the software.

Steps	Detailed Instructions	Comments														
<p>1 Set up six methods for six different fragmentor voltages.</p> <ul style="list-style-type: none"><li>• Change to a SIM experiment.</li><li>• Use 60, 80, 100, 140, 180 and 220 volts as the fragmentor voltages for the six methods.</li><li>• Save the methods as <b>iiiMS2SIMxxx.m</b>, where <b>iii</b> are your initials and <b>xxx</b> is the voltage.</li></ul>	<p>a In the <b>Scan Type</b> dropdown list, click <b>MS2 SIM</b>.</p>  <p>The screenshot shows the 'Properties' window for a method. The 'Scan Type' dropdown menu is open, showing options: MS2 Scan, MS2 SIM (highlighted), MRM, Dynamic MRM, Product Ion, Precursor Ion, Neutral Loss, and Neutral Gain. The 'Time segments' table below shows a single segment with 'MS2 SIM' selected in the 'Scan Type' column.</p> <table border="1"><thead><tr><th>#</th><th>Start Time</th><th>Scan Type</th><th>Div Valve</th><th>Delta EMV (+)</th><th>Delta EMV (-)</th><th>Stored</th></tr></thead><tbody><tr><td>1</td><td>0</td><td>MS2 SIM</td><td>To MS</td><td>0</td><td>0</td><td><input checked="" type="checkbox"/></td></tr></tbody></table>	#	Start Time	Scan Type	Div Valve	Delta EMV (+)	Delta EMV (-)	Stored	1	0	MS2 SIM	To MS	0	0	<input checked="" type="checkbox"/>	
#	Start Time	Scan Type	Div Valve	Delta EMV (+)	Delta EMV (-)	Stored										
1	0	MS2 SIM	To MS	0	0	<input checked="" type="checkbox"/>										

## Exercise 1 – Develop an acquisition method

### Task 3. Find optimum fragmentor voltage for maximum response


Steps	Detailed Instructions	Comments
	<p><b>b</b> In the <b>Acquisition</b> tab, enter the <b>Compound Name</b> and <b>Mass</b> (precursor ion mass) for sulfadimethoxine.</p> <p><b>c</b> Right-click anywhere in the Scan segments section, and click <b>Add Row</b>.</p> <p><b>d</b> Type the <b>Compound Name</b> and the <b>Mass</b> for sulfachloropyridazine.</p> <p><b>e</b> Repeat steps c and d for sulfamethazine and sulfamethizole.</p> <p><b>f</b> <b>Save the method as <i>iii</i>MS2SIM140.m</b>, where <i>iii</i> are your initials.</p> <p><b>g</b> Change the fragmentor voltage to 60, and save the method as <b><i>iii</i>MS2SIM060</b>, where <i>iii</i> are your initials.</p> <p><b>h</b> Repeat <a href="#">step g</a> for voltages 80, 100, 180 and 220, saving the methods as <b><i>iii</i>MS2SIM080</b>, <b><i>iii</i>MS2SIM100</b>, <b><i>iii</i>MS2SIM180</b> and <b><i>iii</i>MS2SIM220</b>, where <i>iii</i> are your initials.</p>	<ul style="list-style-type: none"> <li>With the MS2SIM Scan Type set, a different set of columns appears in the Acquisition window.</li> <li>The Instrument Control and Data Acquisition program creates a SIM experiment for each compound mass, starting with a default fragmentor voltage of 140. See the example below.</li> <li>The Fragmentor column is grayed out if the instrument type is an Agilent 6490.</li> </ul>

Acquisition									
Source	Chromatogram	Instrument	Diagnostics						
Scan segments									
Compound Name	ISTD?	Mass	MS2 Res	Dwell	Fragmentor	Cell Accelerator Voltage	Polarity		
sulfadimethoxine	<input type="checkbox"/>	311	Unit		200	140	7	Positive	
sulfachloropyridazine	<input type="checkbox"/>	285	Unit		200	140	7	Positive	
sulfamethazine	<input type="checkbox"/>	279	Unit		200	140	7	Positive	
sulfamethizole	<input type="checkbox"/>	271	Unit		200	140	7	Positive	

## Exercise 1 – Develop an acquisition method

### Task 3. Find optimum fragmentor voltage for maximum response


Steps	Detailed Instructions	Comments
<p>2 Set up and run the worklist (optional).</p> <ul style="list-style-type: none"><li>Set up six samples with Sample Name SulfaDrugMix to inject 1 ul from vials 1-6 or the ones you choose.</li><li>Specify the data files as <b>iiiSulfaSIM.xxx.d</b>, where <b>iii</b> are your initials and <b>xxx</b> is the voltage.</li></ul>	<p><b>a</b> Click the <b>Worklist</b> icon if necessary to make sure the worklist is visible.</p> <p><b>b</b> Click <b>Worklist &gt; New</b> to start a new worklist. You do not need to save the last worklist.</p> <p><b>c</b> To set up the run, right-click the upper left corner of the worklist, and click <b>Worklist Run Parameters</b>.</p> <p><b>d</b> Type the paths for the method and data files.</p> <p><b>e</b> Type the information for the 60 voltage run.</p> <p><b>f</b> Click <b>Worklist &gt; Add Sample</b>. Another sample is added to the Worklist. Add five samples to the worklist for voltages 80-220.</p> <p><b>g</b> Mark the checkbox to the left of the Sample Name for each of the six samples.</p>	<ul style="list-style-type: none"><li>This step is optional because you can use data files shipped with the system to perform many of the tasks in this exercise.</li></ul>



	Sample Name	Sample Position	Method	Data File	Sample Type	Lev
1	<input checked="" type="checkbox"/> Sample 1	Vial 1	pth:MS2SIM60.m	SulfaSim60.d	Sample	
2	<input checked="" type="checkbox"/> Sample 1	No Injection	pth:MS2SIM80.m	SulfaSim80.d	Sample	
3	<input checked="" type="checkbox"/> Sample 1	No Injection	pth:MS2SIM100.m	SulfaSim100.d	Sample	
4	<input checked="" type="checkbox"/> Sample 1	No Injection	pth:MS2SIM140.m	SulfaSim140.d	Sample	
5	<input checked="" type="checkbox"/> Sample 1	No Injection	pth:MS2SIM180.m	SulfaSim180.d	Sample	
6	<input checked="" type="checkbox"/> Sample 1	No Injection	pth:MS2SIM220.m	SulfaSim220.d	Sample	

- h** Start the worklist.
- Click **Worklist > Run**.
  - Click the  icon in the main toolbar.
  - Click the  icon in the worklist toolbar.
- Note that the program only runs those samples that are marked with a checkmark.
  - You can also run the worklist in locked mode by clicking the  button in the main toolbar.


Steps	Detailed Instructions	Comments
<p><b>3</b> Set up a qualitative method to view the EIC data automatically.</p> <ul style="list-style-type: none"> <li>Open the data file <b>Sulfa_SIM60.d</b> or your own <b>iiiSulfa_SIM60.d</b>, where <b>iii</b> are your initials.</li> <li>In the Method Editor, add in the EICs corresponding to the precursor ion masses of 271, 279, 285, and 311.</li> <li>Save the method as <b>iiiExercise1</b>, where <b>iii</b> are your initials.</li> </ul>	<p><b>a</b> Click <b>File &gt; Open Data File</b>. The system displays the Open Data File dialog box</p> <p><b>b</b> Select either <b>Sulfa_SIM60.d</b> or <b>iiiSulfa_SIM60.d</b>, and click <b>Open</b>.</p> <p><b>c</b> Click <b>Method &gt; Method Editor</b> or <b>View &gt; Method Editor</b>. The system displays the Method Editor window.</p>	<ul style="list-style-type: none"> <li>The Qualitative Analysis program should be open. If not, see <b>“Double-click the Qualitative Analysis icon.”</b> on page 11.</li> </ul>



## Exercise 1 – Develop an acquisition method



### Task 3. Find optimum fragmentor voltage for maximum response

Steps	Detailed Instructions	Comments
	<p><b>d</b> If necessary, click <b>Define Chromatograms</b> in the Chromatogram section of the Method Explorer.</p> <p><b>e</b> To delete the BPC chromatogram, click <b>Delete</b>.</p> <p><b>f</b> Select <b>EIC</b> for the <b>Chromatogram Definition Type</b>.</p> <p><b>g</b> In the MS Chromatogram tab, make sure <b>MS Level</b> is set to <b>All</b> and <b>Scans</b> is set to <b>All scan types</b>.</p> <p><b>h</b> Clear the <b>Do cycle sum</b> check box.</p> <p><b>i</b> Type 271 as the <b>m/z value</b>.</p> <p><b>j</b> Click <b>Add</b>.</p> <p><b>k</b> Repeat steps i and j for the other precursor ions, 279, 285 and 311.</p> <p><b>l</b> Click <b>Method &gt; Save As</b>. The system opens the Save As dialog box</p> <p><b>m</b> Save the method as <i>iiiExercise 1.m</i>.</p> <p><b>n</b> Click <b>Save</b>.</p>	<ul style="list-style-type: none"><li>• The default Method Editor list selection after installation is <b>Integrate (MS)</b>.</li><li>• You can also select Define Chromatograms from the Method Items list in the Method Editor window.</li></ul>




## Exercise 1 – Develop an acquisition method

### Task 3. Find optimum fragmentor voltage for maximum response

Steps	Detailed Instructions	Comments
<p>4 Extract the chromatogram for the data file and view the results.</p> <ul style="list-style-type: none"><li>• Make sure you can see all five chromatograms, the TIC and four EICs.</li></ul>	<p>a Click the <b>Run</b> button on the Method Editor toolbar.</p>  <p>b To see the TIC and four EICs, click the arrow next to the Maximum Number of List Panes icon in the Chromatogram Results toolbar, as shown in the example below.</p> <p>c Select <b>5</b> to view five chromatograms simultaneously. The system displays chromatogram results as shown below.</p>	<ul style="list-style-type: none"><li>• You can also click the <b>Chromatograms &gt; Extract Defined Chromatograms</b> command to extract the defined chromatograms.</li></ul>
		

## Exercise 1 – Develop an acquisition method


### Task 3. Find optimum fragmentor voltage for maximum response

Steps	Detailed Instructions	Comments
<p>5 Extract the remaining ion chromatograms automatically.</p> <ul style="list-style-type: none"><li>• Extract Defined Chromatograms should be the default action for Assign File Open Actions.</li><li>• Open the remaining data files, Sulfa_SIM80.d through Sulfa_SIM220.d.</li><li>• Close the Method Explorer.</li></ul>	<p>a Select <b>File Open Actions</b> from the General section in the Method Explorer.</p> <p>b Make sure that <b>Actions to be run</b> list only contains <b>Extract Defined Chromatograms</b>.</p>	<ul style="list-style-type: none"><li>• The Qualitative Analysis Method Editor lets you define actions to be performed automatically upon opening a data file(s).</li></ul>
		
	<p>c Click <b>File &gt; Open Data File</b>. The system displays the Open Data File dialog box.</p> <p>d Select the data files to be opened, Sulfa_SIM80.d through Sulfa_SIM220.d.</p> <p>e Mark the <b>Run 'File Open' actions from selected method</b> check box. (lower left corner)</p>	

Steps

Detailed Instructions

Comments



f Click **Open**.

The Qualitative Analysis program displays all the EICs for all the data files selected.

g To close the Method Explorer and Method Editor, click the **X** in the upper right corner of each window.

- You can also close the Method Explorer and Method Editor windows by clicking the **View > Method Explorer** command and the **View > Method Editor** command.


## Exercise 1 – Develop an acquisition method


### Task 3. Find optimum fragmentor voltage for maximum response

#### Steps

#### Detailed Instructions

#### Comments




Steps	Detailed Instructions	Comments
<p>6 Select the fragmentor voltage that produces the maximum response for each of the precursor ions.</p> <ul style="list-style-type: none"> <li>Close the data files after you determine the optimum voltage.</li> </ul>	<p>a In the Data Navigator window, highlight the EICs for 271.0 <i>m/z</i>.</p> <p>b Click the <b>Show only the highlighted items</b> icon, . Only the 271 <i>m/z</i> check boxes are now marked.</p> <p>c Look at the relative intensities of each peak to determine which fragmentor voltage setting will be best to use for the 271 precursor.</p>	<ul style="list-style-type: none"> <li>You press the <b>Ctrl</b> key to be able to select multiple objects from the Data Navigator window.</li> <li>You press the <b>Shift</b> key to be able to select a group of objects.</li> <li>A fragmentor voltage of 100 should be sufficient for each precursor ion.</li> <li>You can now determine the product ions that are available for the multiple-reaction monitoring experiments to maximize sensitivity</li> </ul>
		
<p>You can overlay the chromatograms by clicking the Overlaid mode icon in the Chromatogram Results toolbar.</p>		
	<p>d Repeat <a href="#">step a</a> through <a href="#">step c</a> for the other three base peaks or precursor ions.</p> <p>e Click <b>File &gt; Close Data File</b>.</p> <p>f Click <b>Close</b> when the Close Data File dialog box appears.</p>	<ul style="list-style-type: none"> <li>Click the different EICs in the Data Navigator window to change which chromatogram is labeled in the Chromatogram Results window. When the color of the label of the chromatogram matches the color of the chromatogram that has the highest intensity, you use the fragmentor voltage that was used for that file.</li> </ul>

## Exercise 1 – Develop an acquisition method


### Task 4. Determine product ion masses

## Task 4. Determine product ion masses

In this part of the method development, we will use three collision energies to determine the best fragment ions to use for the eventual Multiple Reaction Monitoring (MRMs) acquisition.

Steps	Detailed Instructions	Comments
<p><b>1</b> Set up three product ion acquisition methods and acquire data.</p> <ul style="list-style-type: none"><li>Use the MS parameters in the example below, but change the Fragmentor voltage to the optimum voltage you determined in the previous task.</li><li>Save methods as <i>iiiSulfamix PI_xx.m</i>, where <i>iii</i> are your initials and <i>xx</i> is the collision energy.</li></ul>	<p><b>a</b> Click the <b>000</b> tab in the Method Editor pane.</p> <p><b>b</b> Select <b>Product Ion</b> in the <b>Scan Type</b> combo box to scan each precursor ion for all its product ions.</p> <p><b>c</b> Enter all MS parameters as listed in the example below, making sure the <b>Collision Energy</b> is set to 1.5 and the Fragmentor voltage is set to the optimum voltage determined in Task 3.</p> <p><b>d</b> Save the method as <i>iiiSulfamix PI_15.m</i>.</p> <p><b>e</b> Repeat <b>step c</b> and <b>step d</b> for collision energies of 30 and 45.</p>	<ul style="list-style-type: none"><li>When you change the <b>Scan Type</b> in the <b>Time Segments</b> table, the <b>Scan segments</b> table is reset. If you want to copy the <b>Scan segments</b> to the new <b>Scan segments</b> table, highlight all of the lines in the <b>Scan segments</b> table and then right-click the Scan segments table and click <b>Copy</b>. After you select a new Scan Type, right-click the Scan segments table and click Paste from Clipboard.</li><li>You cannot copy and paste the <b>Scan segments</b> table between all <b>Scan Types</b>.</li></ul>
		
<p><b>2</b> Set up and run the worklist (optional).</p> <ul style="list-style-type: none"><li>Specify the data files as <i>iiiSulfamix PI_xx.d</i>, where <i>iii</i> are your initials and <i>xx</i> is the collision energy.</li></ul>	<p><b>a</b> Click the <b>Worklist</b> tab.</p> <p><b>b</b> Add three samples to the worklist for collision energies 15, 30 and 45.</p> <p><b>c</b> Mark the check box to the left of the Sample Name for each sample you are adding.</p> <p><b>d</b> Click <b>Worklist &gt; Run</b>.</p>	<ul style="list-style-type: none"><li>This step is optional because you can determine the product ion masses from the data files shipped with the system.</li><li>Use the instructions in Step 2 of Task 3 to set up the worklist.</li></ul>


Steps	Detailed Instructions	Comments
<p><b>3</b> Set up a qualitative method to integrate and extract product ion spectra.</p> <ul style="list-style-type: none"> <li>• Use the data files <b>SulfamixPI_xx.d</b>, where <i>xx</i> is the collision energy, or your own data files, <b>iiiSulfamixPI_xx.d</b>.</li> <li>• Open Method Explorer and Method Editor.</li> <li>• Use TICs set up for MS/MS, product ion and each of the precursor ions 271, 279, 285, 311.</li> <li>• Make sure the MS/MS integrator has been selected and the maximum number of peaks has been limited to the largest 100 peaks.</li> <li>• Add the ability to integrate and extract peak spectra to the file actions run upon data opening.</li> <li>• Save the changes to the current method.</li> </ul>	<p><b>a</b> Click the <b>Open Data File</b> icon in the toolbar.</p> <p><b>b</b> Select <b>SulfamixPI_15.d</b>.</p> <p><b>c</b> Make sure that the <b>Run File Open Actions from Specified Method</b> check box is clear, and click <b>Open</b>.</p> <p><b>d</b> Make sure the Method Explorer and the Method Editor windows are displayed; otherwise, click the <b>Method Explorer</b> and then <b>Method Editor</b> icons. </p> <p><b>e</b> In the Chromatogram section in the Method Explorer window, select <b>Define Chromatograms</b>.</p> <p><b>f</b> Delete any existing chromatograms in the <b>Defined Chromatograms</b> list.</p> <p><b>g</b> Select <b>TIC</b> from the <b>Type</b> list in the <b>Define chromatograms</b> section.</p> <p><b>h</b> For <b>MS level</b>, select <b>MS/MS</b>.</p> <p><b>i</b> Mark the <b>Do cycle sum</b> check box.</p> <p><b>j</b> For <b>Scans</b>, select <b>Product ion</b>.</p> <p><b>k</b> For <b>Precursor ion m/z</b>, type 271.</p> <p><b>l</b> Click the <b>Add</b> button.</p> <p><b>m</b> Repeat steps j and k for each ion.</p>	<ul style="list-style-type: none"> <li>• The Qualitative Analysis program should already be open and contain <i>iiexercise 1.m</i> as the method.</li> </ul>



## Exercise 1 – Develop an acquisition method


### Task 4. Determine product ion masses

Steps	Detailed Instructions	Comments
	<ul style="list-style-type: none"><li>n From the Method Explorer in the Chromatogram section, click <b>Integrate (MS/MS)</b>.</li><li>o Select <b>MS/MS</b> as the <b>Integrator selection</b>, if necessary.</li></ul>	<ul style="list-style-type: none"><li>• These data files contain MS/MS data, so you need to modify the parameters in the Integrate (MS/MS) section. If the data file contained only MS data, you would need to modify the parameters in the Integrate (MS) section.</li></ul>




**Figure 3** Integrate (MS/MS) > Integrator Tab

- p Click the **Peak Filters** tab. Make sure that the **Limit (by height) to the largest** check box is marked and set to the value 100 as shown below.



**Figure 4** Integrate (MS/MS) > Peak Filters tab

Steps	Detailed Instructions	Comments
	<p>q Click <b>General</b> in Method Explorer, and then click <b>File Open Actions</b>.</p> <p>r Select <b>Integrate and extract peak spectra</b> from the Available actions list and click  to add this to <b>Actions to be run</b>.</p>	
		
	<p><b>Figure 5</b> General &gt; File Open Actions tab</p>	
	<p>s To apply the changes to the current method, <i>exercise1.m</i>, click the <b>Save Method</b> icon.  You can also click <b>Method &gt; Save</b>.</p>	
<p>4 Run the qualitative method on the current data file.</p>	<ul style="list-style-type: none"> <li>• In the Method Editor toolbar, click the <b>Run</b> button, . When the Assign Actions to Run Opening A Data File section is displayed, the <b>Actions to be run</b> list is executed.</li> </ul>	<ul style="list-style-type: none"> <li>• The program first extracts the product ion chromatograms for each precursor ion in the data file.</li> <li>• Next, it finds the largest peak in the total ion chromatograms, and integrates and extracts peak spectra from each integrated peak.</li> <li>• See <a href="#">Figure 6</a> on page 28.</li> </ul>


## Exercise 1 – Develop an acquisition method

### Task 4. Determine product ion masses

#### Steps

#### Detailed Instructions

#### Comments




**Figure 6** Results for integration and extraction of peak spectra.

- 5 Run the 'File Open' actions on the remaining product ion data files.
- Use either the example files, **Sulfamix PI\_xx.d**, or the data files you acquired in [step 2](#).
- a Click **File > Open Data File**.  
The system displays the Open Data File dialog box.
- b Hold the **Ctrl** key and do one of these:
- Select the two data files **Sulfamix PI\_30.d**, and **Sulfamix PI\_45.d**.
  - Select the data files you acquired in [step 2](#).
- c Mark the **Run 'File Open' actions from selected method** check box in the Open Data File dialog box, and click **Open**.
- After the data files open, the Qual method first extracts the product ion chromatograms for each precursor ion.
  - Next, it finds the largest peak in the total ion chromatograms, and integrates and extracts peak spectra from each integrated peak.

## Exercise 1 – Develop an acquisition method

### Task 4. Determine product ion masses

Steps	Detailed Instructions	Comments
6 Identify product ions. <ul style="list-style-type: none"> <li>View each set of TICs and spectra individually (e.g., 271 m/z first).</li> <li>Close the data files.</li> </ul>	<p><b>a</b> In the Data Navigator, select the TICs and spectra for the 271 m/z precursor ion.</p> <p><b>b</b> Click the <b>Show only the highlighted items</b> icon.</p> <p><b>c</b> Click <b>View &gt; MS Spectrum Peak List 1</b>.</p> <p><b>d</b> Examine the spectra to see which fragment ions are produced at which collision energies.</p> <p><b>e</b> Repeat steps a to d until all the product ions are identified.</p>	<ul style="list-style-type: none"> <li>The m/z 155.7 product ion is the most abundant of any product ion and the highest signal is recorded at 15 V. This means that a good choice for the MRM for sulfamethizole would be 271.0 &gt; 155.7 when the collision energy is around 15 V.</li> <li>The peak may not be labeled if the peak is too wide.</li> </ul>



**f** Click the **Close Data File** icon in the main toolbar, and click **Close** when the dialog box containing the list of data files pops up.

- The product ions appear to be:  
Sulfamethizole-271.0 > 155.7  
Sulfamethazine-279.0 > 185.8  
Sulfachloropyridazine-285.0 > 155.7  
Sulfadimethoxine-311.0 > 155.7

## Exercise 1 – Develop an acquisition method

### Task 5. Find optimum collision energy for MRM acquisition

## Task 5. Find optimum collision energy for MRM acquisition


In this task, you set up MRM acquisition methods for the sulfa drugs for different collision energies. By examining the spectra and comparing peak intensities, you determine the optimal collision energy settings for the compounds.

Steps	Detailed Instructions	Comments
<p><b>1</b> Set up three MRM acquisition methods.</p> <ul style="list-style-type: none"><li>Use all the MS parameters in the example below except for the collision energy value.</li><li>Use collision energies of 10, 15 and 20.</li><li>Save methods as <b>iiiSulfamix MRM_xx.m</b>, where <b>iii</b> are your initials and <b>xx</b> is the collision energy.</li></ul>	<p><b>a</b> Click the <b>Q00</b> tab.</p> <p><b>b</b> Set <b>Scan Type</b> to <b>MRM</b>.</p> <p><b>c</b> Enter all MS parameters shown in the example below except for the collision energy value.</p> <p><b>d</b> In the collision energy column, type 10 for each compound.</p> <p><b>e</b> Save the method as <b>iiiSulfamix MRM_10.m</b>.</p> <p><b>f</b> Repeat <b>step d</b> and <b>step e</b> for collision energies of 15, 20, 25, 30 and 35 saving the methods as <b>iiiSulfamix MRM_xx.m</b>, where <b>iii</b> are your initials and <b>xx</b> is the collision energy.</p>	<ul style="list-style-type: none"><li>Because the largest peaks were produced with a collision energy of 15 in the previous exercise, you will look at only those collision energies to either side of 15.</li></ul>

The screenshot displays the software interface for setting up an MRM acquisition method. On the left, the 'Time segments' table is visible, showing a single segment for MRM. On the right, the 'Acquisition' table lists four compounds with their respective parameters.

Compound Name	ISTD?	Precursor Ion	MS1 Res	Product Ion	MS2 Res	Dwell	Fragmentor	Collision Energy	Cell Accelerator Voltage	Polarity
Sulfadimethoxine	<input type="checkbox"/>	311	Unit	155.7	Unit	50	100	10	7	Positive
Sulfachloropyridazine	<input type="checkbox"/>	285	Unit	155.7	Unit	50	100	10	7	Positive
Sulfamethazine	<input type="checkbox"/>	279	Unit	185.7	Unit	50	100	10	7	Positive
<input checked="" type="checkbox"/> Sulfamethizole	<input type="checkbox"/>	271	Unit	155.8	Unit	50	100	10	7	Positive

<p><b>2</b> Set up and run the worklist (optional).</p> <ul style="list-style-type: none"><li>Specify the data files as <b>iiiSulfamix MRM_xx.d</b>, where <b>iii</b> are your initials and <b>xx</b> is the collision energy.</li></ul>	<p><b>a</b> Click the <b>Worklist</b> tab to make the worklist visible.</p> <p><b>b</b> Add six samples to the worklist for collision energies 10, 15, 20, 25, 30, 35.</p> <p><b>c</b> Mark the checkbox to the left of the Sample Name for each of the three samples.</p> <p><b>d</b> Click <b>Worklist &gt; Run</b>.</p>	<ul style="list-style-type: none"><li>This step is optional because you can use the six example data files in the next step.</li></ul>
--	--	--

Steps	Detailed Instructions	Comments
<p><b>3</b> Compare the compound transition intensities at different collision energies.</p> <ul style="list-style-type: none"> <li>Open the MRM data files: SulfamixMRM_10.d SulfamixMRM_15.d SulfamixMRM_20.d SulfamixMRM_25.d SulfamixMRM_30.d SulfamixMRM_35.d</li> <li>Set the MRM chromatogram extraction parameters as shown at right for all transitions.</li> <li>Disable the TICs for clarity and examine the peak intensities.</li> <li>Compare the intensities of each compound transition obtained at one collision energy with the same compound transition obtained at another collision energy. (Do this in Overlaid Mode with all the MRM chromatograms.)</li> <li>Close the data files but don't save results.</li> <li>Refer to <a href="#">Table 4</a> on page 32 for optimal method settings for each compound.</li> </ul>	<p><b>a</b> Open the <b>Qualitative Analysis</b> program.</p> <p><b>b</b> Clear the <b>Run 'File Open' actions...</b> check box.</p> <p><b>c</b> Open the MRM data files in the Qualitative Analysis program.</p> <p><b>d</b> Right-click the Chromatogram Results window, and click <b>Extract Chromatograms</b> from the shortcut menu.</p> <p><b>e</b> To select all data files, click the last file while holding down the <b>Shift</b> key.</p> <p><b>f</b> Enter the parameters as listed in the example below, and click <b>OK</b>.</p> <p><b>g</b> Clear the TIC check boxes to make the MRM chromatograms easier to view.</p>	<ul style="list-style-type: none"> <li>Why a spectrum for MRM? It's a feature of the program to show spectra even for MRM experiments and can be quite handy for comparing relative intensities of product ions generated from the same precursor.</li> </ul>
		
<p><b>h</b> Click the <b>Overlaid Mode</b> icon, .</p> <p><b>i</b> Compare peak intensities for each compound transition in each data file in the Chromatogram Results window.</p>		<ul style="list-style-type: none"> <li>Compare the colors shown in Chromatogram Results with the color next to the MRM transition name in the Data Navigator.</li> <li>You can also right-click the Chromatogram Results window header and compare the colors of the chromatograms to the colors of the titles in the shortcut menu.</li> </ul>


## Exercise 1 – Develop an acquisition method

### Task 5. Find optimum collision energy for MRM acquisition

#### Steps

#### Detailed Instructions

#### Comments



Unless you decide to acquire MRMs at lower collision energies, you should find that the optimal method settings are as shown in [Table 4](#).

- j Click the **Close Data File** icon in the main toolbar, and click **Close** when the Close Data File dialog box appears.

- You now have all the information you need to do an MRM acquisition experiment of the sulfa drug mixture. Consider doing at least one more run with those settings.

**Table 4** Compounds and Collision Energy

Compounds	MRM Transition	Collision Energy (V)
Sulfamethizole	271.0 > 155.8	10
Sulfamethazine	279.0 > 185.7	15
Sulfachloropyradizine	285.0 > 155.7	10
Sulfadimethoxine	311.0 > 155.7	15

## Exercise 2 – Develop a Dynamic MRM acquisition method from an MRM acquisition data file or an MRM method

The purpose of this exercise is to create a Dynamic MRM method from an acquired MRM data file for sulfamix\_MRM data files with the correct retention times for Dynamic MRM using the Quantitative Analysis program.

For this exercise, you have three main tasks:

- “Task 1. Create a batch file from an existing MRM data file” on page 33
- “Task 2. Print a report in the Quantitative Analysis program” on page 36
- “Task 3. Create a Dynamic MRM method using Update dMRM” on page 37

You can easily create a Dynamic MRM method from an existing MRM method.

- “Task 4. Create a Dynamic MRM method from an MRM method” on page 39

### Task 1. Create a batch file from an existing MRM data file

In this exercise, you create a batch and a method from an existing MRM data file.

Steps	Detailed Instructions	Comments
<b>1</b> Open the Quantitative Analysis program and create a batch file with one sample file, SulfamixMRM_35.d. <ul style="list-style-type: none"><li>• Copy the data file SulfamixMRM_35.d from the installation disk to the \MassHunter\Data\MRM_to_DMRM folder.</li></ul>	<ul style="list-style-type: none"><li><b>a</b> Double-click the <b>QQQ Quantitative Analysis</b> icon.</li><li><b>b</b> Click <b>File &gt; New Batch</b>.</li><li><b>c</b> Navigate to the \MassHunter\Data\MRM_to_DMRM folder.</li><li><b>d</b> Type MRM_to_DMRM in the <b>File Name</b> text box.</li><li><b>e</b> Click <b>Open</b>.</li><li><b>f</b> Click <b>File &gt; Add Samples</b>.</li><li><b>g</b> Select the file SulfamixMRM_35.d.</li><li><b>h</b> Click <b>OK</b>.</li></ul>	<ul style="list-style-type: none"><li>• The file <b>SulfamixMRM_35.d</b> is on the installation disk in the \Support\Data folder. Copy this entire folder to the \MassHunter\Data\MRM_to_DMRM folder.</li></ul>

## Exercise 2 – Develop a Dynamic MRM acquisition method from an MRM acquisition data file or an MRM method


### Task 1. Create a batch file from an existing MRM data file

Steps	Detailed Instructions	Comments
2 Create a method for that batch using MRM data.	<p><b>a</b> Click <b>Method &gt; New &gt; New Method from Acquired MRM data</b>.</p> <p><b>b</b> Select the <b>SulfamixMRM_35.d</b> data file.</p> <p><b>c</b> Click <b>Open</b>.</p>	
<p>3 Set the Concentration Setup, Qualifier Setup, and Calibration Curve Setup.</p> <ul style="list-style-type: none"> <li>• Add calibration level 1 with a concentration of 10000.</li> <li>• Set the <b>Uncertainty</b> to Relative for all qualifiers.</li> <li>• Set the <b>Curve Fit</b> to Linear.</li> <li>• Set the <b>Curve Fit Origin</b> to Include.</li> <li>• Set the <b>Curve Fit Weight</b> to None.</li> </ul>	<p><b>a</b> Select <b>Concentration Setup</b> in the Manual Setup Tasks section in the Method Tasks pane.</p> <p><b>b</b> Select the first compound in the table.</p> <p><b>c</b> Right-click the compound row and click <b>New Calibration Level</b> from the shortcut menu.</p> <p><b>d</b> Enter 1 in the <b>Level</b> column and 10 in the <b>Conc.</b> column.</p> <p><b>e</b> Right-click in the Level box and click <b>Copy Calibration Levels To</b>.</p> <p><b>f</b> Click <b>Select All</b>. Click <b>OK</b>.</p> <p><b>g</b> Select <b>Qualifier Setup</b> in the Manual Setup Tasks section in the Method Tasks pane.</p> <p><b>h</b> Verify that the <b>Uncertainty</b> is Relative.</p> <p><b>i</b> Select <b>Calibration Curve Setup</b> in the Manual Setup Tasks section in the Method Tasks pane.</p> <p><b>j</b> Set <b>Curve Fit</b> to <b>Linear</b> for all compounds.</p> <p><b>k</b> Set <b>CF Origin</b> to <b>Include</b> for all compounds.</p> <p><b>l</b> Set <b>CF Weight</b> to <b>None</b> for all compounds.</p>	<ul style="list-style-type: none"> <li>• Refer to the online Help in the Quantitative Analysis program for additional help on these tasks.</li> </ul>

## Exercise 2 – Develop a Dynamic MRM acquisition method from an MRM acquisition data file or an MRM method

### Task 1. Create a batch file from an existing MRM data file

Steps	Detailed Instructions	Comments
-------	-----------------------	----------



- 4 Verify method and then save the method and apply the method to the batch.
  - a Click **Method > Validate**.
  - b Click OK on the message box. Fix any errors, if necessary.
  - c Click **Method > Save As**.
  - d Enter MRM\_to\_DMRM.
  - e Click the **Save** button.
  - f Click **Method > Exit**.
  - g Click **Yes** to apply the method to the batch.
  
- 5 Analyze and save the batch.
  - a Click **Analyze > Analyze Batch**.
  - b Click **File > Save Batch**.

## Exercise 2 – Develop a Dynamic MRM acquisition method from an MRM acquisition data file or an MRM method

### Task 2. Print a report in the Quantitative Analysis program

## Task 2. Print a report in the Quantitative Analysis program

In this task, you print a report using any template.


You can update a Dynamic MRM method using either a data file or a quantitation report folder, so this task creates the quantitation report folder.

Steps	Detailed Instructions	Comments
1 Print a report using the template MRM_to_DMRM.xltx.	<p><b>a</b> Click <b>File &gt; Save</b>.</p> <p><b>b</b> Click <b>Report &gt; Generate</b>. The system displays the Report dialog box.</p> <p><b>c</b> Select the <b>Template</b> file.</p> <p><b>d</b> Select the <b>Report</b> folder. This folder name will be used in the next task.</p> <p><b>e</b> Click <b>OK</b>.</p>	<ul style="list-style-type: none"><li>• Copy the <b>MRM_to_DMRM.xltx</b> template from the <b>\Support\Data</b> folder on the installation disk.</li><li>• For this report, you do not need to print the report. You need to click <b>Advanced</b> to select a different printer. If you don't want to print this report, click <b>Advanced</b> instead.</li></ul>
2 Check the status of the report using the Queue Viewer program.	<p><b>a</b> Click <b>Report &gt; Queue Viewer</b>.</p> <p><b>b</b> Wait for the report to finish printing.</p> <p><b>c</b> Close the <b>Task Queue Viewer</b> program.</p>	

## Task 3. Create a Dynamic MRM method using Update dMRM

You can create a Dynamic MRM method from an MRM data file or a Quantitative Analysis method. You first set the **Scan Type** to Dynamic MRM, and then you use the Update MRM Method dialog box.

Steps	Detailed Instructions	Comments
1	<p>Open the method <i>iiiSulfamix MRM_10.m</i> and save it to a new name with the format <i>iiiSulfamix dMRM.m</i>, where <i>iii</i> are your initials.</p>	<ul style="list-style-type: none"> <li>In this example, the batch is in the <code>\MassHunter\Data\MRM_to_dMRM</code> folder.</li> </ul>
2	<p>Change the method to a dynamic MRM method with the same compounds. You can either use a data file or the report that was generated in the last task.</p>	<ul style="list-style-type: none"> <li>The Update MRM Method tool automatically sets the Scan type to <b>Dynamic MRM</b>.</li> <li>You can select either a data file that was acquired with a <b>Scan Type</b> of <b>MRM</b> or a Quant Report folder as the input to this dialog box. The Scan segments are created from one of these two input sources.</li> </ul>



You can update the compounds in the Scan segments table by using a QQQ data file or a Quantitative analysis report folder.

## Exercise 2 – Develop a Dynamic MRM acquisition method from an MRM acquisition data file or an MRM method

### Task 3. Create a Dynamic MRM method using Update dMRM

#### Steps

#### Detailed Instructions

#### Comments

The screenshot shows the 'Acquisition' method configuration window. The 'Scan segments' table is populated with the following data:

Compound Group	Compound Name	ISTD?	Precursor Ion	M51 Res	Product Ion	M52 Res	Ret Time (min)	Delta Ret Time	Fragmentor	Collision Energy	Cell Accelerator Voltage	Polarity
	Compound_1	<input type="checkbox"/>	311	Unit	155.7	Unit	2.15	0.81	100	0	0	7 Positive
	Compound_2	<input type="checkbox"/>	285	Unit	155.7	Unit	0.74	0.95	100	0	0	7 Positive
	Compound_3	<input type="checkbox"/>	273	Unit	185.7	Unit	1.16	1.17	100	0	0	7 Positive
	Compound_4	<input type="checkbox"/>	271	Unit	195.8	Unit	0.41	0.37	100	0	0	7 Positive

Dynamic MRM Parameters: Cycle Time [500] ms, Triggered MRM  Triggered Repeats [3]

The compounds from the data file or quantitation report are automatically added to the Scan segments table.

- g Select the original compound in the Scan segments table.
- h Right-click the row and click **Delete Row**.
- i Verify that each row has a **Compound Name**. A blank Compound Name is not allowed.
- j Click **Method > Save**.

## Task 4. Create a Dynamic MRM method from an MRM method

You can create a Dynamic MRM method directly from an MRM method by using the Paste from Clipboard command from the shortcut menu.

Steps	Detailed Instructions	Comments
1	<p>Open the method <i>iii</i>Sulfamix MRM_10.m and save it to a new name with the format <i>iii</i>Sulfamix dMRM_Easy.m, where <i>iii</i> are your initials.</p>	
	<p>a Click <b>File &gt; Open &gt; Method.</b>                      b Select the <i>iii</i>Sulfamix MRM_10.m method.                      c Click <b>OK.</b>                      d Click <b>Method &gt; Save As.</b>                      e Type the new method name with the format <i>iii</i>Sulfamix dMRM2.m.                      f Click the <b>Save</b> button.</p>	
2	<p>Copy all compounds from the Scan segments table in the MRM method.</p>	
	<p>a Click the <b>Acquisition</b> tab in the QQQ tab in the Method Editor.                      b Select all of the rows in the Scan segments table.                      c Right-click the Scan segments table and click <b>Copy.</b></p>	<ul style="list-style-type: none"> <li>To select all of the rows in the Scan segments table, you select the first row in the table. Then, you scroll to the last row in the Scan segments table. Press the <b>Shift</b> key and select the last row in the table.</li> </ul>
3	<p>Change the Scan Type to Dynamic MRM and paste the rows into the new Scan segments table.</p>	
	<p>a Select <b>Dynamic MRM</b> as the Scan Type.                      b Right-click the Scan segments table and click <b>Paste from Clipboard.</b>                      c Select the original compound in the Scan segments table.                      d Right-click and click <b>Delete Row.</b>                      e Click <b>Method &gt; Save.</b></p>	<ul style="list-style-type: none"> <li>To combine multiple Time Segments into one Dynamic MRM Time Segment, you paste the Scan segments into Excel and create one long list. When all of the scan segments have been pasted into Excel, then copy all of the Scan segments in Excel.</li> </ul>

The screenshot shows the 'Acquisition' tab in the software. The 'Scan segments' table is visible with the following data:

Compound Name	ISTD?	Precursor Ion	MS1 Res	Product Ion	MS2 Res	Fragmentor	Collision Energy	Cell Accelerator Voltage	Ret Time (min)	Delta Ret Time	Polarity
sulfamethoxine	<input type="checkbox"/>	350	Unit	200	Unit				0	0	Positive
sulfamethoxazole	<input type="checkbox"/>	311	Unit	155.7	Unit				5	0.81	Positive
sulfachloropyridazin	<input type="checkbox"/>	285	Unit	155.7	Unit				4	0.55	Positive
sulfamethazine	<input type="checkbox"/>	279	Unit	185.7	Unit				6	1.17	Positive
sulfamethizole	<input type="checkbox"/>	271	Unit	155.8	Unit				1	0.37	Positive

A context menu is open over the table, showing options: Insert Row, Append Row, Delete Row, Sort, Import from optimizer..., Update MRM Method..., Edit MRM Method..., Calibrate MRM Method..., Cut, Copy, Paste, and Paste from Clipboard. The 'Copy' option is highlighted.

Below the table, the 'Dynamic MRM Parameters' section shows 'Cycle Time' set to 500 ms.

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 1. Create a Triggered Dynamic MRM method from a Dynamic MRM method manually

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

For this exercise you analyze a mixture of four sulfonamide compounds.

### Task 1. Create a Triggered Dynamic MRM method from a Dynamic MRM method manually

You can create a Triggered Dynamic MRM method directly from a Dynamic MRM method. In a Triggered Dynamic MRM method, you specify some of the transitions to be primary transitions. These transitions are acquired for the entire retention time window. Some of these primary transitions are also marked as triggers. As the data is acquired, the software checks whether or not the abundances of the trigger transitions are higher than the threshold. If the abundances are higher than the thresholds and other additional conditions are met, then the secondary transitions are acquired. These other conditions are described in the *Concepts* guide.

Steps	Detailed Instructions	Comments
1 Open the method <i>iiiSulfamix</i> MRM_10.m and save it to a new name with the format <b><i>iiiSulfamix</i> dMRM_Easy.m</b> , where <i>iii</i> are your initials.	<p><b>a</b> Click <b>File &gt; Open &gt; Method</b>.</p> <p><b>b</b> Select the <i>iiiSulfamix_dMRM2.m</i> method.</p> <p><b>c</b> Click <b>OK</b>.</p> <p><b>d</b> Click <b>Method &gt; Save As</b>.</p> <p><b>e</b> Type the new method name with the format <b><i>iiiSulfamix_TriggeredDMRM.m</i></b>.</p> <p><b>f</b> Click the <b>Save</b> button.</p>	<ul style="list-style-type: none"><li>• A Triggered Dynamic MRM method is a type of Dynamic MRM method. The <b>Scan Type</b> for both methods is <b>Dynamic MRM</b>.</li><li>• The Dynamic MRM method is the template method for the optimization.</li></ul>

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 1. Create a Triggered Dynamic MRM method from a Dynamic MRM method manually

Steps	Detailed Instructions	Comments
2	<p>Change the method to a triggered dynamic MRM method.</p> <ol style="list-style-type: none"> <li><b>a</b> Click the <b>Acquisition</b> tab in the QQQ tab in the Method Editor.</li> <li><b>b</b> Mark the <b>Triggered</b> check box in the Triggered MRM section. This section is only available if the <b>Scan Type</b> is <b>Dynamic MRM</b>.</li> <li><b>c</b> Select whether to automatically mark the highest product ion as the Primary.</li> <li><b>d</b> Enter the value for <b>Repeats</b>.</li> </ol>	<ul style="list-style-type: none"> <li>• Several columns are added to the <b>Scan segments</b> table. These columns only apply to a triggered dynamic MRM method.</li> <li>• The value <b>Repeats</b> is the number of times to acquire each of the secondary transitions when the triggering conditions are met.</li> </ul>
3	<p>Select the transitions that are the <b>Primary</b> transitions.</p> <ol style="list-style-type: none"> <li><b>a</b> For each transition, mark the <b>Primary</b> check box if it is a <b>Primary</b> transition.</li> <li><b>b</b> Verify that you have marked at least one transition as the <b>Primary</b> transition for each <b>Compound Name</b>.</li> </ol>	<ul style="list-style-type: none"> <li>• You can select multiple transitions from each compound to be <b>Primary</b> transitions. If a transition has the same <b>Compound Name</b>, then it is part of the same compound. You must mark at least one transition as a <b>Primary</b> transition for each compound.</li> </ul>
4	<p>Select the transitions that are the <b>Trigger</b> transitions and set the trigger conditions.</p> <ol style="list-style-type: none"> <li><b>a</b> For each compound, mark the <b>Trigger</b> check box if it is a Trigger transition.</li> <li><b>b</b> (optional) Mark a second <b>Trigger</b> transition.</li> <li><b>c</b> (optional) Enter the <b>Threshold</b> value for each <b>Trigger</b> transition.</li> <li><b>d</b> (optional) Enter the <b>Trigger Entrance</b> for each <b>Trigger</b> transition.</li> <li><b>e</b> (optional) Enter the <b>Trigger Delay</b> for each <b>Trigger</b> transition.</li> <li><b>f</b> (optional) Enter the <b>Trigger Window</b> for each <b>Trigger</b> transition.</li> </ol>	<ul style="list-style-type: none"> <li>• For each compound, you can have two Trigger transitions.</li> <li>• If the <b>Trigger</b> transition has an abundance over the <b>Threshold</b>, then that triggering condition is met.</li> <li>• By default, the <b>Trigger Entrance</b>, the <b>Trigger Delay</b> and the <b>Trigger Window</b> are set to 0. If these values are 0, then these triggering conditions are not enabled.</li> <li>• See the <i>Concepts Guide</i> for more information on these trigger conditions.</li> </ul>

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments																	
Scan segments																			
Compound Group	Compound Name	Precursor	Product Ion	Primary	Trigger	Threshold	Ret Time (min)	Delta Ret Time	Fragmentor	Collision Energy	Trigger Entrance	Trigger Delay	Trigger Window	Cell Accelerator	Polarity	MS1 Res	MS2 Res	ISTD?	
1	sulfachloropyridazine	265	197	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	800	0.5	1	380	8	2	0	0	7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfachloropyridazine	265	156	<input checked="" type="checkbox"/>	<input type="checkbox"/>		0.5	1	380	8				7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfachloropyridazine	265	108	<input type="checkbox"/>	<input type="checkbox"/>				380	20				7	Positive	Unit	Unit	<input type="checkbox"/>	
1	sulfadimethoxine	311.1	245.1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1000	1.2	1	380	12	0	1	0	7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfadimethoxine	311.1	173.1	<input checked="" type="checkbox"/>	<input type="checkbox"/>		1.2	1	380	24				7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfadimethoxine	311.1	156	<input type="checkbox"/>	<input type="checkbox"/>				380	16				7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfadimethoxine	311.1	108	<input type="checkbox"/>	<input type="checkbox"/>				380	24				7	Positive	Unit	Unit	<input type="checkbox"/>	
2	sulfamethazine	279.1	186	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		900	0.8	1	380	12	0	0	0.4	7	Positive	Unit	Unit	<input type="checkbox"/>
	sulfamethazine	279.1	155.9	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1000	0.8	1	380	12	0	0	0.4	7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfamethazine	279.1	124.1	<input type="checkbox"/>	<input type="checkbox"/>				380	20				7	Positive	Unit	Unit	<input type="checkbox"/>	
2	sulfamethazine	279.1	108	<input type="checkbox"/>	<input type="checkbox"/>				380	24				7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfamethazole	271	253.4	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1100	0.3	1	380	0	2	1	0.3	7	Positive	Unit	Unit	<input type="checkbox"/>	
	sulfamethazole	271	156	<input checked="" type="checkbox"/>	<input type="checkbox"/>		0.3	1	380	4				7	Positive	Unit	Unit	<input type="checkbox"/>	
1	sulfamethazole	271	108	<input type="checkbox"/>	<input type="checkbox"/>				380	20				7	Positive	Unit	Unit	<input type="checkbox"/>	
Dynamic MRM Parameters																			
Cycle Time	500 ms	Total MRMs = 14	Max Concurrent MRMs = 14	Min/Max Dwell = 32.21 ms/32.21 ms.															
		Primary Only - Total MRMs = 9	Max Concurrent MRMs = 9	Min/Max Dwell = 52.06															
Triggered MRM																			
<input checked="" type="checkbox"/> Triggered																			
Repeats <input type="text" value="3"/>																			

## Task 2. Add/Modify compounds in an existing database

You can also manually add compounds to a database and modify the compounds in the database. In the next task, you create a Triggered Dynamic MRM method from the compounds in the database.

Steps	Detailed Instructions	Comments
1	Review the <b>iiiSulfamix_dMRM2.m</b> , where <b>iii</b> are your initials.	<ul style="list-style-type: none"> <li>A Triggered Dynamic MRM method is a type of Dynamic MRM method. The <b>Scan Type</b> for both methods is <b>Dynamic MRM</b>.</li> </ul>
2	Start the MassHunter Optimizer software.	<ul style="list-style-type: none"> <li>Double-click the <b>Optimizer</b> icon.  .</li> <li>If you are optimizing peptides, use the <b>Optimizer for Peptides</b> program.</li> </ul>

### Exercise 3 – Create a Triggered Dynamic MRM acquisition method


#### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
3 Set parameters on the Optimizer Setup tab.	<ol style="list-style-type: none"><li>Click the <b>Optimizer Setup</b> tab.</li><li>Click the <b>Injection (with or without column)</b> button.</li><li>Set the <b>CE</b> range from 4 to 48.</li><li>Set the <b>Cell Accelerator Voltage</b> to 7.</li><li>Right-click the table and click <b>Add Method</b>.</li><li>Select the <b>iiiSulfamix_dMRM2.m</b> method.</li></ol>	<ul style="list-style-type: none"><li>To create low mass product ions from a precursor ion near 300 <math>m/z</math>, you need fairly high collision energies.</li></ul>
4 Set parameters on the Precursor Ion Selection tab.	<ol style="list-style-type: none"><li>Click the <b>Precursor Ion Selection</b> tab.</li><li>Verify that +H is marked for the <b>Positive ions (with priorities)</b> list.</li></ol>	
5 Set parameters on the Product Ion Selection tab.	<ol style="list-style-type: none"><li>Click the <b>Product Ion Selection</b> tab.</li><li>Click the <b>Mass (m/z)</b> button under Low mass cut-off.</li><li>Enter 60 for the low mass cut-off.</li></ol>	<ul style="list-style-type: none"><li>On the Product Ion Selection, you can automatically add up to 4 product ions per compound (for instance, 2 primaries and 2 secondaries). You want 8 to 10 peaks in the composite spectrum to prove that this is indicative of the compound, so you need to add at least some of the product ions manually.</li></ul>

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
6	<p>Set parameters on the Compound Setup tab and add additional transitions.</p> <ul style="list-style-type: none"> <li>For Precursor ion 311, add the following product ions: 244.8, 229.7, 217.7, 172.9, 107.9, 91.9, 79.9, 64.8</li> <li>For Precursor 285, add the following product ions: 129.9, 107.9, 91.9, 79.8, 64.8</li> <li>For Precursor 279, add the following product ions: 212.8, 155.9, 123.9, 107.9, 91.9, 79.8, 64.9</li> <li>For Precursor 271, add the following product ions: 177.8, 115.9, 107.9, 92, 80, 64.9</li> </ul>	<ul style="list-style-type: none"> <li>For each compound, we are going to add additional transitions.</li> <li>In the Qualitative Analysis program, you examine Product Ion data files which you acquired previously to determine additional transitions to add. See <a href="#">“Task 4. Determine product ion masses”</a> on page 24.</li> <li>You can use the arrow keys to move between rows in the Product table.</li> </ul>
	<p><b>a</b> Click the <b>Compound Setup</b> tab.</p> <p><b>b</b> Click the <b>Import/Export &gt; Import from Acquisition Methods</b> command.</p> <p><b>c</b> Select the <b>iiiSulfamix_dMRM2.m</b> method and click <b>Open</b>.</p> <p><b>d</b> (optional) Right-click the tab and click <b>Expand/Collapse All Rows</b>.</p> <p><b>e</b> Select one of the <b>Product</b> rows for one of the compounds. In this example, select the <b>Product</b> row 155.7 for <b>Precursor</b> 311.</p> <p><b>f</b> Right-click the Product row and click <b>Add Product Ion</b>. In this example, you add 8 product ion rows.</p> <p><b>g</b> Enter the Product in each of the product ion rows that were added. See <a href="#">“To determine product ions in the Qualitative Analysis program:”</a> on page 45.</p> <p><b>h</b> Add product ions for compounds 1, 2 and 3.</p>	



This product ion scan has a precursor mass of 311. You examine the MS spectrum to determine the product ions to add to the Product ion section of the Compound Setup table.

The product ions that are manually added as additional Product ions in the MassHunter Optimizer software are shown in the MS Spectrum Results window. The green boxes were added in this guide to show which product ions were used.

### Exercise 3 – Create a Triggered Dynamic MRM acquisition method


#### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
3 Set parameters on the Optimizer Setup tab.	<ol style="list-style-type: none"><li>Click the <b>Optimizer Setup</b> tab.</li><li>Click the <b>Injection (with or without column)</b> button.</li><li>Set the <b>CE</b> range from 4 to 48.</li><li>Set the <b>Cell Accelerator Voltage</b> to 7.</li><li>Right-click the table and click <b>Add Method</b>.</li><li>Select the <b>iiiSulfamix_dMRM2.m</b> method.</li></ol>	<ul style="list-style-type: none"><li>To create low mass product ions from a precursor ion near 300 <math>m/z</math>, you need fairly high collision energies.</li></ul>
4 Set parameters on the Precursor Ion Selection tab.	<ol style="list-style-type: none"><li>Click the <b>Precursor Ion Selection</b> tab.</li><li>Verify that +H is marked for the <b>Positive ions (with priorities)</b> list.</li></ol>	
5 Set parameters on the Product Ion Selection tab.	<ol style="list-style-type: none"><li>Click the <b>Product Ion Selection</b> tab.</li><li>Click the <b>Mass (m/z)</b> button under Low mass cut-off.</li><li>Enter 60 for the low mass cut-off.</li></ol>	<ul style="list-style-type: none"><li>On the Product Ion Selection, you can automatically add up to 4 product ions per compound (for instance, 2 primaries and 2 secondaries). You want 8 to 10 peaks in the composite spectrum to prove that this is indicative of the compound, so you need to add at least some of the product ions manually.</li></ul>

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
6	<p>Set parameters on the Compound Setup tab and add additional transitions.</p> <ul style="list-style-type: none"> <li>For Precursor ion 311, add the following product ions: 244.8, 229.7, 217.7, 172.9, 107.9, 91.9, 79.9, 64.8</li> <li>For Precursor 285, add the following product ions: 129.9, 107.9, 91.9, 79.8, 64.8</li> <li>For Precursor 279, add the following product ions: 212.8, 155.9, 123.9, 107.9, 91.9, 79.8, 64.9</li> <li>For Precursor 271, add the following product ions: 177.8, 115.9, 107.9, 92, 80, 64.9</li> </ul>	<ul style="list-style-type: none"> <li>For each compound, we are going to add additional transitions.</li> <li>In the Qualitative Analysis program, you examine Product Ion data files which you acquired previously to determine additional transitions to add. See <a href="#">“Task 4. Determine product ion masses”</a> on page 24.</li> <li>You can use the arrow keys to move between rows in the Product table.</li> </ul>
	<p><b>a</b> Click the <b>Compound Setup</b> tab.</p> <p><b>b</b> Click the <b>Import/Export &gt; Import from Acquisition Methods</b> command.</p> <p><b>c</b> Select the <b>iiiSulfamix_dMRM2.m</b> method and click <b>Open</b>.</p> <p><b>d</b> (optional) Right-click the tab and click <b>Expand/Collapse All Rows</b>.</p> <p><b>e</b> Select one of the <b>Product</b> rows for one of the compounds. In this example, select the <b>Product</b> row 155.7 for <b>Precursor</b> 311.</p> <p><b>f</b> Right-click the Product row and click <b>Add Product Ion</b>. In this example, you add 8 product ion rows.</p> <p><b>g</b> Enter the Product in each of the product ion rows that were added. See <a href="#">“To determine product ions in the Qualitative Analysis program:”</a> on page 45.</p> <p><b>h</b> Add product ions for compounds 1, 2 and 3.</p>	




This product ion scan has a precursor mass of 311. You examine the MS spectrum to determine the product ions to add to the Product ion section of the Compound Setup table.

The product ions that are manually added as additional Product ions in the MassHunter Optimizer software are shown in the MS Spectrum Results window. The green boxes were added in this guide to show which product ions were used.

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
<ul style="list-style-type: none"><li>To determine product ions in the Qualitative Analysis program:</li></ul>	<ol style="list-style-type: none"><li>Open the <b>SulfamixPI_15.d</b> from “<b>Task 4. Determine product ion masses</b>” on page 24.</li><li>Click <b>Find &gt; Find Compounds by Targeted MS/MS</b>.</li><li>Close the Compound List window.</li><li>Select a compound in the Data Navigator window. For this example, click <b>Cpd 4</b>.</li><li>Click the Autoscale Y-axis in the MS Spectrum Results toolbar.</li><li>Right-click and drag to zoom in on the MS spectrum.</li></ol>	<ul style="list-style-type: none"><li>If possible, rearrange the windows on the screen so you can see the Optimizer program and the Qualitative Analysis program at the same time.</li></ul>



This product ion scan has a precursor mass of 311. You examine the MS spectrum to determine the product ions to add to the Product ion section of the Compound Setup table.

<p><b>7</b> Set other parameters in the Compound Setup tab and start the optimization.</p> <ul style="list-style-type: none"><li>You cannot perform a multi-compound run.</li><li>You have to mark each row in the table to use.</li></ul>	<ol style="list-style-type: none"><li>Mark the check box in the left column at the top of the table. The check box for every row in the table is marked.</li><li>Clear the <b>Perform multi-compound run</b> check box in the right column.</li><li>Click the <b>Start Optimization</b> button in the Optimizer toolbar.</li></ol>	<ul style="list-style-type: none"><li>You cannot perform a multi-compound run with the number of transitions that were added. If you mark this check box, then the Expected peak width (base) is automatically set to almost 80 seconds wide. If you clear this check box, then the Expected peak width is calculated to be around 9 seconds which is more appropriate.</li></ul>
--	--	---



## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
11 In the Database Browser program, select the transitions.	<p><b>a</b> Mark the <b>Show All Records</b> check box.</p> <p><b>b</b> Click the <b>Select top</b> button under Select Transitions.</p> <p><b>c</b> Type 10 for the ranked transitions.</p> <p><b>d</b> Click the <b>Select Transitions</b> button.</p>	<ul style="list-style-type: none"> <li>All the transitions that you typed in are visible.</li> <li>The tools to allow you to set up <b>Primary transitions</b> and <b>Secondary transitions</b> are available in this program.</li> </ul>

The screenshot shows the 'Select Transitions' dialog box in the Database Browser program. The 'Select top' button is highlighted, and the value '10' is entered in the 'ranked transitions' field. The 'Set primary and trigger flags' section shows 'Set top' set to '2' and 'ranked transitions as primary' checked. The 'Rank transitions by' section shows 'Abundance' selected. Below the dialog is a table of transition data.

Compound Name	Formula	MW	Polarity	Species	Precursor	Product	Frag	CE	CAV	Primary	Trigger	RT	RT Window	Abundance
sulfachloropyridazin			Positive		285	155.7	380	10	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	5723674
sulfachloropyridazin			Positive		285	129.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	315469
sulfachloropyridazin			Positive		285	107.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	5015922
sulfachloropyridazin			Positive		285	91.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	7438199
sulfachloropyridazin			Positive		285	79.8	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	1529348
sulfachloropyridazin			Positive		285	64.8	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	6508111
sulfadimethoxine			Positive		311	155.7	380	15	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	3335529
sulfadimethoxine			Positive		311	244.8	380	12	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	548481
sulfadimethoxine			Positive		311	229.7	380	20	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	135578
sulfadimethoxine			Positive		311	217.7	380	16	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	434965
sulfadimethoxine			Positive		311	172.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	414589
sulfadimethoxine			Positive		311	107.9	380	20	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	2291401
sulfadimethoxine			Positive		311	91.9	380	32	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	3102685
sulfadimethoxine			Positive		311	79.9	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	761344
sulfadimethoxine			Positive		311	64.8	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	2583686
sulfamethazine			Positive		279	185.7	380	11	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	3720936
sulfamethazine			Positive		279	212.8	380	20	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	421592
sulfamethazine			Positive		279	155.9	380	12	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	1377529
sulfamethazine			Positive		279	123.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	2633848
sulfamethazine			Positive		279	107.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	2407737
sulfamethazine			Positive		279	91.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	3643700
sulfamethazine			Positive		279	79.8	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	1242301
sulfamethazine			Positive		279	64.9	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	3793552
sulfamethazole			Positive		271	155.8	380	6	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.37	1	4578854
sulfamethazole			Positive		271	177.8	380	12	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.37	1	89336
sulfamethazole			Positive		271	115.9	380	16	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.37	1	515222
sulfamethazole			Positive		271	107.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.37	1	3108252
sulfamethazole			Positive		271	92	380	28	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.37	1	5192801

12 In the Database Browser program, automatically select the <b>Primary</b> transitions and <b>Trigger</b> transition.	<p><b>a</b> In the <b>Set top ranked transitions as primary</b> box, enter 2.</p> <p><b>b</b> Click the <b>Set Primaries and Trigger</b> button.</p>	<ul style="list-style-type: none"> <li>The software automatically selects the two most abundant transitions as the <b>Primary</b> transitions.</li> <li>The software also selects the most abundant transition as the <b>Trigger</b>.</li> <li>You can manually select a second <b>Trigger</b> transition.</li> </ul>
--	--	---

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

#### Steps

#### Detailed Instructions

#### Comments

Select Transitions:  Select top 10 ranked transitions  Primary transitions  Secondary transitions

Set primary and trigger flags: Set top 2 ranked transitions as primary

Rank transitions by:  Abundance  Response Factor

<input type="checkbox"/>	Compound Name	Formula	MW	Polarity	Species	Precursor	Product	Frag	CE	CAV	Primary	Trigger	RT	RT Window	Abundance
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	155.7	380	10	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	5723674
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	129.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	315469
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	107.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	5015922
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	91.9	380	24	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.65	1	7438199
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	79.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	1529348
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	64.8	380	48	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	6508111
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	155.7	380	15	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2.03	1	3335528
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	244.8	380	12	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	548481
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	223.7	380	20	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	135578
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	217.7	380	16	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	434965
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	172.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	414589
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	107.9	380	20	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	2291401
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	91.9	380	32	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	3102685
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	79.9	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	761344
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	64.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	2583686
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	185.7	380	11	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	3720936
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	212.8	380	20	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	421592
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	155.9	380	12	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	1377529
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	123.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	2633848
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	107.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	2407737
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	91.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	3643700
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	79.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	1242301
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	64.9	380	48	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.98	1	3793552
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	155.8	380	6	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	4578854
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	177.8	380	12	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	89336
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	115.9	380	16	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	515222
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	107.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	3108252
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	92	380	28	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.37	1	5192801

You examine the Primary column and the Trigger column to determine which transitions are selected. You can select one or two Trigger transitions. You can select multiple Primary transitions.

#### 13 Review the Primary transitions and Trigger transitions.

- For sulfachloropyridazine, select 285  $m/z$  -> 156  $m/z$  transition as the Primary and Trigger transition.
- For sulfadimethoxine, select 311  $m/z$  -> 156  $m/z$  transition as the Primary and Trigger transition.
- For sulfamethazine, select 279  $m/z$  -> 186  $m/z$  transition as the Primary and Trigger transition.
- For sulfamethazole, select 271  $m/z$  -> 156  $m/z$  transition as the Primary and Trigger transition.

- Review each compound. Change the Primary and Trigger transitions to the transitions listed in the left column.
- Change the other Primary transitions as shown below.

- The software selected the most abundant transitions which in this example often had a low  $m/z$  for the Product Ion. A very abundant low  $m/z$  ion may be unsuitable as a Primary transition.
- You can select two Primary transitions as triggers for a compound.

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

#### Steps

#### Detailed Instructions

#### Comments

Select Transitions

Select top 10 ranked transitions

Primary transitions

Secondary transitions

Select Transitions

Set primary and trigger flags

Set top 2 ranked transitions as primary

Set Primaries and Trigger

Rank transitions by

Abundance

Response Factor

<input checked="" type="checkbox"/>	Compound Name	Formula	MW	Polarity	Species	Precursor	Product	Frag	CE	CAV	Primary	Trigger	RT	RT Window	Abundance
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	155.7	380	10	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.65	1	5723674
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	129.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	315469
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	107.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	5015322
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	91.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	7438199
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	79.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	1529348
<input checked="" type="checkbox"/>	sulfachloropyridazin			Positive		285	64.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	6608111
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	155.7	380	15	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2.03	1	3335528
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	244.8	380	12	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	548481
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	223.7	380	20	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	135578
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	217.7	380	16	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	434965
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	172.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	414589
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	107.9	380	20	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	2291401
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	91.9	380	32	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	3102685
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	79.9	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	761344
<input checked="" type="checkbox"/>	sulfadimethoxine			Positive		311	64.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	2583686
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	185.7	380	11	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.98	1	3720936
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	212.8	380	20	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	421592
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	155.9	380	12	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	1377529
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	123.9	380	24	7	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	2633848
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	107.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	2407737
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	91.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	3643700
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	79.8	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	1242301
<input checked="" type="checkbox"/>	sulfamethazine			Positive		279	64.9	380	48	7	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	3793552
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	155.8	380	6	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	4578854
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	177.8	380	12	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	89336
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	115.9	380	16	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	515222
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	107.9	380	24	7	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	3108252
<input checked="" type="checkbox"/>	sulfamethazole			Positive		271	92	380	28	7	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.37	1	5192801

**14** Review the **Import List** table on the **Import List** tab.

- a** Click the **Add to Import List** button.
- b** Click the **Import List** tab.
- c** Review the **Import List** table.
- d** Click the **Import** button.

- In this example, you are importing from the database to the **Import List**. Then, you are importing from **Database Browser** to **Optimizer**.

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

#### Steps

#### Detailed Instructions

#### Comments

Compound Name	Formula	Mw	Polarity	Species	Precursor	Product	Frag	CE	Primary	Trigger	RT	RT window	Abundance	RF	Acq Method	Project Name	Use Lists
sulfachloropyridazin			Positive		205	150.7	380	10	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.65	1	5723674		D\MassHunter\M	DefaultProject	
sulfachloropyridazin			Positive		205	129.9	380	24	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	375469		D\MassHunter\M	DefaultProject	
sulfachloropyridazin			Positive		205	107.9	380	24	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.65	1	5015922		D\MassHunter\M	DefaultProject	
sulfachloropyridazin			Positive		205	91.9	380	24	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	7436159		D\MassHunter\M	DefaultProject	
sulfachloropyridazin			Positive		205	79.8	380	48	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	1529248		D\MassHunter\M	DefaultProject	
sulfachloropyridazin			Positive		205	64.0	380	48	<input type="checkbox"/>	<input type="checkbox"/>	0.65	1	6506111		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	155.7	300	15	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2.03	1	3395528		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	244.0	300	12	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	546481		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	229.7	300	20	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	135570		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	217.7	380	16	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	434985		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	172.9	380	24	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	414569		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	107.9	380	20	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2.03	1	2291401		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	91.9	380	32	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	3102695		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	79.9	380	40	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	761344		D\MassHunter\M	DefaultProject	
sulfamethoxazole			Positive		311	64.0	380	48	<input type="checkbox"/>	<input type="checkbox"/>	2.03	1	2582696		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	185.7	300	11	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.98	1	3720936		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	212.0	380	20	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	421592		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	185.9	380	12	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	1377529		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	123.9	380	24	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.98	1	2633849		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	107.9	380	24	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	2407737		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	91.9	380	24	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	3642700		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	79.8	380	48	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	1242301		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		279	64.9	380	48	<input type="checkbox"/>	<input type="checkbox"/>	0.98	1	3792552		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	195.0	300	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	0.37	1	4570954		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	177.0	380	12	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	85336		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	115.9	380	16	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	915222		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	107.9	380	24	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	3106252		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	92	380	28	<input checked="" type="checkbox"/>	<input type="checkbox"/>	0.37	1	5193071		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	80	380	48	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	1196598		D\MassHunter\M	DefaultProject	
sulfanethazole			Positive		271	64.9	380	48	<input type="checkbox"/>	<input type="checkbox"/>	0.37	1	4791359		D\MassHunter\M	DefaultProject	

**15** Review the Compound Setup table in Optimizer. You replace all compounds with the compounds from the Database Browser program.

- Click the **Yes to All** button.
- In the Compound Setup tab in the Optimizer program, review the compounds.

- The compounds in the Optimizer program are overwritten by the compounds that you updated in the Database Browser program.

Replace

Similar record for the compound 'sulfachloropyridazine' is already present in current project.

Do you want to replace it?

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 2. Add/Modify compounds in an existing database

Steps	Detailed Instructions	Comments
<b>16</b> Save the new compound parameters to the database.	<ul style="list-style-type: none"><li>Click the <b>File &gt; Save Compounds</b> command to save all of the changes to the database.</li></ul>	<ul style="list-style-type: none"><li>You cannot see these results by default, but the Primary and Trigger transitions are updated in the project.</li><li>The <b>Primary</b> column, <b>Trigger</b> column, <b>Trigger Entrance Delay</b> column, <b>Trigger Delay</b> column, <b>Trigger Window</b> column and <b>Trigger MRM Threshold</b> column are available in the Compound Setup tab, but they are hidden by default.</li></ul>

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 3. Create a Triggered Dynamic MRM method from an existing database

## Task 3. Create a Triggered Dynamic MRM method from an existing database

You can create a Triggered Dynamic MRM method from a database such as the Pesticides or Forensics/Tox database. These databases can be purchased from Agilent. You can also copy the information from an Excel spreadsheet, but that method is not described in this guide.

Steps	Detailed Instructions	Comments
1 In the Data Acquisition program, you now import the updated compounds from the database. These compounds have optimized collision energies and also Primary and Trigger transitions marked.	<ol style="list-style-type: none"><li>a Switch to the Data Acquisition program.</li><li>b Open the <b>iiiSulfamix_dMRM2.m</b> method.</li><li>c In the <b>QQQ</b> tab, click the <b>Acquisition</b> tab. The Scan segments table contains four rows which are deleted later.</li><li>d Right-click the Scan Segments table and click <b>Import from Database Browser</b>. The Database Browser program opens.</li><li>e Mark the <b>Show All Records</b> check box.</li><li>f Mark all of the transitions for the four sulfa drug compounds. Clear the check boxes next to any unwanted compounds.</li><li>g Click the <b>Add to Import List</b> button.</li><li>h Click the <b>Import List</b> tab.</li><li>i Review the Import List table.</li><li>j Click the <b>Import</b> button.</li><li>k Delete the original compounds from the Scan segments table.</li><li>l Mark the <b>Triggered</b> check box under Triggered MRM.</li></ol>	<ul style="list-style-type: none"><li>• Before you import compounds from Database Browser, the Scan segments table contains at least one row. After importing compounds from the Database Browser, you need to remove any original rows.</li><li>• The Scan segments table always has to have at least one row.</li><li>• The triggering information is loaded from the Database Browser program even if the Triggered check box is clear.</li><li>• See the online Help for the Data Acquisition program and the QQQ Concepts Guide for an explanation of the other triggering conditions: <b>Trigger Entrance, Trigger Delay, and Trigger Window.</b></li></ul>

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 3. Create a Triggered Dynamic MRM method from an existing database

Steps	Detailed Instructions	Comments
-------	-----------------------	----------

Acquisition																
Scan segments																
Compound Group	Compound Name	ISTD?	Precursor Ion	MS1 Res	Product Ion	MS2 Res	Primary	Trigger	Threshold	Ret Time (min)	Delta Ret Time	Fragmentor	Collision Energy	Cell Accelerator Voltage	Polarity	Trigger Enable
	sulfachloropyridazin	<input type="checkbox"/>	285	Unit	155.7	Unit	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	48875	0.65	1	380	10	7	Positive	<input type="checkbox"/>
	sulfachloropyridazin	<input type="checkbox"/>	285	Unit	129.9	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	24	7	Positive	<input type="checkbox"/>
	sulfachloropyridazin	<input type="checkbox"/>	285	Unit	107.9	Unit	<input checked="" type="checkbox"/>	<input type="checkbox"/>		0.65	1	380	24	7	Positive	<input type="checkbox"/>
	sulfachloropyridazin	<input type="checkbox"/>	285	Unit	91.9	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	24	7	Positive	<input type="checkbox"/>
	sulfachloropyridazin	<input type="checkbox"/>	285	Unit	79.8	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	48	7	Positive	<input type="checkbox"/>
	sulfachloropyridazin	<input type="checkbox"/>	285	Unit	64.8	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	48	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	155.7	Unit	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	169666	2.03	1	380	15	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	244.8	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	12	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	229.7	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	20	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	217.7	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	16	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	172.9	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	24	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	107.9	Unit	<input checked="" type="checkbox"/>	<input type="checkbox"/>		2.03	1	380	20	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	91.9	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	32	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	79.9	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	48	7	Positive	<input type="checkbox"/>
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	64.8	Unit	<input type="checkbox"/>	<input type="checkbox"/>				380	48	7	Positive	<input type="checkbox"/>

Dynamic MRM Parameters: Cycle Time 200 ms Total MRMs = 4 Max Concurrent MRMs = 3 Min/Max Dwell = 65.78 ms/193.20 ms

Triggered MRM:  Triggered Repeats 3

- 2 Save the method to a new method name, **iiiSulfas\_TriggerOpt.m**, where **iii** are your initials.
  - a Click the **Method > Save Method** command.
  - b Type **iiiSulfas\_TriggerOpt.m**.
  - c Click the **Save** button.

## Exercise 3 – Create a Triggered Dynamic MRM acquisition method

### Task 3. Create a Triggered Dynamic MRM method from an existing database

#### Steps

3 Review the method in the Dynamic MRM Viewer dialog box.

#### Detailed Instructions

- Right-click the Scan segments table and click **Edit DMRM Method**. The Dynamic MRM Viewer dialog box is opened.
- Type 200 for the **Cycle time**. This value is shown in the Acquisition tab.
- Click between the **Primaries only** button and the **All transitions** button if the **Dynamic MRM Statistics** information is not updating. Then, click the **All transitions** button.

#### Comments

- The compounds in the Optimizer program were overwritten by the compounds that you updated in the Database Browser program.
- You can modify the **Cycle time** and see how the **Minimum Dwell Time** is changed. If the **Minimum Dwell Time** is less than 5 ms, and especially if it is less than 2 ms, then signal-to-noise is poor.
- A **Dwell Time** of 8 ms per transition is fine.

Dynamic MRM Viewer

Compound: [All] Compound Group: [All]

Compound Group	Compound Name	Precursor Ion	Product Ion	RT	RT Window	Primary	Trigger	Threshold	Frag	CE	CAV
	sulfachloropyridazine	285.00	155.70	650	1.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	48875	380	10	
	sulfachloropyridazine	285.00	107.90	650	1.000	<input checked="" type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfachloropyridazine	285.00	129.90	650	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfachloropyridazine	285.00	91.90	650	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfachloropyridazine	285.00	79.90	650	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	48	
	sulfachloropyridazine	285.00	64.90	650	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	48	
	sulfadimethoxine	311.00	155.70	2.030	1.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	169666	380	15	
	sulfadimethoxine	311.00	107.90	2.030	1.000	<input checked="" type="checkbox"/>	<input type="checkbox"/>		380	20	
	sulfadimethoxine	311.00	244.80	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	12	
	sulfadimethoxine	311.00	229.70	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	20	
	sulfadimethoxine	311.00	217.70	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	16	
	sulfadimethoxine	311.00	172.90	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfadimethoxine	311.00	91.90	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	32	
	sulfadimethoxine	311.00	79.90	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	48	
	sulfadimethoxine	311.00	64.80	2.030	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	48	
	sulfamethazine	279.00	165.70	980	1.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	183517	380	11	
	sulfamethazine	279.00	123.90	980	1.000	<input checked="" type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfamethazine	279.00	212.90	980	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	20	
	sulfamethazine	279.00	155.90	980	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	12	
	sulfamethazine	279.00	107.90	980	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfamethazine	279.00	91.90	980	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	24	
	sulfamethazine	279.00	79.90	980	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	48	
	sulfamethazine	279.00	64.90	980	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	48	
	sulfamethazole	271.00	155.80	370	1.000	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	68236	380	6	
	sulfamethazole	271.00	92.00	370	1.000	<input checked="" type="checkbox"/>	<input type="checkbox"/>		380	28	
	sulfamethazole	271.00	177.80	370	1.000	<input checked="" type="checkbox"/>	<input type="checkbox"/>		380	12	
	sulfamethazole	271.00	115.90	370	1.000	<input type="checkbox"/>	<input type="checkbox"/>		380	16	

Dynamic MRM Statistics

Total MRMs	30
Minimum Concurrent MRMs	7
Maximum Concurrent MRMs	21
Minimum Dwell Time	8.61 ms
Maximum Dwell Time	27.66 ms
Minimum Cycle Time	124.26 ms

Parameters

Cycle time: 200 ms

Calculations include:  
 Primaries only  All transitions

Review Tools

Override RT window: 1 min  
 Check minimum data pts: 64 pts

Split Method

Split method

Split by: Minimum Dwell Time

Number of methods: 2  
Max concurrent MRMs: 200  
Min dwell time: 5

Split Method: [ ]

Plot type: Concurrent MRMs  Select transitions on Click

Concurrent MRMs vs Retention Time


Concurrent MRMs

Retention Time (min)

Add Compounds... Save Split Methods... Reset Default Close

**Exercise 3 – Create a Triggered Dynamic MRM acquisition method**  
**Task 3. Create a Triggered Dynamic MRM method from an existing database**

Steps	Detailed Instructions	Comments
4	<p>Review the Trigger Thresholds to verify that they are appropriate.</p> <ol style="list-style-type: none"> <li>a Do an injection to make sure that the Trigger Thresholds are set properly.</li> <li>b Right-click the Scan segments table and click <b>Update DMRM Method</b>.</li> <li>c In the MRM Update Options dialog box, select <b>True</b> for <b>Update Threshold?</b></li> <li>d Enter the value for the <b>Percent of Height</b> for the Trigger Threshold.</li> <li>e Select the data file that you just acquired.</li> <li>f Click <b>OK</b>.</li> </ol>	



## Exercise 4 – Optimize Acquisition parameters

Task 1. Use the Optimizer Software to optimize acquisition parameters

# Exercise 4 – Optimize Acquisition parameters

For this exercise you optimize a mixture of four sulfonamide compounds.

## Task 1. Use the Optimizer Software to optimize acquisition parameters

The Optimizer Software helps you optimize acquisition parameters. Specifically, it automates the selection of the best precursor ions, the optimization of the fragmentor voltage for each precursor ion, selection of the best product ions, and optimization of collision energy values for each transition for a list of compounds you specify.

To do this task, you first need to create the method *iiiSulfamix MRM\_10.m* in “[Task 5. Find optimum collision energy for MRM acquisition](#)” on page 30. You do not need to acquire the data file.

The Fragmentor Voltage for the 6490 is set automatically during Autotune. The Fragmentor voltage for a 6490 is not optimized. The Fragmentor parameters and results will not be displayed for a 6490 instrument.

Steps	Detailed Instructions	Comments
1 Start the MassHunter Optimizer software.	• Double-click the <b>Optimizer</b> icon. 	• If you are optimizing peptides, use the <b>Optimizer for Peptides</b> program.


## Exercise 4 – Optimize Acquisition parameters

### Task 1. Use the Optimizer Software to optimize acquisition parameters

#### Steps

#### Detailed Instructions

#### Comments



## Exercise 4 – Optimize Acquisition parameters


### Task 1. Use the Optimizer Software to optimize acquisition parameters

Steps	Detailed Instructions	Comments
2 Set the optimization parameters.	<p><b>a</b> Click the <b>Optimizer Setup</b> tab.</p> <p><b>b</b> Set the <b>Sample introduction</b> method to <b>Injection</b>.</p> <p><b>c</b> Set the Fragmentor ramp parameters as follows:</p> <ul style="list-style-type: none"><li>• Set the range for ramping the <b>Fragmentor</b> values from 90 to 135.</li><li>• Clear the <b>Fragmentor Fine</b> check box.</li></ul> <p><b>d</b> Set the range for ramping the Collision Energy from 0 to 40 V.</p> <p><b>e</b> Select a Path for data files to store the optimization run data.</p> <p><b>f</b> Right-click the table on the right and select <b>Add Method</b> from the shortcut menu.</p> <p><b>g</b> Click the button on the right side of the Acq Method cell to open the Open Method dialog box.</p> <p><b>h</b> Select the method created in the previous exercise <b>iiiSulfamix MRM_10.m</b> and click <b>OK</b>. The Polarity and Ion Source will be filled in from the values set in the selected method.</p> <p><b>i</b> Check to make sure that the Ion Source from the method matches the physical configuration of your instrument.</p> <p><b>j</b> Repeat <a href="#">step f</a> to <a href="#">step i</a> to select additional methods.</p>	<ul style="list-style-type: none"><li>• Fine optimization refines the coarse ramping values and provides better optimization but takes longer to run.</li><li>• The data can be displayed later with Agilent MassHunter Qualitative Analysis software.</li><li>• The Fragmentor Voltage is not optimized for an Agilent 6490 Triple Quadrupole. It is set automatically when you Autotune. The Fragmentor parameters and results for a 6490 are not shown in the Optimizer program.</li></ul>

## Exercise 4 – Optimize Acquisition parameters

### Task 1. Use the Optimizer Software to optimize acquisition parameters


Steps	Detailed Instructions	Comments
3 Select the precursor ions	<ol style="list-style-type: none"><li>Click the <b>Precursor Ion Selection</b> tab.</li><li>Select the <b>Positive ions</b> +H adduct.</li><li>Select the <b>Charge state</b> of 1.</li><li>Set the search priority of the precursor ions.</li><li>(optional) To exclude certain masses from consideration, click <b>Exclude masses</b> at the bottom of the screen. Enter the <b>m/z Values</b> to exclude separated by commas and/or enter a <b>Minimum abundance value</b> in counts.</li></ol>	<ul style="list-style-type: none"><li>Mark the <b>Use most abundant precursor ion</b> check box to use the most abundant precursor ion.</li><li>Clear the <b>Use most abundant precursor ion</b> check box and use the <b>Up</b> and <b>Down</b> arrow buttons to set the search order (ions at the top of the list are given more priority).</li><li>You can also enter Neutral Losses to exclude (for example H<sub>2</sub>O).</li></ul>



## Exercise 4 – Optimize Acquisition parameters

### Task 1. Use the Optimizer Software to optimize acquisition parameters

Steps	Detailed Instructions	Comments
4 Select the product ions	<ol style="list-style-type: none"><li>Click the <b>Product Ion Selection</b> tab.</li><li>Enter a Low mass cut-off value. Select Mass (<math>m/z</math>) of 60 <math>m/z</math>.</li><li>To exclude certain masses from consideration, click <b>Exclude masses</b> option at the bottom of the screen. Enter the <b><math>m/z</math> Values</b> to exclude separated by commas and/or enter a <b>Minimum abundance value</b> in counts.</li><li>If desired, you can also enter <b>Neutral Losses</b> to exclude, for example <math>H_2O</math>. Enter a formula in the box and click the button to add it to the list.</li></ol>	



## Exercise 4 – Optimize Acquisition parameters

### Task 1. Use the Optimizer Software to optimize acquisition parameters

Steps	Detailed Instructions	Comments
<p>5 Set up a compound list. The formula for the four Sulfa Drugs are:</p> <ul style="list-style-type: none"> <li>• Sulfamethizole <math>C_9H_{10}O_2N_4S_2</math></li> <li>• Sulfamethazine <math>C_{12}H_{14}O_2N_4S</math></li> <li>• Sulfachloropyridazine <math>C_{10}H_9O_2N_4S</math></li> <li>• Sulfadimethoxine <math>C_{12}H_{14}O_4N_4S</math></li> </ul>	<p><b>a</b> Click the <b>Compound Setup</b> tab.</p> <p><b>b</b> Clear the <b>Show results summary</b> check box above the table while you set up the compound list.</p> <p><b>c</b> Right-click the table and select <b>Add Compound</b> from the shortcut menu to add a row to the end of the table.</p> <p><b>d</b> Enter Sulfamethizole as the <b>Compound Name</b>.</p> <p><b>e</b> Enter Sulfa drugs as the group name in the <b>Groups</b> column.</p> <p><b>f</b> Enter <math>C_9H_{10}O_2N_4S_2</math> as the <b>Formula</b> of the compound. The mass is calculated.</p> <p><b>g</b> Enter the <b>Sample Position</b> for the new compound.</p> <p><b>h</b> (optional) Enter an <b>Optimization dwell time</b> value to set longer or shorter cycle times.</p> <p><b>i</b> Repeat the steps above to add the other three sulfa drugs to the table.</p> <p><b>j</b> Mark the <b>Select</b> columns for the compounds (rows) to use for optimization.</p> <p><b>k</b> Save the compound list to the database or to the current project.</p>	<ul style="list-style-type: none"> <li>• Compounds are global to all projects. Compound information such as name, group, formula, and mass in one project will be reflected in the entire database.</li> <li>• If no methods or ions are specified here, then optimization for the compound uses the methods from the Optimizer Setup tab and information from the Precursor Ion Selection and Product Ion Selection tabs to generate the ions.</li> <li>• You can also enter the monoisotopic mass in the <b>Mass</b> column instead of the <b>Formula</b>.</li> </ul>

## Exercise 4 – Optimize Acquisition parameters

### Task 1. Use the Optimizer Software to optimize acquisition parameters

#### Steps

#### Detailed Instructions

#### Comments

**The Fragmentor parameters and results are not displayed for an Agilent 6490 Triple Quadrupole. The Fragmentor voltage for a 6490 is set automatically during Autotune.**

	Compound Name	Groups	Formula	Mass	Sample Position
1	Sulfamethizole	Sulfa drugs	C8H10O2N4S2	270.02	P1-A1
2	Sulfamethazine	Sulfa drugs	C12H14O2N4S	278.08	P1-A1
3	Sulfachloropyridazin	Sulfa drugs	C10H9O2N4SCl	284.01	P1-A1
4	Sulfadimethazine	Sulfa drugs	C12H14O4N4S	310.07	P1-A1

#### 6 Start the optimization process.

- Click the Start Optimization button () on the toolbar
- or
- Click the **Ion Breakdown Profile** button () on the toolbar.

#### 7 Review results.


- a** Click the **Compound Setup** tab.
  - b** Mark the **Show results summary** check box above the table.
  - c** Review the following values for each transition ion in the Compound Table:
    - Fragmentor
    - Collision Energy
  - d** Review the printed optimization report.
- (optional) Use the Agilent MassHunter Workstation Qualitative Analysis program to look at the data.
  - See the online Help for the Optimizer program or the Optimizer Quick Start Guide to learn how to import optimization results to acquisition for MRM time segments.

## Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

The “Source and iFunnel Optimizer” Software helps you optimize acquisition parameters for the source and iFunnel.

To do this task, you first need to create the method *iiiSulfamix\_dMRM2.m* in “[Task 4. Create a Dynamic MRM method from an MRM method](#)” on page 39. You do not need to acquire the data file. When you use this software, a worklist for each of the parameters being optimized is added to the Study Manager program.

Steps	Detailed Instructions	Comments
1 Start the MassHunter Data Acquisition software and load the <i>iiiSulfamix_dMRM2.m</i> method. Save this method to a new name.	<ol style="list-style-type: none"> <li>Start the Data Acquisition program.</li> <li>Make sure that Acquisition appears as the selection in the <b>Context</b> text box. If Tune is the selection, click <b>Acquisition</b> from the <b>Context</b> dropdown menu in the Combo bar.</li> <li>Load the <i>iiiSulfamix_dMRM2.m</i> method.</li> <li>Save this method to the name <i>iiiSulfamix_SourceOpt.m</i>.</li> <li>View the Method Editor window.</li> </ol>	<ul style="list-style-type: none"> <li>The first step is to create a template method. This method is used when you are optimizing the source and iFunnel parameters.</li> </ul>
2 Edit the Source parameters.	<ol style="list-style-type: none"> <li>Click the <b>QQQ</b> tab.</li> <li>Click the <b>Source</b> tab on the QQQ tab.</li> <li>Modify the parameters to the recommended starting parameters for source optimization. These parameters are shown in the following image.</li> </ol>	



## Exercise 4 – Optimize Acquisition parameters

### Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

#### Steps

3 Create template method.

#### Detailed Instructions

- Click the **Acquisition** tab.
- Select a single ion for each compound for the optimization.
- Save the method.

#### Comments

- A transition for each compound is already included in the Dynamic MRM method.

Properties DA Sampler Sampler Pretreatment Binary Pump Column Comp. DAD QQQ

Tune file: ahunes.tune.xml

Stop time:  No limit/As Pump

Ion source: AJS ESI

Time filtering:  Peak width 0.07 min

#	Start Time	Scan Type	Div Valve	Delta EMV (+)	Delta EMV (-)	Stored
1	0	Dynamic MRM	To MS	200	0	<input checked="" type="checkbox"/>

Compound Group	Compound Name	ISTD?	Precursor Ion	MS1 Res	Product Ion	MS2 Res	Ret Time (min)	Delta Ret Time	Fragmentor	Collision Energy	Cell Accelerator Voltage	Polarity
	sulfachloropyridazine	<input type="checkbox"/>	285	Unit	155.7	Unit	0.65	1	380	10	7	Positive
	sulfadimethoxine	<input type="checkbox"/>	311	Unit	155.7	Unit	2.03	1	380	15	7	Positive
	sulfamethazine	<input type="checkbox"/>	279	Unit	185.7	Unit	0.98	1	380	15	7	Positive
	sulfamethizole	<input type="checkbox"/>	271	Unit	155.8	Unit	0.37	1	380	10	7	Positive

Dynamic MRM Parameters: Cycle Time 200 ms Total MRMs = 4 Max Concurrent MRMs = 3 Min/Max Dwell = 65.78 ms/199.20 ms

Triggered MRM:  Triggered Repeats 3

4 Start the MassHunter “Source and iFunnel Optimizer” program.

- Double-click the **Source Optimizer** icon ().

Source and iFunnel Optimizer - SourceOptConfig.opt

File Tool Help

Project parameters

Optimize method: C:\MassHunter\methods\default.m

Project Folder: C:\MassHunter\data

Project Name: MyTest  Append timestamp

Instrument parameters

<input checked="" type="checkbox"/>	Types	PreWait (ms)	Replicate	StepWait (ms)	Start Value	End Value	Step Size
<input checked="" type="checkbox"/>	High Pressure RF	0	1	0	30	210	20
<input checked="" type="checkbox"/>	Low Pressure RF	0	1	0	40	180	20
<input checked="" type="checkbox"/>	Gas Temp	30	1	20	120	230	30
<input checked="" type="checkbox"/>	Gas Flow	30	1	0	11	20	2
<input checked="" type="checkbox"/>	Nebulizer	0	1	0	20	40	5
<input checked="" type="checkbox"/>	Capillary	0	1	0	1500	4500	500

Worklist parameters

Sample Name: MyTestSampleName Sample Position: Val 1

Worklist position of data file used for calibration: 1

Create Methods Submit Close

When this program starts, it automatically selects the default.m method. This method is not set up for an Agilent Jet Stream source, so no Agilent Jet Stream parameters are shown in the Instrument parameters table.


**Steps**

**Detailed Instructions**

**Comments**

5 Select the template method, **iiiSulfamix\_SourceOpt.m**.

- a Click the **Browse** button. The Browse For Folder dialog box is opened.
- b Select the **iiiSulfamix\_SourceOpt.m** method. Click the **OK** button.
- c If the ion source in the method is different than the ion source in the Instrument parameters list, a warning message is opened. Click **OK**.



The Sheath Gas Temp, Sheath Gas Flow, and Nozzle Voltage are all specific to the Agilent Jet Stream. If you do not have an Agilent Jet Stream source, these rows are not included in the table.

The High Pressure RF and Low Pressure RF are only included if the QQQ model is a 6490.

6 Change the order of the rows in the Instrument parameters table to the following:

- Ion Funnel parameters
- Sheath Gas temperature and flow
- Gas temperature and flow
- Nebulizer
- Capillary
- Nozzle Voltage

- a If necessary, select the row that shows the **High Pressure RF** parameter.
- b Drag this row to the top of the table. Both of the Ion Funnel parameters are moved together.
- c Verify that the order of the rows in your table is as indicated.

- The order of the parameters in the Instrument parameters table is the order that the parameters are optimized. You want to optimize the parameters that have the greatest effect on the source optimization first. The Ion Funnel parameters have the greatest effect, so you move those parameters to the top of the list.
- By default, the parameters are in the optimized list.

## Exercise 4 – Optimize Acquisition parameters

Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

### Steps

### Detailed Instructions

### Comments

Instrument parameters

	<input checked="" type="checkbox"/>	Types	PreWait (min)	Replicate	StepWait (min)	StartValue	EndValue	StepSize
▶	<input checked="" type="checkbox"/>	High Pressure RF	0	1	0	70	210	20
	<input checked="" type="checkbox"/>	Low Pressure RF	0	1	0	40	160	20
	<input checked="" type="checkbox"/>	Sheath Gas Temp	30	1	20	200	400	50
	<input checked="" type="checkbox"/>	Sheath Gas Flow	30	1	0	10	12	1
	<input checked="" type="checkbox"/>	Gas Temp	30	1	20	120	230	30
	<input checked="" type="checkbox"/>	Gas Flow	30	1	0	11	20	2
	<input checked="" type="checkbox"/>	Nebulizer	0	1	0	20	40	5
	<input checked="" type="checkbox"/>	Capillary	0	1	0	1500	4500	500
	<input checked="" type="checkbox"/>	Nozzle Voltage	0	1	0	0	2000	500

The High Pressure RF and the Low Pressure RF are always optimized together. If you move one of these rows in the table, the other row is also moved. All possible combinations of the High Pressure RF and Low Pressure RF are tried to find the optimal values. For this example, the High Pressure RF parameter has 8 different parameter settings (70, 90, 110, etc.). The Low Pressure RF parameter has 7 different parameter settings. So, the program automatically creates 56 (8 \* 7) different methods (1 for each parameter combination) just for optimizing the ion Funnel parameters.

7 Review the values for each parameter in the **Instrument parameters** table.

For each row in the table, verify:

- **PreWait** (in minutes).
- **Replicate**.
- **StepWait** (in minutes).
- **StartValue**.
- **EndValue**.
- **StepSize**.

- When the study for each parameter is loaded for the first time but before you run the first run, you wait the **PreWait** number of minutes before starting the run. Some parameters (that are electronic) stabilize almost instantly (in milliseconds), so you do not need to wait. For flows and temperatures, you want to have a **PreWait** before you run the study.
- You also want to wait for temperature parameters in between changing the parameter to a different value, so you also set the **StepWait** (in minutes).

8 Save the Instrument parameters.

- a Click **File > Save As (\*.opt)**.
- b Enter the name for this set of instrument parameters.
- c Click **OK**.

9 Modify the Instrument parameters table to only modify one parameter for this task. This task only optimizes the **Capillary** voltage.

- a Mark the check box next to the **Capillary**.
- b Clear the check boxes next to all of the other parameters.

- For this example, you optimize the **Capillary**. Usually, you optimize the parameters in the order specified in the Instrument parameters table.

Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

Steps	Detailed Instructions	Comments
-------	-----------------------	----------

Instrument parameters

	<input type="checkbox"/>	Types	PreWait (min)	Replicate	StepWait (min)	Start Value	End Value	Step Size
	<input type="checkbox"/>	High Pressure RF	0	1	0	70	210	20
	<input type="checkbox"/>	Low Pressure RF	0	1	0	40	160	20
	<input type="checkbox"/>	Sheath Gas Temp	30	1	20	200	400	50
	<input type="checkbox"/>	Sheath Gas Flow	30	1	0	10	12	1
	<input type="checkbox"/>	Gas Temp	30	1	20	120	230	30
	<input type="checkbox"/>	Gas Flow	30	1	0	11	20	2
	<input type="checkbox"/>	Nebulizer	0	1	0	20	40	5
	<input checked="" type="checkbox"/>	Capillary	0	1	0	1500	4500	500
	<input type="checkbox"/>	Nozzle Voltage	0	1	0	0	2000	500

10 Set the **Project Folder** and the **Project Name**.

- a Select the **Project Folder**. In this example, select `\MassHunter\Data`.
- b Enter a **Project Name**.
- c (optional) Mark the **Append timestamp** check box.

- If you mark the **Append timestamp** check box, then a time stamp is automatically added to the **Project Name** when you click the **Submit** button.

Project parameters

Optimize method

Project Folder

Project Name   Append timestamp

11 Set the **Worklist parameters**.

- a Type the **Sample Name**.
- b Type the **Sample Position**.
- c Type the **Worklist position of data file used for calibration**.

- If you mark the **Append timestamp** check box, then a time stamp is automatically added to the **Project Name** when you click **Submit**.
- For each parameter that is optimized, a batch file is created for Quantitative Analysis. One of the injections is considered 100% of the starting value. The value of **Worklist position of data file used for calibration** states which data file to use. If you enter `1`, then the data file from the first row is used.

Worklist parameters

Sample Name

Sample Position

Worklist position of data file used for calibration

## Exercise 4 – Optimize Acquisition parameters


### Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

Steps	Detailed Instructions	Comments
12 Create the methods and submit the study to the Study Manager.	<ol style="list-style-type: none"><li>Click <b>Create Methods</b>.</li><li>Click <b>Submit</b>.</li></ol>	<ul style="list-style-type: none"><li>When you click <b>Create Methods</b>, a message at the bottom of the main window states how many <b>Methods</b> were created, how many <b>Injections</b> are involved, and the <b>Estimated time</b>. The <b>Estimated time</b> is only an approximation.</li></ul>

Project created: D:\MassHunter\data\S\_Opt\_vcp\_2012912\_1 6 methods 6 injections Estimated time: 21 minutes

**The estimated time includes the Stoptime for the method plus one minute per injection. It does not consider the Posttime specified in the method. Also, it does not include the PreWait nor the StepWait that you entered in the Instrument parameters table.**

13 Review the study (or studies) submitted to Study Manager.	<ol style="list-style-type: none"><li>Open the Study Manager program.</li><li>Select a row in the Pending Studies table.</li><li>Right-click the row and click <b>Edit Worklist From Study</b>.</li><li>Review the worklist in the Edit Worklist dialog box. Click <b>Save</b>.</li></ol>	<ul style="list-style-type: none"><li>A study is submitted for each parameter that you marked in the Instrument parameters table.</li><li>Only one study is created for the <b>High Pressure RF</b> and <b>Low Pressure RF</b> parameters.</li></ul>
--	---	--



**The name of the study is the Instrument parameter that is being optimized. A separate study is added for each parameter that is being optimized.**

**You can examine or edit the worklist for the study. You right-click the line in the Pending Studies table and click Edit Worklist from Study.**

Steps

Detailed Instructions

Comments

	Sample Name	Sample Position	Method	Date File	Sample Type	Level Name	Comment
1	sulfas	Vial 31	S_Opt_Capillary2000 m	S_Opt_Capillary2000_1.d	Calibration	1	
2	sulfas	Vial 31	S_Opt_Capillary2500 m	S_Opt_Capillary2500_1.d	Sample		
3	sulfas	Vial 31	S_Opt_Capillary3000 m	S_Opt_Capillary3000_1.d	Sample		
4	sulfas	Vial 31	S_Opt_Capillary3500 m	S_Opt_Capillary3500_1.d	Sample		
5	sulfas	Vial 31	S_Opt_Capillary4000 m	S_Opt_Capillary4000_1.d	Sample		
6	sulfas	Vial 31	S_Opt_Capillary4500 m	S_Opt_Capillary4500_1.d	Sample		
7	Script SCP_Opt_Quant(D:\MassHunter\data\S_Opt_vcp_2012912_013223(Capillary3))\MPL_Acq_Scripts.exe						

The script that is run at the end of the worklist creates the Quantitative Analysis batch file.

14 Modify the Study Manager parameters to run a standby script when the study completes and then start the Study Manager.

- a Click the **Settings** tab in the Ribbon.
- b Mark the **Enable standby script execution on idle** check box.
- c Click the “...” button to select the script to execute.
- d Select **SCP\_InstrumentStandby** and click the **OK** button.
- e Enter 1 for the **Wait for** time.
- f Click the **Start** button if necessary.

- When the Study Manager is not running a study for the time specified, then the script you select is executed.

Pending Studies		Completed Studies				
Name	Path	Submitter	Plate Assignment	Est. Study Duration	Est. Start Time	Sample Count
▶ Capillary 3	D:\MassHunter\d...	admin	31	00:21:00		7


The name of the study is the Instrument parameter that is being optimized. A separate study is added for each parameter that is being optimized.

You can examine or edit the worklist for the study. You right-click the line in the Pending Studies table and click Edit Worklist from Study.

## Exercise 4 – Optimize Acquisition parameters

### Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

Steps	Detailed Instructions	Comments
15 When the study completes, stop the Study Manager queue and exit from the Study Manager program.	<ol style="list-style-type: none"><li>Click the <b>Stop &gt; Immediately</b> command.</li><li>Close the Study Manager program.</li></ol>	
16 Open the data in the Quantitative Analysis program.	<ol style="list-style-type: none"><li>Start the Quantitative Analysis program.</li><li>Click <b>File &gt; Open Batch</b>.</li><li>Navigate to the location of the study.</li><li>Select the <b>Batch file</b> named <i>Capillary.batch.bin</i> and click <b>Open</b>.</li></ol>	



You specified the Project Folder and the Project Name in the “Source and iFunnel Optimizer” program before you submitted the study.

The batch file is created automatically at the end of the study.

The “system” folder contains all of the methods that were used in this study.

17 Review the Batch Table.

- Switch to **Multiple Compound View**.
- Add the **Area** column to the table.
- For each compound, right-click the **Final Conc.** column and click **Plot this column**.
- Examine the Area column and the Final Conc. graph to determine the best capillary voltage.
- Close the Quantitative Analysis program.

- Refer to the online Help for the Quantitative Analysis program to learn how to do these tasks.
- In this case, all four compounds optimize at the same setting. Often, different compounds have different optimal settings, and you have to compromise.

## Exercise 4 – Optimize Acquisition parameters

### Task 2. Use the “Source and iFunnel Optimizer” program to optimize acquisition parameters

#### Steps

#### Detailed Instructions

#### Comments


Batch Table

Sample: sulfas    Sample Type: <All>    Compound: sulfamethazole    ISTD:

Name	Data File	Type	Level	Acq. Date-Time	sulfamethazole Results				sulfachloropyridazine Results				sulfamethazine Results				sulfadimethoxine Results			
					RT	Final Conc.	Accuracy	Area	RT	Final Conc.	Accuracy	Area	RT	Final Conc.	Accuracy	Area	RT	Final Conc.	Accuracy	Area
sulfas	S_Dpt_Capillary2000_1.d	Cal	1	9/12/2012 1:36 PM	0.365	100.0000	100.0	2395970	0.649	100.0000	100.0	1745306	0.977	100.0000	100.0	6145561	2.025	100.0000	100.0	369776
sulfas	S_Dpt_Capillary2000_1.d	Sample		9/12/2012 1:42 PM	0.365	82.8676		1952336	0.646	85.3296		1489263	0.970	73.8517		4538599	2.021	99.3833		3667020
sulfas	S_Dpt_Capillary3000_1.d	Sample		9/12/2012 1:48 PM	0.365	81.0429		1909346	0.646	82.1921		1434503	0.977	69.8619		4293407	2.025	97.4503		3595698
sulfas	S_Dpt_Capillary3500_1.d	Sample		9/12/2012 1:54 PM	0.365	73.6798		1735779	0.649	73.5986		1284520	0.977	62.5807		3849594	2.025	88.9432		3281804
sulfas	S_Dpt_Capillary4000_1.d	Sample		9/12/2012 2:00 PM	0.365	72.4241		1706291	0.646	71.1248		1241346	0.974	61.3445		3763965	2.025	87.5352		3223852
sulfas	S_Dpt_Capillary4500_1.d	Sample		9/12/2012 2:07 PM	0.365	69.3667		1634259	0.649	67.9969		1186752	0.977	58.4264		3590633	2.025	83.7877		3076487

18 (optional) Review the data files in the Qualitative Analysis program.

- Start the Qualitative Analysis program.
  - Open all of the data files in the study.
  - Click **Find > Find Compounds by MRM**.
  - Select all of the data files and click the **Find** button.
  - Click the **Edit > Auto-Color Mode > Single Color per Data File** menu item.
  - Clear the check boxes next to the **TIC** for each data file.
  - Examine the results in the Chromatogram Results window.
  - Close the Qualitative Analysis program.
- By default, the program selects different colors for different transitions.
  - It is clear that the conditions used for the blue chromatograms are the best, and the blue chromatograms are for the data file with the capillary voltage set to 2000.



[www.agilent.com](http://www.agilent.com)

## In This Book

This exercise helps you use the Agilent 6400 Series Triple Quadrupole LC/MS system. In this guide, you acquire data and then analyze the results using the Qualitative Analysis program to learn how to develop an acquisition method.

© Agilent Technologies, Inc. 2012

Printed in USA  
Revision A, November 2012



G3335-90136



**Agilent Technologies**