



Multiresidue Analysis of Pesticides in Olive Samples Using GC/MS/MS

Application Note

Food Testing and Agriculture

Authors

A. Moreno López, L. Moreno López, and
J.L. Pineda Lucas
Laboratorio Químico Microbiológico
S.A.
Sevilla, Spain
Joan Stevens
Agilent Technologies, Inc.

Abstract

This application note describes an analytical method for the determination of 27 pesticide residues in olives. Olives have a high lipid content of about 80 to 85%, which can adversely affect pesticide recoveries and chromatographic systems. Therefore, a modified QuEChERS method for the extraction and analysis by gas chromatography/triple quadrupole (GC/MS/MS) multiple reaction monitoring method was employed. The method was validated for olives in terms of recovery, repeatability, and reproducibility. The results demonstrated that the method achieved acceptable quantitative recoveries of 70 to 120%, as recommended in SANCO/12571 [1], with RSDs < 20%. Limits of quantification at or below the regulatory maximum residue limits for the pesticides were achieved.

Introduction

Pesticide residue analysis is essential for the health of humans and animals, the import/export trade, and regulated control purposes. Many pesticide classes are used in agriculture, and most pesticides have regulatory guidelines, for example maximum residue levels (MRLs) in food, with analytical determination based on GC/MS or LC/MS/MS [2].

Methods using gas chromatography coupled to mass spectrometry (GC/MS) are based on selected ion monitoring (SIM), which is fairly sensitive, but identification potential and non-target/retrospective analysis capabilities are sacrificed. Determination of GC-amenable pesticides in food samples by using GC/MS/MS has emerged in the last decade as a valuable approach, which allows higher selectivity and sensitivity, while minimizing or even eliminating most chromatographic



Agilent Technologies

interferences [3]. Despite the increased sensitivity and selectivity of GC/MS/MS co-extract matrix components can cause matrix effects that can negatively impact results. An optimized sample preparation method can minimize these effects.

Materials and Methods

A modified EN extraction method was used for the extraction of 27 pesticides from olives. The solvent-modified method, which incorporated a mixture of nonpolar and polar aprotic solvents (ethyl acetate:cyclohexane:acetone) versus the polar aprotic solvent (acetonitrile) described in EN 15662 [4] was used to extract the pesticides from the high-lipid matrix associated with olives. A mixture of nonpolar and polar aprotic solvents can offer a wider polarity range for extractable compounds relative to the use of only a polar aprotic solvent, ACN. The pesticides included organochlorine, organophosphate, and pyrethroid classes.

Consumables and instruments

- Agilent Bond Elut QuEChERS EN Extraction Kit (p/n 5982-5650)
- Agilent Bond Elut QuEChERS Dispersive Kit, EN Method, 15 mL (p/n 5982-5156)
- Agilent SPE Bulk Sorbent, C18 Endcapped (p/n 5982-8082)
- Agilent J&W HP-5ms Ultra Inert GC column, 30 m × 0.25 mm (p/n 19091S-433UI)
- Agilent 7890A GC
- Agilent 7000A Triple Quadrupole GC/MS System
- Agilent 7693 ALS Injector

Solvents and pesticide standards

- Ethyl acetate, residue analysis purity, Baker
- Cyclohexane, residue analysis purity, Sharlau
- Acetone, residue analysis purity, Baker
- Solvent mixture: ethyl acetate:cyclohexane:acetone (1:1:4)
- Pesticide standards, Dr. Ehrenstorfer

GC conditions

Carrier:	Constant pressure, 22.0 psi
Oven temperature:	Initial 70 °C (2 min), 25 °C/min to 150 °C (0 min), 3 °C/min to 200 °C (0 min), 8 °C/min to 280 °C (10 min)
Injector temperature:	250 °C
Septum purge:	On, 3 mL/min
Purge flow split vent:	100 mL/min at 0.75 minutes
Injection:	Splitless, 1.0 µL
RTL compound:	Chlorpyrifos methyl

MS conditions

Ion source:	EPC
Source temperature:	280 °C
Collision gas:	He quench gas on, 2.35 mL/min N ₂ collision gas on, 1.5 mL/min
Transfer line temperature:	280 °C
MS quad 1,2 temperature:	150 °C both
MS1/MS2 resolution:	Wide/wide
MRM settings:	See Table 1

Table 1. Retention time and MRM parameters for 27 pesticides.

Pesticide	RT (locked to chlorpyrifos methyl)
Chlorpropham	11.05
Heptenophos	9.737
Pirimiphos methyl	18.307
Hexachlorobenzene	12.377
Fonofos	13.889
Simazine	12.909
Terbufos	13.796
Terbuthylazine	13.810
Diazinon	14.466
Pirimicarb	15.677
Chlorpyrifos methyl	16.59
Chloropropylate	25.419
Parathion-methyl	16.594
Bifenthrin	28.839
Fenitrothion	18.072
Parathion	19.275
Chlorthal dimethyl	19.433
Chlorfenvinphos	21.557
Pendimethalin	20.991
<i>Alpha</i> -endosulfan	22.637
Procymidone	21.962
<i>Beta</i> -endosulfan	25.3158
Endosulfan sulfate	26.76
Acrinathrin	30.724
Aldrin	18.528
Chlorpyrifos	19.234
4,4'-Dichlorobenzophenone (dicofol)	19.201

Table 2 shows the time segments. MS1 and MS2 resolutions were wide/wide throughout. Dwell time was 10 ms.

Table 2. Time segments.

Compound	Precursor ion	Product ion	Collision energy	Compound	Precursor ion	Product ion	Collision energy
Time segment 2				Time segment 8			
Chlorpropham	213	171	5	4,4'-Dichlorobenzophenone	139	111	15
Chlorpropham	213	127	5	4,4'-Dichlorobenzophenone	139	75	35
Heptenophos	124	89	20	Aldrin	263	193	30
Heptenophos	124	63	35	Aldrin	263	191	30
Time segment 3				Chlorpyrifos	197	169	15
Hexachlorobenzene	284	249	20	Chlorpyrifos	197	107	40
Hexachlorobenzene	284	214	35	Fenitrothion	277	125	15
Time segment 4				Fenitrothion	277	109	20
Fonofos	246	109	15	Parathion	291	109	10
Fonofos	246	81	30	Parathion	291	81	25
Simazine	201	186	5	Pirimiphos methyl	305	290	10
Simazine	201	173	5	Pirimiphos methyl	305	180	5
Terbufos	231	175	10	Time segment 9			
Terbufos	231	129	25	Chlorthal dimethyl	301	223	25
Terbuthylazine	214	132	10	Chlorthal dimethyl	299	221	25
Terbuthylazine	214	104	20	Time segment 10			
Time segment 5				Chlorfenvinphos	267	159	15
Diazinon	179	137	20	Chlorfenvinphos	267	81	30
Diazinon	179	121	40	Pendimethalin	252	162	10
Fonofos	246	109	15	Pendimethalin	252	161	20
Fonofos	246	81	30	Time segment 11			
Terbufos	231	175	10	Endosulfan (alpha isomer)	241	206	10
Terbufos	231	129	25	Endosulfan (alpha isomer)	241	170	20
Terbuthylazine	214	132	10	Procymidone	283	255	10
Terbuthylazine	214	104	20	Procymidone	283	96	10
Time segment 6				Time segment 13			
Pirimicarb	238	166	10	Endosulfan (beta isomer)	241	206	15
Pirimicarb	166	96	15	Endosulfan (beta isomer)	195	159	5
Terbacil	161	144	10	Time segment 14			
Terbacil	161	88	20	Endosulfan (beta isomer)	241	206	15
Time segment 7				Endosulfan (beta isomer)	195	159	5
Chlorpyrifos methyl	286	271	20	Chlorpropylate	251	139	15
Chlorpyrifos methyl	286	93	25	Chlorpropylate	251	111	35
Methyl parathion	263	109	15	Time segment 15			
Methyl parathion	263	79	30	Endosulfan sulfate	272	237	20
				Endosulfan sulfate	272	117	40
				Time segment 16			
				Bifenthrin	181	166	20
				Bifenthrin	181	165	25
				Time segment 17			
				Acrinathrin	289	93	5
				Acrinathrin	208	181	5

Sample preparation

Olives were picked from local trees. The olives were cut into thin slices and a portion of the sample (approximately 1 g) was used for the analysis. The rest was stored frozen. If the olives contained a stone, it was removed and an estimate of the percentage contribution of stone present in the final sample was calculated. This is a legislative requirement for LMR in calculations. This was done by taking a representative portion of the sample (two or three whole pieces) to evaluate the percentage of the stone, as follows. First, weigh the representative samples containing the stone. Then, remove the stones from the samples and weigh them. Finally, calculate the percentage associated with the stones using Equation 1.

$$\% H = \frac{\text{Weight of stone}}{\text{Weight of sample}} \times 100 \quad \text{Eqn.1}$$

Where:

% H = percentage of the stone

Procedure

1. Weigh 10.0 g \pm 0.1 g of chopped sample in a 50-mL centrifuge tube.
2. Add 10 mL of the solvent mixture ethyl acetate:cyclohexane:acetone (1:1:4) then cap the tube and shake for at least 1 minute.
3. Add the Bond Elut QuEChERS extraction salt packet (p/n 5982-5650) to the sample tube; shake for 1 minute.
4. Centrifuge the tube at 4,000 rpm for 5 minutes.
5. Transfer 6 mL of the upper organic layer into a 15-mL Bond Elut dispersive tube (p/n 5982-5156) to which an additional 300 mg of C18EC (p/n 5982-8082) has been added.
6. Shake for at least 1 minute, then centrifuge at 4,000 rpm for 5 minutes.
7. Transfer 1 to 2 mL of the upper organic layer extract to a chromatography vial and seal tightly.

Results and Discussion

All target pesticides were separated and well detected by the HP-5ms Ultra Inert GC column. With the powerful selectivity provided by GC/MS/MS, the MRM chromatograms of the matrix blank did not show any interference peaks with the target analytes. Figure 1 shows the GC/MS/MS chromatogram of a 20 $\mu\text{g}/\text{kg}$ -fortified olive extract processed by the modified EN QuEChERS extraction method. Figure 2 shows the exceptional quantitative analysis of pesticides considered to be more problematic in olives and olive oil analysis, such as terbuthylazine and chlorpyrifos.

Calculations

Once the pesticide concentration in the vial was obtained from the calibration curve, we worked out sample concentrations in mg/kg using Equation 2.

$$C_{\text{Pesticide}} (\text{mg}/\text{kg}) = C_{\text{Pesticide}} (\mu\text{g}/\text{L}) \times \frac{10}{W} \times 1e-3 \times [1 - (\frac{\% H}{100})] \quad \text{Eqn.2}$$

Where:

$C_{\text{Pesticide}}$ = concentration ($\mu\text{g}/\text{L}$) obtained from the calibration curve

W = initial weight of the sample taken for the test in g (10.0 g)

% H = percentage contribution of stones (if necessary)

The final result is expressed in I.S. units, and results given in $\mu\text{g}/\text{kg}$.




Figure 1. Overlay of MRM transitions of 27 pesticides spiked in olives at 20 µg/kg obtained using GC/MS/MS.





Figure 2. A) Terbutylazine at 5 µg/kg g/kg, $R^2 = 0.9999$. B) Chlorpyrifos at 5 µg/kg, $R^2 = 0.9999$.

Linearity and limit of quantification (LOQ)

The linearity calibration range for all the pesticides tested was 5 to 60 µg/kg. Calibration curves using spiked matrix blanks were made at 5, 10, 20, and 60 µg/kg, where 60 µg/kg was at least 120% of the upper limit. The calibration curves were generated by plotting the relative responses of the analytes. The 10 µg/kg quantification limits established for all pesticides were at or below the MRLs of these pesticides in fruits and vegetables. The correlation coefficients (R^2) for all compounds were > 0.9999.

Recovery and reproducibility

Recovery and reproducibility were evaluated by spiking pesticide standards in comminuted olive samples at levels of 5, 10, and 20 µg/kg, which were then prepared by the EN QuEChERS extraction method and the solvent-modified cleanup method. These QC samples were quantitated against the matrix spike calibration curve. The analysis was performed in replicates of five at each level. The recovery and reproducibility (RSD) data are shown in Table 3. It can be seen that the 27 pesticides had excellent recoveries and precision.

Table 3. Recovery %, (RSD%), and reproducibility of pesticides in olives with Agilent Bond Elut QuEChERS EN extraction method and solvent modified dispersive SPE kit for fruits and vegetables with fats and waxes.

Compound	L1	L2	L3
4,4'-Dichlorobenzophenone	93.94 (7.32)	92.95 (7.62)	106.43 (6.28)
Acrinathrin	101.92 (6.66)	95.13 (14.01)	105.25 (2.00)
Aldrin	90.46 (12.86)	90.36 (9.60)	100.00 (10.52)
Bifenthrin	96.50 (9.12)	95.88 (7.87)	104.29 (8.01)
Chlorpyrifos	93.15 (4.43)	86.07 (14.21)	108.10 (5.53)
Chlorthal dimethyl	98.78 (5.78)	96.05 (7.47)	99.96 (16.41)
Chlorfenvinphos	93.10 (6.04)	89.01 (11.80)	106.86 (6.38)
Chloropropylate	91.24 (7.46)	86.68 (7.88)	108.70 (5.41)
Chlorpyrifos methyl	100.40 (11.72)	92.70 (8.63)	110.22 (3.23)
Chlorpropham	98.38 (7.30)	91.00 (15.89)	105.07 (6.86)
Diazinon	94.11 (10.24)	94.68 (14.24)	108.81 (5.98)
Endosulfan (<i>alpha</i> isomer)	89.13 (16.31)	93.97 (11.88)	105.55 (5.64)
Endosulfan <i>beta</i>	83.90 (21.05)	93.06 (10.03)	98.25 (4.41)
Endosulfan sulfate	94.67 (4.76)	93.48 (7.74)	109.58 (4.75)
Fenitrothion	91.48 (8.74)	86.33 (17.91)	109.38 (10.91)
Fonofos	89.55 (5.21)	91.06 (9.75)	104.63 (7.52)
Heptenophos	97.24 (4.02)	93.55 (17.02)	101.78 (12.40)
Hexachlorobenzene	95.25 (3.52)	87.10 (8.51)	105.62 (7.34)
Parathion-methyl	94.00 (15.54)	85.27 (10.50)	104.57 (7.81)
Parathion	86.26 (11.20)	81.57 (10.98)	94.98 (10.91)
Pendimethalin	90.25 (12.40)	81.78 (7.91)	96.00 (4.64)
Pirimicarb	84.52 (5.09)	76.15 (5.43)	89.66 (10.16)
Pirimiphos methyl	93.02 (7.29)	85.34 (13.69)	106.92 (5.54)
Procymidone	90.75 (7.97)	85.04 (3.41)	98.88 (7.29)
Simazine	94.12 (4.59)	90.17 (11.80)	105.81 (8.75)
Terbufos	105.06 (8.45)	98.92 (9.06)	112.62 (6.92)
Terbutylazine	97.31 (9.04)	91.66 (7.59)	105.34 (9.99)

Level	Calibration points (µg/L)	Limit of quantification (µg/kg)
L1	5.0 (at least 70% LQ)	10.0
L2	10	
L3	20	
L4	60 (at least 120% of upper limit)	

Conclusions

The Agilent Bond Elut QuEChERS EN buffered extraction kits and dispersive-SPE kits for fruits and vegetables containing fats and waxes, with a modified-solvent mixture, provided a simple, fast, and effective method for the purification of representative pesticides in olives. The recovery and reproducibility, based on matrix-spiked standards, were acceptable for multiclass pesticide determination in olives. The matrix effects from the olives did not interfere with the quantitation of target compounds. The LOQs of the pesticides were at or below regulated MRLs in foods. As the selected pesticides represented a broad variety of different classes and properties, Bond Elut QuEChERS EN extraction and dispersive-SPE with modified solvent mixtures are excellent choices for other pesticides in similar food matrixes.

References

1. Anon. Guidance document on analytical quality control and validation procedures for pesticide residues analysis in food and feed. SANCO/12571_2013. European Commission (2013).
2. L. Alder, *et al. Mass Spectrom. Rev.* **25**, 838 (2006).
3. F. Hernández, *et al. Anal. Methods* **5**, 5875 (2003).
4. Anon. Foods of plant origin. Determination of pesticide residues using GC-MS and/or LC-MS/MS following acetonitrile extraction/partitioning and clean-up by dispersive SPE. QuEChERS-method. EN 15662. European Standard.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2014
Printed in the USA
June 18, 2014
5991-4800EN



Agilent Technologies