

Lipid Analysis

by Evaporative Light Scattering Detection



Introduction

The composition of the lipid fraction in fats and oils influences the nutritional value of the food and is important for health related issues such as digestion, heart disease, and more.

Triglycerides are a major component of the non-polar lipid fraction. This is usually very complex due to various chain lengths and different degrees of saturation. Phospholipids and sphingolipids have a positive influence on cholesterol levels and are also used for their emulsifying qualities.

Combine the benefits of Evaporative Light Scattering Detection with Alltech's HPLC columns for sensitive, efficient, gradient-compatible separations of both polar and non-polar lipid fractions.



Benefits of Evaporative Light Scattering Detection (ELSD) for Lipid Analysis

- More sensitive than RI and low-wavelength UV
- Gradient compatible
- Eliminates need for derivatization
- Solvent choice not limited by UV-cutoff



The Evaporative Light Scattering Detector (ELSD) is a universal HPLC detector that detects any compound less volatile than the mobile phase to low nanogram levels. It has many advantages over RI and low-wavelength UV detection for lipid analyses.

Since most lipid compounds are not good UV absorbers, they must be monitored at low wavelengths. This can be hindered by the mobile phase solvent's UV cutoff and poor baseline stability. With ELSD, the mobile phase evaporates before detection, so it is blind to the mobile phase's spectroscopic properties and compositional changes.

Although refractive index (RI) detection is universal in nature, its lack of sensitivity and gradient incompatibility make it less than ideal for lipid analysis. Evaporation of the mobile phase before detection makes the ELSD gradient compatible for improved resolution and faster separations. In addition, ELSDs can be as much as 10 to 1000 times more sensitive than RI detectors.

References

[Analysis of Phospho- and Sphingolipids in Dairy Products by a New HPLC Method](#), R. Rombaut, J. V. Camp and K. Dewettinck - Ghent University, Belgium, published in J. Dairy Sci. 88:482-488.

[Triglyceride Analysis of Edible Oils by HPLC-ELSD](#), R. Rombaut, I. Foubert, V. Van Hoed, K. Degraeve, R. Verhé and K. Dewettinck - Ghent University, Belgium, to be published.


How does Evaporative Light Scattering Detection work?

ELSDs detect any compound less volatile than the mobile phase using a simple three-step process:




1. Nebulization

Column effluent passes through a needle and mixes with nitrogen gas to form a dispersion of droplets.



2. Evaporation

Droplets pass through a heated "drift tube" where the mobile phase evaporates, leaving a fine mist of dried sample particles in solvent vapor.



3. Detection

The sample particles pass through a cell and scatter light from a laser beam. The scattered light is detected, generating a signal.



Alltima™ HP C18 Hi-Load Columns for Triglycerides

- Excellent peak symmetry
- Superior retention and loading capacity
- Stable from pH 1-10
- Extremely low bleed



Alltima HP C18 Hi-Load columns use high purity, low metal, pH stable silica for high resolution with no tailing and little to no bleed. The high carbon load, combined with strong organic mobile phases, enhance triglyceride solubility and performance with ELSD or MS detection. Alltima HP C18 Hi-Load columns are ideal for analysis of complex hydrophobic mixtures.



Prevail™ Silica for Polar Lipids

- High capacity and resolution of polar lipids
- Long column life
- Stable from pH 1-8

Prevail silica's high surface area and high purity increase capacity and resolution of polar lipids while maintaining excellent peak shape. Choose 3µm particles for high efficiencies, and a variety of column formats for speed and ELSD and MS compatibility. Prevail Silica is ideal for analysis of polar lipids in dairy products.




Triglycerides



Butter

CHROM
10563



Column: Alltima™ HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:

Time:	0	25	26	27	28
%B:	30	51	70	70	30


Flow Rate: 0.72mL/min

Detector: ELSD

Column Temp: 20°C

Coconut Oil

CHROM
10564



Column: Alltima HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:

Time:	0	25	26	27	28
%B:	30	51	70	70	30

Flow Rate: 0.72mL/min


Detector: ELSD

Column Temp: 20°C

- | | |
|-----------|------------------|
| 1. CCLa | 8. LaLaO |
| 2. C11 | 9. LaMM |
| 3. CLaLa | 10. LaMO |
| 4. CCM | 11. LaMP/MMM |
| 5. LaLaLa | 12. MMP/LaPP/MPO |
| 6. CLaM | 13. C17 |
| 7. LaLaM | |

Cocoa Butter

CHROM
10565



Column: Alltima HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:

Time:	0	25	26	27	28
%B:	30	51	70	70	30

Flow Rate: 0.72mL/min


Detector: ELSD

Column Temp: 20°C

- | | |
|------------|---------|
| 1. C13 | 8. SLS |
| 2. PLO | 9. POS |
| 3. PLP | 10. PPS |
| 4. POO | 11. SOS |
| 5. PLS | 12. PSS |
| 6. POP | 13. SOA |
| 7. SOO/PPP | 14. C19 |

Corn Oil

CHROM
10566



Column: Alltima HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:





Time:	0	25	26	27	28
%B:	30	51	70	70	30

Flow Rate: 0.72mL/min

Detector: ELSD

Column Temp: 20°C

- | | |
|---------|-------------|
| 1. C11 | 8. PLP |
| 2. LLLn | 9. OOO |
| 3. LLL | 10. SLO |
| 4. LLO | 11. POO |
| 5. PLL | 12. POP |
| 6. LOO | 13. SOO/PPP |
| 7. PLO | 14. C17 |






Triglycerides

Palm Oil

CHROM
10571



- | | |
|--------|-------------|
| 1. C13 | 9. POO |
| 2. OLL | 10. PLS |
| 3. PLL | 11. POP |
| 4. LOO | 12. SOO/PPP |
| 5. PLO | 13. POS |
| 6. PLP | 14. PPS |
| 7. MOP | 15. SOS |
| 8. OOO | 16. C19 |

Column: Alltima™ HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:

Time:	0	25	26	27	28
%B:	30	51	70	70	30


Flow Rate: 0.72mL/min

Detector: ELSD

Column Temp: 20°C

Palm Kernel Stearin Fraction

CHROM
10572



- | | |
|-----------|------------------|
| 1. CCLa | 8. LaMM |
| 2. C11 | 9. LaMO |
| 3. CLaLa | 10. LaMP/MMM |
| 4. CMM | 11. MMO/LaPO |
| 5. LaLaLa | 12. MMP/LaPP/MPO |
| 6. LaLaM | 13. C17 |
| 7. LaLaO | |

Column: Alltima HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:

Time:	0	25	26	27	28
%B:	30	51	70	70	30


Flow Rate: 0.72mL/min

Detector: ELSD

Column Temp: 20°C

Palm Kernel Olein Fraction

CHROM
10573



- | | |
|-----------|------------------|
| 1. CCLa | 10. LaMO |
| 2. C11 | 11. LaMP/MMM |
| 3. CLaLa | 12. MMO/LaPO |
| 4. CCM | 13. MMP/LaPP/MPO |
| 5. LaLaLa | 14. PLP |
| 6. CLaM | 15. POO |
| 7. LaLaM | 16. POP |
| 8. LaLaO | 17. C17 |
| 9. LaMM | |

Column: Alltima HP C18 Hi-Load, 3µm, 150 x 3.0mm
(Part No. 87610)

Mobile Phase: A: Acetonitrile
B: Dichloromethane

Gradient:

Time:	0	25	26	27	28
%B:	30	51	70	70	30

Flow Rate: 0.72mL/min

Detector: ELSD

Column Temp: 20°C





Compound Translation Key

Triglycerides

A	Arachidic Acid
C	Caproic Acid
La	Lauric Acid
L	Linoleic Acid
Ln	Linolenic Acid
M	Myristic Acid
O	Oleic Acid
P	Palmitic Acid
S	Stearic Acid

Polar Lipids

PS	Phosphatidylserine
PE	Phosphatidylethanolamine
PI	Phosphatidylinositol
PC	Phosphatidylcholine
PA	Phosphatidic Acid
LPE	lyso-Phosphatidylethanolamine
LPC	lyso-Phosphatidylcholine
GLUCCER	Glucosylceramide
LACCER	Lactosylceramide
SM1	Sphingomyelin 1
SM2	Sphingomyelin 2
U	Unknown



Ordering Information

Triglycerides – Alltima™ HP Columns

Alltima HP HPLC Columns			
PHASE	FORMAT	I.D. x LENGTH	PART No.
C18 Hi-Load, 3µm	Expedite™ MS	2.1 x 10mm	87691
	Expedite MS	2.1 x 20mm	87692
	LC/MS	2.1 x 50mm	87506
	LC/MS	2.1 x 100mm	87687
	LC/MS	2.1 x 150mm	87688
	Solvent-Reducer	3.0 x 150mm	87610
	Expedite MS	4.6 x 10mm	87693
	Expedite MS	4.6 x 20mm	87694
	Analytical	4.6 x 50mm	87827
	Analytical	4.6 x 100mm	87685
	Analytical	4.6 x 150mm	87686
	Rocket™	7.0 x 33mm	87689
	Rocket	7.0 x 53mm	87690
	C18 Hi-Load, 5µm	LC/MS	2.1 x 150mm
LC/MS		2.1 x 250mm	87700
Solvent-Reducer		3.0 x 150mm	87611
Solvent-Reducer		3.0 x 250mm	87612
Analytical		4.6 x 150mm	87697
Analytical		4.6 x 250mm	87698

Alltima HP All-Guard™ Cartridges*			
PHASE	I.D. x LENGTH	QTY	PART No.
C18 Hi-Load	2.1 x 7.5mm	3	87701
	3.0 x 7.5mm	3	87624
	4.6 x 7.5mm	3	87702
All-Guard Cartridge Holder (Includes Direct-Connect™ Column Coupler)		ea	80101

*All-Guard Holder required

Polar Lipids – Prevail™ Columns

Prevail HPLC Columns			
PHASE	FORMAT	I.D. x LENGTH	PART No.
Silica, 3µm	Expedite MS	2.1 x 10mm	43841
	Expedite MS	2.1 x 20mm	43826
	LC/MS	2.1 x 50mm	43868
	LC/MS	2.1 x 100mm	43805
	LC/MS	2.1 x 150mm	99267
	Solvent-Reducer	3.0 x 150mm	99341
	Expedite MS	4.6 x 10mm	43858
	Expedite MS	4.6 x 20mm	43816
	Analytical	4.6 x 50mm	43842
	Analytical	4.6 x 100mm	99269
	Analytical	4.6 x 150mm	99271
	Rocket	7.0 x 33mm	99284
Rocket	7.0 x 53mm	99283	
Silica, 5µm	Analytical	2.1 x 150mm	99273
	Solvent-Reducer	3.0 x 150mm	99339
	Solvent-Reducer	3.0 x 250mm	99340
	Analytical	4.6 x 150mm	99275
	Analytical	4.6 x 250mm	99277

Prevail™ All-Guard Cartridges*			
PHASE	I.D. x LENGTH	QTY	PART No.
Silica	2.1 x 7.5mm	3	99133
	3.0 x 7.5mm	3	99354
	4.6 x 7.5mm	3	99291
All-Guard Cartridge Holder (Includes Direct-Connect Column Coupler)		ea	80101

*All-Guard Holder required

© Copyright 2010 Alltech Associates, Inc. Alltech Associates Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Davison is a business unit of W. R. Grace & Co.-Conn., Global Headquarters 7500 Grace Drive, Columbia, MD 21044.

The information contained herein is based on our testing and experience and is offered for the user's consideration, investigation and verification. Since operating and use conditions vary and since we do not control such conditions, we must DISCLAIM ANY WARRANTY, EXPRESSED OR IMPLIED, with regard to results to be obtained from the use of this product. Test methods are available on request. Alltech reserves the right to change prices and/or specifications without prior notification.

ALL-GUARD™, ALLTECH®, ALLTIMA™, DIRECT-CONNECT™, EXPEDITE™, GROM™, JONES CHROMATOGRAPHY™, MODCOL®, PREVAIL™, ROCKET™, and VYDAC® are trademarks of Alltech Associates, Inc. DAVISIL®, DAVISON®, GRACE®, and GRACE DAVISON® are registered trademarks of W. R. Grace & Co.-Conn. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership.



Grace Davison Discovery Sciences Headquarters

2051 Waukegan Road • Deerfield, Illinois 60015-1899

Tel: 847.948.8600 • Fax: 847.948.1078 • Email: discoverysciences@grace.com • Web: www.discoverysciences.com