




# **Agilent OpenLAB Chromatography Data System (CDS)**

**ChemStation Edition Version C.01.04**

## **Release Notes**



**Agilent Technologies**

## Notices

© Agilent Technologies, Inc. 2010-2011, 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

### Manual Part Number

M8301-90130

### Edition

First Edition, May 2012

Agilent Technologies, Inc.

### Software Revision

This guide is valid for Agilent OpenLAB Chromatography Data System (CDS) ChemStation Edition version C.01.04.

Microsoft ® is a U.S. registered trademark of Microsoft Corporation.

This product may be used as a component of an in vitro diagnostic system if the system is registered with the appropriate authorities and complies with the relevant regulations. Otherwise, it is intended only for general laboratory use.

## Warranty

**The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.**

## Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

## Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

## Safety Notices

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

### WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

## **Introduction**

This document describes the functions, features and changes introduced with this and previous releases of the Agilent OpenLAB Chromatography Data System (CDS) ChemStation Edition product.

References to product documentation regarding known issues and workarounds are also provided.

## Version History

This provides a listing of the major feature modifications made in each release of OpenLAB CDS ChemStation Edition.

### Version C.01.04

#### **OpenLAB Data Store A.01.01 Supported with CDS**

OpenLAB Data Store is a central repository where you can safely archive, store, and share electronic files.

OpenLAB CDS ChemStation Edition and Shared Services support the use of OpenLAB Data Store as the storage location for your scientific data. Use of Data Store storage is supported in networked and distributed configurations of OpenLAB CDS and provides several benefits for your system:

- Central Storage with security for small to medium sized labs
- Handles localized Chinese/Japanese content
- Web based data access for search, share, and review
- 21 CFR 11 Compliant e-signatures and audit trail

#### **Master Installer Updates**

- Master Installer now supports installation on top of an existing C.01.03 system resulting in an automatic upgrade.
- Master Installer now supports the installation of either Data Store or ECM for central data storage system components including servers, AICs, and clients.
- Software installation repair is now supported using the repair utility found in the maintenance section of the Master Installer.
- A new additional software installation tool is now available in the maintenance section of the Master Installer. This currently only supports the addition of OpenLAB Data Store components to an AIC or client upgraded from a previous revision.
- Scripted installation of clients and AICs is now supported. An XML file may be exported at the end of installer parameter definition which can then be used to install an identical configuration on other machines.

#### **OpenLAB Control Panel/Shared Services Updates**

- Enhancements to the Licensing User Interface
- Generic projects for Data Store
- Authentication for Data Store
- Data Store Synchronization

#### **NEW html-page for accessing OpenLAB documentation and manuals**

- OpenLAB CDS Documentation is accessed from the START menu > All Programs > Agilent Technologies > OpenLAB CDS Documentation

## Updates to Intelligent Reporting

### Report Template Editor (RTE)

- User-entered report parameters
- Protection of individual report items in report templates
- Double row table
- Document map
- Export report as DOC and TXT

### Method and Sequence Queue

- Single Sample Queuing
- Enhanced options for single run data file naming

### Method Download Options at Instrument startup

- User options for method download on startup of ChemStation

### Data Analysis Improvements

- Calculation of relative retention time, Resolution, Number of theoretical plates, S/N ratio, and Peak-to-valley ratio according to European Pharmacopeia and Japanese Pharmacopeia and reporting of these values using Intelligent Reporting.
- Peak Performance
- Enhancements to Area Summing

## Instrument Control Drivers

### Updated Agilent LC Drivers

- Improved status dashboard for enhanced control and interaction
- 1290 Quat Pump (G4204A)
- 1220 Infinity LC with DAD (G4294B)
- Flexible Cube (G4227A)(Standalone Driver)
- Universal Interface Box II (G1390B)

### Updated Agilent GC Drivers, 78xx Driver Version 4.02, 68xx Driver Version 6.11

- Support for the 7650 Injector with 78xx GCs (G4567A)
- New menu item to extend run for 68xx GC
- GC Tray user interface enhanced (RC.NET only)
- Method Audit Trail and Method Resolution Audit Trail

## Version C.01.03

### Sequence Planner

The **Queue Planner** is used to set up an ordered set of sequences (either ChemStation sequence templates or Easy Sequences) and pauses. The **Queue Planner** allows you to schedule lengthy routine tasks to run overnight or over a weekend.

### Instrument Control Drivers

#### New G3581A 490 micro GC support

#### Updated Agilent LC Driver

- Support for ISET
- RC.NET drivers for single G1364A/B/C and D Automated Fraction Collectors,
- Low-flow Samplers, pumps and valves (G1377A, G1389A, G1376A, G2226, G1162A, G1163A)
- New LC hardware G5664A Bio-inert AFC, 1290 UVD (G1170A),

#### Updated Agilent GC Drivers, 78xx Driver Version 4.02, 68xx Driver Version 6.11

- Updates to Low Thermal Mass II (LTM-II) Rapid Heating/Cooling System for 7890A GC. This release integrates complete control of the Low Thermal Mass II (LTM II) system into the 7890A GC and Chemstation. The added support of ID tags on the LTM II columns allows for automatic configurations loaded in the GC and Chemstation and the ability to track number of runs, over temperatures, and excursions.

## Version C.01.02

### Agilent OpenLAB ECM connectivity

With the C.01.02 release of OpenLAB CDS ChemStation Edition, ECM is supported as a data storage location. With this feature you get the following enhancements for your data system:

Store and retrieve data directly in the Agilent OpenLAB Enterprise Content Management (ECM) system from the CDS with ECM as a valid storage location and authentication provider.

This provides the following benefits:

- Keyword searches for data, methods and sequences
- Web based data access for sharing and review
- Automated data archival
- Result package version control
- Enhanced security
- Consolidate CDS data with all laboratory content
- Thin client configuration support

### Master Installer Updates

Automated uninstall of all OpenLAB CDS components is now available from the master installer.

### New Topology: Distributed Systems (Networked AIC and Clients)

### Agilent Support information

### Support for LC/MS and CE/MS on Workstation and Networked Workstation

### Updated Agilent GC Drivers, 78xx Driver Version 4.01, 68xx Driver Version 6.10

- Barcode support with the 68xx Software Driver version 6.10 has been updated to transfer the vial from the barcode reader directly to the turret for injection. The intermediate stop back to the tray location has been eliminated.
- Instrument parameter actuals are now available next to the set points in method edit user interface
- 7890A: Digital Signal Switching, Signal Zero, Signal Freeze/Resume available in the Signal Event Table (7890A requires a minimum Firmware revision of A.01.12.)
- Hide and display the instrument actuals view and graphic plot
- 7890A: Shutdown events such as Hydrogen shutdown now included in Software Status (7890A requires a minimum Firmware revision of A.01.12.)
- Column User interface includes Constant Flow, Constant Pressure, Ramped Pressure, and Constant Flow for the Control Mode selections
- G4514A 7693A tray support for 7820A GC
- Start run from GC by pressing the start button now available for 68xx GC

## Version C.01.01

### Intelligent Reporting

OpenLAB CDS provides a new Intelligent Reporting feature for GC and LC ChemStation data. The new intelligent reporting is provided in addition to classic reporting. The new Intelligent Reporting allows you to create state-of-the-art reports in an industry standard reporting format. You can easily create and modify report templates with the new Report Template Editor built into the system that has an easy to use WYSIWYG interface. Report templates are saved in the standardized Report Definition Language (RDL) format that is also used by Microsoft Business Intelligence Studio.

When Intelligent Reporting is activated, a new Review view is available. It allows you to apply any report template to any combination of data files.

### New licensing

With OpenLAB CDS ChemStation Edition, a new licensing strategy is introduced that helps you use your licenses more effectively. All instrument control, driver and Add-on licenses are floating licenses. Any instrument that starts up requests licenses from License Management, and when the instrument is closed, it returns the licenses. The new licensing is based on Flexera Software™ SubScribeNet® technology. It simplifies your software license management and license tracking and allows you to receive new or urgent software and patch download information immediately after release.

### OpenLAB Control Panel / OpenLAB Shared Services

OpenLAB Shared Services allow the management of users, instruments and licenses. OpenLAB Shared Services is managed using the OpenLAB Control Panel.

- **Instrument Management:** OpenLAB Control Panel allows you to set up and configure instruments. Instruments are launched from the OpenLAB Control Panel. In addition, the Control Panel provides a lab-at-a-glance overview of the status of the instruments on a workstation or networked workstation.
- **User Management:** OpenLAB Control Panel allows you to manage users, groups, roles, and privileges. If you manage your users within an ECM system or a Windows domain, you can map those existing users into OpenLAB Shared Services.
- **License Management:** This service includes the administration of licenses that are required for your instrument modules and Add-Ons. You can add or remove licenses and view the status of all licenses.

### New Data Analysis Workflows

OpenLAB ChemStation C.01.01 provides two Data Analysis modes: the Recalculation mode and the Reprocessing mode. The reprocessing mode allows you to reprocess samples in the context of a sequence (e.g. for a bracketed calibration). The recalculation mode allows you to quickly recalculate a sample or a set of samples with a different method. The required functions associated with each mode are available in separate toolbars.



**Result sets:** Containers are now called result sets. Methods that are used in result sets may be modified during acquisition. You can create your own result sets using any collection of existing data. You can use the user-assembled result sets, for example, for cross sample calculations.

**Method Handling:** It is possible to directly load master methods and sequence methods from the ChemStation explorer. The Update Methods dialog box allows you to synchronize master methods and sequence methods.

## Known Issues

For known issues and workarounds in the OpenLAB CDS ChemStation software at the time of release see the Software Status Bulletin file named *83xxAA-SSB.html* located on Disc 2 of the OpenLAB CDS Installation media.

You may also visit [www.Agilent.com](http://www.Agilent.com) for up to date issue information for all Agilent software products:

<http://www.chem.agilent.com/en-US/Support/Downloads/Patches/Pages/default.aspx>