


Parallel GC for Complete RGA Analysis

Application Brief



Chunxiao Wang


A previous application brief [1] has shown that a 7890A GC configured with three parallel channels provides a complete refinery gas analysis (RGA) within six minutes. The configuration for fast RGA in the brief has been updated by adding a fifth valve, which can now be supported by the 7890A GC. The updated configuration is almost the same as the previous one except for the third channel (TCD) for H₂ analysis using N₂ or Ar as carrier gas to improve H₂ detectability and linearity. The updated configuration uses a 10-port valve with a pre-column for backflushing late-eluting components while H₂ is separating on the molsieve column instead of a three-way splitter plus split/splitless inlet.

Refinery gases are mixtures of various gas streams produced in refinery processes. They can be used as a fuel gas, a final product, or a feedstock for further processing. The composition of refinery gas streams is very complex, typically containing hydrocarbons, permanent gases, sulfur compounds, etc. An exact and fast analysis of the components is essential for optimizing refinery processes and controlling product quality.

The Agilent 7890A GC now supports an optional detector (TCD), allowing simultaneous detection across three channels. This provides a complete analysis of permanent gases, including nitrogen, hydrogen, oxygen, carbon monoxide,

Highlights

- One 7890A GC configured with three parallel channels with simultaneous detection provides a comprehensive, fast, and high-resolution analysis of refinery gas in 6 minutes.
- Use of optimized columns allows faster analysis of hydrocarbons and permanent gases using a single oven temperature program without the need for an additional column oven.
- A third TCD channel can be used for improving hydrogen detection and linearity by using nitrogen (or argon) as carrier gas.
- A new, easy-to-use union tubing connector based on capillary flow technology is used to connect valves and capillary columns to improve the chromatographic performance, including peak shape.
- Excellent results are achieved. The lowest detection limit is 50 ppm for all compounds, 500 ppm for hydrogen sulfide.
- ChemStation macro program is supplied for RGA reporting.
- The system can be obtained by ordering option SP1 7890-0322 for the standard fast RGA and 7890-0338 for the fast RGA with Hastelloy valves and nickel tubing for H₂S containing samples on the 7890A.



Column 1 HayeSep Q 80/100 mesh
 Column 2 HayeSep Q 80/100 mesh
 Column 3 Molsieve 5A 60/80 mesh
 Column 4 HayeSep Q 80/100 mesh

Column 5 Molsieve 5A 60/80 mesh
 Column 6 DB-1
 Column 7 HP-PL0T Al₂O₃
 PCM: Electronic pneumatics control (EPC) module

Figure1. RGA valve system.



Agilent Technologies

carbon dioxide, and hydrocarbons to nC6. The total run time is less than 6 minutes. The configuration is suitable for most refinery gas streams such as atmospheric overhead, FCC overhead, fuel gas, and recycle gases.

In this analysis, a single Agilent 7890A GC is configured with three channels, including an FID channel and 2 TCD channels. Light hydrocarbons are determined on the FID channel using an alumina column. One TCD is used with nitrogen or argon carrier gas for improved determination of hydrogen and helium; the other TCD is used with helium carrier for the detection of all other required permanent gases. The configuration is shown in Figure 1. An Agilent union tube connector, based on capillary flow technology, is used to quickly and easily connect the valve and capillary column for improved performance. The system conforms to published methods such as ASTM D1945 [2], D1946 [3], and UOP 539 [4].

Separation resulting from each channel is illustrated in Figure 2. The top chromatogram shows the hydrocarbon analysis. A PLOT AL203 column provides excellent separation of hydrocarbons from C1 to nC5 containing 22 isomers. Components heavier than nC6 are backflushed early in the run as a group (C6+) through a short DB-1 pre-column. The middle chromatogram shows the separation of permanent gases using helium as the carrier gas on the second TCD channel (B TCD). H₂S and COS can be analyzed on the second TCD channel as well, requiring 3 to 4 additional minutes. The bottom chromatogram shows the

separation of hydrogen. Because hydrogen has only a small difference in thermal conductivity compared to helium, it requires an additional TCD with nitrogen or argon as the carrier gas to improve the hydrogen detectability and linearity. All channels operate simultaneously to provide a comprehensive, fast analysis with high resolution of components. A macro program automatically provides the calculation of gas properties. Reports can be generated using formulas specified in the ASTM/GPA and/or ISO standards. Reports in mole%, weight%, volume%, or any combination of the three are available.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.




Figure 2. Refinery gas calibration standards analysis.

Reference

1. Chunxiao Wang, "Parallel GC for Complete RGA Analysis," Agilent application brief, 5989-6103EN, January 19, 2007
2. ASTM D1945-03, "Standard Test Method for Analysis of Natural Gas by Gas Chromatography," ASTM International, 100 Bar Harbor Drive, West Conshohocken, PA 19428 USA.
3. ASTM D1946-90 (2006), "Standard Practice for Analysis of Reformed Gas by Gas Chromatography," ASTM International, 100 Bar Harbor Drive, West Conshohocken, PA 19428 USA.
4. UOP Method 539, "Refinery Gas Analysis by Gas Chromatography," ASTM International, 100 Bar Harbor Drive, West Conshohocken, PA 19428, USA.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2008

Published in the USA
December 12, 2008
5989-7438EN