

# PLATINblue UHPLC/HPLC Plus systems

▶ What limits your LC today?



# Flexibility

More than just UHPLC

## Optimized

Eluent supply ✓  
Eluent mixing ✓  
Injection ✓  
Column ✓  
Detection ✓  
Data processing ✓

The PLATINblue family of liquid chromatography instruments was designed to offer superior results not only in terms of resolution and throughput, but also in terms of reliability, operating costs, and ease-of-use. PLATINblue systems are capable of running under standard HPLC conditions, fast LC conditions, and UHPLC up to 1000 bar.

## Freedom of choice

Depending on your lab's current application focus we offer two versions of PLATINblue premium liquid chromatography systems:


- ▶ **PLATINblue UHPLC**
  - Ultimate performance with sub-2  $\mu\text{m}$  fully porous and core-shell columns
  - Up to 1000 bar
  - Highest sensitivity
  - Ideal for MS
- ▶ **PLATINblue HPLC Plus**
  - Routine HPLC up to 750 bar
  - Upgradeable to UHPLC at any time
  - Most affordable PLATINblue solution

## PLATINblue features

- ▶ Modular system architecture
- ▶ Eluent supply with integrated degasser
- ▶ Optimized system volume for high resolution
- ▶ High speed autosampler
- ▶ Fast detectors, up to 100 Hz data rate
- ▶ PDA available with high sensitivity flow cell
- ▶ User-friendly design and operation
- ▶ Supported by software packages: ChromGate® (incl. OQ/PQ), EZChrom Elite, and Chromeleon®, as well as MS software Xcalibur™ and Analyst™

## Gradient flexibility

- ▶ Binary high pressure gradient for ultra fast composition changes
- ▶ Quaternary low pressure gradient (standard with HPLC Plus)



# Productivity

## High resolution and throughput

Increase your laboratory's productivity with PLATINblue UHPLC/HPLC Plus and achieve:

- ▶ higher resolution and unmatched separation efficiency at equivalent column lengths
- ▶ faster separations, thanks to the use of shorter columns with equal or higher separation performance

### High resolution

- ▶ PLATINblue UHPLC offers ultimate performance with an extended pressure range up to 1000 bar (15000 psi) at flow rates of up to 2 ml/min and of 800 bar (11600 psi) up to 5 ml/min
- ▶ Optimized system volume

### Short cycle time



- ▶ fast gradients with up to 3% composition change rate per second
- ▶ high-performance autosampler for fast injection (15 s or < 60 s with needle wash)
- ▶ optimized communication and software operations for short initialization times

### High throughput

- ▶ process up to 768 samples in 2 x 384-well plates or up to 96 samples in standard 1.8 ml vials

### Fast detection

- ▶ data acquisition rates up to 100 Hz provide higher resolution and peak capacity, as well as a more accurate integration



### Tip

The PLATINblue system includes eluent filters for protection against particulate matter. When working with columns packed with sub-2µm particle sizes, samples should also always be filtered through a 0.2µm membrane.

# Precision

## Sensitivity and reproducibility

High-quality PLATINblue components were designed to work together perfectly to enable accurate measurement with optimal reproducibility – a prerequisite for analysis results you can rely on.

The sharp and narrow peaks produced by fast chromatography in conjunction with low-noise and low-drift PDA detectors ensure sensitive measurements.


### Ultimate detection sensitivity

- ▶ electronically controlled pump drive provides for low pulsation eluent flow and a stable baseline
- ▶ high light intensity over the full UV-Vis range through the combination of a long-life “high brightness” D<sub>2</sub> with a tungsten-halogen lamp
- ▶ optimized flow cells with path lengths up to 50 mm and volumes as low as 2.4 µl allow for most sensitive detection of narrow peaks
- ▶ low noise (5 µAU) and drift specifications
- ▶ excellent linearity of up to 3 AU for precise simultaneous determination of main substances and trace components

### Reproducibility

Precision can only be as good as the weakest link in the chain. All critical components of the PLATINblue systems were optimized for the demanding requirements of fast LC and UHPLC, resulting in systems with excellent reproducibility. The advanced column thermostat T-1 adds to overall system precision by keeping column conditions stable and cooling the effluent before detection.

Excellent reproducibility of overall system




Six successive separations using online mixing of a ternary eluent  
Eluent: 3% CH<sub>3</sub>COOH/MeOH/ACN = 78/20/2  
Flow rate: 1 ml/min, Temp.: 30 °C, Wavelength: 225 nm

#### Retention time

< 0.3% RSD

#### Peak areas

< 1.1% RSD with partial loopfill injection  
(< 0.3% RSD with full loop injection)



### Tip

A very high injection peak can be disadvantageous in fast chromatography since the first sample peaks elute very quickly after the injection. Dissolve the sample in eluent whenever possible to minimize this disturbance.



# Ease of use

Simple to operate and to service

To stay focused on what is really important, the PLATINblue ultra high-performance LC system is equipped with many useful features that make operation more convenient, relieve you from time-consuming routine checks and assist you in your daily work.

**Tip**

The lifespan of the detector lamps can be conserved by having one or both turned off during interruptions in analysis and switched back on automatically when resuming measurements.

## Comfortable operation

- ▶ intuitive touchscreen control\*
- ▶ automatic configuration of communication with PC and software eases setup
- ▶ automatic selection of optimal detector integration time
- ▶ standby mode

## Easy access for quick maintenance


- ▶ solid glass doors provide convenient access to system components
- ▶ removable side panel for easy lamp exchange and maintenance

## Reliable and GLP compliant

- ▶ fast Ethernet communication
- ▶ automatic pump head detection (RFID)
- ▶ column thermostat: RFID recognition for up to six columns simultaneously
- ▶ automatic lamp detection
- ▶ automatic wavelength verification
- ▶ comprehensive GLP reporting functions

## Safety and protection

- ▶ system functions are monitored, minimum pressure shut-off and eluent level control feature in ChromGate® software provide protection against leaks and running dry
- ▶ soft-start flow for protection of columns and system components
- ▶ Leak detection



Simple lamp replacement thanks to detachable side panel

# PLATINblue

## Component overview



▶ **PLATINblue UHPLC**  
HPG Pump: P-1/P-1 A60013/14

Ultra high pressure pump optimized for fast liquid chromatography, capable of pressures up to 1000 bar (15000 psi) and flow rates up to 5 ml/min. The electronically controlled pump drive provides for low pulsation, a highly precise flow rate and a stable baseline.

<b>Flow rate</b>	0.01 – 5.0 ml/min or 0.001 ml increments
<b>Flow rate accuracy</b>	±1 %
<b>Flow rate precision</b>	< 0.1 %
<b>Pressure range</b>	
UHPLC:	1000 bar up to 2 ml/min, 800 bar up to 5 ml/min
<b>Eluent mixing</b>	SmartMix UHPLC mixing technology
<b>binary HPG:</b>	2 x Pump P-1, integrated 2 channel degasser and SmartMix
<b>Features</b>	touchscreen, active pressure and pulsation compensation, RFID automatic pump head detection (GLP), interchangeable pump heads

also available as:

▶ **PLATINblue HPLC Plus**  
HPG Pump P-1/P-1 A60015/16

Gradient pump configuration for pressures up to 750 bar



▶ **PLATINblue HPLC Plus**  
LPG Pump: P-1/M-1 A60015/A60513

High pressure pump optimized for fast and routine HPLC, capable of pressures up to 750 bar (10875 psi). The electronically controlled pump drive provides for low pulsation and a highly precise flow rate. This most affordable PLATINblue solution can be upgraded to UHPLC at any time.

<b>Flow rate</b>	0.01 – 10.0 ml/min 0.001 ml increments
<b>Flow rate accuracy</b>	±1 %
<b>Flow rate precision</b>	< 0.1 %
<b>Pressure range</b>	
HPLC Plus:	750 bar up to 5 ml/min, 400 bar up to 10 ml/min
<b>Eluent mixing</b>	SmartMix HPLC mixing technology
<b>quaternary LPG:</b>	1 x Pump P-1 incl. SmartMix, Modular Eluent Manager M-1, 4 channel degasser and interface

▶ **PLATINblue M-1**  
Modular Eluent Manager A60513

<b>Quaternary LPG</b>	blending of up to 4 eluents
<b>Degasser module</b>	4 channels, Teflon® AF, high efficiency
<b>A/D-D/A interface</b>	24-bit, 4 analog inputs, 4 analog outputs, several digital I/O



▶ **PLATINblue PDA-1**  
Diode Array Detector A62031

This extremely sensitive diode array detector with a data acquisition rate up to 100Hz with dual-lamp configuration has a wavelength range of 190–1000nm. New ultra sensitive fiber optical flow cell provides for highest detection performance. Use spectrum data to determine peak purity or to facilitate the identification of unknown substances.

<b>Lamps</b>	high-brightness D <sub>2</sub> , tungsten-halogen
<b>Wavelength range</b>	190–1000 nm
<b>Wavelength accuracy</b>	< 1 nm
<b>Max. data acquisition rate</b>	100 Hz
<b>Diodes</b>	1024
<b>Pixel pitch</b>	< 1 nm
<b>Channels</b>	max. 6
<b>Noise*</b>	±5 µAU
<b>Drift*</b>	< 300 µAU/h
<b>Linearity</b>	0 – 2 AU
<b>Spectral bandwidth</b>	< 2.5 nm
<b>Special features</b>	touchscreen, exchangeable flow cells
<b>Standard flow cell</b>	sensitive, fiber optical technology, fused silica, 2.4 µl, 10 mm, 100 bar
<b>Optional flow cell</b>	total reflection technology, up to 5 times higher sensitivity than standard flow cell, fiber optics, 10.0 µl, 50 mm, 100 bar

\* according to ASTM




▶ **Fluorescence Detector**  
A59201

With a signal-to-noise ratio of 2000 for the water Raman band, this fluorescence detector is well-suited for trace-level analysis. It features high acquisition speed and a quick response time of just 10 ms both enabling detection of sharpest peaks in ultra-fast analysis.

<b>Lamps</b>	xenon lamp, mercury lamp (wavelength check)
<b>Wavelength range</b>	200–750 nm
<b>Wavelength accuracy</b>	± 2 nm
<b>Wavelength precision</b>	± 0.2 nm
<b>Max. data acquisition rate</b>	50 Hz
<b>Channels</b>	max. 2
<b>Signal/Noise</b>	2000 H <sub>2</sub> O Raman peak
<b>Spectral bandwidth</b>	< 20 nm
<b>Standard flow cell</b>	12 µl, 20 bar

## Tip

The column thermostat PLATINblue T-1 is also available without switching valves.



### ▶ PLATINblue AS-1 Autosampler

A63500

Inject samples from well plates and standard vials accurately and fast with this UHPLC autosampler, which can handle pressures up to 1000 bar (15000 psi). Supports full and partial loopfill injection as well as "microliter-pickup" for flexible selection of injection volumes. The tray cooling option allows for analysis of temperature-sensitive samples.

UHPLC injection	up to 1000 bar (15000 psi)
Sample capacity	max. 768 sample positions with well plates or 96 sample positions with standard autosampler vials
Injection volume	1–100 µl
Standard fitted loop	10 µl
Injection cycle time	15 s, < 60 s with wash
Injection modes	full loop, partial loopfill and "microliter-pickup"
Precision	RSD < 0.3 % full loop injection
Carry over	< 0.05 % with needle wash
Sample tray cooling	4–22 °C



### ▶ PLATINblue T-1 Column Thermostat

A63410/A63412\*

This column thermostat can accommodate up to six HPLC columns and ensures stable temperature conditions in the range of 5 to 80 °C. The eluent is cooled before reaching the detector to minimize noise and drift. It also features RFID column recognition for up to six columns and has two valves for advanced switching tasks.

Temperature range	5 – 80 °C
Temperature stability	0.1 °C
Postcolumn cooling	yes, 15–35 °C
Switching valves	max 2 (available types: 6, 8, 10 port multi position or 6 and 10 port 2-position)
Column compartment	370 mm x 85 mm (L x W)
Column holder	magnetic
Column recognition:	RFID for up to six columns**

\* A63412 version without touchscreen and valves

\*\* Non-KNAUER columns can be fitted with a custom-programmed KNAUER RFID tag. Ask for details.

### ▶ Manual injection valve

A64601

UHPLC injection up to 1000 bar (15000 psi)  
Automatic switching valves and a manual injection valve are also available

# Cleverness

## User-friendly, scalable and compliant

ChromGate® PLATINblue Edition is our easy-to-use software solution for laboratories with high demands for functionality, flexibility and scalability, as well as GLP compliance. Based on EZChrom Elite™ software, ChromGate® PLATINblue Edition software adds extended KNAUER-specific functions.

### Smart product features that save time

- ▶ automatic system configuration for the PLATINblue system via the instrument network
- ▶ navigation with tree structure for quick orientation

Only ONE software for your chromatography laboratory is required to control all KNAUER instruments and a wide range of chromatography instrumentation from other manufacturers. Your lab benefits from reduced training expenses, uniform operability and harmonized data reporting.

### Flexible reporting functions

- ▶ report layouts are easily customized via multiple options in Custom Report feature
- ▶ predefined standard reports for quick output of results

### Scalability and compliance

- ▶ scalable from a single instrument to a multi-user enterprise-wide solution thanks to Client/Server functionality
- ▶ FDA 21 CFR Part 11 compliance

### Easy to use

- ▶ simple method creation via software "wizards"
- ▶ "sequence wizard" provides assistance and saves time when programming a sequence of analyses

## Tip ChromGate® PLATINblue Edition Software

Our full-featured scalable chromatography software includes one client/server licence (additional licences optional). Every complete PLATINblue system also includes a major brand PC with the latest technology, large TFT display and 8-port Ethernet/LAN router.

### Convenient diode-array features

- ▶ clear presentation of 3-D data
- ▶ easy creation of spectrum libraries and fast spectrum search

### NEW

- ▶ PLATINblue systems can also be controlled by Dionex Chromeleon™ software.



# Ultimate

## Front-end for LC/MS

PLATINblue UHPLC in combination with the MSQ Plus™ mass spectrometer offers the most sensitive and compact LC/MS solution for high throughput applications.



### Seamless integration

- ▶ when operating the PLATINblue UHPLC system as part of an LC/MS solution, Xcalibur™ takes over complete LC instrument control and mass spectrometer data handling
- ▶ splitless sample introduction at flow rates up to 2 ml/min for maximum transfer of analytes

### High performance

- ▶ highly sensitive ion optics
- ▶ 60° M-Path™ source for high ruggedness and lowest carryover
- ▶ square quadrupole RF/DC pre-filter lens protects from contamination and enables high-efficiency transmission
- ▶ fast scan speed up to 12000 amu/s for full compatibility with narrow UHPLC peaks

### Work with ease

- ▶ the MSQ Plus™ detector is very robust and easy to handle
- ▶ exchanging the probes between ESI and APCI mode takes only minutes
- ▶ tool-free maintenance of ionization source

### Versatile

- ▶ available with two ionization modes: electrospray ionization (ESI) or atmospheric pressure chemical ionization (APCI)
- ▶ ESI is a very mild method, even for thermally labile species. It is suitable for polar compounds or ionic substances in solution, e.g. anions/cations, peptides, small nucleotides, drugs, or sugars (< 2000 Da)
- ▶ APCI requires higher vaporization temperatures, but is less dependent on solvent type and solvent flow rate. Typical detectable compounds include hydrocarbons, alcohols, aldehydes, ketones, esters – excluding acidic or basic groups (< 1000 Da)



Detail of MSQ Plus with atmospheric pressure chemical ionization (APCI) installed



Detail of MSQ Plus with electrospray ionization (ESI) installed



### ▶ MS detector specifications

<b>Ionization modes</b>	electrospray (ESI) or atmospheric pressure chemical ionization (APCI)
<b>Mass range</b>	17 – 2000 Da
<b>Sensitivity*</b>	
<b>ESI, positive ion</b>	1000:1 RMS, 50 pg erythromycin (10 µl x 5 pg/µl)
<b>ESI, negative ion</b>	500:1 RMS, 20 pg p-nitrophenol (10 µl x 2 pg/µl)
<b>APCI, positive ion</b>	200:1 RMS, 50 pg erythromycin (10 µl x 5 pg/µl)
<b>APCI, negative ion</b>	50:1 RMS, 20 pg p-nitrophenol (10 µl x 2 pg/µl)

\*signal-to-noise, loop injection, 1 ml/min flow rate

### ▶ Xcalibur™ MS Data System

Xcalibur is a very versatile and easy-to-use mass spectrometry data system, offering transparent navigation through the process of instrument setup, sequence setup, and data acquisition. MS data analysis is simplified using any of the three browser environments: Qual Browser, Quan Browser, and Library Browser. Intuitively create reports with XReport software.

The Xcalibur MS data system provides for control and data acquisition of all PLATINblue instruments and the MSQ Plus MS detector, as well as the complete range of Thermo MS detectors.

M-Path, MSQ Plus, and Xcalibur are trademarks of Thermo Fisher Scientific Inc.



# Versatility


Wide range of application areas

## Pharmaceutical

Separation of paracetamol and common synthesis impurities



<b>Separation column</b>	
BlueOrchid C18 1.8µm, 100 x 2mm	
<b>Separation conditions</b>	
Eluent:	A: 1.7 mM NaH <sub>2</sub> PO <sub>4</sub> , pH 3.7 B: ACN
Gradient:	0.0–0.3 min 13% B 0.3–2 min 13–70% B 2–2.5 min 70% B
Flow rate:	0.85 ml/min
Mode:	RP-Mode, gradient
Injection volume:	1 µl
Pressure:	950 bar (13800 psi)
Detection:	MW-1, 254 nm (100Hz, 0.005 s) 10 mm, 2 µl flow cell
Temperature:	50°C




- 1 4-Aminophenol
- 2 Paracetamol
- 3 N-(4-hydroxyphenyl)propanamide
- 4 N-(2-hydroxyphenyl)acetamide
- 5 N-(3-chloro-4-hydroxyphenyl)acetamide
- 6 1-(4-Hydroxyphenyl)ethanoneoxime
- 7 N-phenylacetamide
- 8 4-(Acetylamino)phenyl acetate
- 9 4-Nitrophenol
- 10 Chloroacetanilide
- 11 1-(2-Hydroxyphenyl)ethanone

## Pharmaceutical

Separation of steroids




<b>Separation column</b>	
BlueOrchid C18 1.8µm, 50 x 2mm	
<b>Separation conditions</b>	
Eluent:	A: H <sub>2</sub> O B: ACN
Gradient:	0–1.5 min 35%–75% B 1.5–2 min 75% B
Flow rate:	1.0 ml/min
Mode:	RP-Mode, gradient
Injection volume:	0.5 µl
Pressure:	650 bar (9430 psi)
Detection:	PDA-1, 254 nm (100Hz, 0.005 s) 10 mm, 2 µl flow cell
Temperature:	30°C




- 1 Cortisone
- 2 Corticosterone
- 3 Testosterone
- 4 Deoxycorticosterone
- 5 Norgestrel
- 6 Progesterone

## Pharmaceutical

Separation of Sudan dyes



<b>Separation column</b>	
BlueOrchid C18 A 1.8µm, 50 x 2mm	
<b>Separation conditions</b>	
Eluent:	A: 1.7 mM NaH <sub>2</sub> PO <sub>4</sub> , pH 4 B: ACN
Gradient:	0–0.8 min 8% B 0.8–1.8 min 8%–40% B
Flow rate:	1.4 ml/min
Mode:	RP-Mode, gradient
Injection volume:	1 µl
Pressure:	750 bar (10880 psi)
Detection:	PDA-1, 265 nm (100Hz, 0.005 s) 10 mm, 2 µl flow cell
Temperature:	40°C




- 1 Sulfadiazine
- 2 Sulfamethiazole
- 3 Sulfamethoxypyridazine
- 4 Sulfamethoxazole
- 5 Sulfamethoxine


PLATINblue systems can be flexibly used to increase the throughput and sensitivity of many applications. Benefit from the advantages of high-speed UHPLC applications without having to completely convert your established HPLC methods. The following are only a few examples of the applications possible with PLATINblue.

## Bioanalytical

Separation of nucleosides and bases



<b>Separation column</b>	
BlueOrchid PFP 1.8µm, 100 x 2mm	
<b>Separation conditions</b>	
Eluent:	20mM NH <sub>4</sub> Ac, pH 3.5 10% A/MeOH 90:10 (v/v)
Flow rate:	0.5 ml/min
Mode:	RP-Mode, isocratic
Injection volume:	1 µl
Pressure:	570 bar (8270 psi)
Detection:	PDA-1, 210 nm (100 Hz, 0.001 s) 10 mm, 2 µl flow cell
Temperature:	30°C




- 1 Cytosine
- 2 Uracil
- 3 Cytidine
- 4 Uridine
- 5 Guanine
- 6 Thymine
- 7 Adenine
- 8 Guanosine
- 9 Thymidine
- 10 Adenosine

## Environmental

Determination of pesticides I



<b>Separation column</b>	
BlueOrchid C18 1.8µm, 50 x 2mm	
<b>Separation conditions</b>	
Eluent:	A: H <sub>2</sub> O B: ACN
Gradient:	0–0.5 min 35% B 0.5–1.5 min 35%–100% B 1.5–2.0 min 100% B
Flow rate:	0.6 ml/min
Mode:	RP-Mode, gradient
Injection volume:	1 µl
Pressure:	300 bar (4350 psi)
Detection:	215 nm (50Hz, 0.05 s) 10 mm, 2 µl flow cell
Temperature:	40°C




- 1 Chloridazon
- 2 Metoxuron
- 3 Monuron
- 4 Diuron
- 5 Propazine
- 6 Linuron
- 7 Metolachlor
- 8 Parathion-ethyl
- 9 β-Endosulfan

## Environmental

Determination of 16 EPA PAH



<b>Separation column</b>	
BlueOrchid PAH, 50 x 2mm	
<b>Separation conditions</b>	
Eluent:	A: MeOH/H <sub>2</sub> O 75:25 B: ACN
Flow rate:	1.0 ml/min
Mode:	RP-Mode, gradient
Injection volume:	1 µl EPA standard
Pressure:	< 500 bar (< 7250 psi)
Detection:	254 nm (100Hz, 0.01 s) 10 mm, 2 µl flow cell
Temperature:	25°C



- |                      |                           |
|----------------------|---------------------------|
| 1 Naphthalene        | 11 Benzo(b)fluoranthene   |
| 2 Acenaphthalene     | 12 Benzo(k)fluoranthene   |
| 3 Acenaphthene       | 13 Benzo(a)pyrene         |
| 4 Fluorene           | 14 Dibenzo(a,h)anthracene |
| 5 Phenanthrene       | 15 Benzo(g,h,i)perylene   |
| 6 Anthracene         | 16 Indeno(1,2,3-cd)pyrene |
| 7 Fluoranthene       |                           |
| 8 Pyrene             |                           |
| 9 Benzo(a)anthracene |                           |
| 10 Chrysene          |                           |

## KNAUER UHPLC columns


KNAUER offers three different column product lines, specifically designed for ultra high-performance liquid chromatography applications.

**BlueShell<sup>®</sup>** columns feature particles which consist of a solid core and a porous shell. The 2.6 µm diameter particles enable high speed and high resolution separations without undesirable high backpressure. The reduced depth of the outer porous layer limits the diffusional path of the analytes, leading to minimized mass transfer resistance and minimized peak broadening. BlueShell columns can achieve up to 200.000 plates/meter. Available as a C18, C18A, and HILIC phase.

**BlueOrchid** columns with fully porous 1.8 µm particle are designed to provide high efficiency and stability, making them the ideal choice for demanding applications. They are available in 12 modifications for a wide range of applications.

**Bluespher<sup>®</sup>** columns were developed to feature very similar characteristics and selectivities to the KNAUER Eurospher I and Eurospher II families to make method transfer between UHPLC, HPLC, and even MPLC easy. Based on 2 µm fully porous particles, the manufacturing process employs mono- and multi-functional silanes to achieve a maximum in loading capacity and chemical stability. All application fields of reversed phase chromatography are covered.

For more information, please refer to the detailed column brochures available on [www.knauer.net](http://www.knauer.net)



### Food


Quick determination of tocopherols



**Separation column**  
BlueOrchid Si 1.8µm, 50 x 2mm

#### Separation conditions

Eluent: Heptane/2-butanol 1000:4 (v/v)  
Flow rate: 1.7 ml/min  
Mode: NP-Mode, isocratic  
Injection volume: 1 µl  
Pressure: 350 bar (5080 psi)  
Detection: MW-1, 280 nm (50 Hz, 0.01 s)  
10 mm, 2 µl flow cell  
Temperature: 25 °C



- 1 Trans-retinol
- 2 Alpha-tocopherol
- 3 Beta-tocopherol
- 4 Gamma-tocopherol
- 5 Delta-tocopherol
- 6 Vitamin D<sub>2</sub>

### Food


Determination of water soluble vitamins



**Separation column**  
BlueOrchid C18 A 1.8µm, 100 x 2mm

#### Separation conditions

Eluent: A: 50 mM NaH<sub>2</sub>PO<sub>4</sub> pH 4  
B: MeOH  
Mode: RP-Mode, gradient  
Flow rate: 0.7 ml/min  
Injection volume: 1 µl  
Pressure: 800 bar (11600 psi)  
Detection: PDA-1, 254 nm (50 Hz, 0.01 s)  
10 mm, 2.4 µl flow cell  
Temperature: 30 °C



- 1 Ascorbic acid
- 2 Thiamine
- 3 Nicotinic acid
- 4 Pyridoxal
- 5 Pyridoxin
- 6 Nicotinamid
- 7 Folic acid
- 8 Cyanocobalamin
- 9 Riboflavin

### Food


Separation of Sudan dyes



**Separation column**  
BlueOrchid C18 1.8µm, 50 x 2mm

#### Separation conditions

Eluent: A: H<sub>2</sub>O (+ 0.1 % formic acid)  
B: ACN (+ 0.1 % formic acid)  
Gradient: 0–0.5 min 70% B  
0.5–1.5 min 100% B  
1.5–2.0 min 100% B  
Mode: RP-Mode, gradient  
Flow rate: 1.0 ml/min  
Injection volume: 1 µl  
Pressure: 400 bar (5800 psi)  
Detection: PDA-1, 240 nm (50 Hz, 0.01 s)  
10 mm, 2.4 µl flow cell  
Temperature: 25 °C



- 1 Impurity 1
- 2 4-Phenylazophenol
- 3 α-Naphthyl red
- 4 Impurity 2
- 5 Impurity 3
- 6 Sudan II
- 7 Sudan III
- 8 Sudan IV

### Clinical


Determination of sedatives, e.g. barbiturates



**Separation column**  
BlueOrchid C18 1.8µm, 50 x 2mm

#### Separation conditions

Eluent: MeOH/H<sub>2</sub>O 60:40 (v/v)  
Flow rate: 0.5 ml/min  
Mode: isocratic  
Injection volume: 1 µl  
Pressure: 826 bar (11980 psi)  
Detection: MW-1, 254 nm (80 Hz, 0.005 s)  
Temperature: 40 °C



- 1 barbital
- 2 luminal
- 3 prominal
- 4 revonal

### Clinical


Separation of tt-muconic acid in urine (I)



**Separation column**  
BlueOrchid C18 A 1.8µm, 100 x 2mm  
Phenyl 1.8µm 100 x 2mm (column tandem)

#### Separation conditions

Eluent: A: H<sub>2</sub>O (1% Hac) B: MeOH  
Gradient: 0–7.5 min 5% B–25% B  
7.5–7.6 min 25% B–90% B  
7.6–9.5 min 90% B  
9.5–9.6 min 90% B–5% B  
Flow rate: 0.4 ml/min  
Mode: RP-Mode, gradient  
Injection volume: 10 µl  
Pressure: 590 bar (8560 psi)  
Detection: PDA-1, 259 nm (50 Hz) 10 mm,  
2 µl flow cell;  
Temperature: 40 °C



- 1 tt-Muconic acid

# Support

## Fast and reliable

Excellent support starts at the product development stage. That's why we developed PLATINblue instruments with ultimate reliability in mind. Each component was designed for maximum dependability and durability, using only the best materials available.



**All support requests will be responded** to by local authorized KNAUER support staff. Guaranteed! A detailed report of your system's status can be easily generated by the PLATINblue ServiceTool software for fast problem diagnosis. KNAUER technical support will contact you within 24 hours with a detailed plan for resolution of the problem. If necessary, we will solve the problem at your site – quickly and without complication.

To keep all of the components in your PLATINblue system running with the highest performance, we offer a comprehensive support and maintenance program to ensure your system's reliable operation.



**We offer you a three-year warranty** on every PLATINblue system in conjunction with a service agreement. This entitles you to an all-inclusive full service plan. The reliable and continuous operation of your PLATINblue system is our highest priority.



**So you can get off running right from the start,** we combine the PLATINblue system installation with training for you and your team. A follow-up training session within the following six months is included in the service package at no extra cost.



**We will accompany you** when working out your methods, step by step. Our application experts are available Monday through Friday to assist you with any application questions you might have. We also offer on-site demonstrations and method transfer workshops.

[www.knauer.net](http://www.knauer.net)

## HPLC · SMB · Osmometry

Wissenschaftliche Gerätebau  
Dr. Ing. Herbert Knauer GmbH  
Hegauer Weg 38  
14163 Berlin, Germany

Phone: +49 30 809727-0  
Telefax: +49 30 8015010  
E-Mail: [info@knauer.net](mailto:info@knauer.net)  
Internet: [www.knauer.net](http://www.knauer.net)

