

Flame Ionisation Detector (FID) for Agilent 1260 Infinity-SFC

SFC-FID-System



Distribuido en España por:

sim 
Scientific Instruments Manufacturer GmbH


 **Ingeniería Analítica**
Chromatography & Spectrometry

Tel: (+34) 93.590.28.50 Fax: (+34) 93.675.05.16
www.ingenieria-analitica.com inf@ingenieria-analitica.com


Assembling of FID and Agilent SFC (with restrictor and hexane purging)

SFC Module

FID



Binary Pump




binary pump

channel A: CO₂

channel B: purging (optional, here: hexane)


Development of Electronic Components (I)

EPC/FID-Board





Development of Electronic Components (II)


Communication-Board



OpenLAB CDS ChemStation : FID Module Configuration / FID Parameters



Online Screen: Status of SFC-FID and PAL



ASTM D-5186: Aromatic content of diesel fuels and aviation turbine fuels by SFC-FID

Composition of Performance

Mixture:

75 % n-hexadecane (**C₁₆**)


20% toluene (**T**)

3 % tetrahydronaphthalene (**THN**)

2 % naphthalene (**N**)


Analytical conditions:

system: 1260 LC, SFC-FID
column: YMC-PACK-SIL_06, 250 × 4.5 mm, 6 nm S-, 5 µm or Agilent ZORBAX RX-SIL, 4.6 × 250 mm, 5 µm
column temp.: 25.0 °C
mob. phase: 100 % CO₂
flow : 1.8 ml
SFC (BPR): 25 °C (150 bar)
injection: 0.5 µl
FID : 300 °C
50 ml H₂, 500 ml air,
50 ml N₂ (make-up)




non aromatic mono-aromatic poly-aromatic

Diesel-Sample – Integration and Calculation





content of aromatic components (A): $M (\%) + P (\%) = A (\%)$

ASTM D-5186: Analytical Requirements

	Requirement	SIM SFC-FID
Resolution	$R_{NM}(n-C_{16}/T) \geq 4$ $R_{MD}(THN/N) \geq 2$	$R_{NM}(n-C_{16}/T) = 10.1$ $R_{MD}(THN/N) = 4.1$
Retention Time Reproducibility	$RSD(n-C_{16}) < 0.5$ $RSD(T) < 0.5$	$RSD(n-C_{16}): 0.24 \%$ $RSD(T): 0.28 \%$
Detector Sensitivity	detection of 0.1 mass% T in $n-C_{16}$	

ASTM D-5186: Analytical Requirements

	Requirement	SIM SFC-FID																
Detector Accuracy	relative response factors (RRF): RRF $\pm$ 10 % of theoretical value	example: toluene 																
Detector Linearity Check	aromatic content of diesel fuel and two dilutions	 <table><thead><tr><th>Sample</th><th>measured aromatics (mass%)</th><th>exp. aromatic content (mass%)</th><th>aromatics deviation (mass%)</th></tr></thead><tbody><tr><td>Fuel</td><td>21.8</td><td></td><td></td></tr><tr><td>F-Dilution* (1:1)</td><td>11.1</td><td>10.8</td><td>0.3</td></tr><tr><td>F-Dilution* (1:3)</td><td>5.9</td><td>5.7</td><td>0.2</td></tr></tbody></table>	Sample	measured aromatics (mass%)	exp. aromatic content (mass%)	aromatics deviation (mass%)	Fuel	21.8			F-Dilution* (1:1)	11.1	10.8	0.3	F-Dilution* (1:3)	5.9	5.7	0.2
Sample	measured aromatics (mass%)	exp. aromatic content (mass%)	aromatics deviation (mass%)															
Fuel	21.8																	
F-Dilution* (1:1)	11.1	10.8	0.3															
F-Dilution* (1:3)	5.9	5.7	0.2															

ASTM D-6550 : Olefin Content of Gasoline by SFC




- two columns / two TCCs:
 1. silica column (25 °C)
 2. silver loaded column (80 °C)

➔ two thermost. column compartments necessary
- valve-switching solution


forward-flow to trap the olefins on the Ag⁺ column
back-flush for quantification
- detector:

SIM FID for SFC

ASTM D-6550 : Valve switching solution




Column and valve configuration



- two 6-port valves with 4 valve positions (A-D) for forward-flush and backflush modes to separate and quantify the olefins
- position D allows Test Method **D5186** to be performed without changing the system

ASTM D-6550 : Performance Test Mixture (PTM)



Performance Test Mixture (PTM):

80 % n-hexane (**n-C₁₆**) → **SATURATES**


10 % toluene (**T**) → **AROMATICS**

10 % 3-methyl-2-pentene, cis/trans (**3M2PE**) → **OLEFINS**

→ final method after
system optimization
(switching times!)

ASTM D-6550 : Precision of Olefin Measurement

10 injections of PTM:



requirement:

RT-Repeatability < 0.5 %


measured:

0.04 – 0.13 % (RSD)

→ RT-repeatability is important for quantification:
column switching is performed on time basis


	Hexane		Toluene		Olefin 3M2PE	
	RT	Area	RT	Area	RT	Area
Average	1.8	16161	4.02	2198.6	7.16	2095.5
SD	0.002	130.6	0.002	17.7	0.003	26.0
RSD [%]	0.13	0.81	0.04	0.80	0.04	1.26

ASTM D-6550 : Olefin Calibration



requirement: intercept statistically not different from zero ✓

ASTM D-6550 : Gasoline Sample



repeatability limit: 0.7 mass%
(shall not be exceeded in more than one case in twenty)

gasoline sample:

12.8 mass% olefins (mean value)

19 results, difference $\leq 0.6 \text{ mass\%}$ ✓

SFC-FID Analyzers



Determination of Aromatic Content
in Diesel Fuel According to
ASTM D-5186
Enhancing the Agilent 1260 Infinity Analytical SFC
System with a Flame Ionization Detector

Application Note
Energy & Chemicals

Authors
Martina Noll-Borchers, Timo Höltscher
SIM Scientific Instruments
Manufacturer GmbH
Oberhausen, Germany
Edgar Naegele, Markus Becker
Agilent Technologies, Inc.
Waldbrunn, Germany



ASTM D5186: 5991-5682EN



Determination of Olefin Content in
Gasoline According to ASTM D6550

Application Note

Energy & Chemicals

Authors
Martina Noll-Borchers, Timo Höltscher
SIM Scientific Instruments Manufacturer GmbH
Oberhausen, Germany
SIM Scientific Instruments Manufacturer GmbH
Edgar Naegele, Markus Becker
Agilent Technologies, Inc.
Waldbrunn, Germany



ASTM D6550: Draft

- **consist of**
- 1260 SFC Control Module and Binary Pump
 - Thermostatted Column Compartment
(1x or 2x together with column switching)
 - 1290 Infinity LC-Injector HTC (width 50 cm) incl. DLW Option
 - SIM FID for SFC with OpenLAB Driver (ChemStation)
- ➔ **two applications: ASTM D-5186 + ASTM D-6550** also with the same system



**Thank you for your
attention!**