

# CDSolutions

## APPLICATIONS INFORMATION USING ADVANCED SAMPLE HANDLING TECHNOLOGY

### Quantitation of Low Levels of Methyl Methacrylate in a Copolymer

When a material is made using several different monomers, like a styrene/butadiene rubber, ethylene/propylene copolymer or a latex paint, pyrolysis-GC/MS can be a valuable tool in identifying not only the monomers used, but the relative amounts of each. Whether mixtures, blends, laminates, random or block copolymers, compounds will be produced relating to each of the monomers present, and the peaks they make can be used to study the molecular formula. This is frequently done using copolymers with just a few monomers, each of which represents a significant part of the copolymer. But the same techniques may be used to determine relatively small concentrations of monomers, and even traces of contaminants.

The polymers shown here are primarily styrene, but each contained a small amount of methyl methacrylate - specifically 0.1%, 0.2%, 0.3% and 0.4%. Figure 1 shows a pyrogram of the 0.4% product, showing the styrene monomer, dimer and trimer as well as peaks for alpha-methyl styrene and toluene, also pyrolysis products of polystyrene, and a small, but still measurable peak for the MMA monomer.

Figure 2 contrasts the 0.1% and 0.4% samples. As the amount of MMA increases, the MMA peak becomes larger, relative to other more constant peaks. The ratio of the MMA peak area to the toluene peak area shows a linear relationship to the concentration of MMA, as shown in Figure 3. Using a peak area ratio of two peaks produced by pyrolysis of the sample makes the assay independent of the sample size, so it is not necessary to weigh each sample before pyrolysis. Producing a graph from polymers of known content then makes it simple to determine the amount of MMA contained in an unknown polymer.




Figure 1. Copolymer with trace of MMA.





Figure 2. 0.1% and 0.4% MMA in copolymer.



**Figure 3.** MMA/Toluene peak area ratios vs. MMA concentration.

### Experimental Parameters

All samples were pyrolyzed in quartz tubes, using a CDS Pyroprobe 5200.

#### Pyroprobe

Pyrolysis: 750°C for 15 seconds  
 Interface: 325°C for 4 minutes  
 Carrier flow: 30 ml/min  
 Valve oven: 325°C  
 Transfer line: 325°C

#### GC/MS

Column: 30 m x 0.25 mm 5% phenyl MS  
 Carrier: Helium  
 Split: 50:1  
 Oven program:  
 40°C for 2 minutes  
 10°C/minute to 325°C

FOR MORE INFORMATION  
 CONCERNING THIS APPLICATION,  
 WE RECOMMEND THE  
 FOLLOWING READING:

F. Wang, Composition and Micro-structure Determination of a Latex System by Pyrolysis Gas Chromatography, Anal. Chem. 71 (1999) 4776-4780.

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.

CDS Analytical, LLC has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto [www.cdsanalytical.com](http://www.cdsanalytical.com).