


Seamless method transfer from the Agilent 1200 Series to the Agilent 1220 and 1260 Infinity LC systems

Comparison of system suitability results according to USP method USP/NF 23 for the analysis of Paracetamol (Acetaminophen) and Aspirin in pain relievers

Technical Overview



Introduction

Regulated environments in the pharmaceutical industry use validated methods that must be in compliance with EP or USP regulations. A requirement for all new instrumentation is that validated methods on the new instrument give results within the operating range and acceptance criteria of existing instruments. This eliminates the time consuming and costly need to revalidate an existing method.

The analysis of paracetamol and aspirin in pain relievers was chosen to demonstrate that seamless method transfer is possible (Figure 1). An isocratic method with UV detection was used according to USP/NF 23.

This analysis used the system suitability test, developed on an Agilent 1200 Series LC system to test whether the Agilent 1220 and Agilent 1260 Infinity LC systems fulfilled the acceptance criteria. The chromatographic conditions of the system suitability test transferred to the Agilent 1260 and Agilent 1220 Infinity LC systems without changes. The results obtained were compared for compliance with results typically obtained on the Agilent 1200 Series LC system.



Agilent Technologies




Figure 1
Structures of acetaminophen (paracetamol), aspirin and caffeine.

Experimental

Instrumentation

A standard Agilent 1200 Series LC

system, an Agilent 1220 Infinity LC system and an Agilent 1260 Infinity LC system with the following configurations were used:

| Configuration of the Agilent 1220 Infinity LC system (G4290B) | Configuration of the Agilent 1260 Infinity LC system | Configuration of the Agilent 1200 Series system |
|---|---|--|
| Gradient pump and vacuum degasser | Quaternary pump with integrated vacuum degasser (G1311B) | Quaternary pump (G1354A) with external degasser vacuum degasser (G1322A) |
| Column oven | Column compartment (G1316A) | Column compartment (G1316A) |
| Variable wavelength detector with 10 mm path length cell | Diode array detector (G4212B) with 10 mm path length cell | Diode array detector (G1315D) with 10 mm path length cell |
| Standard Auto sampler, sample cooling not available | Autosampler (G1367E) + sample cooling unit at 20 °C | Standard Autosampler +sample cooling unit at 20 °C |
| Software: ChemStation B.04.02 | Software: ChemStation B.04.02 | Software: ChemStation B.04.02 |

Table 1
Instrument configurations.

Chromatographic conditions according to USP method

| | |
|-------------------|---|
| Column: | Agilent ZORBAX Eclipse Plus, 3 mm × 100 mm, 3.5 μm (internal diameter 35% less than original method, particle size 30% less than original method) |
| Mobile phase: | Water/methanol/acetic acid = 69/28/3 |
| Pump settings: | No gradient (in accordance with EP regulations) |
| Stop time: | 5 min |
| Flow rate: | 1 mL/min, isocratic (50% less than original method) |
| Injection volume: | 10 μL |
| Column temp: | 45 °C |
| Detector: | Agilent 1220 Infinity LC system with 10 mm path length flow cell Agilent 1200 Series system with 10 mm path length Agilent 1260 Infinity LC system with 10 mm path length Peak width 0.05 min (10 Hz) Signal 275 nm |

The original method was changed according to the typical allowed changes for chromatographic parameters (Table 2).

For all experiments the same column, same batch of prepared sample and the same batch of mobile phase was used.

Preparation of samples

The reference solution was prepared as follows (Table 3).

| Chromatographic parameter | Typically allowed changes |
|--|--|
| Mobile phase pH | ± 0.2 units |
| Concentration of salts in buffers | ± 10% |
| Ratio of mobile phase percentages | ± 30% of the minor component, or 0.2% absolute value of that component, whichever is greater. However, a change in any component cannot exceed ± 10% absolute, nor can the final concentration be reduced to zero. |
| Wavelength or UV detector | No change permitted |
| Column length | ± 70% |
| Internal diameter of column | ± 50% |
| Particle size of column packing material | Can be reduced by 50% |
| Flow rate | ± 50% |
| Injection volume | Increased to as much as twice the volume specified, provided no adverse effects. Must be within stated linearity range of the method |
| Column compartment temperature | ± 10 °C |

Table 2
Typically accepted changes for USP methods.

| | Stock solution in mobile phase | 1: 5 diluted in water |
|----------------|--------------------------------|-----------------------|
| Acetaminophen | 5.5 mg/10 mL | 1.1 µg/10 µL |
| Caffeine | 1.3 mg/10 mL | 0.26 µg/10 µL |
| Aspirin | 3.9 mg/10 mL | 0.78 µg/10 µL |
| Benzoic acid | 4 mg/10 mL | 0.8 µg/10 µL |
| Salicylic acid | 4 mg/10 mL | 0.8 µg/10 µL |

Table 3
Sample concentration.

Results and Discussion

System suitability testing verified that the LC system fulfills the acceptance criteria typical for USP methods. The following acceptance criteria were evaluated:

- Precision of areas must be < 2% RSD
- Precision of retention times must be < 0.5% RSD

- Resolution must be > 1.4
- Tailing factor for < 1.3
- Peak symmetry > 0.75

Evaluation of column performance parameters such as symmetry and tailing factors determine the influence of different instruments on these parameters.

Figure 2 shows an overlay of the chromatograms obtained on the three Agilent LC systems.




Figure 2
Overlay of chromatograms.

The difference in retention times is less than 5% from the reference chromatogram of the Agilent 1200 Series LC system. The system suitability test was applied to both new instruments and the results compared to those obtained on the Agilent 1200 Series LC system. In Tables 4–6, the results are combined. Data for retention times, symmetry, USP

tailing and resolution can be obtained from the “Extended Performance” report. This report is based on calibrated peaks. Signal-to-noise data are included in the “Performance and Noise” report. Statistical data are included in the “Sequence Summary” report. All report styles are available from the “Specify Report” menu in the ChemStation.¹

| Compound | Agilent 1200 Series LC system | | Agilent 1220 Infinity LC system | | Agilent 1260 Infinity LC system | |
|----------------|-------------------------------|-----------------|---------------------------------|------------|---------------------------------|------------|
| | Retention time (min) | Resolution (hH) | Retention time (min) | Resolution | Retention time (min) | Resolution |
| Acetaminophen | 0.656 | | 0.684 | | 0.654 | |
| Caffeine | 1.025 | 8.822 | 1.059 | 6.401 | 1.040 | 8.198 |
| Aspirin | 2.480 | 19.535 | 2.541 | 17.419 | 2.569 | 19.604 |
| Benzoic acid | 2.807 | 3.095 | 2.869 | 2.905 | 2.917 | 3.168 |
| Salicylic acid | 3.182 | 3.281 | 3.255 | 3.181 | 3.306 | 3.242 |

Table 4
Comparison of retention times, signal-to-noise and resolution data.

Figure 3 shows an overlay of six consecutive runs performed on the Agilent 1220 Infinity LC system for precision, evaluation of retention times and areas.




Figure 3
Overlay of six consecutive runs using the Agilent 1220 Infinity LC system; evaluation of precision data.

| Compound | Agilent 1200 Series LC System | | Agilent 1220 Infinity LC System | | Agilent 1260 Infinity LC System | |
|----------------|-------------------------------|--------------|---------------------------------|--------------|---------------------------------|--------------|
| | RSD RT (%) | RSD area (%) | RSD RT (%) | RSD area (%) | RSD RT (%) | RSD area (%) |
| Acetaminophen | 0.107 | 0.1012 | 0.043 | 0.0740 | 0.053 | 0.1227 |
| Caffeine | 0.068 | 0.0447 | 0.020 | 0.1291 | 0.036 | 0.0323 |
| Aspirin | 0.032 | 0.0797 | 0.019 | 0.0959 | 0.020 | 0.1260 |
| Benzoic acid | 0.032 | 0.1089 | 0.010 | 0.1034 | 0.018 | 0.2026 |
| Salicylic acid | 0.024 | 0.0303 | 0.025 | 0.0776 | 0.018 | 0.0270 |

Table 5
Comparison of precision data.

| Compound | Agilent 1200 Series LC System | | Agilent 1220 Infinity LC System | | Agilent 1260 Infinity LC System | |
|----------------|-------------------------------|--------------|---------------------------------|--------------|---------------------------------|--------------|
| | Peak symmetry | Peak tailing | Peak symmetry | Peak tailing | Peak symmetry | Peak tailing |
| Acetaminophen | 0.782 | 1.270 | 0.765 | 1.254 | 0.833 | 1.142 |
| Caffeine | 0.841 | 1.136 | 0.818 | 1.168 | 0.891 | 1.094 |
| Aspirin | 0.908 | 1.060 | 0.903 | 1.080 | 0.917 | 1.061 |
| Benzoic acid | 0.866 | 1.122 | 0.905 | 1.082 | 0.857 | 1.137 |
| Salicylic acid | 0.889 | 1.069 | 0.889 | 1.079 | 0.883 | 1.079 |

Table 6
Comparison of peak symmetry and peak tailing.

The results for the Agilent 1220 Infinity LC and the Agilent 1260 Infinity LC systems show that both systems are fully compliant with system suitability acceptance criteria stated for the Agilent 1200 Series LC system.

Conclusion

Seamless method transfer is possible for the analysis of acetaminophen (paracetamol) in pain relievers using an isocratic EP method. The system suitability test, developed on an Agilent 1200 Series LC system for this method, was applied on the new Agilent 1260 Infinity LC and the Agilent 1220 Infinity

LC systems. The obtained results were compared with results typically obtained on the Agilent 1200 Series LC system (Table 6).

Both new instruments fulfilled the acceptance criteria of the system suitability test. The retention time shift for all peaks was < 5% using the same column and the same batch of mobile phase. The precision of retention times were better on the Infinity LC systems than on the original Agilent 1200 Series LC system.

Reference

1. "Understanding Your ChemStation," Agilent Literature, G2070-91124, 2007.

| Parameter and acceptance criteria | Agilent 1200 Series LC System | Agilent 1220 Infinity LC System | Agilent 1260 Infinity LC System |
|--|--------------------------------------|--|--|
| RSD RT < 0.5% RSD | 0.024–0.107% | 0.010–0.043% | 0.018–0.053 % |
| RSD area < 2% RSD | < 0.11% | < 0.13% | < 0.21% |
| Rs > 1.4 Benzoic acid | 3.095 | 2.869 | 3.168 |
| Tailing factors < 1.3 | < 1.28 | < 1.26 | < 1.15 |
| Symmetry factors > 0.75 | > 0.78 | > 0.766 | > 0.853 |

Table 6
Comparison of acceptance criteria.

www.agilent.com/chem/1200

© Agilent Technologies, Inc., 2010
July 1, 2010
Publication Number 5990-6206EN

