

Sensitive Femtogram Determination of Aflatoxins B₁, B₂, G₁ and G₂ in Food Matrices using Triple Quadrupole LC/MS

Application Note

Food Safety

Authors

Yang Chen and Jack Cappozzo
National Center for Food Safety and
Technology
Illinois Institute of Technology
Illinois
USA

Peter JW Stone
Agilent Technologies Inc
Stevens Creek Boulevard
Santa Clara, CA, 95051
USA.

Abstract

A simple and inexpensive sample cleanup procedure based on a dispersive solid phase adsorption approach (C18) is effective in removing background matrix contaminants for reliable determination of aflatoxins in food at the femtogram level by triple quadrupole LC/MS. This application demonstrates fast analysis time (< 6 min) with good chromatographic resolution and separation for all four aflatoxins. Standard curves for each aflatoxin analyte show good linearity (> 0.998) across a wide concentration range (0.1–100 µg/L). Recoveries using the dispersive solid phase adsorption approach were between 85–110% for each aflatoxin for all four spiked food matrices and were comparable to other widely used SPE routines. The limit of detection was determined to be < 0.15 µg/kg and the limit of quantitation < 0.5 µg/kg for all four sample matrices. Precision data was typically below 5% RSD for all analytes.



Agilent Technologies

Introduction

Aflatoxins are a group of mycotoxins produced as metabolites by the fungi *aspergillus flavus* and *aspergillus parasiticus* [1]. They can be found in various foods including grains, nuts, and spices [2]. There are four major naturally occurring aflatoxins: B₁, B₂, G₁ and G₂ (Figure 1). Exposure to them can cause cancer in humans and live stock, therefore reliable and sensitive analytical methods for the determination of aflatoxins are required to safeguard our food supply.

Experimental

These analyses were performed using an Agilent G6460A Triple Quadrupole LC/MS/MS System equipped with Agilent Jet Stream Technology [3] using an Agilent 1200 Series SL LC. The LC system consisted of a binary pump (G1312B), vacuum degasser (G1379B), a low carryover automatic liquid sampler (G1367D), thermostatted column compartment (G1316B) and MassHunter data system.

Aflatoxin standards and foods

Purified aflatoxin standards (B₁, B₂, G₁ and G₂) were obtained from Sigma-Aldrich. Aflatoxin-free corn flour, wheat, peanut and walnut samples obtained from a local grocery store were used for recovery studies.




Figure 1. Structures of aflatoxins B₁, B₂, G₁ and G₂

Instrumentation

Rapid Resolution HPLC Conditions and Configuration:

- Agilent 1200 Series Binary Pump SL (G1312B)
- High Performance WP Sampler SL Plus (G1367D)
Sampler Thermostat (G1330B)
- Thermostatted Column Compartment SL, including 10P/Two-Position switching valve (G1316B with option #057)

Method Conditions:

Column:	Agilent ZORBAX Eclipse Plus C18, 2.1 x 50 mm, 1.8 µm	
Column temperature:	40 °C	
Injection volume:	5 µL	
Autosampler temp:	4 °C	
Needle wash:	Flushport (100% methanol), 5 seconds	
Mobile phase:	A = 10 mM NH ₄ acetate in water B = 100% methanol	
Gradient flow rate:	0.6 mL/min (no split)	
Gradient:	Time (min)	%B
	0	5
	5	100
	6	100
Analysis time:	6 min	
Equilibration time:	1.5 min	
Total run time:	7.5 min	

Mass Spectrometer Source Conditions and Configuration:

Agilent 6460 Triple Quadrupole LC/MS equipped with Agilent Jet Stream Technology.

Ion Source Conditions:

Ion Mode:	ESI/Agilent Jet Stream, Positive ionization
Capillary Voltage:	4000 V
Drying gas (nitrogen):	10 L/min
Drying gas temperature:	325 °C
Nebulizer gas (nitrogen):	50 psi
Sheath Gas temperature:	350 °C
Sheath Gas flow:	11 L/min
Nozzle Voltage:	0 V
Q1 and Q2 Resolution:	0.7 amu [autotune]
Delta EMV:	400V

The Triple Quadrupole MS MRM parameters are listed in Table 1. All fragmentor voltage (frag) settings and respective collision energies (CE) and the most abundant MS/MS product ions per analyte were determined automatically using the Agilent MassHunter Optimizer Software.

Table 1. MRM Transitions for Aflatoxins and Respective Internal Standards

Name	Retention time (min)	Fragmentor voltage (V)	Precursor ion (m/z)	Product ion (m/z)	Collision energy (eV)
Aflatoxin B ₁	4.68	130	313.1	241.1	35
				285.1	20
				269.1	25
Aflatoxin B ₂	4.57	130	315.1	287.1	25
				259.1	25
				243.1	40
				243.1	25
Aflatoxin G ₁	4.40	130	329.1	243.1	25
				311.1	20
				283.1	20
Aflatoxin G ₂	4.26	130	331.1	245.1	30
				285.1	25
				313.1	25
				301.1	20
Isotope B ₁	4.68	130	330.1	255.1	40
				303	25
Isotope B ₂	4.57	130	332.1	273.0	30
				257.1	25
Isotope G ₁	4.40	130	346.1	299.1	25
				299.1	25
Isotope G ₂	4.26	130	348.1	330.1	25
				259.1	30

Sample Preparation and Recovery Studies

Corn flour, ground wheat, peanut and walnut samples (10 g each) were spiked with a mixture of four aflatoxin standards, each at 5 and 25 ng/g. This was then extracted using 40 mL of acetonitrile-water (84:16, v/v) for 30 min with shaking at room temperature. The extract was cleaned up using both C18 powdered adsorbent material (ODS SPE bulk sorbent, Agilent p/n 5982-1182) and MycoSep 226 multifunctional SPE (Romer). Aliquots (0.4 mL) of the cleaned up extracts were diluted with 0.6 mL 10 mM ammonium acetate in water.

The sample was then centrifuged at 14,000 rpm for 3 min prior to LC/MS/MS analysis.

Each food matrix and spike level was conducted in seven replicates to represent and maintain statistical integrity. A schematic of this sample preparation is illustrated in Figure 2.




Figure 2. Schematic matrix sample preparation workflow showing the dispersive solid phase adsorption approach versus a widely used SPE approach.

Results and Discussion

The rapid chromatography conditions as outlined in the experimental section yielded good chromatographic resolution for each aflatoxin analyte and each analysis was completed in six minutes. A typical chromatogram is shown in Figure 3(a), which illustrates 1 ppb concentration level of each aflatoxin together with the corresponding isotopically labeled internal standards at a concentration level of 2.5 ppb (Figure 3(b)). These chromatograms show overlaid extracted ion chromatograms (EICs).

Standard curves for aflatoxins B₁, B₂, G₁ and G₂ all showed a good linearity through the concentration range 0.1 to 100 ppb each with a linear correlation (R^2) of greater than 0.999. Figure 4 illustrates an overlay of each standard curve on the same scale, but without internal standard correction. The use of internal standards effectively adjusted for matrix differences, as shown in Figure 5.




Figure 3. LC/MS/MS chromatogram of aflatoxin B₁, B₂, G₁ and G₂ standards at 1 ppb with corresponding isotopically labelled internal standards at 2.5 ppb.




Figure 4. Overlaid standard curves for aflatoxins B₁, B₂, G₁ and G₂.




Figure 5. Overlaid standard curves for aflatoxins B₁, B₂, G₁ and G₂ with internal standard correction.

In order to determine the limits of detection (LOD) and reporting (LOR), seven separate and mutually exclusive batches of each food matrix were tested via the two sample preparation protocols and analytical methodology outlined previously. The results outlined in this document are derived from the average values across the seven batches of each matrix (N=7.)

Limits of detection were determined using the protocol of chromatographic signal-to-noise ratio of above 3/1 (peak to peak.) Limits of reporting were determined using the protocol of chromatographic signal-to-noise ratio above 10/1 (peak to peak.)

Table 2 details the observed LODs and LORs for each aflatoxin across the series of four food matrices. The limit of detection overall was determined to be < 0.15 µg/kg and the limit of quantitation < 0.5 µg/kg for all four sample matrices and both sample preparation routines.

Tables 3 and 4 summarize the LOD data obtained across the seven batches. Table 3 data is presented with respect to dispersive solid phase adsorption only, and Table 4 data using the Mycosep SPE sample preparation only.

Table 2. Limits of Detection and Reporting Observed for Aflatoxins B₁, B₂, G₁ and G₂ Across Four Food Matrices via Dispersive SPA and Mycosep SPE Sample Preparation Approaches

Food Matrix	Aflatoxin	Mycosep (#226, Romer)	Dispersive C18 ODS SPE Bulk Sorbent, Agilent (p/n 5982-1182)		
		LOD ng/g (S/N>3)	LOR ng/g (S/N>10)	LOD ng/g (S/N>3)	LOR ng/g (S/N>10)
Corn (Ave 7-batches)	B ₁	0.047	0.16	0.060	0.20
	B ₂	0.036	0.12	0.085	0.28
	G ₁	0.08	0.28	0.10	0.35
	G ₂	0.046	0.15	0.033	0.11
Wheat (Ave 7-batches)	B ₁	0.068	0.23	0.012	0.042
	B ₂	0.11	0.36	0.037	0.12
	G ₁	0.14	0.47	0.15	0.50
	G ₂	0.038	0.13	0.11	0.36
Peanut (Ave 7-batches)	B ₁	0.051	0.17	0.056	0.19
	B ₂	0.045	0.15	0.069	0.23
	G ₁	0.07	0.23	0.05	0.15
	G ₂	0.052	0.17	0.14	0.45
Walnut (Ave 7-batches)	B ₁	0.12	0.41	0.093	0.31
	B ₂	0.035	0.12	0.098	0.33
	G ₁	0.03	0.10	0.12	0.40
	G ₂	0.047	0.16	0.04	0.13

Table 3. LOD Results Observed for Aflatoxins B₁, B₂, G₁ and G₂ via the Dispersive SPA

Food Matrix	B ₁ LOD (ng/g)	B ₂ LOD (ng/g)	G ₁ LOD (ng/g)	G ₂ LOD (ng/g)
Corn	0.060	0.085	0.100	0.033
Wheat	0.012	0.037	0.150	0.110
Peanut	0.056	0.069	0.050	0.140
Walnut	0.093	0.098	0.120	0.040
Average	0.055	0.072	0.105	0.080
Mass On-Column (fg)	275	360	525	400

Table 4. LOD Results Observed for Aflatoxins B₁, B₂, G₁ and G₂ via the Mycosep, #226 SPE Sample Preparation Approach

Food Matrix	B ₁ LOD (ng/g)	B ₂ LOD (ng/g)	G ₁ LOD (ng/g)	G ₂ LOD (ng/g)
Corn	0.047	0.036	0.080	0.046
Wheat	0.068	0.110	0.140	0.038
Peanut	0.051	0.045	0.070	0.052
Walnut	0.120	0.035	0.030	0.047
Average	0.072	0.057	0.080	0.046
Mass On-Column (fg)	360	283	400	229

Sample Preparation Approach

Recovery studies were extensively undertaken for both sample cleanup techniques in parallel across the four matrices and across the seven batches for each aflatoxin analyte. Figures 6 (a) and (b) graphically depict the recovery trends

across the four food matrices for the dispersive C18 cleanup approach and the MycoSep SPE cleanup, respectively. As illustrated, the recoveries for both sample cleanup procedures were between 85-110 % for each of the aflatoxins for all four spiked food matrices, with the MycoSep cleanup method only marginally better than the C18 one for walnut samples .




Figure 6(a). Recovery of aflatoxin B₁, B₂, G₁ and G₂ from food matrices using C₁₈ dispersive cleanup.




Figure 6(b). Recovery of aflatoxin B₁, B₂, G₁ and G₂ from food matrices using MycoSep, #226 SPE cleanup.

Aflatoxin analyte recovery data for each separate food matrix is detailed in Tables 5 through 8 and was undertaken at two concentration spiked levels of 5 ng/g and 25 ng/g. Each sample batch tested was split and divided between the

two cleanup procedures outlined in this document following the natural settlement and supernatant transfer step outlined in Figure 2(a).

Table 5. Spiked Corn Samples – Recovery Studies (% Recovery, ± RSD, N=7)

	Aflatoxin	Corn spiked at 5 ng/g C18 cleanup	Corn spiked at 25 ng/g C18 cleanup	Corn spiked at 5 ng/g Mycosep#226	Corn spiked at 25 ng/g Mycosep#226
No Internal Standard	B ₁	101.7 ± 3.7	95.7 ± 3.0	107.8 ± 2.8	105.4 ± 2.8
	B ₂	95 ± 6.3	95.4 ± 1.4	103.0 ± 5.7	105.3 ± 2.0
	G ₁	102.7 ± 7.1	96.8 ± 1.75	110.2 ± 7.9	103.7 ± 3.0
	G ₂	107.9 ± 3.5	97.8 ± 0.88	108.1 ± 5.4	104.3 ± 2.4
Internal Standard	B ₁	102.3 ± 2.9	100.1 ± 2.4	108.2 ± 5.1	97.8 ± 3.3
	B ₂	100.0 ± 7.9	94.0 ± 3.1	101.7 ± 4.7	92.8 ± 3.5
	G ₁	107.3 ± 3.5	97.0 ± 6.0	110.3 ± 3.6	102.5 ± 1.8
	G ₂	101.3 ± 5.6	100.4 ± 3.8	104.9 ± 5.4	97.2 ± 6.3

Table 6. Spiked Wheat Samples – Recovery Studies (% Recovery, ± RSD, N=7)

	Aflatoxin	Wheat spiked at 5 ng/g C18 cleanup	Wheat spiked at 25 ng/g C18 cleanup	Wheat spiked at 5 ng/g Mycosep#226	Wheat spiked at 25 ng/g Mycosep#226
No Internal Standard	B ₁	100.1 ± 4.4	96.6 ± 2.9	113.5 ± 5.9	100.6 ± 1.8
	B ₂	98.2 ± 6.9	96.4 ± 2.6	105.1 ± 4.5	102.1 ± 4.4
	G ₁	100.5 ± 5.5	105.4 ± 3.8	111.5 ± 10.4	106.1 ± 4.0
	G ₂	104.9 ± 3.2	106.7 ± 1.3	108.6 ± 5.2	103.7 ± 2.9
Internal Standard	B ₁	100.9 ± 3.6	109.3 ± 4.7	107.5 ± 4.8	111.7 ± 4.9
	B ₂	85.2 ± 7.7	99.8 ± 2.8	92.4 ± 6.3	101.0 ± 4.0
	G ₁	110.6 ± 7.8	112.8 ± 1.8	117.6 ± 7.7	109.3 ± 5.7
	G ₂	108.4 ± 6.2	108.3 ± 3.9	115.6 ± 7.1	109.8 ± 3.6

Table 7. Spiked Peanut Samples – Recovery Studies (% Recovery, ± RSD, N=7)

	Aflatoxin	Peanut spiked at 5 ng/g C18 cleanup	Peanut spiked at 25 ng/g C18 cleanup	Peanut spiked at 5 ng/g Mycosep#226	Peanut spiked at 25 ng/g Mycosep#226
No Internal Standard	B ₁	96.7 ± 3.4	97.0 ± 4.6	112.0 ± 8.4	104.9 ± 1.7
	B ₂	98.3 ± 4.7	97.4 ± 2.9	108.0 ± 4.6	104.5 ± 2.0
	G ₁	95.0 ± 5.6	95.0 ± 4.9	109.9 ± 2.1	105.7 ± 3.4
	G ₂	100.0 ± 2.3	100.0 ± 2.0	114.7 ± 3.2	106.3 ± 1.1
Internal Standard	B ₁	101.8 ± 3.6	96.1 ± 2.0	100.0 ± 6.8	103.0 ± 3.5
	B ₂	102.5 ± 5.5	100.2 ± 5.0	99.4 ± 4.1	102.9 ± 2.7
	G ₁	105.7 ± 7.3	99.2 ± 2.2	105.2 ± 4.3	101.7 ± 5.2
	G ₂	107.5 ± 10.9	104.9 ± 6.7	109.3 ± 8.7	102.4 ± 3.1

Table 8. Spiked Walnut Samples – Recovery Studies (% Recovery, ± RSD, N=7)

	Aflatoxin	Walnut spiked at 5 ng/g C18 cleanup	Walnut spiked at 25 ng/g C18 cleanup	Walnut spiked at 5 ng/g Mycosep#226	Walnut spiked at 25 ng/g Mycosep#226
No Internal Standard	B ₁	84.9 ± 3.7	85.2 ± 2.2	101.4 ± 3.2	101.0 ± 2.3
	B ₂	91.5 ± 3.9	89.8 ± 2.8	104.2 ± 7.9	106.3 ± 2.9
	G ₁	89.4 ± 4.4	86.7 ± 1.5	103.9 ± 5.9	101.7 ± 4.2
	G ₂	84.0 ± 4.0	83.1 ± 1.3	109.9 ± 3.4	106.3 ± 1.5
Internal Standard	B ₁	106.5 ± 4.9	98.9 ± 4.1	93.8 ± 1.4	100.2 ± 2.9
	B ₂	99 ± 5.4	96.5 ± 3.5	92.4 ± 2.7	98.7 ± 4.4
	G ₁	103.2 ± 5.9	94.9 ± 2.5	102.8 ± 9.0	102.1 ± 3.9
	G ₂	100.2 ± 6.2	97.5 ± 4.6	99.5 ± 6.8	101.2 ± 3.8

Conclusions

An inexpensive and rapid LC/MS/MS method has been developed for the analysis and confirmation of aflatoxins B₁, B₂, G₁ and G₂ in cereals and nuts, with a detection limit of less than 1 ppb. This method is inclusive of sample preparation using a dispersive C18 solid phase adsorption approach (Agilent bulk sorbent p/n 5982-1182.) The performance of this simple sample cleanup procedure was comparable to that of a widely used and generally accepted SPE approach in terms of matrix cleanup and aflatoxin recoveries.

Aflatoxin limits of detection were determined to be less than 0.15 µg/kg and aflatoxin limits of reporting were all less than 0.5 µg/kg for all four sample matrices (corn, wheat, peanut and walnut.)

Standard curves for aflatoxins B₁, B₂, G₁ and G₂ showed a good linearity through the concentration range of 0.1 to 100 ppb with a linear correlation (R^2) of greater than 0.999 for all analytes.

Aflatoxin recoveries were between 85-110% for each of the aflatoxins for all four spiked food matrices using the dispersive C18 solid phase adsorption approach.

References:

1. Microchim Acts 153, 2006, 101-108
2. Rapid Commun. Mass Spectrum, 2009; 23: 3-11
3. "Agilent Jet Stream Thermal Gradient Focussing Technology" Agilent Technologies Publication 5990-3494EN.

For More Information

For more information on Agilent products and services visit our website

www.agilent.com/chem

For more details concerning this application note, please contact Peter JW Stone at Agilent Technologies Inc., 5301 Stevens Creek Blvd, Santa Clara, CA, 95051, USA.

Reference herein to any specific commercial products or non-profit organization, process, or service by trade name, trademark, manufacturer, or other-wise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government and shall not be used for advertising or product endorsement purposes.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011
Printed in the USA
March 9, 2011
5990-6894EN



Agilent Technologies