

Simple and Highly Efficient Drugs of Abuse Testing Methods by SPE and LC/MS/MS

Irina Dioumaeva¹, John Hughes^{*2}, William Hudson¹
¹Agilent Technologies, Lake Forest, CA, USA; ²Agilent Technologies, Pleasanton, CA, USA.

2011 SOFT - TIAFT
 Poster # P175



Introduction

New SAMHSA guidelines (October 2010) allow LC/MS/MS methods to be used for confirmation of initial drug tests by government-certified (NLCP) workplace drug testing labs. SAMHSA requires to analyze 6 classes of drugs and metabolites in urine: 1) 6-acetylmorphine; 2) phencyclidine; 3) benzoylecgonine; 4) a group of 5 amphetamines; 5) morphine and codeine, and 6) 11-nor-9-carboxy- Δ^9 -tetrahydrocannabinol. We present simplified LC/MS/MS confirmation methods, with sample preparation, for all SAMHSA-required categories using the latest instruments and separation tools from Agilent Technologies, including Bond Elut Plexa PCX mixed mode polymeric SPE cartridges and a Poroshell 120 column. Without elaborate sample pretreatment, solid phase extraction on Plexa PCX provides high extraction efficiency and reduced ion suppression due to a combination of a strong cation-exchange functionality, an amide-free polymer surface and a hydrophilic pore gradient. The LC method on a superficially porous column provides excellent separation and short retention times; with max pressure < 400 bar, these methods do not require high pressure UHPLC systems. As required by SAMHSA and NLCP, our MS/MS methods use qualifier ions and are linear across a wide dynamic range.

SAMHSA analytes and cutoffs

Target analyte	Confirmation cutoff, ng/mL	Limit of Detection, ng/mL
Amphetamines		
Amphetamine	250	25
Methamphetamine	250	25
MDA	250	25
MDMA	250	25
MDEA	250	25
Benzoylecgonine (BE)	100	10
11-nor-9-carboxy-Δ^9-THC	15	1.5
Phencyclidine (PCP)	25	2.5
Opiates		
Morphine	2000	200
Codeine	2000	200
6-Acetylmorphine (6-MAM)	10	1

Materials and instrumentation

- Agilent Bond Elut Plexa PCX 30 mg 3 mL cartridges
- Agilent Poroshell 120 EC-C18, 3 x 50 mm, 2.7 μ m column
- Agilent 1260 LC
- Agilent 6460 MS with Jet Stream electrospray source

Experimental

	BE	PCP	6-AM	Ampheta- mines	Opiates	THCA
Sample	1 mL urine spiked with ISTD	1 mL urine spiked with ISTD	1 mL urine spiked with ISTD	0.5 mL urine spiked with ISTD	0.5 mL urine spiked with ISTD	0.5 mL urine spiked with ISTD
Pre-treatment	1 mL 2% formic acid	1 mL 2% formic acid	1 mL 2% formic acid	1 mL 2% formic acid	Add 125 μ L conc. HCl. Incubate 90 min at 95 \pm 5 $^{\circ}$ C, cool and add 2 mL 0.1 M sodium acetate buffer (pH 4.5). Add 250 μ L 7 N KOH. Adjust pH <6. Centrifuge	Add 100 μ L 7 KOH. Incubate 30 min at 60 \pm 5 $^{\circ}$ C, cool and add 125 μ L methanol. Add 100 μ L glacial acetic acid and 1.5 mL 0.2 M sodium acetate buffer (pH 4). Centrifuge
Condition-ing	0.5 mL methanol	0.5 mL methanol	0.5 mL methanol	0.5 mL methanol	0.5 mL methanol	0.5 mL methanol
Wash 1	1 mL 2% formic acid	1 mL 2% formic acid	1 mL 2% formic acid	1 mL 2% formic acid	1 mL 2% formic acid	2 x 2 mL ACN:2% acetic acid (10:90)
Wash 2	1 mL methanol	1 mL methanol	1 mL methanol	1 mL methanol	1 mL methanol	2 mL ACN:2% acetic acid (30:70)
Drying	10 min	10 min	10 min	10 min	10 min	10 min
Rinse						Rinse with 0.2 mL hexane. Pull with low vacuum, dry 3 min
Elution	1 mL MeOH:NH ₄ OH (100:20)	1 mL Ethyl Acetate:MeOH: NH ₄ OH (80:20:5)	1 mL MeOH:NH ₄ OH (100:10)	1 mL ethyl acetate:MeOH:NH ₄ OH (50:50:20)	2 mL MeOH:NH ₄ OH(100:20) Soak with eluent - let eluate drip - apply low vacuum	1.5mL ethyl acetate:IPA (80:20) in 2 steps (0.5+1.0 mL). At each step soak with eluent - let eluate drip - apply low vacuum
Evaporation & Reconstitution	Evaporate to dryness at <40 $^{\circ}$ C. Reconstitute in 1 mL initial mobile phase	Evaporate to dryness at <40 $^{\circ}$ C. Reconstitute in 1 mL initial mobile phase	Evaporate to dryness at <40 $^{\circ}$ C. Reconstitute in 1 mL initial mobile phase	Evaporate to 0.2 mL at 37 $^{\circ}$ C; add 100 μ L 0.025 N HCl in MeOH, vortex; evaporate to dryness. Reconstitute in 0.5 mL initial mobile phase	Evaporate to dryness at <40 $^{\circ}$ C. Reconstitute in 0.5 mL initial mobile phase	Evaporate to dryness at <40 $^{\circ}$ C. Reconstitute in 0.5 mL initial mobile phase


	BE, PCP		6-AM		Amphetamines		Opiates		THCA	
Gradient	0 min	10%B	0 min	10%B	0 min	15%B	0-0.5 min	5%B	0 min	30%B
	0.5 min	10%B	1.5 min	25%B	1.5 min	15%B	1.5 min	25%B	1.0 min	95%B
	2.5 min	70%B	2.0 min	60%B	3.5 min	30%B	2.5 min	55%B	2.5 min	95%B
	2.51 min	90%B	2.1 min	90%B	3.6 min	90%B	2.6 min	90%B		
Mobile phase	A - 0.1% formic acid in water; B - 0.1% formic acid in methanol									A - 5 mM ammonium formate, B - methanol
Injection volume	2 μ L	2 μ L	10 μ L		2 μ L		2 μ L		10 μ L	

	BE, PCP	6-AM	Amphetamines	Opiates	THCA
ES Source	BE, PCP	6-AM	Amphetamines	Opiates	THCA
Ionization mode	positive	positive	positive	positive	negative
Capillary voltage	3000 V	2800 V	4000 V	3000 V	3000 V
Drying gas flow	10 L/min	13 L/min	10 L/min	10 L/min	11 L/min
Drying gas temp.	350 $^{\circ}$ C	350 $^{\circ}$ C	350 $^{\circ}$ C	350 $^{\circ}$ C	350 $^{\circ}$ C
Nebulizer gas	35 psi	35 psi	35 psi	35 psi	35 psi
Sheath gas flow	12 L/min	12 L/min	12 L/min	12 L/min	12 L/min
Sheath gas temp.	400 $^{\circ}$ C	400 $^{\circ}$ C	400 $^{\circ}$ C	400 $^{\circ}$ C	320 $^{\circ}$ C
MS	Delta EMV				
	200 V	200 V	200 V	0	800 V

Experimental

Compound	Parent	Product	Fragmentor	Collision energy
Benzoylecgonine	290.1	168.1	90	15
Benzoylecgonine	290.1	105.1	90	30
Benzoylecgonine-D₈	298.2	171.1	90	15
Benzoylecgonine-D₈	298.2	110.1	90	30
Phencyclidine	244.2	86.1	80	7
Phencyclidine	244.2	159.1	80	7
Phencyclidine-D₅	249.2	86.1	80	7
Phencyclidine-D₅	249.2	164.1	80	7
6-Acetylmorphine	328.2	165.1	140	40
6-Acetylmorphine	328.2	211.1	140	25
6-Acetylmorphine-D₅	334.2	165.1	140	40
6-Acetylmorphine-D₅	334.2	211.1	140	25
Amphetamine	136.1	91.1	64	14
Amphetamine	136.1	119.1	64	4
Amphetamine-D₅	142.1	93.1	66	5
Amphetamine-D₅	142.1	125.1	66	5
MDA	180.1	163.1	92	13
MDA	180.1	105.1	92	17
MDA-D₅	185.1	168.1	68	5
MDA-D₅	185.1	110.1	68	21
MDEA	208.1	163.1	88	8
MDEA	208.1	133.1	88	17
MDEA-D₅	214.2	166.1	90	8
MDEA-D₅	214.2	108.1	90	25
MDMA	194.1	163.1	84	5
MDMA	194.1	135.1	84	17
MDMA-D₅	199.1	165.1	82	4
MDMA-D₅	199.1	107.1	82	25
Methamphetamine	150.1	91.1	80	16
Methamphetamine	150.1	119.1	80	4
Methamphetamine-D₅	159.2	93.1	77	13
Methamphetamine-D₅	159.2	125.1	77	5
Ephedrine-Pseudoephedrine	166.1	133.1	80	21
Phentermine	150.1	133.1	80	6
Phenylpropanolamine	152.1	117.1	80	20
Codeine	300.2	165.1	130	46
Codeine	300.2	215.1	130	23
Codeine-D₅	306.2	165.1	130	44
Codeine-D₅	306.2	218.1	130	23
Morphine	286.1	165.1	130	43
Morphine	286.1	201.1	130	23
Morphine-D₆	292.1	181.1	130	40
Morphine-D₆	292.1	165.1	130	42
Oxycodone	316.2	298.1	130	15
Oxycodone	302.2	284.1	130	17
Hydrocodone	300.2	199.1	130	30
Norcodeine	286.1	225.1	130	20
Hydromorphone	286.1	185.1	130	28
11-nor-9-carboxy-Δ^9-THC	343.2	299.2	135	18
11-nor-9-carboxy-Δ^9-THC	343.2	245.1	135	30
11-nor-9-carboxy-Δ^9-THC-D₅	352.2	308.2	145	18
11-nor-9-carboxy-Δ^9-THC-D₅	352.2	254.2	145	30

Chromatograms, Poroshell 120 EC-C18, 3.0 x 50 mm, 2.7 μ m



Amphetamines and interferences peaks:

- Phenylpropanolamine
- Ephedrine
- Pseudoephedrine
- Amphetamine
- Methamphetamine
- MDA
- MDMA
- MDEA
- Phentermine

Opiates and interferences peaks:

- Morphine
- Oxycodone
- Hydromorphone
- Codeine
- Norcodeine (small blue peak)
- Oxycodone
- Hydrocodone


Method evaluation, n=5

Compound	Process efficiency * %	Extraction recovery * %	Matrix effect * %	Accuracy* %	Precision* (CV), %
Benzoylecgonine	85	86	99	103**	0.8**
Phencyclidine	83	85	98	93	0.5
6-Acetylmorphine	83	83	100	106**	0.6**
Amphetamine	86	94	91	107**	0.6**
Methamphetamine	93	94	99	105	0.5
MDA	91	95	95	92**	1.1**
MDMA	93	97	96	101**	0.5**
MDEA	95	96	98	106**	0.3**
Morphine	83	85	98	108	0.6
Codeine	85	86	99	108	0.7
11-nor-9-carboxy-Δ^9-THC	73	65	113	98	2.2

* Measured according to the method of Matuszewski et al. *Anal. Chem.* **75**: 3019-3030 (2003)
 All values except those denoted by * are obtained for the cutoff level
 ** Values obtained for 40% cutoff

Results and Discussion

Examples of Calibration Curves



Conclusions

- Presented methods meet SAMHSA requirements for LOD, linearity, accuracy and precision, and LC separation from possible interferences. Method accuracy is within 10% of target and precision (CV) is within 2.2%.
- High process efficiencies, extraction recoveries and insignificant matrix effects demonstrate excellent quality of Agilent Bond Elut Plexa PCX.
- Extraction methods are simplified to 3 types of pretreatment (simple formic acid pretreatment for four drug classes, acid hydrolysis for opiates and base hydrolysis for THCA) and to combinations of only three organic eluents.
- All LC separations use the same Agilent Poroshell120 EC-C18 column and utilize low backpressure <400 bar, simple mobile phases and no column heating. The longest separation (amphetamines) is completed within 3.2 minutes.
- MS methods use low sample injection volume (2 – 10 μ L), no sample pre-concentration, and demonstrate excellent S/N ratio for the lowest SAMHSA-required concentrations due to the enhanced sensitivity of Agilent 6460 QQQ with Jet Stream electrospray source.