



Author:

C. Zawodny

Thermal Desorption for PAH Analysis

Application Note

Environment

Polyaromatic hydrocarbons are produced when coal, gas, and oil as well as wood, and other organic materials are combusted. Even the grilling of meat produces some of these compounds. They are present in petroleum products as creosote and asphalt. They are very common environmental contaminants. The few natural sources of PAH's are produced by forest fires and volacanos. Chemically PAH's are a class of semivolatile compounds that are known as fused rings. The common environmental contaminants range in molecular weight from 128 (naph-thlene) to 276 (benzo(g,h,i)perylene). Most of them have mutagenic properties and are suspected carcinogens. Due to the health risks associated with these semivolatile materials, the EPA has developed a number of methods to qualitate and quantitate these compounds. EPA methods such as 8100, 8270, 610, and 625 use GC/MS to qualitate and quantitate semivolatiles present in soils, solids, and liquid wastes.

A quantitative PAH standard mix was purchased at a concentration of $2000\mu g/$ ml in methanol. Dilutions were made to give concentrations of 5,10,15, 35, 100, 300, and 500ng/ μ l respectively. Each concentration was run multiple times. The spiked sampling tube with Tenax was thermally desorbed at 340°C for 10 minutes and analyzed using a GC/MS. Figure 1 is a chromatogram of the standard showing the PAH components. Figures 2 and 3 show selected PAH compounds with area plotted against nanograms and a .99 correlation coefficient.



CDS Autosampler Dynatherm 9300

				4000000	
Value Oven:	20000			4000000	y = 7661x + 119405
Valve Oven:	300°C			3500000	$R^2 = 0.993$
Transfer Line:	325°C		Area	3000000	
Tube Heat: Trap Heat:	340°C 340°C	10 minutes 10 minutes	Alea	2500000	
				2000000	
map nout.	040 0	To minutes		1500000	
				1000000	
				500000	
				0	· · · · · ·
					0 100 200 300

300 400 500 600 nanograms



4500000



Figure 3. Benzo(a)Pyrene Linearity

GC/MS