

## TDTS 20

# Confirming sorbent tube retention volumes and checking for analyte breakthrough

*It is recommended that this Application Note is read in conjunction with Application Note TDTS 5 (Advice on sorbent selection, tube conditioning, tube storage and air sampling).*

### Introduction

Information on how to determine sorbent tube retention volumes, breakthrough volumes and safe sampling volumes is presented in Application Note TDTS 5 and in various international standard methods for thermal desorption, e.g. EN ISO 16017-1, Annexes A and B. However, as sorbent tubes age it is possible for retention volumes to fall, risking breakthrough of the most volatile target compounds. It is also possible for breakthrough volumes to fall under extreme environmental conditions (e.g. high humidity and high temperature), though the impact of such effects varies significantly for different sorbents.

This Application Note describes a simple procedure for checking analyte breakthrough prior to or during field monitoring, using pairs of tubes linked together in series. This guidance also complies with recommendations in standard thermal desorption methods for air monitoring.

### Field validation of quantitative retention during sampling

As a routine check of pumped tube retention during field air monitoring studies, two identical, conditioned sorbent tubes should be linked together in series using an inert stainless steel union with combined PTFE ferrules (part number C-UN010). The sampling end of the back-up tube should be connected to the outlet end of the front (sampling) tube. The sampling pump is attached to the exhaust end of the back-up tube. At least one such tube pair should be deployed in all routine monitoring exercises and, in large scale studies, it is recommended that at least one in ten monitoring locations are sampled using such tube pairs.

After monitoring, the back-up tubes should be analysed in the same analytical sequence as all the sampling tubes and study blanks.

- If the mass of one or more target analytes on the back-up tube is >5% of the mass on the front sampling tube, *breakthrough* of that compound or compounds can be said to have occurred. The amount can be quantified by adding the mass of analyte measured on the back-up tube to that determined on the front sampling tube.

- If the mass of one or more target analytes on the back-up tube is >10% of the mass on the front sampling tube, breakthrough of that compound or compounds is significant, and data for any such compounds should not be viewed as quantitative.

### Extended checks of breakthrough and other tube performance indicators

After a batch of tubes has been used extensively (e.g. 20 field monitoring exercises or thermal cycles) a representative number of tubes from that batch should be more stringently checked for breakthrough. In this case it is recommended that at least six tube pairs are prepared as above and used to sample a standard or real atmosphere under conditions which are as close as possible to the worst-case real-life scenario, *i.e.* highest natural humidity, highest ambient temperature and highest expected VOC concentrations.

The sampling points of all the tube pairs should be placed close together to ensure that, as far as possible, they are all sampling the same atmosphere. The selected sampling location should be well ventilated. Three different pump flow rates between 10 and 200 mL/min should be used for sampling, with duplicates at each level. All the tubes should be analysed in one sequence, and the results evaluated to ensure that compound masses increase relative to the volume of air sampled and are reproducible at each volume. Selective loss of one or more analytes at higher volumes indicates breakthrough. Detection of >5% of one or more analytes on the back-up tubes can be used to confirm breakthrough.

If one or more analyte is observed to increase at a higher rate than the sample volume, this can be an indication of artefact formation or sorbent/analyte breakdown. Such effects can be caused, for example, by the air containing above-ambient concentrations of reactive gases such as ozone, or by sorbent degradation. In the latter case, the tubes would need repacking.

## Appendix – Breakthrough volumes

### Notes on the tables

- Markes has not undertaken any validation work on the following data. It is always recommended that safe sampling volumes are checked using one of the procedures outlined above
- Breakthrough volumes have been extrapolated to relate to a standard ¼" stainless steel TD tube, packed with the typical mass of sorbent
- SSV = Safe sampling volume
- U = Unretained.

### 1. On 200 mg Carbotrap C–Carbograph 2TD at 20 °C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Methane	-164	U	–
Ethane	-89	U	–
Propane	-42	U	–
Butane	-0.5	U	–
Pentane	36	U	–
Hexane	69	0.3	0.15
Heptane	98	1.4	0.7
Octane	125	5.6	2.8
Nonane	151	20	10
Decane	174	100	50
Undecane	196	500	250
Dodecane	216	2600	1300
Tridecane	235	7400	3700
Tetradecane	254	22000	11000
Pentadecane	271	46000	23000
Hexadecane	287	100000	50000
Heptadecane	302	300000	150000
Octadecane	316	800000	400000
<b>Alcohols</b>			
Methanol	65	U	–
Ethanol	78	U	–
Propan-1-ol	97	U	–
Butan-1-ol	118	U	–
Pentan-1-ol	138	0.42	0.21
Hexan-1-ol	158	1.42	0.71
Heptan-1-ol	176	6.6	3.3
Octan-1-ol	195	16	8
Nonan-1-ol	213	40	20
Decan-1-ol	231	90	45
Undecan-1-ol	243	240	120
Dodecan-1-ol	259	800	400
Tridecan-1-ol	–	2000	1000
Tetradecan-1-ol	289	3200	1600
Pentadecan-1-ol	–	7200	3600

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Halogenated compounds</b>			
<i>Monohalogenated compounds</i>			
Chloromethane	-24	U	–
Chloroethene (Vinyl chloride)	-13	U	–
2-Chloroethyl vinyl ether	109	0.1	0.05
Chlorobenzene	132	1	0.5
<i>Dihalogenated compounds</i>			
1,1-Dichloroethane	57	U	–
1,2-Dichloroethane	83	U	–
1,1-Dichloroethene	32	U	–
1,2-Dichloroethene	49	U	–
1,2-Dichloropropane	95	U	–
1,3-Dichloropropene	104–112	U	–
1,2-Dichlorobenzene	181	6	3
1,3-Dichlorobenzene	173	6	3
1,4-Dichlorobenzene	174	6	3
<i>Trihalogenated compounds</i>			
Bromodichloromethane	90	U	–
Dibromochloromethane	119	U	–
Chlorodifluoromethane	-41	U	–
Trichloromethane	61	U	–
Tribromomethane	149	0.1	0.5
1,1,1-Trichloroethane	74	U	–
1,1,2-Trichloroethane	110	U	–
Trichloroethene	87	U	–
<i>Tetrahalogenated compounds</i>			
Dichlorodifluoromethane	-30	U	–
Trichlorofluoromethane	24	U	–
Tetrachloromethane	76	U	–
1,1,2,2-Tetrachloroethane	146	0.2	0.1
Tetrachloroethene	121	0.4	0.02

### 2. On 200 mg Tenax TA at 20 °C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Methane	-164	U	–
Ethane	-89	U	–
Propane	-42	U (0.03)	–
Butane	-0.5	0.2	0.1
Pentane	36	1	0.5
Hexane	69	6.4	3.2
Heptane	98	34	17
Octane	125	160	80
Nonane	151	1400	700
Decane	174	4200	2100
Undecane	196	25000	12500
Dodecane	216	126000	63000

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes (continued)</b>			
Tridecane	235	25000	12500
Tetradecane	254	60000	30000
Pentadecane	271	114000	57000
Hexadecane	287	210000	105000
Heptadecane	302	340000	170000
Octadecane	316	560000	280000
Nonadecane	330	800000	400000
Eicosane	343	1.32 × 10 <sup>6</sup>	660000
n-C <sub>21</sub>	357	2.0 × 10 <sup>6</sup>	1.0 × 10 <sup>6</sup>
n-C <sub>22</sub>	369	3.2 × 10 <sup>6</sup>	1.6 × 10 <sup>6</sup>
n-C <sub>23</sub>	380	5.0 × 10 <sup>6</sup>	2.5 × 10 <sup>6</sup>
n-C <sub>24</sub>	391	8.2 × 10 <sup>6</sup>	4.1 × 10 <sup>6</sup>
n-C <sub>25</sub>	402	1.25 × 10 <sup>8</sup>	6.25 × 10 <sup>6</sup>
n-C <sub>26</sub>	412	2.5 × 10 <sup>7</sup>	1.25 × 10 <sup>6</sup>
n-C <sub>27</sub>	442	5.0 × 10 <sup>7</sup>	2.5 × 10 <sup>6</sup>
n-C <sub>28</sub>	432	1.04 × 10 <sup>8</sup>	5.2 × 10 <sup>6</sup>
n-C <sub>29</sub>	441	2.2 × 10 <sup>8</sup>	1.1 × 10 <sup>8</sup>
n-C <sub>30</sub>	450	5.0 × 10 <sup>8</sup>	2.5 × 10 <sup>8</sup>
<b>Aromatics</b>			
Benzene	80	12.5	6
Toluene	111	76	38
Xylene	138-144	610	305
Ethylbenzene	137	360	180
n-Propylbenzene	159	1700	850
Isopropylbenzene	152	960	480
n-Butylbenzene	183	3600	1800
n-Hexylbenzene	226	87200	43600
Ethyltoluene	162	2000	1000
Styrene	145	600	300
Methylstyrene	167	2400	1200
Nitrobenzene	211	28000	14000
p-Cymene	177	2600	1300
Butylated hydroxy toluene (BHT)	265	2700	1350
Trimethylbenzene	165-176	3600	1800
o-Cresol	191	6600	3300
p-Cresol	202	6200	3100
o-Ethyltoluene	164	2000	1000
m-Ethyltoluene	161	2000	1000
p-Ethyltoluene	162	2000	1000
Naphthalene	218	20000	10000
Biphenyl	256	63200	31600
Phenanthrene	340	112400	56200

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Terpenes and alkaloids</b>			
Limonene	176	2400	1200
Nicotine	247	20000	10000
Geraniol	229	126000	63000
Cotinine	—	280000	140000
<b>Alkenes</b>			
But-1-ene	-6	0.1	0.05
Pent-1-ene	30	0.9	0.45
Pent-2-ene	37	1.2	0.6
Hex-1-ene	63	4.6	2.3
Hex-2-ene	68	6.4	3.2
Hex-3-ene	66	5.6	2.8
Hept-1-ene	94	24	12
Oct-1-ene	121	112	56
Oct-2-ene	123	164	82
Oct-3-ene	121	130	65
Non-1-ene	—	640	320
Dec-1-ene	170	1260	630
Undec-1-ene	—	2820	1410
Dodec-1-ene	214	5600	2800
Tridec-1-ene	—	8400	4200
Tetradec-1-ene	—	12600	6300
Pentadec-1-ene	—	20000	10000
Hexadec-1-ene	—	31600	15800
Heptadec-1-ene	—	50200	25100
Octadec-1-ene	—	80000	40000
Nonadec-1-ene	—	200000	100000
Eicos-1-ene	—	500000	250000
<b>Esters</b>			
Methyl acetate	57	1.6	0.8
Ethyl acetate	71	54	27
Propyl acetate	102	36	18
Isopropyl acetate	90	12	6
n-Butyl acetate	126	170	85
Isobutyl acetate	115	265	132.5
tert-Butyl acetate	98	176	88
Pentyl acetate	149	940	470
Hexyl acetate	172	9000	4500
Methyl acrylate	81	13	6.5
Ethyl acrylate	100	48	24
Methyl methacrylate	100	55	27.5

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Glycol ethers</b>			
Methoxyethanol	125	6	3
Ethoxyethanol	136	10	5
Butoxyethanol	170	70	35
2-(2-Butoxyethoxy)-ethanol	—	50200	25100
Methoxypropanol	118	27	13.5
Methoxyethyl acetate	145	16	8
Ethoxyethyl acetate	156	30	15
Butoxyethyl acetate	92	300	150
<b>Aldehydes</b>			
Acetaldehyde	20	0.1	0.05
Propanal	48	1	0.5
2-Methylpropanal	63	3.4	1.7
Butanal	75	6	3
3-Methylbutanal	91	13	6.5
Pentanal	103	22	11
Hexanal	131	100	50
Heptanal	153	700	350
Octanal	169	2500	1250
Nonanal	191	12600	6300
2-Methylcyclohexane carboxaldehyde	—	500	250
3-Methylcyclohexane carboxaldehyde	—	600	300
4-Methylcyclohexane carboxaldehyde	—	800	400
Furfural	162	600	300
<b>Ketones</b>			
Acetone	56	1.2	0.6
Butan-2-one (Methyl ethyl ketone)	80	6.4	3.2
Pentan-2-one	102	36	18
Pentan-3-one	102	40	20
4-Methylpentan-2-one (Methyl isobutyl ketone)	116	52	26
Hexan-2-one	128	200	100
Heptan-2-one	151	1000	500
Heptan-3-one	147	800	400
Heptan-4-one	144	700	350
Octan-2-one	173	4000	2000
Octan-3-one	167	1160	580
Nonan-2-one	195	10000	5000
Nonan-4-one	187	6400	3200
Cyclohexanone	155	340	170
3,5,5-Trimethylcyclohex-2-en-1-one	214	11200	5600
Benzaldehyde	179	220	110

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alcohols</b>			
Methanol	65	U	—
Ethanol	78	0.4	0.2
Propan-1-ol	97	2.2	1.1
Propan-2-ol	82	1	0.5
2-Methylpropan-2-ol	83	10	5
Butan-1-ol	118	10	5
Butan-2-ol	108	5.6	2.8
2-Methylbutan-1-ol	129	16	8
3-Methylbutan-1-ol	130	32	16
2-Methylbutan-2-ol	101	0.2	0.1
Dimethylaminoethanol	133	36	18
Pentan-1-ol	138	76	38
Pentan-2-ol	119	34	17
Pentan-3-ol	115	9	4.5
Hexan-1-ol	158	360	180
Heptan-1-ol	176	2000	1000
Octan-1-ol	195	2800	1400
Nonan-1-ol	213	17800	8900
Decan-1-ol	231	56200	28100
Undecan-1-ol	243	141600	70800
Dodecan-1-ol	259	280000	140000
Tridecan-1-ol	—	400000	200000
Tetradecan-1-ol	289	640000	320000
Pentadecan-1-ol	—	$1.12 \times 10^6$	560000
Hexadecan-1-ol	—	$1.80 \times 10^6$	900000
Heptadecan-1-ol	—	$3.56 \times 10^6$	$1.78 \times 10^6$
Octadecan-1-ol	—	$7 \times 10^6$	$3.5 \times 10^6$
Nonadecan-1-ol	—	$2 \times 10^7$	$1 \times 10^7$
Eicosan-1-ol	—	$5.62 \times 10^7$	$2.81 \times 10^6$
Docosan-1-ol	—	$4 \times 10^8$	$2 \times 10^8$
Phenol	182	480	240
o-Cresol	191	6600	3300
p-Cresol	201	6200	3100
<b>Glycols</b>			
Ethylene glycol	197	7	3.5
Propane-1,2-diol	186	20	10
Butane-1,3-diol	203	157	78.5
(2-Hydroxyethoxy)-ethan-2-ol (Diethylene glycol)	245	1260	630

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Halogenated compounds</b>			
<i>Monohalogenated compounds</i>			
Chloromethane	-24	U	–
Chloroethane	12	0.8	0.4
Chloroethene (Vinyl chloride)	-13	0.1	0.05
2-Chloroethyl vinyl ether	109	40	20
Chlorobenzene	131	52	26
<i>Dihalogenated compounds</i>			
Dichloromethane	40	0.9	0.45
1,1-Dichloroethane	57	2.2	1.1
1,2-Dichloroethane	84	10.8	5.4
1,1-Dichloroethene	32	0.8	0.4
1,2-Dichloroethene	49	1.8	0.9
1,2-Dichloropropane	95	16	8
1,3-Dichloropropene	104–112	30	15
1,2-Dichlorobenzene	181	660	330
1,3-Dichlorobenzene	173	540	270
1,4-Dichlorobenzene	174	58	290
<i>Trihalogenated compounds</i>			
Bromodichloromethane	90	18	9
Dibromochloromethane	119	63	31.5
Chlorodifluoromethane	-41	U	–
Trichloromethane (Chloroform)	61	3.8	1.9
Tribromomethane (Bromoform)	149	200	100
1,1,1-Trichloroethane	74	2.2	1.1
1,1,2-Trichloroethane	114	68	34
Trichloroethene	87	11	5.5
<i>Tetrahalogenated compounds</i>			
Dichlorodifluoromethane	-30	U	–
Trichlorofluoromethane	24	0.4	0.2
Tetrachloromethane	76	12.4	6.2
Tetrachloromethane	76	4.2	2.1
1,1,1,2-Tetrachloroethane	130	156	78
1,1,1,2-Tetrachloroethane	146	340	170
Tetrachloroethene	121	96	48
<b>Acids</b>			
Acetic acid	118	1.1	0.55
Propanoic acid	141	6.4	3.2
Butanoic acid	164	28	14
Pentanoic acid	186	140	70
Hexanoic acid	205	620	310
<b>Acid anhydrides</b>			
Maleic anhydride	202	176	88

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Amines</b>			
Methylamine	-6	0.1	0.05
Ethylamine	17	0.9	0.45
n-Propylamine	48	3.2	1.6
n-Butylamine	77	13	6.5
sec-Butylamine	83	4.6	2.3
Isobutylamine	68	5	2.5
tert-Butylamine	44	U	–
n-Pentylamine	104	52	26
Isopentylamine	95	26	13
n-Hexylamine	129	200	100
Cyclohexylamine	133	180	90
n-Heptylamine	154	1000	500
n-Octylamine	178	3600	1800
n-Nonylamine	202	9000	4500
n-Decylamine	216	28200	14100
Pyridine	116	16	8
Aniline	183	440	220
Benzylamine	184	1000	500

### 3. On 200 mg Tenax GR at 20 °C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Methane	-164	U	–
Ethane	-89	U	–
Propane	-42	U	–
Butane	-0.5	U (0.08)	–
Pentane	36	0.34	0.17
Hexane	69	1.3	0.65
Heptane	98	6	3
Octane	125	18	9
Nonane	151	50	25
Decane	174	240	120
Dodecane	216	820	410
Tetradecane	254	2400	1200
Hexadecane	287	9400	4700
Octadecane	316	40000	20000
Eicosane	343	170000	85000
<b>Aldehydes</b>			
Propanal	48	0.6	0.3
2-Methylpropanal	63	1.2	0.6
Butanal	75	1.9	0.95
3-Methylbutanal	91	4.4	2.2
Pentanal	103	5	2.5
Hexanal	131	14	7
Heptanal	153	56	28
Octanal	169	120	60
Nonanal	191	160	160

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Ketones</b>			
Acetone	56	0.7	0.35
Butan-2-one (Methyl ethyl ketone)	80	2.4	1.2
Pentan-2-one	102	4.8	2.4
Pentan-3-one	102	5	2.5
Hexan-2-one	128	14	7
2-Methylcyclohexanone	–	24	12
3-Methylcyclohexanone	–	12	6
4-Methylcyclohexanone	–	13	6.5
Heptan-2-one	151	40	20
Heptan-3-one	147	39	19.5
Heptan-4-one	144	32	16
Octan-2-one	173	100	50
Octan-3-one	167	42	21
Nonan-2-one	195	300	150
Nonan-5-one	186	280	140
<b>Alcohols</b>			
Methanol	65	U	–
Ethanol	78	U (0.16)	–
Propan-1-ol	97	0.6	0.16
propan-2-ol	82	0.32	0.3
Butan-1-ol	118	2	0.55
Butan-2-ol	98	1.1	1
Pentan-1-ol	138	4.6	1.6
Pentan-2-ol	119	3.2	2.3
Hexan-1-ol	158	20	10
Heptan-1-ol	180	46	23
Octan-1-ol	195	100	50
Nonan-1-ol	213	170	85
Decan-1-ol	231	400	200
Undecan-1-ol	243	600	300
Dodecan-1-ol	259	1200	600
Tridecan-1-ol	–	2000	1000
Tetradecan-1-ol	289	3600	1800
Pentadecan-1-ol	–	4400	2200

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Halogenated compounds</b>			
<i>Monohalogenated compounds</i>			
Chloromethane	–24	U	–
Chloroethane	12	0.1	0.05
Chloroethene (Vinyl chloride)	–13	U	–
2-Chloroethyl vinyl ether	109	6	3
Chlorobenzene	132	20	10
<i>Dihalogenated compounds</i>			
Dichloromethane	40	0.5	0.25
1,1-Dichloroethane	57	0.8	0.4
1,2-Dichloroethane	83	2.4	1.2
1,1-Dichloroethene	32	0.5	0.25
1,2-Dichloroethene	49	0.6	0.3
1,2-Dichloropropane	95	2.2	1.1
1,3-Dichloropropene	104–112	4	2
1,2-Dichlorobenzene	181	160	80
1,3-Dichlorobenzene	173	120	60
1,4-Dichlorobenzene	174	140	70
<i>Trihalogenated compounds</i>			
Bromodichloromethane	90	2	1
Dibromochloromethane	119	10	5
Chlorodifluoromethane	–41	U	–
Trichloromethane (Chloroform)	61	1.2	0.6
Tribromomethane (Bromoform)	149	40	20
1,1,1-Trichloroethane	74	1.8	0.9
1,1,2-Trichloroethane	110	6	3
Trichloroethene	87	3.4	1.7
<i>Tetrahalogenated compounds</i>			
Dichlorodifluoromethane	–30	U	–
Trichlorofluoromethane	–	0.1	0.05
Tetrachloromethane	76	2	1
1,1,2,2-Tetrachloroethane	146	38	19
Tetrachloroethene	121	7.6	3.8
<b>Amines</b>			
n-Butylamine	77	2.4	1.2
sec-Butylamine	83	1	0.5
Isobutylamine	68	1.3	0.65
tert-Butylamine	44	U	–
n-Pentylamine	104	14	7
Isopentylamine	95	4	2
n-Hexylamine	129	40	20
Cyclohexylamine	133	5	2.5
n-Heptylamine	154	70	35
n-Octylamine	178	170	85
n-Decylamine	216	400	200
Benzylamine	184	70	35

**4. On 200 mg Carbotrap/Carbograph 1TD/Carbopack B at 20 °C**

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Methane	-164	U	—
Ethane	-89	U	—
Propane	-42	U	—
Butane	-0.5	U	—
Pentane	36	1.2	0.6
Hexane	69	16	8
Heptane	98	90	45
Octane	125	1500	750
Nonane	151	14000	7000
Decane	174	200000	100000
Undecane	96	600000	300000
Dodecane	216	1.6 × 10 <sup>6</sup>	800000
<b>Alcohols</b>			
Methanol	65	U	—
Ethanol	78	U	—
Propan-1-ol	97	0.64	0.32
Butan-1-ol	118	6	3
Pentan-1-ol	138	20	10
Hexan-1-ol	158	60	30
Heptan-1-ol	176	180	90
Octan-1-ol	180	600	300
Nonan-1-ol	213	3000	1500
Decan-1-ol	231	6000	3000
<b>Halogenated compounds</b>			
<i>Monohalogenated compounds</i>			
Chloromethane	-24	U	—
Chloroethene (Vinyl chloride)	-13	U	—
2-Chloroethyl vinyl ether	109	2.6	1.3
Chlorobenzene	132	20	10
<i>Dihalogenated compounds</i>			
Dichloromethane	40	0.04	0.02
1,1-Dichloroethane	57	0.2	0.1
1,2-Dichloroethane	83	0.3	0.15
1,1-Dichloroethene	32	0.1	0.05
1,2-Dichloroethene	49	0.2	0.1
1,2-Dichloropropane	95	1.2	0.6
1,3-Dichloropropene	104-112	1.3	0.65
1,2-Dichlorobenzene	181	102	51
1,3-Dichlorobenzene	173	102	51
1,4-Dichlorobenzene	174	102	51

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Halogenated compounds (continued)</b>			
<i>Trihalogenated compounds</i>			
Dibromochloromethane	119	1.5	0.75
Trichloromethane (Chloroform)	61	0.2	0.1
Tribromomethane (Bromoform)	149	3.8	1.9
1,1,1-Trichloroethane	74	0.6	0.3
1,1,2-Trichloroethane	110	1.6	0.8
Trichloroethene	87	1.6	0.8
<i>Tetrahalogenated compounds</i>			
Dichlorodifluoromethane	-30	U	—

**5. On 300 mg Carbopack X**

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
Carbon disulfide <sup>a</sup>	46	6.5	3
Buta-1,3-diene	-4	>25	>17
Benzene	80	10800	5400

<sup>a</sup> At 30 °C, 80% relative humidity.

**6. On 300 mg Chromosorb 106**

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Pentane	35	11.2	56
Hexane	69	60	30
Heptane	98	325	162
Octane	125	2076	1038
Nonane	151	14000	7000
Decane	174	74000	37000
<b>Aromatics</b>			
Benzene	80	53	26
Toluene	111	165	82
Xylene	138-144	1554	777
Ethylbenzene	137	730	365
Trimethylbenzene	165-176	5650	2825
<b>Esters</b>			
Methyl acetate	58	5.2	2.6
Ethyl acetate	71	39	19.5
Propyl acetate	102	297	148.5
Isopropyl acetate	90	147	73.5
n-Butyl acetate	126	1460	730
Isobutyl acetate	115	880	440
tert-Butyl acetate	98	327	163.5

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Glycol ethers</b>			
Methoxyethanol	125	9.6	4.8
Ethoxyethanol	136	150	75
Methoxyethyl acetate	145	1720	860
Ethoxyethyl acetate	156	8100	4050
<b>Ketones</b>			
Acetone	56	2.9	1.45
Butan-2-one (Methyl ethyl ketone)	80	21	10.5
4-Methylpentan-2-one (Methyl isobutyl ketone)	118	490	245
<b>Alcohols</b>			
Ethanol	78	2.4	1.2
Propan-1-ol	97	17	8.5
Propan-2-ol	82	9	4.5
Butan-1-ol	118	96	48
Butan-2-ol	108	60	30
<b>Halogenated compounds</b>			
Tetrachloromethane	76	44	22
1,2-Dichloroethane	84	34	17
1,1,1-Trichloroethane	74	17	8.5

### 7. On 500 mg PoraPak N

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Pentane	35	8.2	4.1
Hexane	69	32	16
Heptane	98	90	45
<b>Aromatics</b>			
Benzene	80	52	26
Pyridine	116	390	195
<b>Ketones</b>			
Butan-2-one (Methyl ethyl ketone)	80	95	47.5
<b>Alcohols</b>			
Ethanol	78	7.5	3.75
Propan-1-ol	97	40	20
Butan-1-ol	118	10	5
Butan-2-ol	108	5.6	2.8
Octan-1-ol	180	2800	1400
Phenol	182	480	240
<b>Acids</b>			
Acetic acid	116	97	48.5
<b>Nitriles</b>			
Acetonitrile	82	7	3.5
Acrylonitrile	82	7	3.5
Propionitrile	97	23	11.5

### 8. On 300 mg UniCarb

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
Butane	-0.5	1640	820
Pentane	35	63000	31500
Hexane	69	3.9 × 10 <sup>6</sup>	1.95 × 10 <sup>6</sup>
Benzene	80	1 × 10 <sup>6</sup>	500000
Dichloromethane	40	395	197
1,1,1-Trichloroethane	74	17600	8800
Methanol	65	264	132
Ethanol	78	6900	3450
Carbon disulfide <sup>a</sup>	46	26	13

<sup>a</sup> At 30 °C, 80% relative humidity.

### 9. On 600 mg Carbosieve SIII at 20 °C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Alkanes</b>			
Methane	-164	U (0.05)	—
Ethane	-89	0.7	—
Propane	-42	5.1	2.55
Butane	-0.5	37.8	18.9
Pentane	36	360	180
Hexane	69	3000	1500
Heptane	98	6000	3000

### 10. On 500 mg Carboxen 569 at 20 °C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Aromatics</b>			
Benzene	80	42.5	21.25
Toluene	111	1350	675
Ethylbenzene	137	1250	625
Isopropylbenzene	152	3750	1875
Trimethylbenzene	165-176	3750	1875
Xylene	138-144	5500	2750
Benzaldehyde	179	6000	3000
<i>p</i> -Cymene	177	10000	5000
Propylbenzene	159	11000	5500
<b>Terpenes</b>			
Limonene	176	8000	4000
Terpinene	182	10000	5000



Compound	b.p. (°C)	Retention volume (L)	SSV (L)
<b>Aldehydes</b>			
Propanal	48	1.2	0.6
Butanal	75	2.6	1.3
2-Methylpropanal	63	4.25	2.125
Pentanal	103	35	17.5
Hexanal	131	600	300
Heptanal	153	10500	5250
Octanal	169	30000	15000
Nonanal	191	90000	45000
<b>Ketones</b>			
Acetone	56	3.5	1.75
Butan-2-one (Methyl ethyl ketone)	80	2	1
Pentan-2-one	102	10.5	5.25
Pentan-3-one	102	20	10
Hexan-2-one	128	400	200
2-Methylcyclohexanone	—	65000	32500
3-Methylcyclohexanone	—	400	200
4-Methylcyclohexanone	—	400	200
Heptan-2-one	151	60000	30000
Heptan-3-one	147	60000	30000
Heptan-4-one	144	60000	30000
Octan-3-one	167	35000	17500
<b>Halogenated compounds</b>			
Tribromomethane (Bromoform)	149	1000	500
1,1,2,2-Tetrachloroethane	146	2000	1000
<b>Amines</b>			
n-Butylamine	77	5.5	2.75
sec-Butylamine	83	5.2	2.6
Isobutylamine	68	4.1	2.05
tert-Butylamine	44	1.6	0.8
n-Pentylamine	104	95	47.5
Isopentylamine	95	10.5	5.25
n-Hexylamine	129	3 × 10 <sup>6</sup>	1.5 × 10 <sup>6</sup>
Benzylamine	184	60000	30000
<b>Other</b>			
Carbon disulfide <sup>a</sup>	46	72	36

<sup>a</sup> At 30°C, 80% relative humidity.

### 11. On 350 mg Carbograph 5TD at 20°C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
Buta-1,3-diene	-4	~10	~7

### 12. On 600 mg Molecular sieve 13X at 20°C

Compound	b.p. (°C)	Retention volume (L)	SSV (L)
Buta-1,3-diene	-4	>100	>70

### 13. Breakthrough volumes for water

It is important to use at least twice the volumes below when purging water from the sorbent bed, to ensure complete elution. It is also important that the water is kept in the gas phase.

Sorbent	Temp. (°C)	Breakthrough volume	
		Manufacturers' data (mL/g) <sup>a</sup>	Extrapolated to typical mass of sorbent in 1/4" stainless steel tube (mL)
Carbotrap C	20	25	5
Tenax TA	20	65	13
	40	—	7 <sup>a</sup>
	60	—	3.6 <sup>a</sup>
Tenax GR	20	47	9.4
Carbotrap	20	39	7.8
Carboxen 569	20	60	12
Carbosieve SIII	20	320	64

<sup>a</sup> For 200 mg Tenax TA.

### Trademarks

UniCarb™ is a trademark of Markes International Ltd, UK.

Carbograph™ is a trademark of LARA s.r.l., Italy.

Carbopack™ is a trademark of Supelco Inc., USA.

Carbosieve™ is a trademark of Supelco Inc., USA.

Carbotrap® is a registered trademark of Sigma-Aldrich Co. LLC, USA.

Carboxen™ is a trademark of Supelco Inc., USA.

Chromosorb® is a registered trademark of Manville Corporation, USA.

PoraPak™ is a trademark of Waters Associates Inc., USA.

Version 8  
June 2012