



## USEPA Method 8260C using the Tekmar Atomx Automated VOC Sample Prep System and a Thermo Focus/DSQ™ II GCMS

### Application Note

By Roger Bardsley

#### **Abstract**

The United States Environmental Protection Agency (USEPA) developed Method 8260 to determine the concentration of Volatile Organic Compounds (VOCs) in water and soil matrices. USEPA Methods 5030 and 5035<sup>1</sup> are used to prepare water, soil and waste samples for analysis by USEPA Method 8260. USEPA Method 8260 was recently updated from 8260B to 8260C<sup>2</sup>. In this study water and soil sample analysis will be performed using USEPA Method 8260C<sup>3</sup> guidelines.

A calibration curve and Method Detection Limits (MDLs) for a target compound list will be performed on the Atomx, an Automated VOC Sample Prep System in conjunction with a Thermo Focus/DSQ II Gas Chromatography-Mass Spectrometry (GC/MS). A 5mL purge volume will be utilized for the water study, while an in-vial purge will be utilized for the soil samples. All conditions outlined in Method 8260C will be followed.

#### **Introduction**

Teledyne Tekmar developed the Atomx, an Automated VOC Sample Prep System that integrates both a purge and trap concentrator with an 80-position autosampler to allow for sampling of multiple matrices such as water, soil and methanol extraction with a single instrument. This "all-in-one" instrument allows for increased throughput and efficiency through the features that it provides.

The Atomx utilizes the #9 U-shape trap and three standard addition vessels. The Atomx, in conjunction with a Thermo Focus GC with the optional Purge and Trap adapter inlet and a Thermo DSQ II MS is an exceptional instrument to meet the new requirements for the detection of the target analytes by USEPA Method 8260C.

In this study, a calibration curve and MDLs were preformed for both water and soil matrices. The water samples employed a calibration curve from 0.5ppb to 200ppb, while the soil samples used a calibration curve from 1.0ppb to 200ppb. A 5mL purge volume was used for the water samples, and the soil samples utilized 5g of sample with 10mL of reagent water added by the Atomx. The conditions that are outlined in USEPA Method 8260C were followed for both water and soil matrices.



Image 1: Atomx, Automated VOC Sample Prep System



## Experimental-Instrument Conditions

The instrument conditions for both the waters and soils application developments and results are listed in Tables 1 through 4.

| GC Parameters |  | MS Parameters         |                     |
|---------------|--|-----------------------|---------------------|
| GC            | Thermo Scientific Focus GC   | MSD                   | Thermo DSQ II       |
| Column        | Restek Rtx®-VMS 20m x 0.18mmID x 1um   | Source                | 230°C               |
| Oven Program  | 40°C for 4 min; 16°C/min to 100°C for 0 min; 25°C /min to 200°C for 3 min, 14.75 min runtime | MS Transfer Line Temp | 230°C               |
| Inlet         | 220°C  | Solvent Delay         | 0.5 min             |
| Column Flow   | 1.0mL/min  | Scan Range            | 35.0 m/z to 270 m/z |
| Gas           | Helium   | Scans                 | 5.5586 scans/sec    |
| Split         | 50mL/min   | Chrom Filter          | 2.5 Sec             |
| Inlet         | Thermo optional Purge and Trap adapter   |                       |                     |

Tables 1 and 2: Thermo Scientific GC and MS Parameters

| Atomx Water Parameters   |           |                                 |                 |
|--------------------------|-----------|---------------------------------|-----------------|
| Variable                 | Value     | Variable                        | Value           |
| Valve Oven Temp          | 140°C     | Dry Purge Flow                  | 100mL/min       |
| Transfer Line Temp       | 140°C     | Dry Purge Temp                  | 20°C            |
| Sample Mount Temp        | 90°C      | Methanol Needle Rinse           | Off             |
| Water Heater Temp        | 90°C      | Methanol Needle Rinse Volume    | 3.0mL           |
| Sample Vial Temp         | 20°C      | Water Needle Rinse Volume       | 7.0mL           |
| Sample Equilibrate Time  | 0.00 min  | Sweep Needle Time               | 0.25min         |
| Soil Valve Temp          | 100°C     | Desorb Preheat Time             | 245°C           |
| Standby Flow             | 10mL/min  | GC Start Signal                 | Start of Desorb |
| Purge Ready Temp         | 40°C      | Desorb Time                     | 2.00 min        |
| Condensate Ready Temp    | 45°C      | Drain Flow                      | 300mL/min       |
| Presweep Time            | 0.25 min  | Desorb Temp                     | 250°C           |
| Prime Sample Fill Volume | 3.0mL     | Methanol Glass Rinse            | On              |
| Sample Volume            | 5.0mL     | Number of Methanol Glass Rinses | 1               |
| Sweep Sample Time        | 0.25 min  | Methanol Glass Rinse Volume     | 3.0mL           |
| Sweep Sample Flow        | 100mL/min | Number of Bake Rinses           | 1               |
| Sparge Vessel Heater     | Off       | Water Bake Rinse Volume         | 7.0mL           |
| Sparge Vessel Temp       | 20°C      | Bake Rinse Sweep Time           | 0.25 min        |
| Prepurge Time            | 0.00 min  | Bake Rinse Sweep Flow           | 100mL/min       |
| Prepurge Flow            | 0mL/min   | Bake Rinse Drain Time           | 0.40 min        |
| Purge Time               | 11.00 min | Bake Time                       | 2.00 min        |
| Purge Flow               | 40mL/min  | Bake Flow                       | 200mL/min       |
| Purge Temp               | 20°C      | Bake Temp                       | 280°C           |
| Condensate Purge Temp    | 20°C      | Condensate Bake Temp            | 200°C           |
| Dry Purge Time           | 0.50 min  |                                 |                 |

Table 3: Atomx Water Parameters (parameters highlighted in yellow were not used.)



| Atomx Soil Parameters |           |                              |                 |
|-----------------------|-----------|------------------------------|-----------------|
| Variable              | Value     | Variable                     | Value           |
| Valve oven Temp       | 150°C     | Purge Time                   | 11.00 min       |
| Transfer Line Temp    | 150°C     | Purge Flow                   | 40mL/min        |
| Sample Mount Temp     | 90°C      | Purge Temp                   | 20°C            |
| Water Heater Temp     | 90°C      | Condensate Purge Temp        | 20°C            |
| Sample Vial Temp      | 40°C      | Dry Purge Time               | 2.00 min        |
| Prepurge Time         | 0.00 min  | Dry Purge Flow               | 100mL/min       |
| Prepurge Flow         | 0mL/min   | Dry Purge Temp               | 20°C            |
| Preheat Mix Speed     | Slow      | Methanol Needle Rinse        | Off             |
| Sample Preheat Time   | 2.00 min  | Methanol Needle Rinse Volume | 3.0mL           |
| Soil Valve Temp       | 120°C     | Water Needle Rinse Volume    | 10.0mL          |
| Standby Flow          | 20mL/min  | Sweep Needle Time            | 0.50 min        |
| Purge Ready Temp      | 40°C      | Desorb Preheat Time          | 245°C           |
| Condensate Ready Temp | 45°C      | GC Start Signal              | Start of Desorb |
| Presweep Time         | 0.25 min  | Desorb Time                  | 2.00 min        |
| Water Volume          | 10mL      | Drain Flow                   | 300mL/min       |
| Sweep Water Time      | 0.25 min  | Desorb Temp                  | 250°C           |
| Sweep Water Flow      | 100mL/min | Bake Time                    | 2.00 min        |
| Sparge Vessel Heater  | Off       | Bake Flow                    | 400mL/min       |
| Sparge Vessel Temp    | 20°C      | Bake Temp                    | 280°C           |
| Purge Mix Speed       | Medium    | Condensate Bake Temp         | 200°C           |

Table 4: Atomx Soil Parameters (parameters highlighted in yellow were not used.)

### Calibration and Sample Preparation

A stock standard solution was prepared at 50ppm ( $\mu\text{g/mL}$ ) utilizing 6 Restek Stock standards to provide 94 compounds related to USEPA method 8260C and other state environmental methods. The Restek standards include those listed in Table 5. The amounts listed in Table 5 were transferred to a 10mL volumetric flask containing P&T methanol. The flask was diluted to volume with P&T methanol.

| Cat#  | Name                          | Concentration                  | Amount            | Vol. | Final Conc. |
|-------|-------------------------------|--------------------------------|-------------------|------|-------------|
| 30633 | 8260B MegaMix®                | 2000 $\mu\text{g/mL}$          | 250 $\mu\text{L}$ | 10mL | 50ppm       |
| 30489 | 8260B Acetate Mix             | 2000 $\mu\text{g/mL}$          | 250 $\mu\text{L}$ | 10mL | 50ppm       |
| 30465 | California Oxygenates Mix     | 2000 – 10,000 $\mu\text{g/mL}$ | 250 $\mu\text{L}$ | 10mL | 50ppm       |
| 30042 | 502.2 Calibration Mix (Gases) | 2000 $\mu\text{g/mL}$          | 250 $\mu\text{L}$ | 10mL | 50ppm       |
| 30265 | 2-Chloroethyl Vinyl Ether     | 2000 $\mu\text{g/mL}$          | 250 $\mu\text{L}$ | 10mL | 50ppm       |
| 30006 | VOA Calibration Mix (Ketones) | 5000 $\mu\text{g/mL}$          | 100 $\mu\text{L}$ | 10mL | 50ppm       |

Table 5: Stock Standard Solution Preparation

Working standards were prepared ranging from 0.5ppb to 200ppb for the water analysis and 1ppb to 200ppb for the soil analysis by diluting the stock standard to a final volume in reagent water. This method of standard preparation creates varying amounts of methanol in the samples to ensure that the higher methanol concentration does not interfere with the mass spectrometer.

A 0.5ppb solution for the water analysis and a 1.0ppb solution for the soil analysis were used to calculate the minimum detectable limits of the compounds.

A 125ppm Internal Standard (IS) and Surrogate Standard (SS) solution was prepared by diluting 500 $\mu$ L of both the Restek 8260 Internal Standard Mix and 8260 Surrogate Mix to 10mL with P&T methanol. This was placed in the standard vessel on the Atomx. 2 $\mu$ L of the IS/SS solution was automatically added to each sample by the Atomx system.

### **Minimum Detection Limits**

The water analysis standard curve and MDL samples were prepared by filling a standard 40mL VOA vial with the appropriate standard. Seven samples at 0.5ppb were used for the calculation of the MDL for the water assay.

The soil analysis standard curve and MDL samples were prepared by transfer 5mL of the appropriate standard of sample to a 40mL VOA vial, Teflon stir bar added and the vial capped. Seven samples at 1.0ppb were used for the MDL calculation.

The samples were placed on the Atomx automated VOC instrument which purged either 5mL of the water sample in the standard purge and trap vessel or added 10mL of water to each soil sample along with IS/SS solution to perform an in vial purge and trap of the soil samples. The trapped VOC's were rapidly desorbed from the trap and transferred to the Focus/DSQ II GC/MS for analysis.

### **Results and Discussion**

Thermo Scientific EnviroLab Forms™ 3.0 GC/MS analysis and reporting software was used to process the calibration standards and MDL sample data. The relative response factors (RF) of all target analytes were evaluated for average RF and calibration %RSD. The minimum detection limit was calculated for the water and the soil samples. These results are presented in Table 6.

Example chromatograms of a 50ppb water standard and a 50ppb soil standard are presented in Figures 1 and 2.

| Compound                       | Water       |       |         |                  | Soil        |       |         |                  |
|--------------------------------|-------------|-------|---------|------------------|-------------|-------|---------|------------------|
|                                | Spike Level | MDL   | Avg. RF | Calibration %RSD | Spike Level | MDL   | Avg. RF | Calibration %RSD |
| Dichlorodifluoromethane        | 0.5         | 0.169 | 0.633   | 9.49             | 1           | 0.329 | 0.436   | 11.24            |
| Chloromethane                  | 0.5         | 0.107 | 2.348   | 7.60             | 1           | 0.121 | 1.975   | 9.24             |
| Vinyl Chloride                 | 0.5         | 0.061 | 0.750   | 5.57             | 1           | 0.142 | 0.684   | 9.44             |
| Bromomethane                   | 0.5         | 0.110 | 0.334   | 8.89             | 1           | 0.395 | 0.293   | 7.97             |
| Chloroethane                   | 0.5         | 0.228 | 0.236   | 5.86             | 1           | 0.187 | 0.205   | 11.87            |
| Trichloromonofluoromethane     | 0.5         | 0.150 | 0.432   | 11.57            | 1           | 0.086 | 0.414   | 9.14             |
| Diethyl Ether                  | 0.5         | 0.112 | 0.284   | 8.11             | 1           | 0.196 | 0.278   | 3.68             |
| 1,1-Dichloroethene             | 0.5         | 0.125 | 0.311   | 4.26             | 1           | 0.169 | 0.288   | 7.05             |
| Carbon disulfide               | 0.5         | 0.296 | 1.251   | 7.98             | 1           | 0.120 | 1.389   | 16.63            |
| 1,1,2-Trichlorotrifluoroethane | 0.5         | 0.149 | 0.169   | 8.19             | 1           | 0.104 | 0.170   | 10.13            |
| Iodomethane*                   | 0.5         | 0.054 | 0.410   | 13.58            | 1           | 0.035 |         | 0.9997           |
| Allyl chloride                 | 0.5         | 0.102 | 0.264   | 6.77             | 1           | 0.321 | 0.240   | 4.86             |
| Methylene Chloride             | 0.5         | 0.107 | 0.299   | 6.21             | 1           | 0.181 | 0.298   | 7.00             |
| Acetone*                       | 0.5         |       |         | 0.9987           | 1           |       |         | 0.9994           |



| Compound                         | Water       |        |         |                  | Soil        |       |         |                  |
|----------------------------------|-------------|--------|---------|------------------|-------------|-------|---------|------------------|
|                                  | Spike Level | MDL    | Avg. RF | Calibration %RSD | Spike Level | MDL   | Avg. RF | Calibration %RSD |
| trans-1,2-Dichloroethene         | 0.5         | 0.070  | 0.352   | 7.82             | 1           | 0.186 | 0.347   | 7.15             |
| Methyl Acetate                   | 0.5         | 0.450  | 0.831   | 5.23             | 1           | 0.154 | 0.570   | 1.77             |
| Methyl-tert-butyl Ether (MTBE)   | 0.5         | 0.65   | 1.298   | 7.03             | 1           | 0.136 | 1.244   | 3.59             |
| tert-Butyl Alcohol (TBA)         | 2.5         | 0.441  | 0.060   | 11.50            | 5           | 1.237 | 0.049   | 7.84             |
| Acetonitrile                     | 0.5         | 0.222  | 0.685   | 5.98             | 1           | 0.132 | 0.654   | 5.32             |
| Diisopropyl Ether (DIPE)         | 0.5         | 0.0435 | 4.083   | 5.16             | 1           | 0.096 | 3.925   | 4.45             |
| Chloroprene                      | 0.5         | 0.099  | 0.794   | 4.61             | 1           | 0.156 | 0.671   | 7.55             |
| 1,1-Dichloroethane               | 0.5         | 0.074  | 1.297   | 4.37             | 1           | 0.114 | 1.154   | 6.10             |
| Acrylonitrile                    | 0.5         | 0.124  | 0.264   | 5.97             | 1           | 0.150 | 0.168   | 2.84             |
| Ethyl-tert-butyl Ether (ETBE)    | 0.5         | 0.047  | 2.165   | 4.7              | 1           | 0.124 | 2.036   | 4.84             |
| Vinyl Acetate                    | 0.5         | 0.051  | 3.434   | 5.94             | 1           | 0.195 | 2.508   | 16.29            |
| cis-1,2-Dichloroethene           | 0.5         | 0.125  | 0.752   | 3.16             | 1           | 0.130 | 0.696   | 4.73             |
| 2,2-Dichloropropane              | 0.5         | 0.081  | 1.093   | 3.35             | 1           | 0.195 | 0.855   | 6.41             |
| Bromochloromethane               | 0.5         | 0.051  | 0.419   | 5.58             | 1           | 0.145 | 0.359   | 8.02             |
| Chloroform                       | 0.5         | 0.085  | 1.044   | 4.31             | 1           | 0.105 | 0.934   | 4.71             |
| Carbon Tetrachloride             | 0.5         | 0.093  | 0.799   | 6.25             | 1           | 0.161 | 0.636   | 6.52             |
| 1,1,1-Trichloroethane            | 0.5         | 0.076  | 0.894   | 4.06             | 1           | 0.109 | 0.744   | 7.05             |
| Dibromofluoromethane (Surrogate) | 25          | 2.264  | 0.605   | 4.27             | 25          | 1.683 | 0.624   | 2.19             |
| Methyl Acrylate                  | 0.5         | 0.098  | 0.672   | 6.84             | 1           | 0.168 | 0.509   | 11.74            |
| Tetrahydrofuran                  | 0.5         | 0.099  | 0.139   | 6.51             | 1           | 0.207 | 0.100   | 4.52             |
| Ethyl Acetate                    | 0.5         | 0.356  | 0.036   | 6.82             | 1           | 1.312 | 0.027   | 7.47             |
| 1,1-Dichloropropene              | 0.5         | 0.143  | 1.029   | 4.29             | 1           | 0.103 | 0.857   | 8.47             |
| 2-Butanone (MEK)                 | 0.5         | 0.097  | 0.149   | 8.52             | 1           | 0.449 | 0.108   | 16.21            |
| Benzene                          | 0.5         | 0.089  | 3.516   | 3.65             | 1           | 0.124 | 3.085   | 5.79             |
| Methacrylonitrile                | 0.5         | 0.132  | 1.026   | 7.00             | 1           | 0.199 | 0.727   | 2.69             |
| Pentafluorobenzene (IS)          | 25          |        |         | 2.48             | 25          |       |         | 3.43             |
| tert-Amyl Methyl Ether (TAME)    | 0.5         | 0.078  | 2.540   | 4.40             | 1           | 0.137 | 2.281   | 4.09             |
| 1,2-Dichloroethane               | 0.5         | 0.075  | 0.888   | 3.63             | 1           | 0.123 | 0.770   | 4.89             |
| Isopropyl Acetate                | 0.5         | 0.057  | 0.185   | 6.76             | 1           | 0.192 | 0.148   | 2.76             |
| Trichloroethylene                | 0.5         | 0.112  | 0.690   | 5.04             | 1           | 0.159 | 0.587   | 7.84             |
| Dibromomethane                   | 0.5         | 0.049  | 0.468   | 5.91             | 1           | 0.127 | 0.365   | 6.91             |
| 1,2-Dichloropropane              | 0.5         | 0.078  | 0.747   | 4.60             | 1           | 0.131 | 0.683   | 3.41             |
| Bromodichloromethane             | 0.5         | 0.075  | 0.767   | 2.99             | 1           | 0.129 | 0.699   | 3.72             |
| Methyl Methacrylate              | 0.5         | 0.085  | 0.680   | 4.86             | 1           | 0.157 | 0.523   | 4.44             |
| n-Propyl acetate                 | 0.5         | 0.128  | 2.491   | 6.23             | 1           | 0.133 | 1.968   | 6.97             |
| 2-Chloroethyl Vinyl Ether        | 0.5         | 0.100  | 0.297   | 7.37             | 1           | 0.102 | 0.208   | 7.77             |
| cis-1,3-Dichloropropene          | 0.5         | 0.055  | 0.728   | 5.75             | 1           | 0.140 | 0.596   | 6.17             |
| Toluene-d8 (Surrogate)           | 25          | 2.049  | 1.312   | 2.63             | 25          | 1.958 | 1.321   | 2.25             |
| Toluene                          | 0.5         | 0.079  | 0.964   | 3.10             | 1           | 0.113 | 0.801   | 8.18             |
| 2-Nitropropane*                  | 0.5         | 1.137  |         | 0.9970           | 1           | 2.652 |         | 0.9972           |
| Tetrachloroethylene              | 0.5         | 0.140  | 0.323   | 9.35             | 1           | 0.347 | 0.201   | 8.41             |
| 4-Methyl-2-pentanone             | 0.5         | 0.172  | 0.139   | 8.12             | 1           | 0.204 | 0.106   | 5.30             |
| trans-1,3-Dichloropropene        | 0.5         | 0.095  | 0.650   | 6.26             | 1           | 0.117 | 0.515   | 5.29             |
| 1,1,2-Trichloroethane            | 0.5         | 0.054  | 0.274   | 5.40             | 1           | 0.131 | 0.214   | 5.20             |
| Ethyl Methacrylate               | 0.5         | 0.068  | 0.666   | 7.35             | 1           | 0.108 | 0.514   | 5.20             |
| Dibromochloromethane             | 0.5         | 0.031  | 0.443   | 8.39             | 1           | 0.126 | 0.362   | 7.43             |
| 1,3-Dichloropropane              | 0.5         | 0.066  | 0.748   | 5.22             | 1           | 0.125 | 0.592   | 4.63             |
| 1,2-Dibromoethane (EDB)          | 0.5         | 0.047  | 0.413   | 5.67             | 1           | 0.114 | 0.303   | 6.86             |



| Compound                       | Water       |       |         |                  | Soil        |       |         |                  |
|--------------------------------|-------------|-------|---------|------------------|-------------|-------|---------|------------------|
|                                | Spike Level | MDL   | Avg. RF | Calibration %RSD | Spike Level | MDL   | Avg. RF | Calibration %RSD |
| Butyl Acetate                  | 0.5         | 0.078 | 1.485   | 9.58             | 1           | 0.093 | 1.104   | 6.63             |
| 2-Hexanone                     | 0.5         | 0.065 | 0.727   | 9.94             | 1           | 0.339 | 0.524   | 9.56             |
| Chlorobenzene-d5 (IS)          | 25          |       |         | 5.50             | 25          |       |         | 1.91             |
| Chlorobenzene                  | 0.5         | 0.069 | 1.131   | 4.04             | 1           | 0.142 | 0.927   | 5.53             |
| Ethylbenzene                   | 0.5         | 0.071 | 1.594   | 3.45             | 1           | 0.179 | 1.296   | 6.99             |
| 1,1,1,2-Tetrachloroethane      | 0.5         | 0.064 | 0.400   | 4.13             | 1           | 0.088 | 0.341   | 5.42             |
| m-,p-Xylene                    | 1           | 0.121 | 0.748   | 3.73             | 2           | 0.331 | 0.582   | 8.04             |
| o-Xylene                       | 0.5         | 0.062 | 0.759   | 3.94             | 1           | 0.119 | 0.599   | 5.40             |
| Bromoform                      | 0.5         | 0.101 | 0.205   | 7.88             | 1           | 0.150 | 0.156   | 6.56             |
| Styrene                        | 0.5         | 0.062 | 1.139   | 5.79             | 1           | 0.172 | 0.936   | 6.06             |
| Isopropylbenzene (Cumene)      | 0.5         | 0.073 | 1.679   | 3.92             | 1           | 0.206 | 1.288   | 7.69             |
| n-Amyl Acetate                 | 0.5         | 0.084 | 1.535   | 10.28            | 1           | 0.116 | 1.146   | 5.56             |
| Bromofluorobenzene (Surrogate) | 25          | 1.550 | 0.487   | 3.11             | 25          | 1.889 | 0.483   | 2.10             |
| Bromobenzene                   | 0.5         | 0.112 | 1.010   | 6.86             | 1           | 0.206 | 0.844   | 6.73             |
| cis-1,4-Dichloro-2-butene      | 0.5         | 0.150 | 0.371   | 5.34             | 1           | 0.403 | 0.595   | 6.57             |
| n-Propylbenzene                | 0.5         | 0.093 | 5.319   | 7.94             | 1           | 0.233 | 4.092   | 8.77             |
| 1,1,2,2-Tetrachloroethane      | 0.5         | 0.067 | 1.344   | 5.13             | 1           | 0.188 | 1.022   | 7.24             |
| 2-Chlorotoluene                | 0.5         | 0.079 | 3.484   | 4.29             | 1           | 0.196 | 2.669   | 6.57             |
| 1,2,3-Trichloropropane         | 0.5         | 0.116 | 1.552   | 5.35             | 1           | 0.237 | 1.098   | 5.52             |
| 1,3,5-Trimethylbenzene         | 0.5         | 0.101 | 4.779   | 4.96             | 1           | 0.206 | 3.546   | 7.02             |
| trans-1,4-dichloro-2-Butene    | 0.5         | 0.148 | 0.358   | 12.26            | 1           | 0.376 | 0.234   | 10.06            |
| 4-Chlorotoluene                | 0.5         | 0.134 | 3.547   | 6.25             | 1           | 0.226 | 2.732   | 6.80             |
| tert-Butylbenzene              | 0.5         | 0.095 | 4.404   | 4.64             | 1           | 0.189 | 3.613   | 6.21             |
| Pentachloroethane              | 0.5         | 0.287 | 0.292   | 24.61            | 1           | 0.583 | 0.400   | 10.10            |
| 1,2,4-Trimethylbenzene         | 0.5         | 0.090 | 4.794   | 5.52             | 1           | 0.178 | 3.504   | 6.63             |
| sec-Butylbenzene               | 0.5         | 0.092 | 5.315   | 7.28             | 1           | 0.225 | 4.065   | 8.20             |
| p-Isopropyltoluene (p-Cymene)  | 0.5         | 0.094 | 4.563   | 8.75             | 1           | 0.187 | 3.311   | 7.67             |
| 1,3-Dichlorobenzene            | 0.5         | 0.131 | 1.733   | 5.68             | 1           | 0.234 | 1.307   | 7.78             |
| 1,4-Dichlorobenzene-d4 (IS)    | 25          |       |         | 7.78             | 25          |       |         | 7.96             |
| 1,4-Dichlorobenzene            | 0.5         | 0.161 | 1.749   | 5.97             | 1           | 0.248 | 1.308   | 9.83             |
| n-Butylbenzene                 | 0.5         | 0.157 | 3.252   | 11.20            | 1           | 0.234 | 2.390   | 9.63             |
| 1,2-Dichlorobenzene            | 0.5         | 0.074 | 1.697   | 5.10             | 1           | 0.176 | 1.269   | 4.49             |
| 1,2-Dibromo-3-chloropropane    | 0.5         | 0.116 | 0.410   | 5.24             | 1           | 0.232 | 0.257   | 3.70             |
| Hexachlorobutadiene            | 0.5         | 0.166 | 0.412   | 4.66             | 1           | 0.384 | 0.304   | 11.82            |
| 1,2,4-Trichlorobenzene         | 0.5         | 0.153 | 1.128   | 11.38            | 1           | 0.297 | 0.686   | 9.51             |
| Naphthalene                    | 0.5         | 0.081 | 5.418   | 9.62             | 1           | 0.334 | 3.621   | 4.61             |
| 1,2,3-Trichlorobenzene         | 0.5         | 0.094 | 1.143   | 11.00            | 1           | 0.341 | 0.730   | 6.15             |

Table 6: Minimum Detection Limit, Method Minimum Response Factor, Method Response Factor, and Calibration %RSD for the Water and Soil Analysis

\* The %RSD value was greater than the 20% allowed by the method. The alternate linear regression curve passed the method requirement of greater than 0.99.




Figure 1: 50ppb Water Chromatogram




Figure 2: 50ppb Soil Chromatogram

## Conclusions

The Atomx, purge and trap concentrator with a multi-matrix autosampler in conjunction with Thermo Scientific Focus/DSQ II GC/MS was shown to be an excellent system for detection for the analytes by USEPA Method 8260C. In this study all analytes have met and pass the requirements set forth by this method.

## References

1. USEPA Method 5030 Purge-And-Trap For Aqueous Samples Revision 3, May 2003
2. USEPA Method 5035 Closed-System Purge-And Trap and Extractions For Volatile Organics In Soil and Waste Samples Revision 1, July 2002
3. USEPA Method 8260C Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS) Revision 3, August 2006