

# Agilent 1200 Infinity Series Multisamplers



User Manual



**Agilent Technologies**

# Notices

© Agilent Technologies, Inc. 2014

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

## Manual Part Number

G7167-90001 Rev. B

## Edition

12/2014

Printed in Germany

Agilent Technologies  
Hewlett-Packard-Strasse 8  
76337 Waldbronn

## Warranty

**The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.**

## Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

## Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

## Safety Notices

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

### WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

# In This Guide

This manual covers the following modules:

- Agilent 1290 Infinity II Multisampler (G7167B)
- Agilent 1260 Infinity Multisampler (G7167A)

## 1 Introduction

This chapter gives an introduction to the Multisampler.

## 2 Site Requirements and Specifications

This chapter provides information on environmental requirements, physical and performance specifications.

## 3 Using the Module

This chapter explains the essential operational parameters of the module.

## 4 Preparing the module

This chapter explains the operational parameters of the module.

## 5 Optimizing Performance

This chapter gives hints on how to optimize the performance or use additional devices.

## 6 Troubleshooting and Diagnostics

This chapter gives an overview about the troubleshooting and diagnostic features and the different user interfaces.

## 7 Error Information

This chapter describes the meaning of error messages, and provides information on probable causes and suggested actions how to recover from error conditions.

### **8 Test Functions and Calibration**

This chapter describes the built in test functions.

### **9 Maintenance**

This chapter describes the maintenance of the Multisampler

### **10 Parts for Maintenance and Upgrade or Options**

This chapter provides information on parts material required for the module.

### **11 Identifying Cables**

This chapter provides information on cables used with the modules.

### **12 Hardware Information**

This chapter describes the module in more detail on hardware and electronics.

### **13 LAN Configuration**

This chapter provides information on connecting the detector to the Agilent ChemStation PC.

### **14 Appendix**

This chapter provides addition information on safety, legal and web.


# Contents

<b>1</b>	<b>Introduction</b>	<b>9</b>
	Product Description (G7167B)	10
	Product Description (G7167A)	11
	Features (G7167B)	12
	Features (G7167A)	13
	Overview of the Module	14
	Multisampler Principle	16
	System Overview	22
<b>2</b>	<b>Site Requirements and Specifications</b>	<b>29</b>
	Site Requirements	30
	Physical Specifications	34
	Performance Specifications	35
	Physical Specifications of the Sample Cooler	40
<b>3</b>	<b>Using the Module</b>	<b>43</b>
	Magnets	44
	Turn on/off	45
	Status indicators	46
	Drawer Status Indicator	47
	Insert vial trays/wellplates	48
	Remove vial trays/wellplates	49
	Installing the Sample Cooler	50
	Transporting the Multisampler with a Sample Cooler Installed	58
<b>4</b>	<b>Preparing the module</b>	<b>61</b>
	Leak and Waste Handling	62
	Preparing the Multisampler	64
	Solvent Information	66
	Recommended Mats and Vials	71
	Capillary Color Coding Guide	73
	Installing Capillaries	74

Setting up the Autosampler with Agilent Open Lab ChemStation	77
<b>5 Optimizing Performance</b>	<b>87</b>
Delay Volume and Extra-Column Volume	88
How to Configure the Optimum Delay Volume	89
How to Achieve Higher Injection Volumes	94
How to Achieve High Throughput	97
How to Achieve Higher Resolution	98
How to Achieve Higher Sensitivity	101
How to Achieve Lowest Carry Over	102
<b>6 Troubleshooting and Diagnostics</b>	<b>105</b>
User Interfaces	106
Agilent Lab Advisor Software	107
<b>7 Error Information</b>	<b>109</b>
What Are Error Messages	110
General Error Messages	111
<b>8 Test Functions and Calibration</b>	<b>117</b>
Introduction	118
System Pressure Test	119
Auto Referencing	122
Maintenance Positions	124
Injector Steps	127
<b>9 Maintenance</b>	<b>129</b>
Introduction to Maintenance	131
Warnings and Cautions	133
Overview of Maintenance	135
Clean the Module	136
Removal and Installation of the Front Door	137
Remove the Needle Assembly	140
Install the Needle Assembly	144
Exchange the Needle Seat	148
Replace the Rotor Seal	151

Remove the Metering Seal	157
Install the Metering Seal	162
Replace the Peristaltic Pump Cartridge	168
Replace the Flushhead Seal	172
Replace the Injection Valve	176
Removing the Sample Loop-Flex	180
Installing the Sample Loop-Flex	184
Replace the Dummy Drawer	187
Remove the Sample Cooler	196
Install the Sample Cooler	197
Replace the Module Firmware	200
<b>10 Parts for Maintenance and Upgrade or Options</b>	<b>201</b>
Overview of Maintenance Parts	202
Sampler Main Assemblies	203
Hotel Drawer	205
Analytical Head Assembly 40 $\mu$ L	206
Analytical Head Assembly 100 $\mu$ L	207
Analytical Head Assembly 900 $\mu$ L	208
Flush Head Assembly 500 $\mu$ L	209
Injection Valve Assembly	211
Injection Valve with Actuator	213
Needle Port Assembly	214
Door Assy	215
Accessory Kit	216
Tubing Kit Sampler Standard	217
Sample Cooler	218
<b>11 Identifying Cables</b>	<b>219</b>
Cable Overview	220
Analog Cables	222
Remote Cables	224
CAN/LAN Cables	228
Agilent Module to PC	229
USB	230

<b>12 Hardware Information</b>	<b>231</b>
Firmware Description	232
Electrical Connections	235
Interfaces	238
Setting the 6-bit Configuration Switch	246
Instrument Layout	248
Early Maintenance Feedback	249
<b>13 LAN Configuration</b>	<b>251</b>
Setting up the module in a LAN environment	252
Connecting the module via LAN	253
<b>14 Appendix</b>	<b>255</b>
General Safety Information	256
Waste Electrical and Electronic Equipment Directive	262
Refrigerant	263
Radio Interference	265
Sound Emission	266
Solvent Information	267
Agilent Technologies on Internet	268



# 1 Introduction

Product Description (G7167B)	10
Product Description (G7167A)	11
Features (G7167B)	12
Features (G7167A)	13
Overview of the Module	14
Multisampler Principle	16
System Overview	22
Leak and Waste Handling	22
Leak and Waste Handling in a Mixed Configuration	26

This chapter gives an introduction to the Multisampler.



## Product Description (G7167B)

The Agilent 1290 Infinity II Multisampler can handle both vials and microtiter plates with ease and efficiency up to 1300 bar system pressure, optimized on chromatographic performance.

In fact, this compact module has the capacity to house up to 6144 samples, all inside the Agilent stack footprint and the robotics to smoothly inject each into the chromatograph in turn.

With Agilent's unique dual-needle design, cycle time is just 5 seconds.

With the multi-wash capability, you can reduce carryover to less than 9 parts per million.




Figure 1 Overview of the Multisampler

## Product Description (G7167A)

The Agilent 1260 Infinity Multisampler can handle both vials and microtiter plates with ease and efficiency up to 600 bar system pressure, optimized on high flexibility.

This compact module can house up to 6144 samples, all inside the Agilent stack footprint and the robotics to inject each into the chromatograph in turn.

With Agilent's dual-needle design, cycle time is about 5 seconds.

With the multi-wash capability, you can reduce carryover to less than 9 parts per million.




Figure 2 Overview of the Multisampler

## Features (G7167B)

- *Unmatched flexibility* – You choose how you want to introduce samples for injection, whether you prefer vials, microtiter plates, or any combination of formats. Sample drawers are available in three heights, and you can mix shallow drawers with deeper ones to accommodate different sample sizes.
- *High capacity* – Using shallow well-plate drawers, the 1290 Infinity II Multisampler takes a maximum load of 16 microtiter plates and up to 6144 samples—the most of any single system.
- *Seamless automation* – Internal robotics move microtiter plates and other sample containers from the sample hotel to the central workspace for sample processing steps and injections.
- *Dual-needle injection* – By running samples alternately through one or the other injection path, you can reduce cycle times to mere seconds, virtually eliminating conventional wait times—whether for large volume loadings or flushing procedures.
- *Scalable injection volumes* – The Agilent unique dual-needle integ setup also enhances flexibility by providing two differently optimized injectors in a single instrument. You can, for example, optimize one path for large volume injections and the other for low delay volumes.
- *Ultralow carryover* – 1290 Infinity II Multisampler is designed for low carryover, but you can take clean to a whole new level with our multi-wash capability, cleansing all relevant injection parts between runs. This sophisticated, integrated feature flushes the injection needle outside with three solvents, and uses seat backflush procedures to reduce carryover to less than 10 ppm.
- *Efficient temperature control* – For temperature-sensitive samples, simply add Agilent's new highly efficient compressor-based cooling system. It allows you to maintain perfect temperature control on all vials and plates inserted into the 1290 Infinity Multisampler.
- *Instant information* – Lights on each drawer tell you all you need to know about loading status, current activity, and accessibility.

## Features (G7167A)

- *Unmatched flexibility* – You can choose how you want to introduce samples for injection, whether you prefer vials, microtiter plates, or any combination of formats. Sample drawers are available in three heights, and you can mix shallow drawers with deeper ones to accommodate different sample sizes.
- *High capacity* – Using shallow well-plate drawers, the 1260 Infinity Multisampler takes a maximum load of 16 microtiter plates and up to 6144 samples—the most of any single system.
- *Seamless automation* – Internal robotics move microtiter plates and other sample containers from the sample hotel to the central workspace for sample processing steps and injections.
- *Dual-needle injection* – By running samples alternately through one or the other injection path, you can reduce cycle times to mere seconds, virtually eliminating conventional wait times—whether for large volume loadings or flushing procedures.
- *Scalable injection volumes* – The Agilent dual-needle setup enhances flexibility by providing two differently optimized injectors in a single instrument. You can, for example, optimize one path for large volume injections and the other for low delay volumes.
- *Ultralow carryover* – The 1260 Infinity Multisampler has a low carryover, and a multi-wash capability, cleansing all relevant injection parts between runs. This integrated feature flushes the injection needle outside with three solvents, and uses seat backflush procedures to reduce carryover to less than 9 ppm.
- *Efficient temperature control* – For temperature-sensitive samples, add Agilent’s compressor-based cooling system. It maintains temperature control on all vials and plates inserted into the 1260 Infinity Multisampler.
- *Instant information* – Lights on each drawer tell you about loading status, current activity, and accessibility.

## Overview of the Module

The Multisampler transport mechanism uses a Cartesian robot. The X-Y drive together with the Z drive optimize the grabbing and positioning for the sample trays and the needle handling inside of the Multisampler. The sample coupler moves the sample container from the sample hotel which stores all the samples and place it on the central workspace. Then the needle coupler of the Z drive take over and grab the needle assembly from the needle station and perform the analytical procedures inside of the Multisampler. Due the uncoupled needle design the robot can do other liquid handling jobs during the analysis.

The multisampler employ an active vial/plate pusher mechanism to hold down the vial or the plate while the needle is drawn back from the sample vessel (a must in the case a septum is used). This active vial/plate pusher employs a sensor to detect the presence of a plate and to ensure accurate movement regardless of plate used. All axes of the transport mechanism are driven by very fast BLCD motors. Optical encoders ensure the correct operation of the movement.

The standard metering device provides injection volumes from 0.1 – 20  $\mu\text{L}$ . A 0.1 – 40  $\mu\text{L}$  injection volume metering device is installed in the G7167A/B, with a 20  $\mu\text{L}$ , low restriction loop capillary restricting the injection volume. The entire flowpath including the metering device is always flushed by the mobile phase after injection for minimum internal carry-over.

To reduce the carry over you have two different possibilities. First the external needle wash. In the Standard configuration the needle flush station is equipped with a peristaltic pump to wash the outside of the needle. This reduces the already low carry-over for very sensitive analysis. The bottle containing the mobile phase for the wash procedure will be located in the solvent bottle cabinet. Produced waste during this operation is channeled safely away through a waste drain. In the Multi-Wash configuration the external needle wash will be done by a micro piezo pump combined with a solvent selection valve, where you can select between three different solvents. If this is not sufficient to reduce the carry over there is an additional and perfect way to achieve the lowest carry over in the Multi-Wash configuration by using the intergrated flush

pump. This high pressure pump can also select between three different solvents and it is capable of reducing the carry over to a minimum by using the seat backflushing. The flush pump outlet capillary is connected to the port 4 of the Multisampler's injection valve, which normally holds the waste line. If the Multisampler is in bypass mode the flush pump connects to the needle seat and can flush backwards through the needle seat into the waste line attached to the needle seat outlet port.

The six-port (only 5 ports are used) injection valve unit is driven by a high-speed hybrid stepper motor. During the sampling sequence, the valve unit bypasses the Multisampler, and connects flow from the pump to the column directly. During injection and analysis, the valve unit directs the flow through the Multisampler which ensures that all of the sample is injected onto the column, and that the metering unit and needle are always free of sample residue before the next sampling sequence begins.

The Cooling Control of the vial/plate temperature in the Multisampler is achieved using an additional Agilent Sample Cooler module. The sample cooler is a micro compressor based refrigerator. A fan draws air from the central workstation above the sample container of the Multisampler. It is then blown through the fins of the cooling module. There it is cooled according the temperature setting. The cooled air enters the Sampler Hotel through a recess underneath the special designed base plate. The air is then distributed evenly through the Sample Hotel ensuring effective temperature control, regardless of how many sample containers are in the drawer. In cooling mode condensation is generated on the cooled side of the Sample Cooler. This condensed water is safely guided into a waste bottle for condensed water located underneath the working bench.

## Multisampler Principle

The movements of the Multisampler components during the sampling sequence are monitored continuously by the Multisampler processor. The processor defines specific time windows and mechanical ranges for each movement. If a specific step of the sampling sequence is not completed successfully, an error message is generated. Solvent is bypassed from the Multisampler by the injection valve during the sampling sequence. After the required sample container was automatically loaded from the sample hotel and placed on the central workspace. The Needle assembly moves via robot to the desired sample position and is lowered into the sample liquid in the sample to allow the metering device to draw up the desired volume by moving its plunger back a certain distance. The needle assembly is then raised again and moved to the needle park station onto the seat to close the sample loop. Sample is applied to the column when the injection valve returns to the mainpass position at the end of the sampling sequence.

The standard sampling sequence occurs in the following order:

- 1 The robot loads the required sample container on the central workspace
- 2 The injection valve switches to the bypass position.
- 3 The plunger of the metering device moves to the initialization position.
- 4 The robot couples into the needle assembly from the needle parkstation.
- 5 The robot unlocks the needle assembly and moves up.
- 6 The coupled needle assembly/robot moves to the desired sample vial (or well plate) position on the central workstation.
- 7 The needle lowers into the sample vial (or well plate).
- 8 The metering device draws the preset sample volume.
- 9 The needle lifts out of the sample vial (or well plate).
- 10 The coupled needle assembly/robot is then moved to the park station onto the seat to close the sample loop.
- 11 The needle assembly is locked into the park station and moves down.
- 12 The injection cycle is completed when the injection valve switches to the mainpass position.

13 The robot moves the sample container back into the sample hotel if the sampling sequence is done. If needle wash is required it will be done between step 9 and 10.

**NOTE**

For the needle seat backflush the Multisampler must be in bypass mode.

If an additional needle seat backflush is required this step must also be done between step 9 and 10.

## Injection Sequence

Before the start of the injection sequence, and during an analysis, the injection valve is in the mainpass position. In this position, the mobile phase flows through the Multisampler metering device, sample loop, and needle, ensuring all parts in contact with sample are flushed during the run, thus minimizing carry-over.




Figure 3 Valve in mainpass, flow through


When the sample sequence begins, the valve unit switches to the bypass position. Solvent from the pump enters the valve unit at port 1, and flows directly to the column through port 6.

The standard injection starts with draw sample from vial/wellplate from the central workstation. In order to do this the needle assembly moves via robot to the desired sample position and is lowered into the sample liquid

## 1 Introduction

### Multisampler Principle


in the sample to allow the metering device to draw up the desired volume by moving its plunger back a certain distance. The needle assembly is then raised again and moved to the needle park station onto the seat to close the sample loop. In case of an injector program several steps are interspersed at this point.



**Figure 4** Valve in bypass, drawing sample

#### *Flush the Needle*

Before injection and to reduce the carry-over for very sensitive analysis, the outside of the needle can be washed in a flush port located behind the injector port. As soon as the needle is on the flush port a wash pump delivers some solvent during a defined time to clean the outside of the needle. At the end of this process the needle assembly returns to the needle port.




**Figure 5** Valve in bypass, washing needle

*Inject- and- Run*

The final step is the inject- and run-step. The six-port valve is switched to the mainpass position, and directs the flow back through the sample loop, which now contains a certain amount of sample. The solvent flow transports the sample onto the column, and separation begins. This is the beginning of a run within an analysis. In this stage, all major performance-influencing hardware is flushed internally by the solvent flow. For standard applications no additional flushing procedure is required.

# 1 Introduction


## Multisampler Principle



**Figure 6** Valve in mainpass, sample injected

### *Needle seat back flush*

After the injection to reduce the carry-over for very sensitive analysis, the needle seat can be flushed by an integrated flush pump with up to 3 different solvents which may have different properties and solvent strengths. As soon as injection valve is in bypass mode the flush pump delivers some solvent during a defined time to clean the needle seat. The back flushing solvent will be guided into the waste line attached on the needle wash port. At the end of this process the injection valve switches back into the mainpass position ready for the next injection. The last rinsing step should always include the mobile phase as solvent to get the initial conditions again.




**Figure 7** Valve in bypass, needle backflush (Multiwash)

## System Overview

### Leak and Waste Handling

The 1290 Infinity II Series has been designed for safe leak and waste handling. It is important that all security concepts are understood and instructions are carefully followed.



**Figure 8** Leak and waste handling concept (overview - typical stack configuration as an example)

The solvent cabinet is designed to store a maximum volume of 6 L solvent. The maximum volume for an individual bottle stored in the solvent cabinet should not exceed 4 L. For details, see the usage guideline for the Agilent 1200 Infinity Series Solvent Cabinets (a printed copy of the guideline has been shipped with the solvent cabinet, electronic copies are available on the Internet).

All leak plane outlets are situated in a consistent position so that all Infinity and Infinity II modules can be stacked on top of each other. Waste tubes are guided through a channel on the right hand side of the instrument, keeping the front access clear from tubes.

The leak plane provides leak management by catching all internal liquid leaks, guiding them to the leak sensor for leak detection, and passing them on to the next module below, if the leak sensor fails. The leak sensor in the leak plane stops the running system as soon as the leak detection level is reached.

Solvent and condensate is guided through the waste channel into the waste container:

- from the detector's flow cell outlet
- from the Multisampler needle wash port
- from the Sample Cooler (condensate)
- from the Seal Wash Sensor
- from the pump's Purge Valve or Multipurpose Valve

The waste tube connected to the leak pan outlet on each of the bottom instruments guides the solvent to a suitable waste container.


## Waste Concept

- 1 Agilent recommends using the 6 L waste can with 1 Stay Safe cap GL45 with 4 ports (5043-1221) for optimal and safe waste disposal. If you decide to use your own waste solution, make sure that the tubes don't immerse in the liquid.




## Leak and Waste Handling in a Mixed Configuration

### Leak and Waste Handling in a Mixed Configuration



**Figure 9** Leak and waste handling with multisampler in a mixed configuration as an example




**Figure 10** Leak and waste handling with multisampler in a mixed configuration as an example (two stack configuration)

**NOTE**

Do not place the multisampler directly on the bench. If a sample cooler is installed, install Base plate (G1328-44121) underneath.

**1 Introduction**  
System Overview



## 2 Site Requirements and Specifications

Site Requirements	30
Physical Specifications	34
Performance Specifications	35
Performance Specifications (G7167B)	37
Performance Specifications (G7167A)	37
Physical Specifications of the Sample Cooler	40

This chapter provides information on environmental requirements, physical and performance specifications.



## Site Requirements

A suitable environment is important to ensure optimal performance of the instrument.

### Power Considerations

The module power supply has wide ranging capability. It accepts any line voltage in the range described in [Table 1](#) on page 34. Consequently there is no voltage selector in the rear of the module. There are also no externally accessible fuses, because automatic electronic fuses are implemented in the power supply.

#### **WARNING**

**Hazard of electrical shock or damage of your instrumentation can result, if the devices are connected to a line voltage higher than specified.**

→ Connect your instrument to the specified line voltage only.

---

#### **WARNING**

**The module is partially energized when switched off, as long as the power cord is plugged in.**

**Repair work at the module can lead to personal injuries, e.g. electrical shock, when the cover is opened and the module is connected to power.**

→ Always unplug the power cable before opening the cover.

→ Do not connect the power cable to the instrument while the covers are removed.

---

**CAUTION**

Inaccessible power plug.

In case of emergency it must be possible to disconnect the instrument from the power line at any time.

- Make sure the power connector of the instrument can be easily reached and unplugged.
  - Provide sufficient space behind the power socket of the instrument to unplug the cable.
- 

## Power Cords

Different power cords are offered as options with the module. The female end of all power cords is identical. It plugs into the power-input socket at the rear. The male end of each power cord is different and designed to match the wall socket of a particular country or region.

**WARNING**

**Absence of ground connection or use of unspecified power cord**

**The absence of ground connection or the use of unspecified power cord can lead to electric shock or short circuit.**

- Never operate your instrumentation from a power outlet that has no ground connection.
  - Never use a power cord other than the Agilent Technologies power cord designed for your region.
- 

**WARNING**

**Use of unsupplied cables**

**Using cables not supplied by Agilent Technologies can lead to damage of the electronic components or personal injury.**

- Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.
-

**WARNING**

**Unintended use of supplied power cords**

Using power cords for unintended purposes can lead to personal injury or damage of electronic equipment.

- Never use the power cords that Agilent Technologies supplies with this instrument for any other equipment.
- 

## Bench Space

The module dimensions and weight (see [Table 1](#) on page 34) allow you to place the module on almost any desk or laboratory bench. It needs an additional 2.5 cm (1.0 inches) of space on either side and approximately 8 cm (3.1 inches) in the rear for air circulation and electric connections.

If the bench shall carry a complete HPLC system, make sure that the bench is designed to bear the weight of all modules.

The module should be operated in a horizontal position, especially if a sample cooler is installed. Check position with a bulb.

## Condensation

**CAUTION**

Condensation within the module

Condensation can damage the system electronics.

- Do not store, ship or use your module under conditions where temperature fluctuations could cause condensation within the module.
  - If your module was shipped in cold weather, leave it in its box and allow it to warm slowly to room temperature to avoid condensation.
-

## Warranty Provisions and Reference

### Warranty provisions

Agilent standard warranty does not cover:

- Improper use or operation outside the product specification
- Defects due to adjustments or maintenance by customer or unauthorized parties
- Visible abuse, negligence or shipping damage
- Unsupported configurations
- Improper site preparations
- Use of non Agilent parts

### Reference

- The removal of nameplates leads to loss of Warranty claim.
- Any unauthorized destruction of seals will invalidate the Warranty claim.

## Physical Specifications

**Table 1** Physical Specifications

Type	Specification	Comments
Weight	22 kg (48.5 lbs)	w/o sample cooler
Dimensions (height × width × depth)	320 x 396 x 468 mm (12.6 x 15.6 x 18.4 inches)	
Line voltage	100 – 240 V~, ± 10 %	Wide-ranging capability
Line frequency	50 or 60 Hz, ± 5 %	
Power consumption	180 VA, 180 W	
Ambient operating temperature	4 - 40 °C (39 - 104 °F)	
Ambient non-operating temperature	-40 – 70 °C (-40 – 158 °F)	
Humidity	< 95 % r.h. at 40 °C (104 °F)	Non-condensing
Operating altitude	Up to 3000 m (9842 ft)	
Non-operating altitude	Up to 4600 m (15092 ft)	For storing the module
Safety standards: IEC, EN, CSA, UL	Installation category II, Pollution degree 2	For indoor use only.

## Performance Specifications

**Table 2** Agilent 1290 Infinity II Multisampler (G7167B) Performance Specifications

Type	Specification	Comment
Injection range for <i>Single-needle</i> instruments	0.1 – 100 µL optional: 20 µL, 40 µL in 0.1 µL increments (default)	Up to 1300 bar using 100 µL analytical head
	0.1 – 500 µL or 900 µL in 0.1 µL increments (default)	Pressure range up to 400 bar due to 900 µL analytical head
	0.1 – 120 µL in 0.1 µL increments with 1290 Infinity large volume injection kit (hardware modification required) G4216-68711	Pressure range up to 1300 bar Multi-draw modus (Injection into needle-seat capillary)
	0.1 – 500 µL, 1500 µL in 0.1 µL increments with 100 µL upgrade kit (hardware modification required) G7167-68711	
Injection range for <i>Dual-needle</i> instruments	0 – 20 µL optional: 40 µL in 0.1 µL increments (default)	Up to 1300 bar using 40 µL analytical head
	0 – 100 µL optional: 500 µL in 0.1 µL increments (default)	Up to 1300 bar using 100 µL analytical head + Multi-load
Precision	<0.15 % RSD or SD <10 nL, whatever is greater	Measured caffeine
Pressure range	Up to 1300 bar (G7167B)	Max pressure for basic instrument
Sample viscosity range	0.2 – 5 cp	
Sample capacity	1H Drawer up to 8 drawers and 16 positions Shallow well plates (MTP)	Max. 6144/1536 samples (384MTP/96)
	2H Drawer up to 4 drawers and 8 positions MTP, deep well plates, vials, Eppendorf	3072 samples, 432 vials (2 mL)
	3H Drawer up to 2 drawers and 4 positions MTP, deep well plates, vials up to 6 mL, Eppendorf	1536 samples, 60 vials (6 mL), 384 vials (1 mL), 216 vials (2 mL)
Injection cycle time	<10 s using following standard conditions: Default draw speed: 100 µL/min Default eject speed: 400 µL/min Injection volume: 1 µL	Using standard Single-needle setup  Time between 2 injections is not mechanically limited, time delay depends on communication speed of software, OS or network connections
	< 5 s using a Dual-needle setup	

## 2 Site Requirements and Specifications

### Performance Specifications

**Table 2** Agilent 1290 Infinity II Multisampler (G7167B) Performance Specifications

Type	Specification	Comment
Carry Over	<0.003 % (30 ppm) Multisampler Standard <0.0009 % (9 ppm) Multisampler Multiwash	Using the following conditions: <ul style="list-style-type: none"> <li>• Column: Agilent ZORBAX SB-C18, 2.1x 50 mm1.8 mm (827700-902)</li> <li>• Mobile Phase: <ul style="list-style-type: none"> <li>• A: 0.1 % TFA in water</li> <li>• B: 0.1 % TFA in Acetonitrile</li> </ul> </li> <li>• Isocratic : % B=35 %</li> <li>• Flow rate: 0.5 mL/min</li> <li>• Temperature: 25 °C</li> <li>• Wavelength: 257 nm</li> <li>• Sample: 1200 ng/μL Chlorhexidine (dissolved with mobile phase A), 1 μL injected and measured on G4212A DAD</li> <li>• Wash solution: H2O with 0.1 % TFA (5 s)</li> </ul>
Multiwash	Outer needle wash and seat backflush for carryover reduction with up to 3 different solvents	
Control and data evaluation	Agilent Open Lab CDS MassHunter QQQ MassHunter QTOF Lab Advisor ICF for 3rd party SW control LC and CE Drivers	A.02.01 or above (A.02.02 supports Sample Entry UI) B.07.00 SP1 <sup>1</sup> or above B.05.01 SP3 <sup>1</sup> or above B.02.05 or above A.02.01 or above A.02.10 or above
Local Control	Agilent Instant Pilot (G4208A)	B.02.17 or above (currently not supported/official release 2015)
Communications	Controller-area network (CAN), Local Area Network (LAN) ERI: ready, start, stop and shut-down signals	
Safety and maintenance	Extensive support for troubleshooting and maintenance is provided by the Instant Pilot, Agilent Lab Advisor, and the Chromatography Data System. Safety-related features are leak detection, safe leak handling, leak output signal for shutdown of pumping system, and low voltages in major maintenance areas.	

**Table 2** Agilent 1290 Infinity II Multisampler (G7167B) Performance Specifications

Type	Specification	Comment
GLP features	Early maintenance feedback (EMF) for continuous tracking of instrument usage with user-settable limits and feedback messages. Electronic records of maintenance and errors.	
Housing	All materials recyclable.	
Metering device	Metering device in high pressure flow path	

<sup>1</sup> only for the integration in an Infinity I LC setup

## Performance Specifications (G7167A)

**Table 3** Agilent 1260 Infinity Multisampler (G7167A) Performance Specifications

Type	Specification	Comment
Injection range for <i>Single-needle</i> instruments	0.1 – 100 µL optional: 20 µL, 40 µL in 0.1 µL increments (default) 0.1 – 500 µL or 900 µL in 0.1 µL increments (default) 0.1 – 120 µL in 0.1 µL increments with 1290 Infinity large volume injection kit (hardware modification required) G4216-68711 0.1 – 500 µL, 1500 µL in 0.1 µL increments with 100 µL upgrade kit (hardware modification required) G7167-68711	Up to 600 bar using 100 µL analytical head Pressure range up to 400 bar due to 900 µL analytical head Pressure range up to 600 bar Multi-draw modus (Injection into needle-seat capillary)
Injection range for <i>Dual-needle</i> instruments	0 – 20 µL optional: 40 µL in 0.1 µL increments (default) 0 – 100 µL optional: 500 µL in 0.1 µL increments (default)	Up to 600 bar using 40 µL analytical head Up to 600 bar using 100 µL analytical head + Multi-load
Precision	<0.15 % RSD or SD < 10 nL, whatever is greater	Measured caffeine
Pressure range	Up to 600 bar (G7167A)	Max pressure for basic instrument
Sample viscosity range	0.2 – 5 cp	

## 2 Site Requirements and Specifications

### Performance Specifications

**Table 3** Agilent 1260 Infinity Multisampler (G7167A) Performance Specifications

Type	Specification	Comment
Sample capacity	1H Drawer up to 8 drawers and 16 positions Shallow well plates (MTP)	Max. 6144/1536 samples (384MTP/96)
	2H Drawer up to 4 drawers and 8 positions MTP, deep well plates, vials, Eppendorf	3072 samples, 432 vials (2 mL)
	3H Drawer up to 2 drawers and 4 positions MTP, deep well plates, vials up to 6 mL, Eppendorf	1536 samples, 60 vials (6 mL), 384 vials (1 mL), 216 vials (2 mL)
Injection cycle time	<10 s using following standard conditions: Default draw speed: 100 µL/min Default eject speed: 400 µL/min Injection volume: 1 µL	Using standard Single-needle setup  Time between 2 injections is not mechanically limited, time delay depends on communication speed of software, OS or network connections
	< 5 s using a Dual-needle setup	
Carry Over	<0.003 % (30 ppm) Multisampler Standard <0.0009 % (9 ppm) Multisampler Multiwash	Using the following conditions: <ul style="list-style-type: none"> <li>• Column: Agilent ZORBAX SB-C18, 2.1x 50 mm/1.8 mm (827700-902)</li> <li>• Mobile Phase: <ul style="list-style-type: none"> <li>• A: 0.1 % TFA in water</li> <li>• B: 0.1 % TFA in Acetonitrile</li> </ul> </li> <li>• Isocratic : % B=35 %</li> <li>• Flow rate: 0.5 mL/min</li> <li>• Temperature: 25 °C</li> <li>• Wavelength: 257 nm</li> <li>• Sample: 1200 ng/µL Chlorhexidine (dissolved with mobile phase A), 1 µL injected and measured on G4212A DAD</li> <li>• Wash solution: H2O with 0.1 % TFA (5 s)</li> </ul>
Multiwash	Outer needle wash and seat backflush for carryover reduction with up to 3 different solvents	
Control and data evaluation	Agilent Open Lab CDS	A.02.01 or above (A.02.02 supports Sample Entry UI)
	MassHunter QQQ MassHunter QTOF Lab Advisor ICF for 3rd party SW control LC and CE Drivers	B.07.00 SP1 <sup>1</sup> or above B.05.01 SP3 <sup>1</sup> or above B.02.05 or above A.02.01 or above A.02.10 or above

**Table 3** Agilent 1260 Infinity Multisampler (G7167A) Performance Specifications

Type	Specification	Comment
Local Control	Agilent Instant Pilot (G4208A)	B.02.17 or above (currently not supported/official release 2015)
Communications	Controller-area network (CAN), Local Area Network (LAN) ERI: ready, start, stop and shut-down signals	
Safety and maintenance	Extensive support for troubleshooting and maintenance is provided by the Instant Pilot, Agilent Lab Advisor, and the Chromatography Data System. Safety-related features are leak detection, safe leak handling, leak output signal for shutdown of pumping system, and low voltages in major maintenance areas.	
GLP features	Early maintenance feedback (EMF) for continuous tracking of instrument usage with user-settable limits and feedback messages. Electronic records of maintenance and errors.	
Housing	All materials recyclable.	
Metering device	Metering device in high pressure flow path	

<sup>1</sup> only for the integration in an Infinity I LC setup

## Physical Specifications of the Sample Cooler

Cooling unit is designed as vapor-compression refrigeration system. Contains fluorinated greenhouse gas (refrigerant) according to the Kyoto protocol. For specifications of refrigerant, charge capacity, carbon dioxide equivalent (CDE), and global warming potential (GWP) see instrument label.

**Table 4** Physical Specification of the Sample Cooler

Type	Specification	Comments
Weight	< 6 kg	
Dimensions (height × width × depth)	205 mm X 340 mm X 370 mm	
Refrigerant gas	HFKW-134a (0.042 kg)	Ozone depletion potential (ODP) = 0
Line voltage	24 VDC (nominal)	
Current	10 A max.	
Ambient operating temperature	4 – 40 ° C (39.2 – 131 ° F)	
Ambient non-operating temperature	-40 – 70 ° C (-20 – 158 ° F)	
Humidity	4 – 35 ° C, 20 – 80 % RH	Non-condensing
Operating altitude	Up to 3000 m (9842 ft)	
Non-operating altitude	Up to 4600 m (15091 ft)	
Safety standards: IEC, CSA, UL	Installation category II, Pollution degree 2	For indoor use only.

**CAUTION**

General hazards and improper disposal

Improper disposal of the media and components used pollutes the environment.


- The breakdown of the sample cooler unit must be carried out by specialist refrigeration company.
- All media must be disposed of in accordance with national and local regulations.
- Please contact your local Agilent Service Center in regard to safe environmental disposal of the appliance.

**Table 5** Performance Specifications Agilent 1290 Sample Cooler

Type	Specifications
Operating principle	High performance, low-energy consumption micro-compressor based cooler with ozone-friendly R134A coolant (42 g), user-upgradable.
Temperature range	from 4 °C to ambient
Temperature settable	from 4 – 40 °C in 1 ° increments
Temperature accuracy (<25 °C, <50 % r.H.)	2 °C to 6 °C at a setpoint of 4 °C
Temperature Stability	± 2 °C at a setpoint of 4 °C

## 2 Site Requirements and Specifications

### Physical Specifications of the Sample Cooler



## 3 Using the Module


Magnets	44
Turn on/off	45
Status indicators	46
Drawer Status Indicator	47
Insert vial trays/wellplates	48
Remove vial trays/wellplates	49
Installing the Sample Cooler	50
Damaged Packaging	50
Install the Sample Cooler	51
Transporting the Multisampler with a Sample Cooler Installed	58

This chapter explains the essential operational parameters of the module.






# Magnets

1




# Turn on/off

<p>1</p> 	<p>2</p>  <p>Power switch (1) On (2) Off</p>
<p>3</p> 	

## Status indicators

- 1 The module status indicator indicates one of six possible module conditions:



### Status indicators

1. Idle
2. Run mode
3. Not-ready. Waiting for a specific pre-run condition to be reached or completed.
4. Error mode - interrupts the analysis and requires attention (for example a leak or defective internal components).
5. Resident mode (blinking) - for example during update of main firmware.
6. Bootloader mode (fast blinking). Try to re-boot the module or try a cold-start. Then try a firmware update.


## Drawer Status Indicator

The module status indicator indicates one of three possible module conditions:

- When the status indicator is *OFF* no sample containers are loaded.
- When the upper, lower or both semi circle status indicators are *ON*, indicates the rear or front position of the drawer or both positions are loaded with a sample containers.
- When semi circle indicators are *blinking* the robot interacts with a drawer.


### NOTE

During blinking of the drawer status indicator. Do not try to open the drawer at this point.




3 Using the Module  
Insert vial trays/wellplates

## Insert vial trays/wellplates



## Remove vial trays/wellplates



## Installing the Sample Cooler

### Damaged Packaging

If the delivery packaging shows signs of external damage, please call your Agilent Technologies sales and service office immediately. Inform your service representative that the instrument may have been damaged during shipment.

#### CAUTION

"Defective on arrival" problems

If there are signs of damage, please do not attempt to install the module. Inspection by Agilent is required to evaluate if the instrument is in good condition or damaged.

- Notify your Agilent sales and service office about the damage.
  - An Agilent service representative will inspect the instrument at your site and initiate appropriate actions.
-

## Install the Sample Cooler

Parts required	p/n	Description
		Multisampler
	G7167-60005	Sample cooler
		Power cord

**Hardware required** Other cables see below and “[Cable Overview](#)” on page 220

**Software required** OpenLabCDS and/or Instant Pilot G4208A with the appropriate revisions, see [Table 2](#) on page 35.

### NOTE


If the sample cooler is disconnected from the power supply, you should wait for at least five minutes before replugging and switching on the compressor again.

### NOTE

Even under average humidity conditions, a significant amount of condensed water gathers every day. A suitable container must be provided and emptied regularly in order to avoid overflow.

**1** Place the Autosampler on the bench.

**2** Ensure that the power switch on the front of the module is OFF (switch stands out).



### 3 Using the Module

#### Installing the Sample Cooler

3 Ensure that the power cable is removed from the instrument.



4 Open the 4 screws on the rear of the module.



5 Remove the sample cooler mainframe cover.



6 Slide in the sample cooler the halfway.



### WARNING

Module is partially energized when switched off, as long as the power cord is plugged in. Repair work at the module can lead to personal injuries, e.g. shock hazard, when the cover is opened and the module is connected to power.

- Make sure that it is always possible to access the power plug.
- Do not use the sample cooler if it is not operating correctly or has been damaged. Disconnect it from the power supply and call your local service center.
- Remove the power cable from the module before opening the cover.
- Do not connect the power cable to the module while the covers are removed.
- If the sample cooler is disconnected from the power supply, you should wait for at least five minutes before switching on the compressor.

### CAUTION


Damaged electronics

- To avoid damages of the electronics of the module make sure the power cords are unplugged before disconnecting or reconnecting the multisampler to the sample cooler cables.

### 3 Using the Module

#### Installing the Sample Cooler

7 Connect power cable and signal/data cable.




#### CAUTION

Damage to the cables

- Do not bend or pinch the cables.
- Fit in the sample cooler perfectly.

8 Slide in the whole unit.



- 9 Tighten the 4 screws which holds the sample cooler unit in place.



### CAUTION

**Routing of the condensation tubing**  
Proper routing of the condensation tubing is critical for correct condensate drainage.

- Do not place the multisampler directly on the bench. Install Base plate (G1328-44121) underneath.

### CAUTION

**Damage through condensation**  
If the condensation tube is located in liquid the condensed water cannot flow out of the tube and the outlet is blocked. Any further condensation will then remain in the instrument. This may damage the instruments electronics.

- Make sure the condensation tube is always above the liquid level in the vessel.
- Make sure the waste container is not sealed.
- Ensure the drain tube has no loops and leads directly into the waste bottle.

### WARNING

**Heavy weight**  
The module is heavy (>22 kg (>46 lbs)).

- Carry the module at least with 2 people.
- Avoid back strain or injury by following all precautions for lifting heavy objects.
- Ensure that the load is as close to your body as possible.
- Ensure that you can cope with the weight of your load.

- 10 Lift the multisampler with the sample cooler installed into the LC stack.

- 11 Use a bubble level to check the leveling of the Multisampler.


### NOTE

The sample cooler should operated in a proper horizontal position.

### 3 Using the Module

#### Installing the Sample Cooler


**12** Install the condensate tubing on the drain outlet of the sample cooler and guide the condensation tubing out to the central waste.



#### NOTE

Do not kink the drain tubing (use an elbow connector) and avoid siphoning effects of the drainage.

**13** Ensure the power switch on the front of the module is OFF (switch stands out).



#### CAUTION

Damage to the sample cooler

- Wait at least 30 min before switching on the compressor of the sample cooler.
- This allows the refrigerant and system lubrication to reach equilibrium.

- 14 Connect the CAN interface cables to other modules in the system (see section *Recommended Stack Configurations* in the technical note *Use of Multisampler in Mixed Configurations* or the service manual).
- 15 If required, connect additional interface and control cables to the autosampler (see section *Recommended Stack Configurations* in the technical note *Use of Multisampler in Mixed Configurations* or the service manual). Refer to the documentation of the Agilent 1200 Infinity Series Instant Pilot or ChemStation for LC for more information.

**NOTE**


In an Agilent 1290 Infinity or 1260 Infinity system, the individual modules are connected by a CAN cable. The Agilent 1200 Infinity Series Instant Pilot can be connected to the CAN bus at any of the modules in the system. If an Agilent detector is part of the system, usually the LAN connection should be at the detector. For more information about connecting the instant pilot or control software refer to the respective user manual.

For connecting the Agilent 1290 Infinity equipment to non-Agilent 1290 Infinity equipment, see "[Cable Overview](#)" on page 220).

- 16 Connect the power cable to the power connector at the rear of the module.



17



Power switch

- (1) On
- (2) Off

### 3 Using the Module

#### Transporting the Multisampler with a Sample Cooler Installed

## Transporting the Multisampler with a Sample Cooler Installed

### NOTE

There are magnets in the front area of the multisampler, see “Magnets” on page 44.

---

### NOTE

When moving the multisampler around the laboratory, make sure that any condensed water inside the thermostat is removed.

- Tilt the module to the back, so that the water inside the thermostat can safely flow into the leak funnel.
  - Otherwise no special precautions are needed for the modules.
- 

### WARNING

#### Heavy weight

The module is heavy (>22 kg (>46 lbs)).

- Carry the module at least with 2 people.
  - Avoid back strain or injury by following all precautions for lifting heavy objects.
  - Ensure that the load is as close to your body as possible.
  - Ensure that you can cope with the weight of your load.
- 

### CAUTION

#### Mechanical damage of the module

If the transport assembly is not parked and not protected by the transport foam, the module could be damaged due to excessive shock of the shipping container during transport.


- Always park the transport assembly before shipment.
  - Store the installation foam in a save place, to use it for later transport of the module.
-

If the multisampler with a sample cooler needs to be shipped to another location via carrier, ensure:

- The two modules are shipped in separate boxes.
- The Sample handler of the multisampler is parked, see **Park Robot** in Agilent Lab Advisor online help for more information.
- The sample containers (vial trays) are removed from the sample hotel.
- The condensed water inside of the sample cooler is removed.

### **3 Using the Module**

#### Transporting the Multisampler with a Sample Cooler Installed



## 4 Preparing the module

Leak and Waste Handling	62
Preparing the Multisampler	64
Solvent Information	66
Recommended Mats and Vials	71
Capillary Color Coding Guide	73
Installing Capillaries	74
Setting up the Autosampler with Agilent Open Lab ChemStation	77
Control Settings	81
Method Parameter Settings	82
Module Configuration View	85

This chapter explains the operational parameters of the module.



## Leak and Waste Handling

### WARNING

**Toxic, flammable and hazardous solvents, samples and reagents**

The handling of solvents, samples and reagents can hold health and safety risks.

- When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
- The volume of substances should be reduced to the minimum required for the analysis.
- Do not operate the instrument in an explosive atmosphere.
- Never exceed the maximal permissible volume of solvents (6 L) in the solvent cabinet.
- Do not use bottles that exceed the maximum permissible volume as specified in the usage guideline for the Agilent 1200 Infinity Series Solvent Cabinets.
- Arrange the bottles as specified in the usage guideline for the solvent cabinet.
- A printed copy of the guideline has been shipped with the solvent cabinet, electronic copies are available on the Internet.
- Ground the waste container.
- The residual free volume in the appropriate waste container must be large enough to collect the waste liquid.
- Check the filling level of the waste container regularly.
- To achieve maximal safety, check the correct installation regularly.
- Do not use solvents with an auto-ignition temperature below 200 °C (392 °F).

### NOTE

**Recommendations for Solvent Cabinet**

For details, see the usage guideline for the Agilent 1200 Infinity Series Solvent Cabinets.

For correct installation of your system contact your Agilent service representative.

## Preparing the Multisampler

For best performance of the multisampler

- When using the multisampler in a system with a vacuum degassing unit, shortly degas your samples before using them in the multisampler.
- Filter samples before use in a 1200 Infinity Series system. Use High pressure filter kit (5067-4638) for inline filtering.
- When using buffer solutions, flush the system with water before switching it off.
- Check the multisampler plungers for scratches, grooves and dents when changing the piston seal. Damaged plungers cause micro leaks and will decrease the lifetime of the seal.
- Solvent Information - Observe recommendations on the use of solvents, see “[Solvent Information](#)” on page 66.
- Priming and Purging the System - When the solvents have been exchanged or the system has been turned off for a certain time (for example, overnight) oxygen will re-diffuse into the solvent channel. Therefore priming and purging of the system is required before starting an application.

**Table 6** Choice of Priming Solvents for Different Purposes

Activity	Solvent	Comments
After an installation	Isopropanol	Best solvent to flush air out of the system
When switching between reverse phase and normal phase (both times)	Isopropanol	Best solvent to flush air out of the system
After an installation	Ethanol or methanol	Alternative to isopropanol (second choice) if no isopropanol is available

**Table 6** Choice of Priming Solvents for Different Purposes

Activity	Solvent	Comments
To clean the system when using buffers	Bidistilled water	Best solvent to re-dissolve buffer crystals
After a solvent change	Bidistilled water	Best solvent to re-dissolve buffer crystals

## Solvent Information

Observe the following recommendations on the use of solvents.

- Follow recommendations for avoiding the growth of algae, see pump manuals.
- Small particles can permanently block capillaries and valves. Therefore, always filter solvents through 0.4  $\mu\text{m}$  filters.
- Avoid or minimize the use of solvents that may corrode parts in the flow path. Consider specifications for the pH range given for different materials like flow cells, valve materials etc. and recommendations in subsequent sections.

## Recommended Wash Solvents

- water
- ethanol
- methanol
- water/acid (especially for basic compounds)
- water/base (especially for acidic compounds)
- water/acetonitrile

## Material Information

Materials in the flow path are carefully selected based on Agilent's experiences in developing highest quality instruments for HPLC analysis over several decades. These materials exhibit excellent robustness under typical HPLC conditions. For any special conditions, please consult the material information section or contact Agilent.

## Disclaimer

Subsequent data were collected from external resources and are meant as a reference. Agilent cannot guarantee the correctness and completeness of such information. Data is based on compatibility libraries, which are not specific for estimating the long-term life time under specific but highly variable conditions of UHPLC systems, solvents, solvent mixtures and samples. Information can also not be generalized due to catalytic effects of impurities like metal ions, complexing agents, oxygen etc. Apart from pure chemical corrosion, other effects like electro corrosion, electrostatic charging (especially for non-conductive organic solvents), swelling of polymer parts etc. need to be considered. Most data available refers to room temperature (typically 20 – 25 °C, 68 – 77 °F). If corrosion is possible, it usually accelerates at higher temperatures. If in doubt, please consult technical literature on chemical compatibility of materials.

## PEEK

PEEK (Polyether-Ether Ketones) combines excellent properties regarding biocompatibility, chemical resistance, mechanical and thermal stability. PEEK is therefore the material of choice for UHPLC and biochemical instrumentation.

It is stable in a pH range between 1 – 12, and inert to many common solvents.

There is still a number of known incompatibilities with chemicals such as chloroform, methylene chloride, THF, DMSO, strong acids (nitric acid > 10 %, sulphuric acid > 10 %, sulfonic acids, trichloroacetic acid), halogenes or aqueous halogene solutions, phenol and derivatives (cresols, salicylic acid etc.).

## Polyimide

Agilent uses semi-crystalline polyimide for rotor seals in valves and needle seats in autosamplers. One supplier of polyimide is DuPont, which brands polyimide as Vespel, which is also used by Agilent.

Polyimide is stable in a pH range between 1 and 10 and in most organic solvents. It is incompatible with concentrated mineral acids (e.g. sulphuric acid), glacial acetic acid, DMSO and THF. It is also degraded by nucleophilic substances like ammonia (e.g. ammonium salts in basic conditions) or acetates.

## Polyethylene (PE)

Agilent uses UHMW (ultra-high molecular weight)-PE/PTFE blends for yellow piston and wash seals, which are used in 1290 Infinity pumps and for normal phase applications in 1260 Infinity pumps.

Polyethylene has a good stability for most common inorganic solvents including acids and bases in a pH range of 1 to 12.5. It is compatible to many organic solvents used in chromatographic systems like methanol, acetonitrile and isopropanol. It has limited stability with aliphatic, aromatic and halogenated hydrocarbons, THF, phenol and derivatives, concentrated acids and bases. For normal phase applications, the maximum pressure should be limited to 200 bar.


## Tantalum (Ta)

Tantalum is inert to most common HPLC solvents and almost all acids except fluoric acid and acids with free sulfur trioxide. It can be corroded by strong bases (e.g. hydroxide solutions > 10 %, diethylamine). It is not recommended for the use with fluoric acid and fluorides.

## Stainless Steel (ST)

Stainless steel is inert against many common solvents. It is stable in the presence of acids and bases in a pH range of 1 to 12.5. It can be corroded by acids below pH 2.3. It can also corrode in following solvents:

- Solutions of alkali halides, their respective acids (for example, lithium iodide, potassium chloride, and so on) and aqueous solutions of halogens.
- High concentrations of inorganic acids like nitric acid, sulfuric acid and organic solvents especially at higher temperatures (replace, if your chromatography method allows, by phosphoric acid or phosphate buffer which are less corrosive against stainless steel).
- Halogenated solvents or mixtures which form radicals and/or acids, for example:



This reaction, in which stainless steel probably acts as a catalyst, occurs quickly with dried chloroform if the drying process removes the stabilizing alcohol.

- Chromatographic grade ethers, which can contain peroxides (for example, THF, dioxane, di-isopropylether). Such ethers should be filtered through dry aluminium oxide which adsorbs the peroxides.
- Solutions of organic acids (acetic acid, formic acid, and so on) in organic solvents. For example, a 1 % solution of acetic acid in methanol will attack steel.
- Solutions containing strong complexing agents (for example, EDTA, ethylene diamine tetra-acetic acid).
- Mixtures of carbon tetrachloride with 2-propanol or THF.

### **Diamond-Like Carbon (DLC)**

Diamond-Like Carbon is inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

### **Fused silica and Quartz (SiO<sub>2</sub>)**

Fused silica is used in 1290 Infinity Flow Cells and capillaries. Quartz is used for classical flow cell windows. It is inert against all common solvents and acids except hydrofluoric acid and acidic solvents containing fluorides. It is corroded by strong bases and should not be used above pH 12 at room temperature. The corrosion of flow cell windows can negatively affect measurement results. For a pH greater than 12, the use of flow cells with sapphire windows is recommended.

### **Gold**

Gold is inert to all common HPLC solvents, acids and bases within the specified pH range. It can be corroded by complexing cyanides and concentrated acids like aqua regia.

### **Zirconium Oxide (ZrO<sub>2</sub>)**

Zirconium Oxide is inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

### **Platinum/Iridium**

Platinum/Iridium is inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

### **Fluorinated polymers (PTFE, PFA, FEP, FFKM)**

Fluorinated polymers like PTFE (polytetrafluorethylene), PFA (perfluoroalkoxy) and FEP (fluorinated ethylene propylene) are inert to almost all common acids, bases, and solvents. FFKM is perfluorinated rubber, which is also resistant to most chemicals. As an elastomer, it may swell in some organic solvents like halogenated hydrocarbons.

TFE/PDD copolymer tubings, which are used in all Agilent degassers except 1322A, are not compatible with fluorinated solvents like Freon, Fluorinert, or Vertrel. They have limited life time in the presence of Hexafluoroisopropanol (HFIP). To ensure the longest possible life with HFIP, it is best to dedicate a particular chamber to this solvent, not to switch solvents, and not to let dry out the chamber. For optimizing the life of the pressure sensor, do not leave HFIP in the chamber when the unit is off.

### **Sapphire, Ruby and Al<sub>2</sub>O<sub>3</sub>-based ceramics**

Sapphire, ruby and ceramics based on aluminum oxide Al<sub>2</sub>O<sub>3</sub> are inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

## Recommended Mats and Vials

**Table 7** Recommended plates and closing mat

Description (Part Number)	Rows	Columns	Plate height	Volume (μL)	Package
384Agilent (5042-1388)	16	24	14.4	80	30
384Corning (No Agilent PN)	16	24	14.4	80	
384Nunc (No Agilent PN)	16	24	14.4	80	
96 well plate 0.5 ml, PP (pack of 10) (5042-1386)	8	12	14.3	500	10
96 well plate 0.5 ml, PP (pack of 120) (5042-1385)					120
96Agilent conical (5042-8502)	8	12	17.3	150	25
96CappedAgilent (5065-4402)	8	12	47.1	300	1
96Corning (No Agilent PN)	8	12	14.3	300	
96CorningV (No Agilent PN)	8	12	14.3	300	
96DeepAgilent31mm (5042-6454)	8	12	31.5	1000	50
96DeepNunc31mm (No Agilent PN)	8	12	31.5	1000	
96DeepRitter41mm (No Agilent PN)	8	12	41.2	800	
96Greiner (No Agilent PN)	8	12	14.3	300	
96GreinerV (No Agilent PN)	8	12	14.3	250	
96Nunc (No Agilent PN)	8	12	14.3	400	
Closing mat for all 96 Agilent plates (5042-1389)	8	12			50

**Recommended  
Vial Plates**

p/n	Description
G2255-68700	Vial plate for 54 x 2 mL vials (6/pk)

## 4 Preparing the module

### Recommended Mats and Vials

p/n	Description
5022-6539	Vial plate for 15 x 6 mL vials (1/pk)
5022-6538	Vial plate for 27 Eppendorf tubes (1/pk)

#### NOTE

Agilent Technologies recommends to use preslit septa.

---

#### NOTE

Bottom sensing is a feature to detect the depth of vials or plates via the software. If the bottom sensing feature is used, the bottom of the plates and vials must resist the needle. Make sure that the material supports this feature.

---

# Capillary Color Coding Guide

Type		Material		Fitting Left/Fitting Right	
Key	Description	Key	Description	Key	Description
Capillary	Connection capillaries	ST	Stainless steel	W	Swagelok + 0.8 mm Port id
Loop	Loop capillaries	Ti	Titanium	S	Swagelok + 1.6 mm Port id
Seat	Autosampler needle seats	PK	PEEK	M	Metric M4 + 0.8 mm Port id
Tube	Tubing	FS/PK	PEEK-coated fused silica*	E	Metric M3 + 1.6 mm Port id
Heat exchanger	Heat exchanger	PK/ST	Stainless steel-coated PEEK**	U	Swagelok union
		PTFE	PTFE	L	Long
		FS	Fused silica	X	Extra long
				H	Long head
				G	Small head SW 4 mm
				N	Small head SW 5 mm
				F	Fingertight
				V	1200 bar
				B	Bio
				P	PEEK

\*Fused silica in contact with solvent  
 \*\*PEEK in contact with solvent

The **type** gives some indication on the primary function, like a loop or a connection capillary.  
 The **material** indicates which raw material is used.  
 The **fitting** left/right indicate which fitting is used on both ends of the capillary.

### At-a-glance color-coding keys

The color of your capillary will help you quickly identify the capillary id – see the chart to the right for reference.


Internal Diameter in mm	Color code
0.015	Orange
0.025	Yellow
0.05	Beige
0.075	Black
0.1	Purple
0.12	Red
0.17	Green
0.20/0.25	Blue
0.3	Grey
0.50	Bone White

**Tip:** As you move to smaller-volume, high efficiency columns, you'll want to use narrow id tubing, as opposed to the wider id tubing used for conventional HPLC instruments.

Figure 11 Syntax for capillary description

## Installing Capillaries

For correct installation of capillary connections of the multisampler it's important to choose the correct SL/SX fittings, see “[Capillary Color Coding Guide](#)” on page 73.



**Figure 12** Capillary connections for the multisampler

Note that the SL fittings are backward compatible to the IDEX valves.


Parts required	p/n	Description
	5067-4650	Capillary ST 0.12 mm x 150 mm SL/SX
	5067-4651	Capillary ST 0.12 mm x 280 mm SL/SX
	5067-4720	Capillary ST 0.17 mm x 150 mm SL/SX
	5067-4722	Capillary ST 0.17 mm x 280 mm SL/SX
	5065-4454	Fitting screw long 10/pk Quantity depends on configuration of the module (number of connections to the multisampler).

The capillaries mentioned above are examples only.


- 1 Select a nut that is long enough for the fitting you'll be using.




- 2 Slide the nut over the end of the tubing.



- 3 Carefully slide the ferrule components on after the nut and then finger-tighten the assembly while ensuring that the tubing is completely seated in the bottom of the end fitting.



- 4 Use a column or injection valve to gently tighten the fitting which forces the ferrule to seat onto the tubing.



**NOTE**

Don't overtighten. Overtightening will shorten the lifetime of the fitting.

## 4 Preparing the module

### Example of a perfect fitting

- 5 Once you believe you have the fitting complete, loosen the nut, and inspect the ferrule for the correct position on the tubing.


### Example of a perfect fitting



### Examples of incorrect fittings

Ferrule cannot seat properly

Mixing chamber



If dimension X is too long, leaks will occur

If dimension X is too short, a dead volume or mixing chamber will occur

Figure 13 Examples of incorrect fittings

## Setting up the Autosampler with Agilent Open Lab ChemStation

The setup of the Multisampler is shown with the Agilent OpenLab ChemStation C.01.06. Depending on the controller (e.g. Agilent Instant Pilot, OpenLab EZChrom, Masshunter) the screens look different.

### NOTE

This section describes the autosampler settings only. For information on the Agilent OpenLab ChemStation or other 1290 Infinity modules refer to the corresponding documentation.





Figure 14 ChemStation Method and Run Control

After successful load of the OpenLab ChemStation, you should see the module as an active item in the graphical user interface (GUI).

## 4 Preparing the module

### Setting up the Autosampler with Agilent Open Lab ChemStation

**Table 8** The Autosampler User Interface



Within the Multisampler user interface, there are active areas. If you move the mouse cursor across the icons (tray, EMF button), the cursor will change and you may click on the icon to

- 1 Turn on/off the autosampler
- 2 Configure the sample hotel
- 3 Get the status of the EMF (Early Maintenance Feature)
- 4 Cooling Temperatur


Current instrument information on:

- **Injection volume**
- **Sample location**

## 4 Preparing the module

### Setting up the Autosampler with Agilent Open Lab ChemStation

**Table 8** The Autosampler User Interface



A right-click into the Active Area will open a menu to

- Show the **Control** User Interface (special module settings)
- Show the **Method** User interface (same as via menu **Instrument > Set up Instrument Method > Setup G7167B**)
- **Injector Program**

When you activate a pretreatment/injector program, it replaces the standard injection cycle.

- **Identify Device**
- **Home All**
- **Switch on Tray Illumination**
- **Auto Clean**
- **Prime**
- **Drawer Configuration**

Changing the load capacity of the Sample Hotel

- **Right Capillaries**
- **Reference Vial Bar**
- **Assign Wellplates**

Wellplate Configuration (same as click on the Tray icon)

**Table 8** The Autosampler User Interface

	<p><b>Module Status</b> shows Run / Ready / Error state and “Not Ready text” or “Error text”</p> <ul style="list-style-type: none"> <li>• <b>Error</b> (Red)</li> <li>• <b>Not ready</b> (yellow)</li> <li>• <b>Ready</b> (green)</li> <li>• <b>Pre run, Post run</b> (purple)</li> <li>• <b>Run</b> (blue)</li> <li>• <b>Idle</b> (green)</li> <li>• <b>Offline</b> (dark gray)</li> <li>• <b>Standby</b> (light gray)</li> </ul>
	<p><b>EMF Status</b> shows Run / Ready / Error state and “Not Ready text” or “Error text”</p> <ul style="list-style-type: none"> <li>• Offline (gray)</li> <li>• Ok</li> <li>• No Maintenance required (green)</li> <li>• EMF warning. Maintenance might be required (yellow)</li> <li>• EMF warning. Maintenance required (red)</li> </ul>


## Control Settings

These settings are available via right click on the Active Area of the ALS GUI.

## 4 Preparing the module

Setting up the Autosampler with Agilent Open Lab ChemStation

Table 9 Control settings



**Missing Vessel:** The handling of missing vessels can be configured.

**Illumination:** Switch on/off the internal illumination  
**Linked Pump:** To configure which pump delivers flow to the Autosampler

**Prime Flush Pump:** Priming the Needle wash flush pump


**Work space:** **Clear the workspace > Immediately, At the End of the Analysis and Never**

## Method Parameter Settings

These settings are available via **Menu > Instrument > Set up Instrument Method Multisampler** or via right click on the Active area.

### NOTE

The signal window in the lower part is not shown when opening the parameter settings via right mouse on the Multisampler user interface.






Figure 15 Method parameter settings

**NOTE**

For additional help and support. Highlight the desired cell and press the **F1** key. A help screen will open with additional information and documentation about the topic.

## 4 Preparing the module


### Setting up the Autosampler with Agilent Open Lab ChemStation



#### Injection Mode/ Needle Wash

The settable **Injection volume** is depending on what kind of configuration is installed. Default configuration 0.1 – 20 µL.

It is possible to select between using the **Standard Wash** or **Standard Wash off**. Using needle wash is one option to obtain minimum carry-over.



The Injection cleaning section allows you to select between the **Standard Wash** option and the **Multi-wash** option. With the **Standard Wash** (default configuration) you can choose between two modes the Flush port or Wash Vial. If the **Multi-Wash** option is installed (additional hardware is required) you can use **Needle Wash** and **Seat Back Flush** together to obtain the lowest carry-over.



#### Multi-wash (Multisampler Injection Cleaning)

The Multi-wash table allows you to specify up to four steps that will be used to clean the system.

The Start Cond. step is not always executed. Therefore it is recommend to check the box to ensure that, at the end of the cleaning procedure, the flow path of the sampler is filled with the starting solvent conditions for the next sample.

For each cleaning step, Click the Solvent down arrow and select the solvent to use (S1, S2, S3) or switch the step Off.

Specify a duration (in seconds) in the Time [s] field.

Mark the check boxes for **Seat Back Flush** and/or **Needle Wash** to include these actions. If both are selected, they are carried out simultaneously; if neither is selected, the step is ignored (equivalent to selecting **Off**). Add a comment in the **Comment** column, if necessary.



#### Stoptime/Posttime

A Multisampler **Stoptime** can be set. For equilibration of the Multisampler a **Posttime** can be set.

## Module Configuration View

The settings are available via menu **Instrument > Instrument Configuration > Multisampler Configuration**.

Multisampler Configuration Instrument 1

Communication:

Device name: Multisampler  
 Type ID: G7167B  
 Serial number: DEA0000001  
 Firmware revision: B.06.60  
 Connection settings...

Options:

Left needle installed      Alternating needle usage not available

**Left needle:**  
 Metering: G4267-60042: 40 µL Analytical Head  
 Loop: G4267-60301: Sample Loop-Flex 20 µL left  
 Seat: G4267-87012: Seat assembly 0.12 mm 1290 In6  
 Max. Injection Volume: 20 µL  
 (Multi-draw disabled)

**Right needle:**  
 Metering: G4267-60042: 40 µL Analytical Head  
 Loop: G4267-60300: Sample Loop-Flex 20 µL right  
 Seat: G4267-87012: Seat assembly 0.12 mm 1290 In6  
 Max. Injection Volume: 20 µL  
 (Multi-draw disabled)

Cooler installed       Multi-wash installed

Reference vial array:  
 G4267-40071: 5 Vial bar

Define Sample Containers...


- **Device name:** based on the module.
- **Type ID:** based on the module (product number). Some modules may allow changing the type based on hardware/firmware. This results in a change of features and functions.
- **Serial number:** based on the module.
- **Firmware revision:** based on the module.
- **Options:** lists installed options.

### NOTE

Changes in the sampler configuration can only be done in the online view of the CDS system, see [Table 8](#) on page 79.

## 4 Preparing the module

### Setting up the Autosampler with Agilent Open Lab ChemStation



## 5 Optimizing Performance

Delay Volume and Extra-Column Volume 88

    Delay Volume 88

How to Configure the Optimum Delay Volume 89

How to Achieve Higher Injection Volumes 94

How to Achieve High Throughput 97

How to Achieve Higher Resolution 98

How to Achieve Higher Sensitivity 101

How to Achieve Lowest Carry Over 102

This chapter gives hints on how to optimize the performance or use additional devices.



## Delay Volume and Extra-Column Volume

The *delay volume* is defined as the system volume between the point of mixing in the pump and the top of the column.

The *extra-column volume* is defined as the volume between the injection point and the detection point, excluding the volume in the column.

### Delay Volume

In gradient separations, this volume causes a delay between the mixture changing in the pump and that change reaching the column. The delay depends on the flow rate and the delay volume of the system. In effect, this means that in every HPLC system there is an additional isocratic segment in the gradient profile at the start of every run. Usually the gradient profile is reported in terms of the mixture settings at the pump and the delay volume is not quoted even though this will have an effect on the chromatography. This effect becomes more significant at low flow rates and small column volumes and can have a large impact on the transferability of gradient methods. It is important, therefore, for fast gradient separations to have small delay volumes, especially with narrow bore columns (e.g., 2.1 mm i.d.) as often used with mass spectrometric detection.


## How to Configure the Optimum Delay Volume

For very fast gradients over 0.5 min the delay volume of the system can be easily reduced without changing the physical configuration of the system. The change is achieved by changing the behavior of the multisampler. The 180  $\mu\text{L}$  delay volume of the autosampler is due to the flow path from the injection valve through the metering device, needle, needle seat and connecting capillaries back to the injection valve (see [Table 10](#) on page 90). To make an injection the valve switches from mainpass to bypass so that the metering device can draw the sample into the needle capillary. The injection is made when the valve switches back to mainpass and the sample is flushed onto the column. The valve remains in this position during analysis so that the autosampler is continually flushed and hence the gradient has to flow through this delay volume to reach the column. This can be eliminated by switching the injection valve from mainpass to bypass after the injection has been made and the injected sample has been flushed onto the column. In practice this can be done a few seconds after injection and is activated by selecting the **Automatic Delay Volume Reduction (ADVR)** function in the autosampler setup menu. The Flush-out Factor (typically 5 times injection volume) ensures that enough time is allowed to flush the sample out of the injector before switching to bypass. For instance a 1  $\mu\text{L}$  injection under standard conditions effectively reduces the system delay volume by approximately 160  $\mu\text{L}$ .


## 5 Optimizing Performance

### How to Configure the Optimum Delay Volume

**Table 10** Schematic of injection steps in 1290 Infinity II Multisampler




**Figure 16** Valve in mainpass, flow through



**Figure 17** Valve in bypass, drawing sample


**Table 10** Schematic of injection steps in 1290 Infinity II Multisampler




**Figure 18** Valve in bypass, washing needle

5 Optimizing Performance  
How to Configure the Optimum Delay Volume

**Table 10** Schematic of injection steps in 1290 Infinity II Multisampler



**Figure 19** Valve in mainpass, sample injected



**Figure 20** Valve in bypass, needle backflush (Multiwash)

When using ADVR it should be noted that the gradient has already started at the pump at the instant of injection. The question should be asked whether the gradient has already reached the autosampler, in which case a small step in the gradient will result. This happens when the delay volume is less than the flush-out volume and is not necessarily a problem but may be a factor to be considered in a method transfer. With a flush-out factor of 5 and an injection volume of 10  $\mu\text{l}$ , the autosampler will allow 50  $\mu\text{l}$  to pass through before switching to bypass which, with a delay volume of 50  $\mu\text{l}$ , means the gradient just reached the injection valve. Smaller injection volumes will have no effect but for larger injection volumes this will introduce a small step in the gradient. The flow rate in use will also have an impact on the decision to use ADVR or not. At 0.2 ml/min the delay time saved is 21 seconds while at 1.0 ml/min it is 4 seconds.

The ADVR function is unlikely to be suitable for applications involving compounds which are known to cause carry-over problems. The best solution to reduce the delay volume is to install the 40  $\mu\text{L}$  Analytical Head and the 20  $\mu\text{L}$  Loop. To get the best results it is also recommended to order the Low dispersion heat exchanger and the micro flow cell for UV. This will reduce the the delay volume by 120  $\mu\text{L}$ .

## How to Achieve Higher Injection Volumes

The standard configuration of the Multisampler can inject a maximum volume of 20  $\mu\text{L}$  with the standard loop capillary. To increase the injection volume the Multidraw upgrade kit (G4216-68711) can be installed. With this kit you can add a maximum of 80  $\mu\text{L}$  to the injection volume of your injector. The total volume for the standard Multisampler is then 100  $\mu\text{L}$  or 120  $\mu\text{L}$  depending on the loop size with 40  $\mu\text{L}$  analytical head installed.

For higher injection volume you can choose between further options. This requires additional hardware modifications. One way to increase the injection volume is to change the analytical head volume. There are a 100  $\mu\text{L}$  and 900  $\mu\text{L}$  analytical heads available. Additionally you can install the Multidraw kit (G7167-68711). With the kit you can add a maximum of 400  $\mu\text{L}$  or 1400  $\mu\text{L}$  to the injection volume of your injector. The total volume is then 500  $\mu\text{L}$  or 1500  $\mu\text{L}$  for the Multisampler with a 100  $\mu\text{L}$  analytical head setup. Note the delay volume of your Multisampler is extended when using the extended seat capillaries from the multi-draw kit. When calculating the delay volume of the Multisampler you have to double the volume of the extended capillaries. The system delay volume due to the Multisampler will increase accordingly.

Whenever a method is scaled down from a larger column to a smaller column it is important that the method translation makes an allowance for reducing the injection volume in proportion to the volume of the column to maintain the performance of the method. This is to keep the volume of the injection at the same percentage volume with respect to the column. This is particularly important if the injection solvent is stronger (more eluotropic) than the starting mobile phase and any increase will affect the separation particularly for early running peaks (low retention factor). In some cases it is the cause of peak distortion and the general rule is to keep the injection solvent the same or weaker than the starting gradient composition. This has a bearing on whether, or by how much, the injection volume can be increased and the user should check for signs of increased dispersion (wider or more skewed peaks and reduced peak resolution) in trying to increase the injection size. If an injection is made in a weak solvent then the volume can probably be increased further because the effect will be to concentrate the analyte on the head of the column at the

start of the gradient. Conversely if the injection is in a stronger solvent than the starting mobile phase then increased injection volume will spread the band of analyte down the column ahead of the gradient resulting in peak dispersion and loss of resolution.

Perhaps the main consideration in determining injection volume is the diameter of the column as this will have a big impact on peak dispersion. Peak heights can be higher on a narrow column than with a larger injection on a wider column because there is less peak dispersion. With 2.1 mm i.d. columns typical injection volumes might range up to 5 to 10  $\mu\text{L}$  but it is very dependent on the chemistry of the analyte and mobile phase as discussed above. In a gradient separation injection volumes of about 5 % of the column volume might be achieved whilst maintaining good resolution and peak dispersion. One way to achieve larger injections is to use a trapping column selected by a switching valve to capture and concentrate the injection before switching it, i.e. injecting it, onto an analytical column, see [Figure 21](#) on page 96. The valve can be conveniently located in the Multicolumn Thermostat.

## 5 Optimizing Performance

### How to Achieve Higher Injection Volumes




Figure 21 Sample enrichment

## How to Achieve High Throughput

The injection can be optimized for speed remembering that drawing the sample too fast can reduce the reproducibility. Marginal gains are to be made here as the sample volumes used tend towards the smaller end of the range in any case. A significant portion of the injection time is the time taken with the needle movements to and from the vial and into the flush port. These manipulations can be performed while the previous separation is running. This is known as *overlapped injection* and it can be easily turned on from the Multisampler setup screen in the control software. The Multisampler can be told to switch the flow through the Multisampler to bypass after the injection has been made and then after, for example, 3 minutes into a 4 minutes run to start the process of aspirating the next sample and preparing for injection. This can typically save 0.5 to 1 minute per injection.

## How to Achieve Higher Resolution

Increased resolution in a separation will improve the qualitative and quantitative data analysis, allow more peaks to be separated or offer further scope for speeding up the separation. This section explains how resolution can be increased by examining the following points:

- Optimize selectivity
- Smaller particle-size packing
- Longer Columns
- Shallower gradients, faster flow

Resolution between two peaks is described by the resolution equation:

$$R_s = \frac{1}{4} \sqrt{N} \frac{(\alpha - 1)(k_2 + 1)}{\alpha k_2}$$

where

- $R_s$ =resolution,
- $N$ =plate count (measure of column efficiency),
- $\alpha$ =selectivity (between two peaks),
- $k_2$ =retention factor of second peak (formerly called capacity factor).

The term that has the most significant effect on resolution is the selectivity,  $\alpha$ , and practically varying this term involves changing the type of stationary phase (C18, C8, phenyl, nitrile etc.), the mobile phase and temperature to maximize the selectivity differences between the solutes to be separated. This is a substantial piece of work which is best done with an automated method development system which allows a wide range of conditions on different columns and mobile phases to be assessed in an ordered scouting protocol. This section considers how to get higher resolution with any chosen stationary and mobile phases. If an automated method development system was used in the decision on phases it is likely that short columns were used for fast analysis in each step of the scouting.

The resolution equation shows that the next most significant term is the plate count or efficiency,  $N$ , and this can be optimized in a number of

ways. N is inversely proportional to the particle size and directly proportional to the length of a column and so smaller particle size and a longer column will give a higher plate number. The pressure rises with the inverse square of the particle size and proportionally with the length of the column. This is the reason that the 1290 Infinity LC system was designed to go to 1200 bar so that it can run sub-two-micron particles and column length can be increased to 100 mm or 150 mm. There are even examples of 100 mm and 150 mm columns linked to give 250 mm length. Resolution increases with the square root of N so doubling the length of the column will increase resolution by a factor of 1.4. What is achievable depends on the viscosity of the mobile phase as this relates directly to the pressure. Methanol mixtures will generate more back pressure than acetonitrile mixtures. Acetonitrile is often preferred because peak shapes are better and narrower in addition to the lower viscosity but methanol generally yields better selectivity (certainly for small molecules less than about 500 Da). The viscosity can be reduced by increasing the temperature but it should be remembered that this can change the selectivity of the separation. Experiment will show if this leads to increase or decrease in selectivity. As flow and pressure are increased it should be remembered that frictional heating inside the column will increase and that can lead to slightly increased dispersion and possibly a small selectivity change both of which could be seen as a reduction in resolution. The latter case might be offset by reducing the temperature of the thermostat by a few degrees and again experiment will reveal the answer.

The van Deemter curve shows that the optimum flow rate through an STM column is higher than for larger particles and is fairly flat as the flow rate increases. Typical, close to optimum, flow rates for STM columns are: 2 ml/min for 4.6 mm i.d.; and 0.4 ml/min for 2.1 mm i.d. columns.

In isocratic separations, increasing the retention factor, k, results in better resolution because the solute is retained longer. In gradient separations the retention is described by  $k^*$  in the following equation:

$$k^* = \frac{t_G}{\Delta\%B} \cdot \frac{F}{V_m} \cdot \frac{100}{S}$$

where:

- $k^*$  = mean k value,
- $t_G$  = time length of gradient (or segment of gradient) (min),

## 5 Optimizing Performance

### How to Achieve Higher Resolution

- $F$  = flow (ml/min),
- $V_m$  = column delay volume,
- $\Delta\%B$  = change in fraction of solvent B during the gradient,
- $S$  = constant (ca. 4-5 for small molecules).

This shows that  $k$  and hence resolution can be increased by having a shallower gradient (2 to 5 %/min change is a guideline), higher flow rate and a smaller volume column. This equation also shows how to speed up an existing gradient – if the flow is doubled but the gradient time is halved,  $k^*$  remains constant and the separation looks the same but happens in half the time. Recently published research has shown how a shorter STM column (at temperatures above 40 °C) can generate higher peak capacity than a longer STM column by virtue of running it faster. (Refer to *Petersson et al., J.Sep.Sci, 31, 2346-2357, 2008, Maximizing peak capacity and separation speed in liquid chromatography*).

## How to Achieve Higher Sensitivity

The sensitivity of a separation method is linked to the choice of stationary and mobile phases as good separation with narrow peaks and a stable baseline with minimal noise are desirable. The choice of instrument configuration will have an effect and a major impact is the setup of the detector. This section considers how sensitivity is affected by:

- Pump mixer volume
- Narrower columns
- Detector flow cell
- Detector parameters

In addition, the discussion on detector parameters also mentions the related topics of selectivity and linearity.

### Columns

Sensitivity is specified as a signal-to-noise ratio (S/N) and hence the need to maximize peak height and minimize baseline noise. Any reduction in peak dispersion will help to maintain peak height and so extra-column volume should be minimized by use of short, narrow internal diameter, connection capillaries and correctly installed fittings. Using smaller inner diameter columns should result in higher peak height and is therefore ideal for applications with limited sample amounts. If the same sample amount can be injected on a smaller i.d. column, then the dilution due to column diameter will be less and the sensitivity will increase. For example, decreasing the column i.d. from 4.6 mm to 2.1 mm results in a theoretical gain in peak height of 4.7 times due to the decreased dilution in the column. For a mass spectrometer detector, the lower flow rates of narrow columns can result in higher ionization efficiencies and therefore higher sensitivity.

## How to Achieve Lowest Carry Over

Carryover is measured when residual peaks from a previous active-containing injection appear in a subsequent blank solvent injection. There will be carry over between active injections which may lead to erroneous results. The level of carryover is reported as the area of the peak in the blank solution expressed as a percentage of the area in the previous active injection. The Multisampler is optimized for lowest carryover by careful design of the flow path and use of materials in which sample adsorption is minimized. A carryover figure of 0.001 % should be achievable even when a triple quadrupole mass spectrometer is the detector. Operating settings of the Multisampler allow the user to set appropriate parameters to minimize carryover in any application involving compounds liable to stick in the system. The following functions of the Multisampler can be used to minimize carryover:

- Internal needle wash
- External needle wash
- Needle seat backflush
- Injection valve cleaning

The flow path, including the inside of the needle, is continuously flushed in normal operation, providing good elimination of carryover for most situations. Automated delay volume reduction (ADVR) will reduce the delay volume but will also reduce the flushing of the Standard Multisampler and should not be used with analytes where carryover might be a problem.

The outside of the needle can be washed using a wash vial in a specific location or the needle can be washed using the flush port. If a wash vial in a tray location specified by the user is chosen then this vial should have no septum and should contain a solvent suitable for washing the sample from the needle. The septum is not used to avoid wiping contamination off the needle on the downstream only to re-apply it on the upstroke. The needle can be dipped into the vial multiple times. This will be effective in removing a small degree of carryover but for more effective washing of the outside of the needle use the flushport.

The flush port is located above and behind the needle seat and in the standard hardware configuration a peristaltic pump delivers the wash solvent. It has a volume of 0.68 mL and the peristaltic pump delivers 5 mL/min, which means the flush port volume is completely refilled with fresh solvent in 7 s.

If the flush port is selected, the user can set how long the outside of the needle is to be washed with fresh solvent. This can last two or three seconds in routine situations where carryover is less of a problem and 10 – 20 s for more complete washing.

It is recommended that washing the outside of the needle in the flush port should be standard procedure to avoid contaminating the needle seat. If the needle seat becomes contaminated it will have to be back-flushed. In the standard setup it must be done by manually changing the flow connections to clean it.

In this standard configuration the task can be done automated by using the Flexible Cube module. If you have installed the Multisampler with Multi-Wash option the flushport will be primed with a micro piezo pump. This pump can choose between 3 different wash solvents.


The flush port and its solvent delivery pump and tubing should be regularly flushed to ensure the lowest carryover. For example, before using the system each day, prime the flush pump for three minutes with appropriate solvent.

When other measures have failed to eliminate carryover it might be that analyte is sticking inside the injector valve. With auto clean feature in the CDS system the injector valve can be set to make additional switching movements to clean out the flow path in the valve if problems occur here with carryover. If the problem compounds need a high percentage of organic phase for elution, it is recommended to switch the injection valve at the high percentage of organic phase after the last peak has eluted. It is also recommended to switch the injection valve again after the initial conditions for the mobile phase have stabilized. This ensures that the bypass groove in the rotor seal of the valve contains the gradient start conditions, which is especially important for flow rates below 0.5 mL/min. For samples where the outside of the needle cannot be cleaned sufficiently with water or alcohol from the flush pump use wash vials with an appropriate solvent. With an injector program several wash vials can be used for cleaning.

## 5 Optimizing Performance

### How to Achieve Lowest Carry Over

The optimum carry-over performance of the Multisampler is achieved after a run-in period of new instruments or after the exchange of consumable parts (like needle, needle seat and valve parts). During injections in this period, surfaces of these parts adjust to each other. After this period, we recommend back-flushing the needle seat in order to get the sealing areas between needle and needle seat clean. Regular Preventive Maintenance service is recommended as the carry-over performance of the Autosampler depends on the integrity of these consumable parts.



## 6 Troubleshooting and Diagnostics

User Interfaces 106

Agilent Lab Advisor Software 107

This chapter gives an overview about the troubleshooting and diagnostic features and the different user interfaces.



## User Interfaces

- Depending on the user interface, the available tests and the screens/reports may vary.
- Preferred tool should be Agilent Lab Advisor Software, see “[Agilent Lab Advisor Software](#)” on page 107.
- The Agilent OpenLab ChemStation C.01.03 and above do not include any maintenance/test functions.
- Screenshots used within these procedures are based on the Agilent Lab Advisor Software.

## Agilent Lab Advisor Software

The Agilent Lab Advisor Software is a standalone product that can be used with or without chromatographic data system. Agilent Lab Advisor helps to manage the lab for high-quality chromatographic results by providing a detailed system overview of all connected analytical instruments with instrument status, Early Maintenance Feedback counters (EMF), instrument configuration information, and diagnostic tests. By the push of a button, a detailed diagnostic report can be generated. Upon request, the user can send this report to Agilent for a significantly improved troubleshooting and repair process.

The Agilent Lab Advisor software is available in two versions:

- Lab Advisor Basic
- Lab Advisor Advanced

Lab Advisor Basic is included with every Agilent 1200 Infinity Series and Infinity II Series pump.


The Lab Advisor Advanced features can be unlocked by purchasing a license key, and include real-time monitoring of instrument actuals, all various instrument signals, and state machines. In addition, all diagnostic test results, calibration results, and acquired signal data can be uploaded to a shared network folder. The Review Client included in Lab Advisor Advanced allows to load and examine the uploaded data no matter on which instrument it was generated. This makes Data Sharing an ideal tool for internal support groups and users who want to track the instrument history of their analytical systems.

The optional Agilent Maintenance Wizard Add-on provides an easy-to-use, step-by-step multimedia guide for performing preventive maintenance on Agilent 1200 Infinity and Infinity II Series instruments.

The tests and diagnostic features that are provided by the Agilent Lab Advisor software may differ from the descriptions in this manual. For details, refer to the Agilent Lab Advisor software help files.

## 6 Troubleshooting and Diagnostics

Agilent Lab Advisor Software



## 7 Error Information

What Are Error Messages	110
General Error Messages	111
Timeout	111
Shutdown	111
Remote Timeout	112
Lost CAN Partner	113
Leak Sensor Short	113
Leak Sensor Open	114
Compensation Sensor Open	114
Compensation Sensor Short	115
Fan Failed	115
Leak	116

This chapter describes the meaning of error messages, and provides information on probable causes and suggested actions how to recover from error conditions.



## What Are Error Messages

Error messages are displayed in the user interface when an electronic, mechanical, or hydraulic (flow path) failure occurs which requires attention before the analysis can be continued (for example, repair, or exchange of consumables is necessary). In the event of such a failure, the red status indicator at the front of the module is switched on, and an entry is written into the module logbook.

If an error occurs outside a method run, other modules will not be informed about this error. If it occurs within a method run, all connected modules will get a notification, all LEDs get red and the run will be stopped. Depending on the module type, this stop is implemented differently. For example, for a pump the flow will be stopped for safety reasons. For a detector, the lamp will stay on in order to avoid equilibration time. Depending on the error type, the next run can only be started, if the error has been resolved, for example liquid from a leak has been dried. Errors for presumably single time events can be recovered by switching on the system in the user interface.

Special handling is done in case of a leak. As a leak is a potential safety issue and may have occurred at a different module from where it has been observed, a leak always causes a shutdown of all modules, even outside a method run.

In all cases, error propagation is done via the CAN bus or via an APG remote cable (see documentation for the APG interface).

# General Error Messages

General error messages are generic to all Agilent series HPLC modules and may show up on other modules as well.

## Timeout

**Error ID: 0062**

The timeout threshold was exceeded.

**Probable cause**

- 1 The analysis was completed successfully, and the timeout function switched off the module as requested.
- 2 A not-ready condition was present during a sequence or multiple-injection run for a period longer than the timeout threshold.

**Suggested actions**

- Check the logbook for the occurrence and source of a not-ready condition. Restart the analysis where required.
- Check the logbook for the occurrence and source of a not-ready condition. Restart the analysis where required.

## Shutdown

**Error ID: 0063**

An external instrument has generated a shutdown signal on the remote line.

The module continually monitors the remote input connectors for status signals. A LOW signal input on pin 4 of the remote connector generates the error message.

Probable cause	Suggested actions
1 Leak detected in another module with a CAN connection to the system.	Fix the leak in the external instrument before restarting the module.
2 Leak detected in an external instrument with a remote connection to the system.	Fix the leak in the external instrument before restarting the module.
3 Shut-down in an external instrument with a remote connection to the system.	Check external instruments for a shut-down condition.
4 The degasser failed to generate sufficient vacuum for solvent degassing.	Check the vacuum degasser for an error condition. Refer to the <i>Service Manual</i> for the degasser or the 1260 pump that has the degasser built-in.

## Remote Timeout

### Error ID: 0070

A not-ready condition is still present on the remote input. When an analysis is started, the system expects all not-ready conditions (for example, a not-ready condition during detector balance) to switch to run conditions within one minute of starting the analysis. If a not-ready condition is still present on the remote line after one minute the error message is generated.

Probable cause	Suggested actions
1 Not-ready condition in one of the instruments connected to the remote line.	Ensure the instrument showing the not-ready condition is installed correctly, and is set up correctly for analysis.
2 Defective remote cable.	Exchange the remote cable.
3 Defective components in the instrument showing the not-ready condition.	Check the instrument for defects (refer to the instrument's documentation).

## Lost CAN Partner

**Error ID: 0071**

During an analysis, the internal synchronization or communication between one or more of the modules in the system has failed.

The system processors continually monitor the system configuration. If one or more of the modules is no longer recognized as being connected to the system, the error message is generated.

**Probable cause**

- 1 CAN cable disconnected.
- 2 Defective CAN cable.
- 3 Defective main board in another module.

**Suggested actions**

- Ensure all the CAN cables are connected correctly.
  - Ensure all CAN cables are installed correctly.
- Exchange the CAN cable.
- Switch off the system. Restart the system, and determine which module or modules are not recognized by the system.

## Leak Sensor Short

**Error ID: 0082**

The leak sensor in the module has failed (short circuit).

The current through the leak sensor is dependent on temperature. A leak is detected when solvent cools the leak sensor, causing the leak sensor current to change within defined limits. If the current increases above the upper limit, the error message is generated.

**Probable cause**

- 1 Defective leak sensor.
- 2 Leak sensor incorrectly routed, being pinched by a metal component.

**Suggested actions**

- Please contact your Agilent service representative.
- Please contact your Agilent service representative.

## Leak Sensor Open

**Error ID: 0083**

The leak sensor in the module has failed (open circuit).

The current through the leak sensor is dependent on temperature. A leak is detected when solvent cools the leak sensor, causing the leak-sensor current to change within defined limits. If the current falls outside the lower limit, the error message is generated.

Probable cause	Suggested actions
1 Leak sensor not connected to the Power Switch board.	Please contact your Agilent service representative.
2 Defective leak sensor.	Please contact your Agilent service representative.
3 Leak sensor incorrectly routed, being pinched by a metal component.	Please contact your Agilent service representative.

## Compensation Sensor Open

**Error ID: 0081**

The ambient-compensation sensor (NTC) on the power switch board in the module has failed (open circuit).

The resistance across the temperature compensation sensor (NTC) on the power switch board is dependent on ambient temperature. The change in resistance is used by the leak circuit to compensate for ambient temperature changes. If the resistance across the sensor increases above the upper limit, the error message is generated.

Probable cause	Suggested actions
1 Loose connection between the power switch board and the main board	Please contact your Agilent service representative.
2 Defective power switch board	Please contact your Agilent service representative.

## Compensation Sensor Short

**Error ID: 0080**

The ambient-compensation sensor (NTC) on the power switch board in the module has failed (open circuit).

The resistance across the temperature compensation sensor (NTC) on the power switch board is dependent on ambient temperature. The change in resistance is used by the leak circuit to compensate for ambient temperature changes. If the resistance across the sensor falls below the lower limit, the error message is generated.

Probable cause	Suggested actions
1 Defective power switch board	Please contact your Agilent service representative.
2 Loose connection between the power switch board and the main board	Please contact your Agilent service representative.

## Fan Failed

**Error ID: 0068**

The cooling fan in the module has failed.

The hall sensor on the fan shaft is used by the main board to monitor the fan speed. If the fan speed falls below a certain limit for a certain length of time, the error message is generated.

Depending on the module, assemblies (e.g. the lamp in the detector) are turned off to assure that the module does not overheat inside.

Probable cause	Suggested actions
1 Fan cable disconnected.	Please contact your Agilent service representative.

Probable cause	Suggested actions
2 Defective fan.	Please contact your Agilent service representative.
3 Defective main board.	Please contact your Agilent service representative.


## Leak

### Error ID: 0064

A leak was detected in the module.

The signals from the two temperature sensors (leak sensor and board-mounted temperature-compensation sensor) are used by the leak algorithm to determine whether a leak is present. When a leak occurs, the leak sensor is cooled by the solvent. This changes the resistance of the leak sensor which is sensed by the leak-sensor circuit on the main board.

Probable cause	Suggested actions
1 Loose fittings.	Ensure all fittings are tight.
2 Broken capillary.	Exchange defective capillaries.



## 8 Test Functions and Calibration

Introduction	118
System Pressure Test	119
System Pressure Test Evaluation	121
Auto Referencing	122
Maintenance Positions	124
Change Needle Assembly	124
Change Sample Loop Capillary	125
Arm Position	126
Change Metering Device	126
Injector Steps	127

This chapter describes the built in test functions.



## Introduction

All tests are described based on the Agilent Lab Advisor Software B.02.06 or above. Other user interfaces may not provide any test or just a few. For details on the use of the interface refer to the interface documentation.

**Table 11** Interfaces and available test functions

Interface	Comment	Available Function
Agilent Lab Advisor	All tests are available	<ul style="list-style-type: none"> <li>• System Pressure test</li> <li>• Maintenance</li> <li>• Drawer Detection/Auto Referencing</li> </ul>
	Adding of pressure to chromatographic signals possible	
Agilent ChemStation	No tests available	<ul style="list-style-type: none"> <li>• Drawer Detection/Auto Referencing</li> <li>• Temperature mainboard</li> <li>• Pressure/Pressure ripple</li> </ul>
	Adding of pressure to chromatographic signals possible	

For details on the use of the interface refer to the interface documentation.

## System Pressure Test

The test determines the leak rate of the system between pump outlet valves and a blank nut. The blank nut can be positioned at different locations in the system before the flow cell, to determine and verify the leak rate of individual modules and components. The test allows for setting the pressure at which the test is performed. The leak rate of high pressure parts is not always a linear function and therefore it is recommended to perform the test at a pressure that corresponds to the normal operating pressure of the system.

- When**
- In case of a suspected leak
  - To verify successful execution of maintenance

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	5067-6127	Blank Nut SL

## 8 Test Functions and Calibration

### System Pressure Test

- 1 Run the System pressure test with the Agilent Lab Advisor (for further information see Online-Help of user interface).




Figure 22 System Pressure Test – Result




Figure 23 System Pressure Test – Dynamic pressure input

## System Pressure Test Evaluation

### Test Failed

#### Probable cause

- 1 Damaged blank nut (poorly shaped from over tightening)
- 2 Pump leakages
- 3 Loose or leaky fittings
- 4 Autosampler leakages
- 5 Themostatted Column Compartment valve leakages

#### Suggested actions

- Before investigating any other possible sources of failure make sure that the blank nut you are using is in a good condition and properly tightened.
- Perform the Pump Head Leak test.
- Tighten the fittings or replace capillaries.
- Perform the Autosampler Leak test.
- Replace the TCC valve rotor seal.

### NOTE

Notice the difference between *error* in the test and a *failed* result! An *error* is caused by an abnormal termination during the operation of the test, whereas a *failed* result indicates that the test result were not within the specified limits.

---

## Auto Referencing


The multisampler auto referencing uses predefined positions on the base plate and the sample hotel to calibrate the positioning of the needle parkstation and the sample hotel. The auto referencing is required to compensate deviations in positioning the needle assembly and the sample tray. The auto referencing is required after disassembling the system or when you exchange the sample handler, the sample hotel, the needle parkstation, the needle assembly or one of the main boards. This function is implemented in the drawer detection and in the needle exchange routine.

**When** After disassembling the module or an exchange of the needle assembly.

**Preparations**

- Workspace of the multisampler is empty
- All drawers are closed properly

- 1 In the Agilent Lab Advisor software select **Service & Diagnostics** in the system screen **Maintenance Positions > Change Needle**, select **YES** click **Start** and wait until the needle assembly is in maintenance position.



- 2 If no needle assembly have to be changed click **Next**. The Agilent Lab Advisor software will perform an auto referencing and reset the EMF counters.

- 3 Click the **Back** button to leave the **Service & Diagnosis** menu.



## Maintenance Positions

Some maintenance procedures require the needle assembly, the sample loop flex, the metering device and the needle seat to be moved to specific positions to enable easy access to components. The maintenance functions move these assemblies into the appropriate maintenance position. In the Agilent Lab Advisor Software the maintenance positions can be selected in the **Service & Diagnostics** view.

**When** Performing maintenance on the module

- 1 Run the Maintenance Positions in the **Service & Diagnostics** View in the Agilent Lab Advisor (for further information see Online-Help of user interface).




Figure 24 Maintenance Positions

## Change Needle Assembly

The Sample handler is positioning the needle assembly so that there is easy access for changing needle assembly or needle seat. The position is far to the left of the needle parkstation, and the current to the motors are off, so that the Z-drive of the robot can be moved while servicing the module.




Figure 25 Change Needle Assembly

## Change Sample Loop Capillary

The **Change Loop** command positions the Z-drive of the robotarm far to the left of the needle parkstation to enable easy exchange of the sample loop cartridge.




Figure 26 Change Sample Loop Capillary

## Arm Position

The home position of the multisampler ensures a better access to the workspace. When transporting the module it is highly recommended to use the **Instrument Control > Park Position** command, in order to place the Sample Handler in a position for safe transport.



Figure 27 Park Position Button

## Change Metering Device

When removing the metering device is necessary (by exchanging the metering seal for instance), the metering drive needs to be moved to a position at the far back, in order to prevent seal and/or piston damage.





Figure 28 Change Metering Device

# Injector Steps

Each movement of the sampling sequence can be done under manual control. This is useful during troubleshooting, where close observation of each of the sampling steps is required to confirm a specific failure mode or verify successful completion of a repair. Each injector step command actually consists of a series of individual commands that move the multisampler components to predefined positions, enabling the specific step to be done.

- 1 Run the **Injector Steps** in the **Service & Diagnostics** View in the Agilent Lab Advisor (for further information see Online-Help of user interface).




## 8 Test Functions and Calibration

### Injector Steps

- 2 Select the individual step command like needle selection and needle position (for further information see Online-Help of user interface).

The screenshot displays the 'Injector Steps' control panel. It includes sections for Tray Selection, Needle Selection, Needle Position, Draw Parameters, Valve, and Device Status. The Device Status section shows a 'Not Ready' error with a 'Clear Error' button. Below the settings is a table with two columns: 'Action' and 'Result'.

Action	Result
User interaction	Clear Error
Device command accepted	Ok



## 9 Maintenance

Introduction to Maintenance	131
Warnings and Cautions	133
Overview of Maintenance	135
Clean the Module	136
Removal and Installation of the Front Door	137
Remove the Needle Assembly	140
Install the Needle Assembly	144
Exchange the Needle Seat	148
Replace the Rotor Seal	151
Remove the Metering Seal	157
Install the Metering Seal	162
Replace the Peristaltic Pump Cartridge	168
Replace the Flushhead Seal	172
Replace the Injection Valve	176
Removing the Sample Loop-Flex	180
Installing the Sample Loop-Flex	184
Replace the Dummy Drawer	187
Optional Configurations	187
Installing and Replacing of Drawers (Upgrade Drawer Kit)	189
Configuration of the Hotel Drawers	192
Remove the Sample Cooler	196
Install the Sample Cooler	197
Replace the Module Firmware	200



## 9 Maintenance

### Injector Steps

This chapter describes the maintenance of the Multisampler

## Introduction to Maintenance

Figure 29 on page 131 shows the main user accessible assemblies of the multisampler. These parts can be accessed from the front (simple repairs) and don't require to remove the multisampler from the system stack.




Figure 29 Main user accessible assemblies (standard)

## 9 Maintenance

### Introduction to Maintenance




Figure 30 Main user accessible assemblies (multiwash)




Figure 31 Overview of drawer, sample tray and sample container

## Warnings and Cautions

### WARNING

#### Personal injury or damage to the product

Agilent is not responsible for any damages caused, in whole or in part, by improper use of the products, unauthorized alterations, adjustments or modifications to the products, failure to comply with procedures in Agilent product user guides, or use of the products in violation of applicable laws, rules or regulations.

- Use your Agilent products only in the manner described in the Agilent product user guides.
- 

### WARNING

#### Electrical shock

Repair work at the module can lead to personal injuries, e.g. shock hazard, when the cover is opened.

- Do not remove the cover of the module.
  - Only certified persons are authorized to carry out repairs inside the module.
- 

### WARNING

#### Sharp metal edges

Sharp-edged parts of the equipment may cause injuries.

- To prevent personal injury, be careful when getting in contact with sharp metal areas.
-

**WARNING**

**Toxic, flammable and hazardous solvents, samples and reagents**

**The handling of solvents, samples and reagents can hold health and safety risks.**

- When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
  - The volume of substances should be reduced to the minimum required for the analysis.
  - Do not operate the instrument in an explosive atmosphere.
- 

**CAUTION**

**Safety standards for external equipment**

- If you connect external equipment to the instrument, make sure that you only use accessory units tested and approved according to the safety standards appropriate for the type of external equipment.
-

## Overview of Maintenance

It is necessary to perform periodic inspection of this instrument to ensure its safe use. It is possible to have these periodic inspections performed by Agilent service representatives on a contractual basis. For information regarding the maintenance inspection contract, contact your Agilent representative.

The following pages describe the maintenance (simple repairs) of the module that can be carried out without opening the main cover.

**Table 12** Overview of maintenance

Procedure	Typical interval (minimum)	Notes
Change needle/needle seat	60000 needle into seat	
Change peristaltic pump cartridge	3000 min on time	
Change rotor seal	30000 injections	

## Clean the Module

To keep the module case clean, use a soft cloth slightly dampened with water, or a solution of water and mild detergent.

### **WARNING**

**Liquid dripping into the electronic compartment of your module can cause shock hazard and damage the module**

- Do not use an excessively damp cloth during cleaning.
  - Drain all solvent lines before opening any connections in the flow path.
-

## Removal and Installation of the Front Door

**When** If the front door is defective or the hinge are damaged.

**Tools required** **Description**  
Flat screwdriver

**Parts required** **p/n** **Description**  
5067-5415 Door Assy

**Preparations** Finish any pending acquisition job and return any plate on the workspace back to the hotel.

### NOTE

For detailed information on position of the magnets, refer to [“Magnets”](#) on page 44

### CAUTION

Magnetic fields

Magnets produce a far-reaching, strong magnetic field.

You can damage for example televisions, laptops, computer harddisks, credit cards, magnetic cards may be damaged as well.

→ Keep magnets at least 25 mm away from devices and objects that could be damaged by strong magnetic fields.

### WARNING

Heart pacemakers


Magnets may derogate the function of heart pacemakers and implanted defibrillators. Heart pacemakers may switch to the test modus which may cause malaise. Defibrillators may malfunction.

→ Bearers of heart pacemakers or implanted defibrillators must stay off at least 55 mm from the magnets.



## 9 Maintenance

### Removal and Installation of the Front Door


1 Open the front door.



2 Press the release buttons and pull the front door out.



- 3 For the Installation of the front door. Insert the hinges into their guides and move the door in until the release buttons click into their final position.



## Remove the Needle Assembly

<b>When</b>	When the limit in the needle into seat counter in the EMF is exceeded or when needle shows indications of damage, blockage or leaks.		
<b>Tools required</b>	<b>p/n</b>	<b>Description</b>	
	8710-0510	Wrench open 1/4 — 5/16 inch	
<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	G4267-87201	Needle Assembly
OR	1	G4267-87210	Needle Assembly (slotted) for high injection volumes
<b>Preparations</b>	In order to avoid leaks, stop the pump running and remove the tubings from the solvent bottles. If available close the shutoff valves.		

### WARNING

#### Risk of injury by uncovered needle

**An uncovered needle is a risk of harm to the operator.**

- Do not open the safety lock of the needle assembly
- Be careful working at the z-robot.
- Wear safety goggles, when removing the needle assembly.

### WARNING

#### Toxic, flammable and hazardous solvents, samples and reagents

**The handling of solvents, samples and reagents can hold health and safety risks.**

- When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.

### NOTE


It is recommended to always exchange the needle assembly and the needle seat at the same time to prevent premature leakage.

- 1 In the Instant Pilot start the maintenance mode and select **Change needle/seat** function.

OR

In the Agilent Lab Advisor software select **Service & Diagnostics** in the system screen (**Tools**) **Maintenance Positions > Change Needle/Loop**, click **Start** and wait until the needle assembly is in maintenance position.

- 2 Open the front door.



**WARNING**


Sharp needle  
Uncovered needles may cause injuries

→ Do not unlock the safety lock of needle assembly.


## 9 Maintenance

### Remove the Needle Assembly

- 3 Remove the needle assembly by slightly pulling the needle cartridge.



- 4 Z-Robot (Z-arm coupler) without the needle assembly.




#### CAUTION

**Damage of the loop**  
The loop shape may be damaged if the loop is stretched or bent too far.


- Avoid to change the loop shape.
- Do not pull or bend the loop too far.

Remove the Needle Assembly


5 The needle assembly is still connected to the loop capillary.




6 Remove the loop plastic adapter.



7 Use a 1/4 inch wrench to loosen the fitting of the loop capillary.



8 Remove the needle assembly.



## Install the Needle Assembly

<b>When</b>	When the limit in the needle into seat counter in the EMF is exceeded or when needle shows indications of damage, blockage or leaks.		
<b>Tools required</b>	<b>p/n</b>	<b>Description</b>	
	8710-0510	Wrench open 1/4 — 5/16 inch	
<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	G4267-87201	Needle Assembly
OR	1	G4267-87210	Needle Assembly (slotted) for high injection volumes
<b>Preparations</b>	In order to avoid leaks, stop the pump running and remove the tubings from the solvent bottles. If available close the shutoff valves.		

### WARNING

#### Risk of injury by uncovered needle


**An uncovered needle is a risk of harm to the operator.**

- Do not open the safety lock of the needle assembly
- Be careful working at the z-robot.
- Wear safety goggles, when removing the needle assembly.

### NOTE

It is recommended to always exchange the needle assembly and the needle seat at the same time to prevent premature leakage.

- 1 Install the loop capillary on top of the needle cartridge (1.) and tighten the fitting hand tight (2.).



### CAUTION


#### Blockages

- Do not overtighten the fitting. A quarter turn should be sufficient.
- 2 Use a 1/4 inch wrench to tighten the fitting of the loop capillary.

## 9 Maintenance

### Install the Needle Assembly



3 Install loop plastic adapter.



4 Pinch and reinsert the needle assembly and the connected loop capillary into the z-arm coupler.

#### NOTE

Check the tension of the loop capillary. This must be forced and guided to the hydraulic box to prevent it from being caught by the Z-drive.




5 In the Instant Pilot close **Change needle /seat**.

OR

In the Agilent Lab Advisor software **Change needle/loop**  
> **End**, click **End** and wait until the needle assembly is in  
the needle park station.

6 Close the front door.



7 Perform a pressure test.

## Exchange the Needle Seat

**When** When seat is visibly damaged, blocked or leaks.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	Wrench open 1/4 — 5/16 inch Flat head screwdriver

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	G4267-87012	High Pressure Needle Seat, 0.12 mm (PEEK)

**Preparations** In order to avoid leaks, stop the pump running and remove the tubings from the solvent bottles. If available close the shutoff valves.

### WARNING


#### Risk of injury by uncovered needle

**An uncovered needle is a risk of harm to the operator.**


- Do not open the safety lock of the needle assembly
- Be careful working at the z-robot.
- Wear safety goggles, when removing the needle assembly.

- 1 In the Instant Pilot start the maintenance mode and select **Change needle/seat** function.  
  
OR  
  
In the Agilent Lab Advisor software select **Service & Diagnostics** in the system screen **Maintenance Positions > Change Needle**, click **Start** and wait until the needle assembly is in maintenance position.


- 2 Open the front door.




3 Disconnect the seat capillary from the Injection valve.




4 With a flat head screw driver carefully lift out the needle seat (1.). Then remove the complete leak tube needle seat assembly from the holder (2.).



5 Insert the new Needle seat (1.). Press it firmly in position (2.).



6 Reconnect the seat capillary to the injection valve.



## 9 Maintenance


### Exchange the Needle Seat

7 In the Instant Pilot close **Change needle /seat**.

OR

In the Agilent Lab Advisor software **Change needle** click **End** and wait until the needle assembly is in the needle park position.

8 Close the front door.



9 Perform a pressure test.

## Replace the Rotor Seal

**When** When poor injection volume reproducibility or when injection valve is leaking.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	Wrench open 1/4 — 5/16 inch
	8710-2394	Hex key 9/64 inch 15 cm long T-handle
		Cleaning tissue and appropriate solvent like isopropanol or methanol

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	5068-0198	Rotor Seal for VICI Injection Valve (PAEK)

### CAUTION

Reduced life time of the injection valve


Component cleanliness is crucial for the life time of the injection valve.

→ Replace the rotor seal in a clean environment.

## 9 Maintenance

### Replace the Rotor Seal


- 1 Open the front door.




- 2 Remove all capillaries from the injection valve with a 1/4 inch wrench.

#### NOTE

Remember the correct plumbing.



- 3 Use a 9/64 inch hex driver to unscrew the two socket screws which hold the stator head in place.




**CAUTION**

**Damage to the stator head**

The polished sealing surface of the stator head contains six ports that access handling can easily damage.

- Avoid touching the polished surface of the stator head.
- Never place the polished surface on a hard surface.


- 4 Carefully remove the stator head. To ensure that the sealing surface of the stator head is not damaged, place it on its outer face.



## 9 Maintenance

### Replace the Rotor Seal

5 Remove the rotor seal.



#### NOTE

Remove the rotor seal with a small tool, gently pry the rotor seal away from the drive.

Examine the rotor sealing surface for scratches and nicks.



- If scratches are visible the rotor seal must be replaced.
- If no scratches are visible clean all the parts with an appropriate solvent, taking care that no surfaces get scratched.

#### CAUTION

Damage to the rotor seal and cross-port leaks

- Before you replace the rotor seal, clean the stator.
- Inspect the stator head and swab it with the appropriate solvent. If more stringent cleaning is required, use a sonicator. Inspect the remaining valve components for contamination. Clean them as necessary.
- If the stator head is scratched, replace the valve.


6 Install new rotor seal.




#### NOTE

Make sure that the rotor sealing surface with its engraved flow passages is facing out. The pattern is asymmetrical to prevent improper placement.

- 7 Reinstall the stator head. The index pins on the drive and the stator head must engage in the corresponding holes. Insert the two socket head screws.



- 8 Using a 9/64 in. L-Hex wrench, tighten each screw gently until you feel resistance (approximately fingertight). Tighten each screw by 1/8 turn, and then tighten each screw again, until the stator is secured to the driver.




**NOTE**


Do not over-tighten the screws. The screws hold the assembly together and do not affect the sealing force. The sealing force is automatically set as the screws close the stator head against the valve body.

**9 Maintenance**  
Replace the Rotor Seal

**9** Reconnect all capillaries to the proper injection valve ports with a 1/4 inch wrench



**10** Close the front door.



**11** Perform a pressure test.

## Remove the Metering Seal


**When** When poor injection volume reproducibility or when metering device / analytical head is leaking.

Tools required	p/n	Description
	8710-0510	Wrench open 1/4 — 5/16 inch
	8710-2392	4 mm Hex key
	01018-23702	Insert tool

Parts required	#	p/n	Description
	1	0905-1717	Metering seal 40 µL for 40 µL analytical head
	1	0905-1719	Metering seal 100 µL for 100 µL analytical head
	1	5067-5920	Piston ceramic If previous piston is scratched


**1** In the Instant Pilot start the maintenance mode and select **Change metering device** function.  
OR  
In the Agilent Lab Advisor software select **Service & Diagnostics** in the system screen (**Tools**) > **Maintenance Positions** > **Change Metering Device**, click start and wait until the metering device is in maintenance position.




## 9 Maintenance

### Remove the Metering Seal


3 Disconnect all capillaries from the metering device.




4 To release the bayonet lock, push (1.) and rotate (2.) the analytical head a quarter left. Then you can pull and detach the analytical head assembly from the actuator (3.).




5 Remove the metering device.




6 Take the metering device. Push against the rear side of the metering device and rotate a quarter left to release the bayonet lock.



7 Now you can separate the analytical head and head body.




8 Remove the piston out of the head body.



## 9 Maintenance

### Remove the Metering Seal

9 Inspect the piston for cleanliness and scratches.




- If dirty:  
Clean the piston with an appropriate solvent.
- If scratched:  
Replace the piston by a new one.

10 Take the analytical head and remove the three screws on the rear side, which holds the support ring in place. Check the support ring for any damages.



**11** Carefully remove the metering seal using the steel side of the insert tool. Clean the chamber with an appropriate solvent and ensure that all particulate matter is removed.



## Install the Metering Seal


**When** After removing the metering seal.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	Wrench open 1/4 — 5/16 inch
	8710-2392	4 mm Hex key
	01018-23702	Insert tool
		Cleaning tissue and appropriate solvent like isopropanol or methanol


<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	0905-1717	Metering seal 40 $\mu$ L for 40 $\mu$ L analytical head
	1	0905-1719	Metering seal 100 $\mu$ L for 100 $\mu$ L analytical head
	1	5067-5920	Piston ceramic If previous piston is scratched

**Preparations** Removing the metering seal, see [“Remove the Metering Seal”](#) on page 157

1 Install the new metering seal using the plastic side of the insert tool. Press it firmly into position. Avoid any offset angle as it might deform the seal.



2 Reassemble the support ring.




## 9 Maintenance

### Install the Metering Seal


- 3 Make sure to comply to the following order of actions:
- a Tighten the three screws fingerthight, then
  - b Tighten the screws a little at a time to keep the support ring surface *parallel* (important!) to the surface of the analytical head.




- 4 Use the twist and lock bayonet mechanisms to reassemble the analytical head assembly. Push the two parts together to couple the head body with the analytical head. Once the pin reaches the bottom of the slot, one or both parts are rotated so that the pin slides along the horizontal arm of the L until it reaches the *serif*. The spring then pushes the male connector up into the *serif* to keep the pin locked into place.



- 5 Press the piston carefully into the housing of the head body and the seal.



- 6 Reinstall the complete analytical head with the actuator housing





**NOTE**

For proper installation, check the correct position of the tag.


## 9 Maintenance

### Install the Metering Seal


- 7 Fix the analytical head by pushing (1.) and rotating (2.) via twist and lock bayonet mechanism.



- 8 Reconnect the capillaries.



9 Close the front door.



In the Instant Pilot exit the maintenance mode and select **Change metering device** function.

OR

In Agilent Lab Advisor software system screen exit **Service & Diagnostics (Tools) > Maintenance Positions > Change Metering Device** click **End** and wait until the metering device is in **Home** position.

10 Perform a pressure test.

## 9 Maintenance

### Replace the Peristaltic Pump Cartridge

# Replace the Peristaltic Pump Cartridge

**When** Tubing blocked or broken

**Tools required** Description

Parts required	#	p/n	Description
	1	5065-4445	Peristaltic pump with Pharmed tubing (default)
OR	1	5042-8507	Peristaltic pump cartridge, silicone tubing
OR	1	5042-9952	Peristaltic pump with Chemsure tubing

**Preparations** Remove the inlet filter of the solvent bottle which guides the solvent to the peristaltic pump to avoid syphoning effects.

### WARNING

**When opening capillary or tube fittings solvents may leak out.**

**The handling of toxic and hazardous solvents and reagents can hold health risks.**


→ Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

### NOTE


The peristaltic pump cartridge is a replaceable unit. The tubing inside the pump is not replaceable.

## Replace the Peristaltic Pump Cartridge


1 Open the front door.




2 Press the two clips on the front of the peristaltic pump cartridge.



3 Pull the cartridge forward off the motor shaft.




4 Disconnect the tubing coupler leading to the wash port and the tubing coupler coming from the solvent bottle.




## 9 Maintenance

### Replace the Peristaltic Pump Cartridge

- 5 Push the new cartridge onto the motor shaft until the clips click into place.




- 6 Connect the wash port tubing to the upper tubing of the new cartridge (use sand paper to get a good grip on the tubing).




## Replace the Peristaltic Pump Cartridge

- 7 Connect the inlet filter of the solvent bottle again. Use the syringe to draw enough solvent for completely filling of the peristaltic pump tubing before continuing to prime the peristaltic pump.



- 8 Close the front door.



## Replace the Flushhead Seal

**When** Flush head is leaking


<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	Wrench open 1/4 — 5/16 inch
	8710-2392	Hex key 4 mm15 cm long T-handle
	01018-23702	Insert tool
OR	G4226-43800	Seal insert tool

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	5067-5918	Seal 500 µL


- Preparations**
- Cleaning tissue
  - Appropriate solvent like isopropanol or methanol

- 1** In the Instant Pilot start the maintenance mode and select **Change metering device** function.
- OR
- In the Agilent Lab Advisor software select **Service & Diagnostics** in the system screen (**Tools**) > **Maintenance Positions** > **Change Metering Device**, click start and wait until the metering device is in maintenance position.


- 2** Open the front door.




3 Remove capillaries and valves from the flush head.




4 Press and turn the Flush Head a quarter left (bayonet fitting) and detach the metering device from the actuator.



5 Pull the flush head away from the hydraulic box



6 Press against the rear side of flush head and turn a quarter left (bayonet fitting) and separate the flush head, head body and the piston.




**NOTE**

Be careful not to break the piston.


## 9 Maintenance

### Replace the Flushhead Seal


7 Remove the piston from the head body.




8 Carefully remove the metering seal from the tip of the piston.




9 Reassemble the flush head and the head body (without piston).




10 Carefully insert the piston with the new metering seal into the flush head assembly.




11 Reinstall the flush head to the actuator housing.




12 Fix the flush head.



13 Connect the capillaries.



14 Close the front door.



## 9 Maintenance

### Replace the Injection Valve

# Replace the Injection Valve


**When** Add new injection valve or replace defective injection valve.

**Tools required** **Description**  
Wrench 9/64


<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	5067-4232	2pos/6port Injection Valve (VICI) 1300 bar (G7167B)
OR	1	5067-4230	Injection Valve Idex 600 bar (G7167A)

**Preparations** Switch off the power of the Multisampler

**1** Disconnect the capillaries.




**2** Turn the spanner nut counter clockwise until the injection valve head detaches from the hydraulic box (Do not use wrenches on the spanner nut).




## Replace the Injection Valve

- 3 Remove the spanner nut from the injection valve head.




- 4 Take the replacement injection valve head and insert it into the open actuator slot of the hydraulic box. Rotate until the unions at the base of the replacement injection valve head and the valve actuator engage



OR

If the outside pin does not fit into the outside groove, you have to turn the valve head until you feel that the two pins snap into the grooves. Now you should feel additional resistance from the valve drive while continue turning the valve head until the pin fits into the groove.



**NOTE**



Check the orientation of the rear side.

Verify the correct position of the Valve TAG.


## 9 Maintenance

### Replace the Injection Valve


- 5 Continue to rotate until the clocking pin in the injection valve head align with the notch in the housing and press the replacement injection valve head into the actuator.



- 6 Replace the spanner nut (1.) and tighten clockwise (2.) (Hand tighten only, do not use wrenches on the spanner nut).



7 Reconnect the capillaries



## Removing the Sample Loop-Flex

**When** If the sample loop flex is defective or damaged.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	Wrench open 1/4 — 5/16 inch

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	G4267-60300	Sample Loop Flex 20 µL, right (red coded)
	G4267-60400	Sample Loop Flex 40 µL, right (green coded)
	G4267-60500	Sample Loop Flex 100 µL, right (blue coded)

**Preparations** Finish any pending acquisition job and return any plate on the workspace back to the hotel.

### WARNING

#### Risk of injury by uncovered needle

**An uncovered needle is a risk of harm to the operator.**


- Do not open the safety lock of the needle assembly
- Be careful working at the z-robot.
- Wear safety goggles, when removing the needle assembly.

**1** In the Instant Pilot start the maintenance mode and select **Change needle/seat** function.


OR

In the Agilent Lab Advisor software select **Service & Diagnostics** in the system screen (**Tools**) **Maintenance Positions > Change Needle/Loop**, click **Start** and wait until the needle assembly is in maintenance position.

**2** Open the front door.



- 3 The needle assembly is still connected to the loop capillary. Use a 1/4 inch wrench to loosen the fitting of the loop capillary connected to the analytical head.



**CAUTION**

**Damage of the loop**  
The loop shape may be damaged if the loop is stretched or bent too far.


- Avoid to change the loop shape.
- Do not pull or bend the loop too far.

**WARNING**

**Sharp needle**  
Uncovered needles may cause injuries

- Do not unlock the safety lock of needle assembly.



- 4 Remove the needle assembly by slightly pulling the needle cartridge.




## 9 Maintenance

### Removing the Sample Loop-Flex


- 5 Remove the cartridge out of its proper position. By gently tilting and pulling it out of the work space of the multisampler.




- 6 Remove the loop plastic adapter.



7 Use a 1/4 inch wrench to loosen the fitting of the loop capillary.



8 Remove the needle assembly.



## Installing the Sample Loop-Flex

**When** If the sample loop flex is defective or damaged.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	Wrench open 1/4 — 5/16 inch

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	G4267-60300	Sample Loop Flex 20 $\mu$ L, right (red coded)
	G4267-60400	Sample Loop Flex 40 $\mu$ L, right (green coded)
	G4267-60500	Sample Loop Flex 100 $\mu$ L, right (blue coded)

**Preparations** Finish any pending acquisition job and return any plate on the workspace back to the hotel.


### **WARNING**

#### **Risk of injury by uncovered needle**

**An uncovered needle is a risk of harm to the operator.**

- Do not open the safety lock of the needle assembly
  - Be careful working at the z-robot.
  - Wear safety goggles, when removing the needle assembly.
-

- 1 Install the loop capillary on top of the needle cartridge (1.) and tighten the fitting hand tight (2.).




**CAUTION**


**Blockages**

- Do not overtighten the fitting. A quarter turn should be sufficient.
- 2 Then use a 1/4 inch wrench to tighten the fitting of the loop capillary.

- 3 Install loop plastic adapter.




- 4 Click the sample loop cartridge in the designated location and keep the right orientation.




## 9 Maintenance

### Installing the Sample Loop-Flex

- 5 Install the shorter capillary of the sample loop cartridge to the analytical head.




- 6 Pinch and reinsert the needle assembly and the connected sample loop capillary into the z-arm coupler.



#### NOTE

Check the tension of the loop capillary. This must be forced and guided to the hydraulic box to prevent it from being caught by the Z-drive.

- 7 Close the front door.



In the Instant Pilot close **Change needle /seat**.

OR

In Agilent Lab Advisor software **Change needle/loop**. Click **NEXT** and wait until the needle is in the needle park station.

Click **Back** to leave the Maintenance window.



#### NOTE

If you need an autoreferencing step included you must choose the change needle procedure

# Replace the Dummy Drawer

## Optional Configurations

Table 13 Overview on optional configurations (examples for uniform types)

	1H	2H	3H	Dummy-Drawer
 <p><b>Delivery Status</b></p>	-	G7167-60020 1x	-	G4267-60024 3x
 <p><b>Up to 8 single height drawers</b> 16 positions Shallow wellplates and MTP Max Sample capacity 1536 / 6144 samples (96 Shallow Wellplates / 384 MTP)</p>	G7167-60021 8x	-	-	-
 <p><b>Up to 4 Dual Height drawers</b> 8 positions Vials (2 mL), deep well plates, MTP, Eppendorf Max Sample capacity 432 / 3072 samples (2 mL Vials/ 384 MTP)</p>	-	G7167-60020 4x	-	-
 <p><b>Up to 2 Drawers Triple Height</b> 4 positions (2H or 2*1H option left over) Vials (6 ml), deep well plates, MTP, Eppendorf Max Sample capacity 60 / 216/ 1536 samples (6 mL Vials/ 2 mL Vials/ 384 MTP)</p>	-	G7167-60020 1x	G7167-60022 2x	-

## 9 Maintenance

### Replace the Dummy Drawer

#### NOTE

Mixed configurations are possible (for example 1x3H- with 1x2H- and 3x1H-drawer).

All positions in the Sample Hotel must be filled either with dummies or drawers. The drawers must be installed from bottom to top.

---

## Installing and Replacing of Drawers (Upgrade Drawer Kit)


<b>Tools required</b>	<b>Description</b>	
	Screwdriver	
<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	G7167-60020	Drawer 2H
	G7167-60021	Drawer 1H
	G7167-60022	Drawer 3H

**NOTE**

Before you start the new drawer installation you have to remove the lower drawer (2H drawer = default configuration) from the Sample Hotel.

**NOTE**


For the best cooling performance avoid the installation of a 3H drawer in the lowest position.




## 9 Maintenance

### Replace the Dummy Drawer


- 3 Unlatch the drawer: Use a screwdriver to press the clamping lever lightly to the left.




- 4 Remove the drawer from the rail guide.



The drawer is now out of the hotel.




- 5 Grab in the recession below the dummy drawer front panel (1.) and lift the left side (2.).



- 6 Remove the dummy drawer.

**NOTE**


At this stage remove all other dummies that will be replaced by hotel drawers.



## 9 Maintenance

### Replace the Dummy Drawer


- 7 Place the new drawer horizontally into the sample hotel. Check that the drawer matches the middle bracket of the sample hotel.



- 8 Push until the complete drawer locks in place.

#### NOTE

Take care that the clamping lever locks.



#### NOTE

Always fill sample hotel completely (no empty drawer slots). Otherwise the drawers can't be configured in the software.

- 9 Configure the hotel drawers in the controller software (see the Online Help of the software for details).

## Configuration of the Hotel Drawers


The configuration of your drawers is necessary to detect the new drawer configuration for your CDS system. When a wrong configuration is detected there will be a mismatch in your CDS system and you are not able to use the new drawers. The new drawer configuration is active and stored after you have done the Drawer Configuration.

## Configure the Hotel Drawers in the Control Software

**Software required** OpenLAB (A.02.01 or above)  
 LC driver (A.02.10 or above)

- Preparations**
- Stop the acquisition run.
  - Remove the sample containers (trays and well plates) from workspace.
  - Complete the drawer installation.
  - Remove the sample containers (trays and well plates) from the drawers.
  - Verify that all sample trays (palettes) are installed in their drawers.
  - All open drawers and dummies have to be closed and installed properly.

- 1 Start OpenLAB CDS ChemStation Edition.
- 2 Right-click on the **Multisampler** GUI.




- 3 Select **Modify > Drawer Configuration** in the GUI screen.

**NOTE**

For correct detection, it is necessary to remove all sample containers (for example 54 vial tray or well plates).

## 9 Maintenance

### Replace the Dummy Drawer




- 4 Follow the Setup or Change configuration screen.
- 5 System is ready after the robot has done Auto Referencing.

## Configure the Hotel Drawers in Lab Advisor

**Software required** Lab Advisor (B.02.05 or above)

- Preparations**
- Stop the acquisition run.
  - Remove the sample containers (trays and well plates) from workspace.
  - Complete the drawer installation.
  - Remove the sample containers (trays and well plates) from the drawers.
  - Verify that all sample trays (palettes) are installed in their drawers.
  - All open drawers and dummies have to be closed and installed properly.

- 1 Start the Lab Advisor Software.
- 2 Connect the instrument and select **Instrument Control** in the system screen.
- 3 Switch In the **Configuration** menu of the Multisampler. Select **Detect Drawers** in the **Hotel Configuration**.



- 4 Follow the Detect Hotel Configuration screen to detect the physically available drawers.

**NOTE**

For correct detection, it is necessary to remove all sample containers (for example 54 vial tray or well plates).

- 5 System is ready after the robot has done Auto Referencing.

## Remove the Sample Cooler

### WARNING

#### Heavy weight

The module is heavy (>22 kg (>46 lbs)).

- Carry the module at least with 2 people.
- Avoid back strain or injury by following all precautions for lifting heavy objects.
- Ensure that the load is as close to your body as possible.
- Ensure that you can cope with the weight of your load.


### CAUTION

#### Routing of the condensation tubing

Proper routing of the condensation tubing is critical for correct condensate drainage.

- Do not place the multisampler directly on the bench. Install Base plate (G1328-44121) underneath.

1



Power switch

- (1) On
- (2) Off

#### Next Steps:

- 2 Remove the power cable from the module.
- 3 Open the four screws on cooler cover.
- 4 Slide the sample cooler the half way out.
- 5 Remove power and the signal cable.
- 6 Slide the cooler completely out.
- 7 Place the sample cooler on the bench.

## Install the Sample Cooler

**When** If the cooler is damaged or defective.

**Tools required** **Description**  
Screwdriver, Pozidriv #1 PT3

**Parts required**

p/n	Description
G4267-60005	Sample Cooler

<p><b>1</b> Slide in halfway</p> 	<p><b>2</b></p>  <p><b>CAUTION</b></p> <p>Damage to the cables</p> <ul style="list-style-type: none"><li>→ Do not bend or pinch the cables.</li><li>→ Fit in the cooler perfectly.</li></ul>
---	---

## 9 Maintenance

### Install the Sample Cooler


3



4



5 Install the condensate tubing correctly.




#### CAUTION

Damage to the sample cooler

→ Wait at least 30 min before switching on the compressor of the sample cooler.

6





## Replace the Module Firmware

<b>When</b>	The installation of newer firmware might be necessary <ul style="list-style-type: none"><li>• if a newer version solves problems of older versions or</li><li>• to keep all systems on the same (validated) revision.</li></ul> The installation of older firmware might be necessary <ul style="list-style-type: none"><li>• to keep all systems on the same (validated) revision or</li><li>• if a new module with newer firmware is added to a system or</li><li>• if third party control software requires a special version.</li></ul>
-------------	---

<b>Tools required</b>	<b>Description</b>
	Agilent Lab Advisor software
OR	Instant Pilot G4208A (only if supported by module)

<b>Parts required</b>	<b>#</b>	<b>Description</b>
	1	Firmware, tools and documentation from Agilent web site


**Preparations** Read update documentation provided with the Firmware Update Tool.

**To upgrade/downgrade the module's firmware carry out the following steps:**

- 1 Download the required module firmware, the latest LAN/USB FW Update Tool and the documentation from the Agilent web.  
[http://www.chem.agilent.com/\\_layouts/agilent/downloadFirmware.aspx?wid=69761](http://www.chem.agilent.com/_layouts/agilent/downloadFirmware.aspx?wid=69761)
- 2 For loading the firmware into the module follow the instructions in the documentation.

*Module Specific Information*

There is no specific information for this module.



## 10 Parts for Maintenance and Upgrade or Options

Overview of Maintenance Parts	202
Sampler Main Assemblies	203
Hotel Drawer	205
Analytical Head Assembly 40 $\mu$ L	206
Analytical Head Assembly 100 $\mu$ L	207
Analytical Head Assembly 900 $\mu$ L	208
Flush Head Assembly 500 $\mu$ L	209
Injection Valve Assembly	211
Injection Valve with Actuator	213
Needle Port Assembly	214
Door Assy	215
Accessory Kit	216
Tubing Kit Sampler Standard	217
Sample Cooler	218

This chapter provides information on parts material required for the module.



## 10 Parts for Maintenance and Upgrade or Options

### Overview of Maintenance Parts

## Overview of Maintenance Parts

<b>p/n</b>	<b>Description</b>
0905-1717	Metering seal 40 µL
0905-1719	Metering seal 100 µL
5067-5918	Metering seal 500 µL
5068-0007	Injection valve rotor seal
G4267-87201	Needle Assembly
G4267-87012	High Pressure Needle Seat, 0.12 mm (PEEK)
G4267-60300	Sample Loop Flex 20 µL, right (red coded)
G4267-60400	Sample Loop Flex 40 µL, right (green coded)
G4267-60500	Sample Loop Flex 100 µL, right (blue coded)
5065-4445	Peristaltic pump with Pharmed tubing

## Sampler Main Assemblies

p/n	Description
G4267-60002	Hotel Assembly
G7167-60020	Drawer 2H
G7167-60021	Drawer 1H
G7167-60022	Drawer 3H
G4267-60001	Sample Handler Assembly
5067-4232	2pos/6port Injection Valve (VICI)
G4267-87201	Needle Assembly
G4267-60044	Needle Port Assembly Station
G4267-60045	Flush head, 500 $\mu$ L
G4267-87012	High Pressure Needle Seat, 0.12 mm (PEEK)
G4267-60300	Sample Loop Flex 20 $\mu$ L, right (red coded)
G4267-60400	Sample Loop Flex 40 $\mu$ L, right (green coded)
G4267-60500	Sample Loop Flex 100 $\mu$ L, right (blue coded)
G4267-60050	Baseplate Assembly

## 10 Parts for Maintenance and Upgrade or Options

### Sampler Main Assemblies




Figure 32 Main user accessible assemblies (standard)




Figure 33 Main user accessible assemblies (multiwash)

# Hotel Drawer

Item	p/n	Description
1	G7167-60021	Drawer 1H (2 p/k (including 2*G4267-60206 Sample Container))
2	G7167-60020	Drawer 2H (including 2*G4267-60205 Sample Container)
3	G7167-60022	Drawer 3H (2 p/k (including 2*G4267-60205 Sample Container))
	G4267-60024	Dummy Drawer (not shown)




Figure 34 Hotel drawer

## Analytical Head Assembly 40 $\mu$ L

Item	p/n	Description
	G4267-60042	Analytical Head, 40 $\mu$ L
1	G4267-60423	Head Assembly, 40 $\mu$ L
2	0905-1717	Metering seal 40 $\mu$ L
3	G4267-60422	Seal Support Assembly, 40 $\mu$ L
4	0515-4384	Screw
5	G4267-60432	Spring Adapter Assembly
6	5067-5920	Piston ceramic




Figure 35 Analytical head assembly, 40  $\mu$ L

## Analytical Head Assembly 100 $\mu$ L

Item	p/n	Description
	G4267-60043	Analytical Head, 100 $\mu$ L
1	G4267-60433	Head Assembly, 100 $\mu$ L
2	0905-1719	Metering seal 100 $\mu$ L
3	G4267-60434	Seal Support Assembly, 100 $\mu$ L
4	0515-1052	Screw 2.5 mm hex
5	G4267-60432	Spring Adapter Assembly
6	5067-5678	Piston 1290 Infinity Pumps, ceramic




Figure 36 Analytical head assembly, 100  $\mu$ L

## Analytical Head Assembly 900 $\mu$ L

Item	p/n	Description
	G4267-60046	Analytical head, 900 $\mu$ L, 400 bar
1	G4267-60461	Head Assembly, 900 $\mu$ L
2	0905-1294	Metering seal, 900 $\mu$ L
3	G4267-60463	Seal Support Assembly, 900 $\mu$ L
4	0515-4864	Screw for heat exchangers
5	G4267-60432	Spring Adapter Assembly
6	G4267-60462	Piston Assembly, 900 $\mu$ L





Figure 37 Analytical head assembly, 900  $\mu$ L

## Flush Head Assembly 500 $\mu$ L

Item	p/n	Description
	G4267-60049	Flush head, 500 $\mu$ L
1	G4267-60491	Flush Head Assembly, 500 $\mu$ L
2	5023-2473	Sealing Plate 500 $\mu$ L
3	G4267-60482	Cylinder Assembly, 500 $\mu$ L
4	5067-5918	Seal 500 $\mu$ L
5		Screw
6	1410-1881	Bearing-Sleeve 8 mm-ID 10 mm-OD 10 mm-LG PI
7	G4267-60432	Spring Adapter Assembly
8	5067-5919	Piston Assembly 500 $\mu$ L
9	G4267-60451	Pump Valve IN
10	G4267-60452	Pump Valve Out

**10** Parts for Maintenance and Upgrade or Options  
Flush Head Assembly 500  $\mu$ L



**Figure 38** Flush head assembly, 500  $\mu$ L


## Injection Valve Assembly

Item	p/n	Description
	5067-4232	2pos/6port Injection Valve (VICI)
1	5068-0019	Stator screws
2	5068-0197	Stator head
3	5068-0198	Rotor Seal for VICI Injection Valve (PAEK)
	5500-1159	Capillary ST 0.17x100 SX/S-2.3 Metering Device to Injection Valve
	5067-4650	Capillary ST 0.12 mm x 150 mm SL/SX Pump to sampler
	5500-1157	Capillary, ST, 0.12 mmx500 mm Sampler to column compartment
	5067-6127	Blank Nut SL

### NOTE

For the VICI Valve SL/SX fittings are mandatory.

**10** Parts for Maintenance and Upgrade or Options  
Injection Valve Assembly



**Figure 39** Injection valve assembly (VICI)

## Injection Valve with Actuator

Item	p/n	Description
1	5067-4232	2pos/6port Injection Valve (VICI)
2	5043-0291	Lock Nut
3	5188-8030	Tag Reader
4	5067-4162	Direct-Actuator-50 Assembly





Figure 40 Injection valve with actuator

## 10 Parts for Maintenance and Upgrade or Options

### Needle Port Assembly

# Needle Port Assembly

Item	p/n	Description
1	G4267-60044	Needle Port Assembly Station
2	G4267-40045	Needle port Adapter



**Figure 41** Needle port assembly

# Door Assy

Item	p/n	Description
	5067-5415	Door Assy
1	5021-1879	Permanent Magnet
2	G4267-68713	Pressure Spring
3	5067-5412 (2x)	Hinge Universal




Figure 42 Door assy

## Accessory Kit

Item	p/n	Description
	G4267-68705	Accessory Kit
1	G4220-60007	Bottle Head Assembly
2	5063-6527	Tubing assembly, i.d. 6 mm, o.d. 9 mm, 1.2 m (to waste)
3	5500-1157	Capillary, ST, 0.12 mmx500 mm
4	5043-1013	Tubing Clip
5	5181-1519	CAN cable, Agilent module to module, 1 m




Figure 43 Accessory kit (standard)

## Tubing Kit Sampler Standard

Item	p/n	Description
	G4267-60061	Tubing-Kit-Sampler-Standard contains:
1	5042-9974	Tubing Flex (1.5 m)
2	5500-1155	Tube Connector, 90 degree, ID 6.4
3	0890-1760	Tubing Flexible 1 ea / 1 meter
4	5042-6422	Tubing connector, 1 mm o.d.
5	0100-1708	Nut 1/8 PPS
6	0100-1700	FERRULE-AY-18IN
7	0100-1846	UNION-TEFZEL





Figure 44 Tubing kit sampler standard

## 10 Parts for Maintenance and Upgrade or Options


### Sample Cooler

# Sample Cooler

Item	p/n	Description
1	G4267-60005	Sample Cooler
	G4267-81015	Cable Power Sample Cooler not shown
	G4267-81014	Cable-Ribbon Sample Cooler not shown
	2110-1519	Fuse 3.50 A125 V not shown



**Figure 45** Sample cooler



## 11 Identifying Cables

Cable Overview	220
Analog Cables	222
Remote Cables	224
CAN/LAN Cables	228
Agilent Module to PC	229
USB	230

This chapter provides information on cables used with the modules.



# Cable Overview

**NOTE**

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

---

**Analog cables**

p/n	Description
35900-60750	Agilent 35900A A/D converter
01046-60105	Analog cable (BNC to general purpose, spade lugs)

**Remote cables**

p/n	Description
5188-8029	ERI to general purpose
5188-8044	Remote Cable ERI – ERI
5188-8045	Remote Cable APG – ERI
5061-3378	Remote Cable to 35900 A/D converter
01046-60201	Agilent module to general purpose

**CAN cables**

p/n	Description
5181-1516	CAN cable, Agilent module to module, 0.5 m
5181-1519	CAN cable, Agilent module to module, 1 m

**LAN cables**

p/n	Description
5023-0203	Cross-over network cable, shielded, 3 m (for point to point connection)
5023-0202	Twisted pair network cable, shielded, 7 m (for point to point connection)

**RS-232 cables  
(not for FUSION  
board)**

p/n	Description
G1530-60600	RS-232 cable, 2 m
RS232-61601	RS-232 cable, 2.5 m Instrument to PC, 9-to-9 pin (female). This cable has special pin-out, and is not compatible with connecting printers and plotters. It's also called "Null Modem Cable" with full handshaking where the wiring is made between pins 1-1, 2-3, 3-2, 4-6, 5-5, 6-4, 7-8, 8-7, 9-9.
5181-1561	RS-232 cable, 8 m

**USB cables**

p/n	Description
5188-8050	USB A M-USB Mini B 3 m (PC-Module)
5188-8049	USB A F-USB Mini B M OTG (Module to Flash Drive)

## Analog Cables



One end of these cables provides a BNC connector to be connected to Agilent modules. The other end depends on the instrument to which connection is being made.


### Agilent Module to 35900 A/D converters

p/n 35900-60750	35900	Pin Agilent module	Signal Name
	1		Not connected
	2	Shield	Analog -
	3	Center	Analog +

### Agilent Module to BNC Connector

p/n 8120-1840	Pin BNC	Pin Agilent module	Signal Name
	Shield	Shield	Analog -
	Center	Center	Analog +

### Agilent Module to General Purpose

p/n 01046-60105	Pin	Pin Agilent module	Signal Name
	1		Not connected
	2	Black	Analog -
	3	Red	Analog +

## Remote Cables

### ERI (Enhanced Remote Interface)

5188-8029 ERI to general purpose

p/n 5188-8029	pin	Enhanced Remote	Classic Remote
<p>D-Sub female 15way user's view to connector</p> <p>I08 I07 I06 I05 I04 I03 I02 I01</p> <p>8 15 9</p> <p>+24V +24V PGND PGND +5V DGND 1Wire DATA</p>	1	I01	START REQUEST
	2	I02	STOP
	3	I03	READY
	4	I04	POWER ON
	5	I05	NOT USED
	6	I06	SHUT DOWN
	7	I07	START
	8	I08	PREPARE
	9	1wire DATA	
	10	DGND	
	11	+5V ERI out	
	12	PGND	
	13	PGND	
	14	+24V ERI out	
	15	+24V ERI out	

5188-8044 ERI to ERI (Connector D\_Subminiature 15 pin)

**Table 14** 5188-8044 ERI to ERI

p/n 5188-8044	Pin (ERI)	Signal	Pin (ERI)
	10	GND	10
	1	Start Request	1
	2	Stop	2
	3	Ready	3
	5	Power on	5
	4	Future	4
	6	Shut Down	6
	7	Start	7
	8	Prepare	8
	Ground Connection	Cable Shielding	NC

5188-8045 ERI to APG (Connector D\_Subminiature 15 pin (ERI), Connector D\_Subminiature 9 pin (APG))

p/n 5188-8045	Pin (ERI)	Signal	Pin (APG)
	10	GND	1
	1	Start Request	9
	2	Stop	8
	3	Ready	7
	5	Power on	6
	4	Future	5
	6	Shut Down	4
	7	Start	3
	8	Prepare	2
	Ground Connection	Cable Shielding	NC


## 11 Identifying Cables

### Remote Cables



One end of these cables provides a Agilent Technologies APG (Analytical Products Group) remote connector to be connected to Agilent modules. The other end depends on the instrument to be connected to.


### Agilent Module to Agilent 35900 A/D Converters

p/n 5061-3378	Pin 35900 A/D	Pin Agilent module	Signal Name	Active (TTL)
	1 - White	1 - White	Digital ground	
	2 - Brown	2 - Brown	Prepare run	Low
	3 - Gray	3 - Gray	Start	Low
	4 - Blue	4 - Blue	Shut down	Low
	5 - Pink	5 - Pink	Not connected	
	6 - Yellow	6 - Yellow	Power on	High
	7 - Red	7 - Red	Ready	High
	8 - Green	8 - Green	Stop	Low
	9 - Black	9 - Black	Start request	Low

## Agilent Module to General Purpose

p/n 01046-60201	Wire Color	Pin Agilent module	Signal Name	Active (TTL)
	White	1	Digital ground	
	Brown	2	Prepare run	Low
	Gray	3	Start	Low
	Blue	4	Shut down	Low
	Pink	5	Not connected	
	Yellow	6	Power on	High
	Red	7	Ready	High
	Green	8	Stop	Low
	Black	9	Start request	Low

## CAN/LAN Cables



Both ends of this cable provide a modular plug to be connected to Agilent modules CAN or LAN connectors.

### CAN Cables

p/n	Description
5181-1516	CAN cable, Agilent module to module, 0.5 m
5181-1519	CAN cable, Agilent module to module, 1 m

### LAN Cables

p/n	Description
5023-0203	Cross-over network cable, shielded, 3 m (for point to point connection)
5023-0202	Twisted pair network cable, shielded, 7 m (for point to point connection)

## Agilent Module to PC

p/n	Description
G1530-60600	RS-232 cable, 2 m
RS232-61601	RS-232 cable, 2.5 m Instrument to PC, 9-to-9 pin (female). This cable has special pin-out, and is not compatible with connecting printers and plotters. It's also called "Null Modem Cable" with full handshaking where the wiring is made between pins 1-1, 2-3, 3-2, 4-6, 5-5, 6-4, 7-8, 8-7, 9-9.
5181-1561	RS-232 cable, 8 m


## 11 Identifying Cables

### USB

# USB

To connect a USB Flash Drive use a USB OTG cable with Mini-B plug and A socket.

<b>p/n</b>	<b>Description</b>
5188-8050	USB A M-USB Mini B 3 m (PC-Module)
5188-8049	USB A F-USB Mini B M OTG (Module to Flash Drive)



## 12 Hardware Information

Firmware Description	232
Electrical Connections	235
Rear view of the module	236
Information on Instrument Serial Number	236
Interfaces	238
Overview Interfaces	239
ERI (Enhanced Remote Interface)	243
Setting the 6-bit Configuration Switch	246
Instrument Layout	248
Early Maintenance Feedback	249

This chapter describes the module in more detail on hardware and electronics.



## Firmware Description

The firmware of the instrument consists of two independent sections:

- a non-instrument specific section, called *resident system*
- an instrument specific section, called *main system*

### Resident System

This resident section of the firmware is identical for all Agilent 1100/1200/1220/1260/1290 series modules. Its properties are:

- the complete communication capabilities (CAN, LAN and RS-232C)
- memory management
- ability to update the firmware of the 'main system'

### Main System

Its properties are:

- the complete communication capabilities (CAN, LAN and RS-232C)
- memory management
- ability to update the firmware of the 'resident system'

In addition the main system comprises the instrument functions that are divided into common functions like

- run synchronization through APG remote,
- error handling,
- diagnostic functions,
- or module specific functions like
  - internal events such as lamp control, filter movements,
  - raw data collection and conversion to absorbance.

### Firmware Updates

Firmware updates can be done using the following tools (latest version should be used):

- Agilent Lab Advisor software with files on the hard disk (\*)
- Firmware Update Tool with local files on the hard disk (\*)
- Instant Pilot (G4208A) with files on a USB Flash Disk

(\*) Required tools, firmware and documentation are available from the Agilent web:

[http://www.chem.agilent.com/\\_layouts/agilent/downloadFirmware.aspx?whid=69761](http://www.chem.agilent.com/_layouts/agilent/downloadFirmware.aspx?whid=69761)

The file naming conventions are:

PPPP\_RVVV\_XXX.dlb, where

PPPP is the product number, for example, 1315B for the G1315B DAD,

R the firmware revision, for example, A for G1315B or B for the G1315C DAD,

VVV is the revision number, for example 650 is revision 6.50,

XXX is the build number of the firmware.

For instructions on firmware updates refer to section *Replacing Firmware* in chapter *"Maintenance"* or use the documentation provided with the *Firmware Update Tools*.

**NOTE**

Update of main system can be done in the resident system only. Update of the resident system can be done in the main system only.

Main and resident firmware must be from the same set.




Figure 46 Firmware Update Mechanism

**NOTE**

Some modules are limited in downgrading due to their main board version or their initial firmware revision. For example, a G1315C DAD SL cannot be downgraded below firmware revision B.01.02 or to a A.xx.xx.

Some modules can be re-branded (e.g. G1314C to G1314B) to allow operation in specific control software environments. In this case the feature set of the target type are use and the feature set of the original are lost. After re-branding (e.g. from G1314B to G1314C), the original feature set is available again.

All these specific informations are described in the documentation provided with the firmware update tools.

The firmware update tools, firmware and documentation are available from the Agilent web.

- [http://www.chem.agilent.com/\\_layouts/agilent/downloadFirmware.aspx?whid=69761](http://www.chem.agilent.com/_layouts/agilent/downloadFirmware.aspx?whid=69761)

## Electrical Connections

- The CAN bus is a serial bus with high-speed data transfer. The two connectors for the CAN bus are used for internal module data transfer and synchronization.
- One analog output provides signals for integrators or data handling systems.
- The ERI/REMOTE connector may be used in combination with other analytical instruments from Agilent Technologies if you want to use features such as start, stop, common shutdown, prepare, and so on.
- With the appropriate software, the LAN connector may be used to control the module from a computer through a LAN connection. This connector is activated and can be configured with the configuration switch.
- With the appropriate software, the USB connector may be used to control the module from a computer through a USB connection.
- The power input socket accepts a line voltage of 100 – 240 VAC  $\pm$  10 % with a line frequency of 50 or 60 Hz. Maximum power consumption varies by module. There is no voltage selector on your module because the power supply has wide-ranging capability. There are no externally accessible fuses because automatic electronic fuses are implemented in the power supply.

**NOTE**

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

---

## Rear view of the module




Figure 47 Rear view of multisampler - electrical connections and label

## Information on Instrument Serial Number

### Serial Number Information 1200 Series and 1290 Infinity

The serial number information on the instrument labels provide the following information:

CCYWWSSSSS

Format

CC

country of manufacturing

- DE = Germany
- JP = Japan
- CN = China

YWW	year and week of last major manufacturing change, e.g. 820 could be week 20 of 1998 or 2008
SSSS	real serial number

## Serial Number Information 1260 Infinity

The serial number information on the instrument labels provide the following information:

CCXZZ00000	Format
CC	Country of manufacturing <ul style="list-style-type: none"> <li>• DE = Germany</li> <li>• JP = Japan</li> <li>• CN = China</li> </ul>
X	Alphabetic character A-Z (used by manufacturing)
ZZ	Alpha-numeric code 0-9, A-Z, where each combination unambiguously denotes a module (there can be more than one code for the same module)
00000	Serial number

## Interfaces

The Agilent 1200 Infinity Series II modules provide the following interfaces:

**Table 15** Agilent 1200 Infinity II Series Interfaces

Module	CAN	USB	LAN (on-board)	RS-232	Analog	APG (A) / ERI (E)	Special
<b>Pumps</b>							
G7104A Flexible Pump	2	No	Yes	Yes	1	A	
G7120A High Speed Pump	2	No	Yes	Yes	1	A	
<b>Samplers</b>							
G7129A/B ALS	2	Yes	Yes	No	No	E	
G7167A/B Multisampler	2	Yes	Yes	No	No	E	
<b>Detectors</b>							
G7114A/B VWD	2	Yes	Yes	No	1	E	
G7117A/B DAD	2	Yes	Yes	No	1	E	
G7115A/B DAD	2	Yes	Yes	No	1	E	
<b>Others</b>							
G7116B MCT	2	No	No	No	No	No	Requires a HOST module via CAN

**NOTE**

The detector (DAD/MWD/FLD/VWD/RID) is the preferred access point for control via LAN. The inter-module communication is done via CAN.

- CAN connectors as interface to other modules
- LAN connector as interface to the control software
- RS-232C as interface to a computer

- USB (Universal Series Bus) as interface to a computer
- REMOTE connector as interface to other Agilent products
- Analog output connector(s) for signal output

## Overview Interfaces

### CAN

The CAN is inter-module communication interface. It is a 2-wire serial bus system supporting high speed data communication and real-time requirement.

### LAN

The modules have either an interface slot for an LAN card (e.g. Agilent G1369B/C LAN Interface) or they have an on-board LAN interface (e.g. detectors G1315C/D DAD and G1365C/D MWD). This interface allows the control of the module/system via a PC with the appropriate control software. Some modules have neither on-board LAN nor an interface slot for a LAN card (e.g. G1170A Valve Drive or G4227A Flex Cube). These are hosted modules and require a Host module with firmware B.06.40 or later or with additional G1369C LAN Card.

#### NOTE

If an Agilent detector (DAD/MWD/FLD/VWD/RID) is in the system, the LAN should be connected to the DAD/MWD/FLD/VWD/RID (due to higher data load). If no Agilent detector is part of the system, the LAN interface should be installed in the pump or autosampler.

### RS-232C (Serial)

The RS-232C connector is used to control the module from a computer through RS-232C connection, using the appropriate software. This connector can be configured with the configuration switch module at the rear of the module. Refer to *Communication Settings for RS-232C*.

## 12 Hardware Information

### Interfaces

#### NOTE

There is no configuration possible on main boards with on-board LAN. These are pre-configured for

- 19200 baud,
- 8 data bit with no parity and
- one start bit and one stop bit are always used (not selectable).

The RS-232C is designed as DCE (data communication equipment) with a 9-pin male SUB-D type connector. The pins are defined as:

**Table 16** RS-232C Connection Table

Pin	Direction	Function
1	In	DCD
2	In	RxD
3	Out	TxD
4	Out	DTR
5		Ground
6	In	DSR
7	Out	RTS
8	In	CTS
9	In	RI




Figure 48 RS-232 Cable

## Analog Signal Output

The analog signal output can be distributed to a recording device. For details refer to the description of the module's main board.

## APG Remote

The APG Remote connector may be used in combination with other analytical instruments from Agilent Technologies if you want to use features as common shut down, prepare, and so on.

Remote control allows easy connection between single instruments or systems to ensure coordinated analysis with simple coupling requirements.

The subminiature D connector is used. The module provides one remote connector which is inputs/outputs (wired- or technique).

To provide maximum safety within a distributed analysis system, one line is dedicated to **SHUT DOWN** the system's critical parts in case any module detects a serious problem. To detect whether all participating modules are switched on or properly powered, one line is defined to summarize the **POWER ON** state of all connected modules. Control of analysis is maintained by signal readiness **READY** for next analysis, followed by **START** of run and optional **STOP** of run triggered on the respective lines. In addition **PREPARE** and **START REQUEST** may be issued. The signal levels are defined as:

- standard TTL levels (0 V is logic true, + 5.0 V is false),

## 12 Hardware Information

### Interfaces

- fan-out is 10,
- input load is 2.2 kOhm against + 5.0 V, and
- output are open collector type, inputs/outputs (wired- or technique).

#### NOTE

All common TTL circuits operate with a 5 V power supply. A TTL signal is defined as "low" or L when between 0 V and 0.8 V and "high" or H when between 2.0 V and 5.0 V (with respect to the ground terminal).

**Table 17** Remote Signal Distribution

Pin	Signal	Description
1	DGND	Digital ground
2	PREPARE	(L) Request to prepare for analysis (for example, calibration, detector lamp on). Receiver is any module performing pre-analysis activities.
3	START	(L) Request to start run / timetable. Receiver is any module performing run-time controlled activities.
4	SHUT DOWN	(L) System has serious problem (for example, leak: stops pump). Receiver is any module capable to reduce safety risk.
5		Not used
6	POWER ON	(H) All modules connected to system are switched on. Receiver is any module relying on operation of others.
7	READY	(H) System is ready for next analysis. Receiver is any sequence controller.
8	STOP	(L) Request to reach system ready state as soon as possible (for example, stop run, abort or finish and stop injection). Receiver is any module performing run-time controlled activities.
9	START REQUEST	(L) Request to start injection cycle (for example, by start key on any module). Receiver is the autosampler.

### Special Interfaces


There is no special interface for this module.

## ERI (Enhanced Remote Interface)

ERI replaces the AGP Remote Interface that is used in the HP 1090/1040/1050/1100 HPLC systems and Agilent 1100/1200/1200 Infinity HPLC modules. All new 1200 Infinity II products using the FUSION core electronics use ERI. This interface is already used in the Agilent Universal Interface Box 2 (UIB2)

### ERI Description

The ERI interface contains eight individual programmable input/output pins. In addition, it provides 24 V power and 5 V power and a serial data line to detect and recognize further add-ons that could be connected to this interface. This way the interface can support various additional devices like sensors, triggers (in and out) and small controllers, etc.



**Figure 49** Location of the ERI interface (example shows a G7114A/B VWD)

	Pin	Enhanced Remote
<p>D-Sub female 15way user's view to connector</p> <p>IO1 IO2 IO3 IO4 IO5 IO6 IO7 IO8 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 1WEprom DGND +5V PGND +24V +24V</p>	1	IO 1 (START REQUEST)
	2	IO 2 (STOP)
	3	IO 3 (READY)
	4	IO 4 (POWER ON)
	5	IO 5 (NOT USED)
	6	IO 6 (SHUT DOWN)
	7	IO 7 (START)
	8	IO 8 (PREPARE)
	9	1 wire DATA
	10	DGND
	11	+5 V ERI out
	12	PGND
	13	PGND
	14	+24 V ERI out
	15	+24 V ERI out

## IO (Input/Output) Lines

- Eight generic bi-directional channels (input or output).
- Same as the APG Remote.
- Devices like valves, relays, ADCs, DACs, controllers can be supported/controlled.

## 1-Wire Data (Future Use)

This serial line can be used to read out an EPROM or write into an EPROM of a connected ERI-device. The firmware can detect the connected type of device automatically and update information in the device (if required).

### 5V Distribution (Future Use)

- Available directly after turn on of the hosting module (assures that certain base functionality of the device can be detected by firmware).
- For digital circuits or similar.
- Provided 500 mA maximum.
- Short-circuit proof with automatic switch off (by firmware).

### 24V Distribution (Future Use)

- Available by firmware command (defined turn on/off).
- For devices that need higher power
  - Class 0: 0.5 A maximum (12 W)
  - Class 1: 1.0 A maximum (24 W)
  - Class 2: 2.0 A maximum (48 W)
- Class depends on hosting module's internal power overhead.
- If a connected device requires more power the firmware detects this (overcurrent detection) and provides the information to the user interface.
- Fuse used for safety protection (on board).
- Short circuit will be detected through hardware.

## Setting the 6-bit Configuration Switch

The 6-bit configuration switch is located at the rear of the module with FUSION electronics. Switch settings provide configuration parameters for LAN and instrument specific initialization procedures.

All modules with FUSION electronics:

- Default is ALL switches DOWN (best settings).
  - Default IP address for LAN 192.168.254.11
- For specific LAN modes switches 4-5 must be set as required.
- For boot resident/cold start modes switches 1+2 or 6 must be UP.




Figure 50 Location of Configuration switch (example shows a G7114A/B VWD)

Table 18 6-bit Configuration Switch

	Mode	Function/Setting				
	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6
COM <sup>1</sup>	0	n.a. <sup>2</sup>	n.a.	LAN Init Mode		n.a.
Use Default IP Address <sup>3</sup>		0	0	0	0	0
Use Stored IP Address		0	0	0	1	0
Use DHCP to request IP Address <sup>4</sup>		0	0	1	0	0

**Table 18** 6-bit Configuration Switch

Test	Mode	Function/Setting				
	1	System	n.a.	n.a.	n.a.	ColdStart
Boot Main System / Keep Data		0	0	0	0	0
Boot Resident System / Keep Data		1	0	0	0	0
Boot Main System / Revert to Default Data		0	0	0	0	1
Boot Resident System / Revert to Default Data		1	0	0	0	1

- <sup>1</sup> When selecting mode COM, settings are stored to non-volatile memory. When selecting mode TEST, COM settings are taken from non-volatile memory.
- <sup>2</sup> not assigned - Always keep these switches on position '0' (off)
- <sup>3</sup> Default IP Address is 192.168.254.11
- <sup>4</sup> Host Name will be the MAC address.

## Instrument Layout

The industrial design of the module incorporates several innovative features. It uses Agilent's E-PAC concept for the packaging of electronics and mechanical assemblies. This concept is based upon the use of expanded polypropylene (EPP) layers of foam plastic spacers in which the mechanical and electronic boards components of the module are placed. This pack is then housed in a metal inner cabinet which is enclosed by a plastic external cabinet. The advantages of this packaging technology are:

- virtual elimination of fixing screws, bolts or ties, reducing the number of components and increasing the speed of assembly/disassembly,
- the plastic layers have air channels molded into them so that cooling air can be guided exactly to the required locations,
- the plastic layers help cushion the electronic and mechanical parts from physical shock, and
- the metal inner cabinet shields the internal electronics from electromagnetic interference and also helps to reduce or eliminate radio frequency emissions from the instrument itself.

## Early Maintenance Feedback

Maintenance requires the exchange of components which are subject to wear or stress. Ideally, the frequency at which components are exchanged should be based on the intensity of usage of the module and the analytical conditions, and not on a predefined time interval. The early maintenance feedback (EMF) feature monitors the usage of specific components in the instrument, and provides feedback when the user-selectable limits have been exceeded. The visual feedback in the user interface provides an indication that maintenance procedures should be scheduled.

### EMF Counters

**EMF counters** increment with use and can be assigned a maximum limit which provides visual feedback in the user interface when the limit is exceeded. Some counters can be reset to zero after the required maintenance procedure.

### Using the EMF Counters


The user-settable **EMF** limits for the **EMF Counters** enable the early maintenance feedback to be adapted to specific user requirements. The useful maintenance cycle is dependent on the requirements for use. Therefore, the definition of the maximum limits need to be determined based on the specific operating conditions of the instrument.

### Setting the EMF Limits

The setting of the **EMF** limits must be optimized over one or two maintenance cycles. Initially the default **EMF** limits should be set. When instrument performance indicates maintenance is necessary, take note of the values displayed by the **EMF counters**. Enter these values (or values slightly less than the displayed values) as **EMF** limits, and then reset the **EMF counters** to zero. The next time the **EMF counters** exceed the new **EMF** limits, the **EMF** flag will be displayed, providing a reminder that maintenance needs to be scheduled.

## 12 Hardware Information

### Early Maintenance Feedback



## 13 LAN Configuration

Setting up the module in a LAN environment [252](#)

Connecting the module via LAN [253](#)

This chapter provides information on connecting the detector to the Agilent ChemStation PC.



## 13 LAN Configuration

Setting up the module in a LAN environment

### Setting up the module in a LAN environment


It is not recommended to connect a system via an autosampler. The detector is producing the most data in the stack, followed by the pump, and it is therefore highly recommended to use either of these modules for the LAN connection.

## Connecting the module via LAN

If the module is being operated as a standalone module or if a connection via LAN is required regardless of above mentioned recommendation, a G1369B/C LAN card has to be used. For installation and configuration, see the G1369B/C documentation.

## 13 LAN Configuration

Connecting the module via LAN



## 14 Appendix

General Safety Information	256
General Safety Information	256
Safety Standards	256
General	256
Before Applying Power	257
Ground the Instrument	257
Do Not Operate in an Explosive Atmosphere	258
Do Not Remove the Instrument Cover	258
Do Not Modify the Instrument	258
In Case of Damage	258
Solvents	259
Symbols	260
Waste Electrical and Electronic Equipment Directive	262
Refrigerant	263
Radio Interference	265
Sound Emission	266
Solvent Information	267
Agilent Technologies on Internet	268

This chapter provides addition information on safety, legal and web.



## General Safety Information

### General Safety Information

The following general safety precautions must be observed during all phases of operation, service, and repair of this instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended use of the instrument. Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

#### **WARNING**

Ensure the proper usage of the equipment.

The protection provided by the equipment may be impaired.

→ The operator of this instrument is advised to use the equipment in a manner as specified in this manual.

---

### Safety Standards

This is a Safety Class I instrument (provided with terminal for protective earthing) and has been manufactured and tested according to international safety standards.

### General

Do not use this product in any manner not specified by the manufacturer. The protective features of this product may be impaired if it is used in a manner not specified in the operation instructions.

## Before Applying Power

### WARNING

Wrong voltage range, frequency or cabling

Personal injury or damage to the instrument

- Verify that the voltage range and frequency of your power distribution matches to the power specification of the individual instrument.
  - Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.
  - Make all connections to the unit before applying power.
- 

### NOTE

Note the instrument's external markings described under “Symbols” on page 260.

---

## Ground the Instrument

### WARNING

Missing electrical ground

Electrical shock

- If your product is provided with a grounding type power plug, the instrument chassis and cover must be connected to an electrical ground to minimize shock hazard.
  - The ground pin must be firmly connected to an electrical ground (safety ground) terminal at the power outlet. Any interruption of the protective (grounding) conductor or disconnection of the protective earth terminal will cause a potential shock hazard that could result in personal injury.
-

## Do Not Operate in an Explosive Atmosphere

### WARNING

Presence of flammable gases or fumes

Explosion hazard

→ Do not operate the instrument in the presence of flammable gases or fumes.

---

## Do Not Remove the Instrument Cover

### WARNING

Instrument covers removed

Electrical shock

→ Do Not Remove the Instrument Cover

→ Only Agilent authorized personnel are allowed to remove instrument covers. Always disconnect the power cables and any external circuits before removing the instrument cover.

---

## Do Not Modify the Instrument

Do not install substitute parts or perform any unauthorized modification to the product. Return the product to an Agilent Sales and Service Office for service and repair to ensure that safety features are maintained.

## In Case of Damage

### WARNING

Damage to the module

Personal injury (for example electrical shock, intoxication)

→ Instruments that appear damaged or defective should be made inoperative and secured against unintended operation until they can be repaired by qualified service personnel.

---

## Solvents

### WARNING

**Toxic, flammable and hazardous solvents, samples and reagents**

**The handling of solvents, samples and reagents can hold health and safety risks.**

- When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
  - The volume of substances should be reduced to the minimum required for the analysis.
  - Do not operate the instrument in an explosive atmosphere.
  - Never exceed the maximal permissible volume of solvents (6 L) in the solvent cabinet.
  - Do not use bottles that exceed the maximum permissible volume as specified in the usage guideline for the Agilent 1200 Infinity Series Solvent Cabinets.
  - Arrange the bottles as specified in the usage guideline for the solvent cabinet.
  - A printed copy of the guideline has been shipped with the solvent cabinet, electronic copies are available on the Internet.
  - Ground the waste container.
  - The residual free volume in the appropriate waste container must be large enough to collect the waste liquid.
  - Check the filling level of the waste container regularly.
  - To achieve maximal safety, check the correct installation regularly.
  - Do not use solvents with an auto-ignition temperature below 200 °C (392 °F).
-

## Symbols

Table 19 Symbols

	The apparatus is marked with this symbol when the user should refer to the instruction manual in order to protect risk of harm to the operator and to protect the apparatus against damage.
	Indicates dangerous voltages.
	Indicates a protected ground terminal.
	The apparatus is marked with this symbol when hot surfaces are available and the user should not touch it when heated up.
	Cooling unit is designed as vapor-compression refrigeration system. Contains fluorinated greenhouse gas (refrigerant) according to the Kyoto protocol. For specifications of refrigerant, charge capacity, carbon dioxide equivalent (CDE), and global warming potential (GWP) see instrument label.
	Confirms that a manufactured product complies with all applicable European Community directives. The European Declaration of Conformity is available at: <a href="http://regulations.corporate.agilent.com/DoC/search.htm">http://regulations.corporate.agilent.com/DoC/search.htm</a>
	Manufacturing date.
	Power symbol indicates On/Off. The apparatus is not completely disconnected from the mains supply when the power switch is in the Off position

Table 19 Symbols

	<p>Pacemaker Magnets could affect the functioning of pacemakers and implanted heart defibrillators. A pacemaker could switch into test mode and cause illness. A heart defibrillator may stop working. If you wear these devices keep at least 55 mm distance to magnets. Warn others who wear these devices from getting too close to magnets.</p>
	<p>Magnetic field Magnets produce a far-reaching, strong magnetic field. They could damage TVs and laptops, computer hard drives, credit and ATM cards, data storage media, mechanical watches, hearing aids and speakers. Keep magnets at least 25 mm away from devices and objects that could be damaged by strong magnetic fields.</p>
	<p>Indicates a pinching or crushing hazard</p>
	<p>Indicates a piercing or cutting hazard.</p>

**WARNING**

**A WARNING**

alerts you to situations that could cause physical injury or death.

- Do not proceed beyond a warning until you have fully understood and met the indicated conditions.

**CAUTION**

**A CAUTION**

alerts you to situations that could cause loss of data, or damage of equipment.

- Do not proceed beyond a caution until you have fully understood and met the indicated conditions.

## Waste Electrical and Electronic Equipment Directive

### Abstract

The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC), adopted by EU Commission on 13 February 2003, is introducing producer responsibility on all electric and electronic appliances starting with 13 August 2005.

#### NOTE

This product complies with the WEEE Directive (2002/96/EC) marking requirements. The affixed label indicates that you must not discard this electrical/electronic product in domestic household waste.

Product Category:

With reference to the equipment types in the WEEE Directive Annex I, this product is classed as a Monitoring and Control Instrumentation product.

---



#### NOTE

Do not dispose of in domestic household waste

To return unwanted products, contact your local Agilent office, or see [www.agilent.com](http://www.agilent.com) for more information.

---

## Refrigerant

The refrigerant HFKW-134a is used only in the Agilent Infinity II Sample Cooler.

**Table 20** Physical properties of refrigerant HFKW-134a

Molecular weight	102
Critical temperature	101.1 °C
Critical pressure	40.6 bar
Boiling point	-26.5 °C

**WARNING**

Refrigerant



Refrigerant HFKW-134a is known as a safe refrigerant, however accidents can occur if it is handled incorrectly. For this reason, the following instructions must be observed:

- Avoid contact with liquid refrigerant HFKW-134a. At atmospheric pressure HFKW-134a evaporates at approximately -26 °C and causes frost bite.
  - After skin contact, rinse the affected area with water.
  - After eye contact, rinse the eye(s) with plenty of water for at least 15 minutes and consult a doctor.
  - HFKW-134a must not be allowed to escape in enclosed areas. Although HFKW-134a is not toxic, there is a danger of suffocation as gaseous refrigerant is heavier than air.
  - Please observe the following first aid instructions. After inhalation, move the affected person to fresh air, keep him warm and allow him to rest. If necessary, he should be supplied with oxygen. If he has stopped breathing or is breathing erratically, he should be given artificial respiration. In the case of cardiac arrest, carry out heart massage. Send for a doctor immediately.
  - Moreover, it must be noted that HFKW-134a must always be extracted from the system and collected. It must never be discharged into the atmosphere on environmental grounds (greenhouse effect).
- 

**CAUTION**

General hazards and improper disposal

Improper disposal of the media and components used pollutes the environment.

- The breakdown of the sample cooler unit must be carried out by specialist refrigeration company.
  - All media must be disposed of in accordance with national and local regulations.
  - Please contact your local Agilent Service Center in regard to safe environmental disposal of the appliance.
-

## Radio Interference

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

### Test and Measurement

If test and measurement equipment is operated with equipment unshielded cables and/or used for measurements on open set-ups, the user has to assure that under operating conditions the radio interference limits are still met within the premises.

## Sound Emission

### Manufacturer's Declaration

This statement is provided to comply with the requirements of the German Sound Emission Directive of 18 January 1991.


This product has a sound pressure emission (at the operator position) < 70 dB.

- Sound Pressure  $L_p < 70$  dB (A)
- At Operator Position
- Normal Operation
- According to ISO 7779:1988/EN 27779/1991 (Type Test)

## Solvent Information

Observe the following recommendations on the use of solvents.

- Brown glass ware can avoid growth of algae.
- Avoid the use of the following steel-corrosive solvents:
  - Solutions of alkali halides and their respective acids (for example, lithium iodide, potassium chloride, and so on),
  - High concentrations of inorganic acids like sulfuric acid and nitric acid, especially at higher temperatures (if your chromatography method allows, replace by phosphoric acid or phosphate buffer which are less corrosive against stainless steel),
  - Halogenated solvents or mixtures which form radicals and/or acids, for example:



This reaction, in which stainless steel probably acts as a catalyst, occurs quickly with dried chloroform if the drying process removes the stabilizing alcohol,

- Chromatographic grade ethers, which can contain peroxides (for example, THF, dioxane, di-isopropyl ether) such ethers should be filtered through dry aluminium oxide which adsorbs the peroxides,
- Solvents containing strong complexing agents (e.g. EDTA),
- Mixtures of carbon tetrachloride with 2-propanol or THF.

## Agilent Technologies on Internet

For the latest information on products and services visit our worldwide web site on the Internet at:

<http://www.agilent.com>

# Index

## A

- accessory kit
  - standard 216
- Agilent Lab Advisor software 107
- Agilent Lab Advisor 107
- Agilent
  - on internet 268
- algae 267
- ambient non-operating temperature 34
- ambient operating temperature 34
- analog signal 241
- analytical head assembly
  - 100 207
  - 40 206
- apg remote 241
- assemblies
  - main 203
- assembly
  - door 215
  - needle port 214

## B

- bench space 32

## C

- cable
  - CAN 228
  - LAN 228
  - overview 220
  - RS-232 229
- cables
  - analog 222
  - remote 224
- CAN

- cable 228
- cleaning 136
- coding
  - color 73
  - guide 73
- color
  - coding 73
  - guide 73
- compensation sensor open 114
- compensation sensor short 115
- condensation 32
- configuration
  - module 85
  - options 187
- control settings 81

## D

- defect on arrival 50
- delay volume
  - description 88
- dimensions 34
- door
  - assembly 215
- doors
  - installation 137
  - removal 137
- drawer upgrade
  - installation 189
- drawer
  - status indicator 47

## E

- electrical connections
  - descriptions of 235

- electronic waste 262
- EMF
  - early maintenance feedback 249
- error messages
  - compensation sensor open 114
  - compensation sensor short 115
  - fan failed 115
  - leak sensor open 114
  - leak sensor short 113
  - leak 116
  - lost CAN partner 113
  - remote timeout 112
  - shutdown 111
  - system pressure test 121
  - timeout 111
- extra-column volume 88

## F

- fan failed 115
- firmware
  - description 232
  - main system 232
  - resident system 232
  - update tool 233
  - updates 232, 200
  - upgrade/downgrade 200
- flush head 209
- frequency range 34

## G

- general error messages 111
- guide
  - coding 73
  - color 73

## Index

### H

- hotel drawer
  - configure in control software 193
  - configure in Lab Advisor 195
- hotel drawer 205
- humidity 34

### I

- indicator
  - drawer status 47
- injection valve
  - actuator 213
- installation
  - bench space 32
  - sample cooler 51
  - upgrade drawer 189
- installing
  - sample loop-flex 184
- instrument layout 248
- interfaces 106
  - Infinity II 238
- internet 268

### L

- LAN
  - cable 228
  - configuration 251
  - environment 252
  - setup 252
- leak sensor open 114
- leak sensor short 113
- leak 116
- line frequency 34
- line voltage 34
- lost CAN partner 113

### M

- maintenance
  - feedback 249
  - overview 202, 135
  - replacing firmware 200
- mats 71
- message
  - remote timeout 112
- metering seal
  - install 162
  - remove 157
- method parameter settings 82
- module firmware
  - replace 200
- module
  - configuration 85
  - overview 14
- multisampler
  - transport 58

### N

- needle assembly
  - installation 144
  - remove 140
- needle port
  - assembly 214
- needle seat
  - exchange 148
- non-operating altitude 34
- non-operating temperature 34

### O

- operating Altitude 34
- operating temperature 34
- optimization
  - achieving higher resolution 98
  - achieving higher sensitivity 101
  - column use 101
- overview

- cable 220
- module 14

### P

- packaging
  - damaged 50
- physical specifications 34
- power consideration 30
- power consumption 34
- power cords 31
- principle 16
- product description 11

### R

- remote
  - cables 224
- removing
  - sample loop-flex 180
- repairs
  - replacing firmware 200
- resolution
  - Optimization 98
- rotor seal
  - replace 151,
- RS-232C
  - cable 229

### S

- safety class I 256
- safety
  - general information 256
  - standards 34
  - symbols 260
- sample cooler 218
  - installation 51
  - install 197
  - remove 196
- sample loop-flex
  - install 184

- removal 180
- sampler
  - main assemblies 203
  - tubing kit 217
- sensitivity
  - optimization 101
- serial number
  - information 236, 237
- setup
  - LAN 252
- shutdown 111
- site requirements
  - power cords 31
- solvents 267
- sound emission 266
- special interfaces 242
- specification
  - physical 34
- System pressure test
  - error message 121

## T

- temperature sensor 116
- timeout 111
- transport 58
- troubleshooting
  - error messages 110
- tubing kit
  - sampler 217

## V

- vials 71
- voltage range 34

## W

- waste
  - electrical and electronic equipment 262
- WEEE directive 262

## In This Book

This manual contains technical reference information about the Agilent 1290 Infinity II Multisampler G7167B and the Agilent 1260 Infinity Multisampler (G7167A).

The manual describes the following:

- Introduction,
- Site requirements and specifications,
- Using the module,
- Preparing the module,
- Optimizing performance,
- Troubleshooting and diagnostics,
- Error information,
- Test functions,
- Maintenance,
- Parts,
- Hardware information,
- LAN configuration,
- Safety and related information.

© Agilent Technologies 2014

Printed in Germany  
12/2014

\* G7167-90001 \*  
\* G7167-90001 \*

G7167-90001  
Rev. B