

Cracking Products of Oleic Acid and Olive Oil

Unmodified vegetable oil can be used as fuel in adapted diesel engines, but because vegetable oil is more viscous than diesel fuel, it can cause problems in such engines over time. To create a fuel with a lower viscosity, oil can be modified. Thermal cracking, or controlled pyrolysis of oil decreases its molecular weight and is thereby effective in converting oil to a more useable fuel. Here, microscale pyrolysis is used to look at cracking products of supermarket bought olive oil, and a free fatty acid found abundantly in olive oil, oleic acid.

Each sample was heated inside a quartz tube using the platinum coil of the autosampler, to a setpoint 750°C for 15 seconds. The resulting pyrogram of oleic acid and olive oil are shown in Figures 1 and 2. Because oils and fatty acids are semivolatiles, it can be difficult to get them to pyrolyze before they evaporate. Oleic acid is particularly stable because it has only one double bond. The large peak at 30 minutes in Figure 1 is oleic acid that desorbed from the coil before it had a chance to pyrolyze.

Despite oleic acid's stability, we obtained abundant pyrolysis products; many of them were not well identified with the search engine. The few listed in Tables 1 and 2 had a library quality match of 80% or above. Oleic acid produced a few aromatics like benzene and toluene. Most products in the olive oil and oleic acid are long chain alkenes, alkanes, and alkynes.

Octane found in Oleic Acid (shown in Figure 3) may be from the cleavage of the weaker single bond next to oleic acid's double bond. Octane is found in olive oil too, and could be the directly related to its high oleic acid content.




Figure 1: Oleic Acid 750°C for 15 seconds.




Figure 2: Olive Oil 750°C for 15 seconds.




Figure 3: m/z 43 extracted to show Octane in Oleic Acid.

Table 1: Oleic Acid pyrolysis products, by retention time.

| | | | | | | | |
|------|-----------------------------|-------|---|-------|-------------------------|-------|---------------------------|
| 1.76 | 1,3-Butadiene | 5.47 | Cyclohexane, ethenyl- | 12.91 | Cyclodecene | 18.53 | Z-1,6-Tridecadiene |
| 2.08 | 1,3-Cyclopentadiene | 6.96 | 1-Nonene | 13.38 | 1-Dodecene | 19.43 | 1-Octadecyne |
| 2.31 | 1-Hexene | 7.72 | Cyclopropane, 1,1'-(1,2-ethanediyl)bis- | 14.05 | 5-Decyne | 20.17 | 1-Hexadecyne |
| 2.82 | 2,4-Hexadiene | 8.23 | Cyclopentene, 1-butyl- | 14.72 | 2,4-Dodecadiene, (E,Z)- | 20.59 | cis-7-Tetradecen-1-ol |
| 3.01 | 1,3-Cyclohexadiene | 9.18 | 1-Decene | 15.20 | Cyclodecene | 21.71 | 8-Heptadecene |
| 3.72 | Cyclopentane, ethenyl- | 10.83 | 1,3-Nonadiene, (E)- | 15.81 | 2,4-Decadienal, (E,E)- | 22.38 | 9-Octadecyne |
| 4.35 | 1,3,5-Hexatriene, 2-methyl- | 11.35 | 1-Undecene | 16.52 | 1,12-Tridecadiene | 26.12 | 9-Octadecenal, (Z)- |
| 4.68 | 1,7-Octadiene | 12.05 | 1,3-Nonadiene, (E)- | 17.10 | 1-Tetradecene | 28.44 | 9-Octadecenoic acid, (E)- |

Table 2: Olive Oil pyrolysis products, by retention time.

| | | | | | | | |
|------|--------------------------|-------|----------------------------|-------|--------------------------|-------|----------------------------------|
| 1.84 | 2-Propenal | 8.23 | Cyclopentene, 1-butyl- | 17.23 | Tetradecane | 24.81 | 9-Octadecenal, (Z)- |
| 2.02 | 1-Buten-3-yne, 2-methyl- | 9.37 | Decane | 18.13 | Cyclotetradecane | 26.87 | 9,17-Octadecadienal, (Z)- |
| 2.26 | 1-Hexene | 9.93 | 1,3-Octadiene | 18.54 | Z-1,6-Tridecadiene | 27.23 | 9,12-Octadecadienoic acid (Z,Z)- |
| 2.42 | 1,4-Hexadiene, (Z)- | 11.65 | 5-Undecene | 18.67 | Cyclododecane | 27.68 | Oleic Acid |
| 2.77 | 1,3-Hexadiene,c&t | 12.08 | 1,3-Nonadiene, (E)- | 18.80 | 2-Tetradecene, (E)- | 27.84 | Z,E-2,13-Octadecadien-1-ol |
| 2.97 | 1,3-Cyclohexadiene | 12.91 | trans-Bicyclo[5.1.0]octane | 19.83 | n-Nonylcyclohexane | 28.04 | 2-Methyl-Z,Z-3,13-octadecadienol |
| 3.18 | 1-Heptene | 13.39 | 1-Dodecene | 20.14 | Cyclododecene, (E)- | 28.14 | 11-Hexadecen-1-ol, acetate, (Z)- |
| 3.39 | 2-Heptene | 13.55 | Dodecane | 20.29 | Cyclohexadecane | 28.26 | 11-Hexadecen-1-ol, acetate, (Z)- |
| 4.82 | 1-Octene | 14.08 | 5-Decyne | 20.51 | Hexadecane | 29.26 | E,Z-2,13-Octadecadien-1-ol |
| 5.12 | 2-Octene, (Z)- | 15.21 | Cyclodecene, (Z)- | 21.74 | 8-Heptadecene | 29.68 | Oleic Acid |
| 6.95 | 1-Nonene | 15.45 | Tridecane | 22.28 | Cyclododecene, 1-methyl- | 31.13 | 9-Octadecenal, (Z)- |
| | | 16.99 | n-Decanoic acid | 23.14 | 5-Octadecene, (E)- | 32.43 | 13-Docosen-1-ol, (Z)- |

Equipment

These samples were analyzed using a CDS Model 5250 Pyroprobe, interfaced to an Agilent 6890/5975 gas chromatograph/mass spectrometer.

Model 5250 Conditions

Valve Oven: 225°C
Transfer Line: 300°C
Temperature: 750°C
Time: 15 seconds

GC Conditions

Carrier: Helium
Column: 5% phenyl, methyl silicone (30m X 0.25mm)
Detector: Agilent 5975 MSD
Scan: 25 to 550 amu
GC Program:
Initial: 40°C for 2 minutes
Ramp: 8°C/min.
Final: 300°C for 5 minutes

FOR MORE INFORMATION
CONCERNING THIS APPLICATION, WE
RECOMMEND THE
FOLLOWING READING:

D. G. Lima et al., Diesel-like fuel obtained by pyrolysis of vegetable oils., J. Anal. Appl. Pyrolysis 71 (2004) 987-996.

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.

CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto www.cdsanalytical.com.