

Creating Linear and Serial Dilution Methods on an Agilent 7696A WorkBench

User's Guide

This user guide describes procedures for creating linear and serial dilutions for a four level calibration curve complete with internal standards. It can be used as a template for creating similar calibration curves with an Agilent 7696A WorkBench.

Linear dilutions: In this example, we start with a 2,000 ppb stock and make a series of four linear dilutions such that the resulting concentrations are: 1,000 ppb, 500 ppb, 250 ppb, and 125 ppb.




Figure 1. Linear dilutions.

Serial dilutions: In this example, we start with a 2,000 ppb stock and make a series of four serial dilutions such that the resulting concentrations are: 1,000 ppb, 500 ppb, 250 ppb, and 125 ppb.

Each of the four calibration standards have the same volume of internal standard added.




Figure 2. Serial dilutions.

Resource Layout Editor

Before building the method for WorkBench, the resources must be allocated using the Resource Layout (Figure 3).




Figure 3. Resource layout.

- Empty vials – calibration standards are made in these
- Methanol - diluent
- ISTD – internal standard
- 2,000 ppb standard stock – starting stock solution

Methanol, ISTD, and 2,000 ppb standard stock can be allocated by volume (the user selects the amount of volume available in the vial[s]) or by use (the user selects how many times the resource can be accessed). Most chemical resources are allocated as “by volume”, however “by use” is helpful when the resource evaporates or degrades quickly after the vial has been pierced.



Create Prep Method - Linear Dilution

Using the Resources allocated above the Agilent 7696A WorkBench method is created stepwise (shown below and in Figure 4).

1. Add 400 μL of methanol to an empty vial, which is renamed '1,000 ppb std'.
2. Add 100 μL of the resource 'ISTD' to the newly created vial '1,000 ppb std'.
3. Add 500 μL of the resource '2,000 ppb Standard Stock' to '1,000 ppb std'.
4. Mix '1,000 ppb std' (1,000 ppb std is now complete).
5. Add 650 μL of methanol to an 'empty vial' which is renamed '500 ppb std'.
6. Add 100 μL of the resource 'ISTD' to the newly created vial '500 ppb std'.
7. Add 250 μL of the resource '2,000 ppb Standard Stock' to '500 ppb std'.
8. Mix '500 ppb std' (500 ppb std is now complete).
9. Add 775 μL of methanol to an empty vial, which is renamed '250 ppb std'.
10. Add 100 μL of the resource 'ISTD' to the newly created vial '250 ppb std'.
11. Add 125 μL of the resource '2,000 ppb Standard Stock' to '250 ppb std'.
12. Mix '250 ppb std' (250 ppb std is now complete).
13. Add 830 μL of methanol to an empty vial, which is renamed '125 ppb std'.
14. Add 7.5 μL of methanol to the newly created vial '125 ppb std'.
15. Add 100 μL of the resource 'ISTD' to '125 ppb std'.
16. Add 62.5 μL of the resource '2,000 ppb Standard Stock' to '125 ppb std'.
17. Mix '125 ppb std' (125 ppb std is now complete).
18. Flag the vials created: 1,000 ppb std, 500 ppb std, 250 ppb std, and 125 ppb std as 'result vials'.

Setup Method
Agilent 7696A Sample Prep Method | Agilent 7696A Configuration

Import | Export

Process in Batch Mode

Version 3.1.36.0

Actions: Add, Mix, Heat, Wait, Flag as result, Move, Wash

Program:

1. Add
2. Add
3. Add
4. Mix
5. Add
6. Add
7. Add
8. Mix
9. Add
10. Add
11. Add
12. Mix
13. Add
14. Add
15. Add
16. Add
17. Mix
18. Flag as result

Steps:

1. Add 400 μL of Methanol (CH₃OH) to 1000ppb std at Back Tower
2. Add 100 μL of ISTD to 1000ppb std at Back Tower
3. Add 500 μL of 2000ppb Standard Stock to 1000ppb std at Back Tower
4. Mix 1000ppb std at 1000 RPM for 0 min 10 sec
5. Add 650 μL of Methanol (CH₃OH) to 500ppb std at Back Tower
6. Add 100 μL of ISTD to 500ppb std at Back Tower
7. Add 250 μL of 2000ppb Standard Stock to 500ppb std at Back Tower
8. Mix 500ppb std at 1000 RPM for 0 min 10 sec
9. Add 775 μL of Methanol (CH₃OH) to 250ppb std at Back Tower
10. Add 100 μL of ISTD to 250ppb std at Back Tower
11. Add 125 μL of 2000ppb Standard Stock to 250ppb std at Back Tower
12. Mix 250ppb std at 1000 RPM for 0 min 10 sec
13. Add 830 μL of Methanol (CH₃OH) to 125ppb std at Back Tower
14. Add 7.5 μL of Methanol (CH₃OH) to 125ppb std at Front Tower
15. Add 100 μL of ISTD to 125ppb std at Back Tower
16. Add 62.5 μL of 2000ppb Standard Stock to 125ppb std at Front Tower
17. Mix 125ppb std at 1000 RPM for 0 min 10 sec
18. Flag 1000ppb std, ... as Result.

Available Resources Tracked By Use

Resource Name	Resource Type	Uses/Vial	Vial Range
empty vial	Empty Container	1	91-94

Available Resources Tracked By Volume

Resource Name	Resource Type	Usable Volume/V
Methanol (CH ₃ OH)	Chemical Resource	1700 μL
ISTD	Chemical Resource	1700 μL
2000ppb Standard Stock	Chemical Resource	1500 μL
Front Solvent A	Turret Location	2000 μL
Front Solvent B	Turret Location	2000 μL

OK | Apply | Cancel | Help

Figure 4. Create prep method, linear dilution.

Create Prep Method - Serial Dilution

Using the Resources allocated above the Agilent 7696A WorkBench method is again created stepwise (shown below and in Figure 5).

1. Add 400 μL of methanol to an empty vial, which is renamed '1,000 ppb std'.
2. Add 100 μL of the resource 'ISTD' to the newly created vial '1,000 ppb std'.
3. Add 500 μL of the resource '2,000 ppb Standard Stock' to '1,000 ppb std'.
4. Mix '1,000 ppb std' (1,000 ppb std is now complete).
5. Add 400 μL of methanol to an empty vial, which is renamed '500 ppb std'.
6. Add 100 μL of the resource 'ISTD' to the newly created vial '500 ppb std'.
7. Add 500 μL of the resource '1,000 ppb std' to '500 ppb std'.
8. Mix '500 ppb std' (500 ppb std is now complete).
9. Add 400 μL of methanol to an empty vial, which is renamed '250 ppb std'.
10. Add 100 μL of the resource 'ISTD' to the newly created vial '250 ppb std'.
11. Add 500 μL of the resource '500 ppb std' to '250 ppb std'.
12. Mix '250 ppb std' (250 ppb std is now complete).
13. Add 400 μL of methanol to an empty vial, which is renamed '125 ppb std'.
14. Add 100 μL of the resource 'ISTD' to the newly created vial '125 ppb std'.
15. Add 500 μL of the resource '250 ppb std' to '125 ppb std'.
16. Mix '125 ppb std' (125 ppb std is now complete).
17. Flag the vials created: 1,000 ppb std, 500 ppb std, 250 ppb std, and 125 ppb std as 'result vials'.

Setup Method
Agilent 7696A Sample Prep Method | Agilent 7696A Configuration

Import | Export

Process in Batch Mode

Version 3.1.36.0

Actions: Add, Mix, Heat, Wait, Flag as result, Move, Wash

Program

Steps

1. Add 400 μL of Methanol (CH₃OH) to 1000ppb std at Back Tower
2. Add 100 μL of ISTD to 1000ppb std at Back Tower
3. Add 500 μL of 2000ppb Standard Stock to 1000ppb std at Back Tower
4. Mix 1000ppb std at 1000 RPM for 0 min 5 sec
5. Add 400 μL of Methanol (CH₃OH) to 500ppb std at Back Tower
6. Add 100 μL of ISTD to 500ppb std at Back Tower
7. Add 500 μL of 1000ppb std to 500ppb std at Back Tower
8. Mix 500ppb std at 1000 RPM for 0 min 5 sec
9. Add 400 μL of Methanol (CH₃OH) to 250ppb std at Back Tower
10. Add 100 μL of ISTD to 250ppb std at Back Tower
11. Add 500 μL of 500ppb std to 250ppb std at Back Tower
12. Mix 250ppb std at 1000 RPM for 0 min 5 sec
13. Add 400 μL of Methanol (CH₃OH) to 125ppb std at Back Tower
14. Add 100 μL of ISTD to 125ppb std at Back Tower
15. Add 500 μL of 250ppb std to 125ppb std at Back Tower
16. Mix 125ppb std at 1000 RPM for 0 min 5 sec
17. Flag 1000ppb std, ... as Result.

Available Resources Tracked By Use

Resource Name	Resource Type	Uses/Vial	Vial Range
empty vial	Empty Container	1	91-94

Available Resources Tracked By Volume

Resource Name	Resource Type	Usable Volume/ μL
Methanol (CH ₃ OH)	Chemical Resource	1700 μL
ISTD	Chemical Resource	1700 μL
2000ppb Standard Stock	Chemical Resource	1500 μL
Front Solvent A	Turret Location	2000 μL
Front Solvent B	Turret Location	none

OK | Apply | Cancel | Help

Figure 5. Create prep method, serial dilution.

Easy Sequence

To run the method created to generate the four calibration standards (Figure 6):

1. Select **Easy Sequence**.
2. Select your saved method.
3. Under **Starting Vial Location** enter any vial position not used in the Resource Layout (for this example we chose vial 1).
4. Under **Number of Samples** enter **1**.

You must enter values here even though there are no 'samples' when running a method to create calibration standards is used.

5. Under **Sample Name** enter a name (optional).
6. Select **Fill Samples** and start the sequence as you would normally.

This sequence will be added to your sequence queue.




Figure 6. Easy sequence.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2012
Printed in the USA
January 17, 2012
5990-9615EN