

Detecting Differences in PEEK Preparation

Pyrolysis, a good technique for competitive analysis, quickly yields information on subtle differences between polymers. In this application note, we spot small differences in two poly ether ether ketone (PEEK) samples.

PEEK is high-strength, radiation-resistant plastic. Thermally stable and highly resistant to chemicals, it is used in machine parts, nuclear power-plant equipment, automobile parts, aerospace components, cable insulation, and medical implants.

Figures 1 and 2 are pyrograms of 2 different PEEK samples. Although largely similar, (they both contain typical pyrolysis products of PEEK such as benzene, phenol, biphenyl, fluorene, and diphenyl methanone), we can see tiny differences between the them. The first PEEK sample has siloxanes (peak #s 2 and 4).

The second sample instead has difluorobenzophenone, dichlorobenzophenone, and diphenyl sulfone (DPS). These products could be residual starting components or solvent. One way of making PEEK involves reacting benzenediol's di-potassium salt with difluorobenzophenone using DPS as a solvent.

Using pyrolysis can make work in the lab easier by providing a quick, simple way for competitive analysis. With almost no sample preparation, you get a complete profile of the unknown polymer, and clues as to how it was made.




Figure 1




Figure 2

Peak Identification

1. Benzene
2. Hexamethyl Cyclotrioxane
3. Toluene
4. Octamethyl Cyclotetrasiloxane
5. Phenol
6. 1,1'-Biphenyl
7. Diphenyl ether
8. Difluorobenzophenone
9. Fluorene
10. Diphenyl Methanone
11. Dichlorobenzophenone
12. Diphenyl Sulfone

Both polymers were pyrolyzed at 750°C for 10 seconds using the CDS Model 5200 interfaced via a transfer line to an Agilent 6890/5973 GC/MS.

Model 5200 Conditions

Pyrolysis to GC mode (no trapping)

Interface zone: 300°C

Transfer Line: 290°C

Pyrolysis probe: Coil

Pyrolysis time: 10 seconds

Valve oven: 250°C

Gas Chromatograph Conditions

Carrier: Helium

Column: 35% Phenyl, methyl silicone

Detector: MSD

GC Program:

Initial: 40°C for 2 minutes

Ramp: 7°C/min.

Final: 300°C for 10 minutes

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.



CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto [www. cdsanalytical.com](http://www.cdsanalytical.com).