

Dual Channel Simulated Distillation of Carbon and Sulfur with the Agilent 7890A GC and 355 Sulfur Chemiluminescence Detector

Application Note

Hydrocarbon Processing

Authors

ChunXiao Wang
Agilent Technologies (Shanghai) Co., Ltd.
412 Ying Lun Road
Waigaoqiao Free Trade Zone
Shanghai 200131
China

Roger Firor and Paul Tripp
Agilent Technologies, Inc.
2850 Centerville Road
Wilmington DE 19808
USA

Abstract

Two-channel simulated distillation by gas chromatography (GC) for both hydrocarbons and sulfur is described. The method utilizes a 7890A GC configured with a high-temperature programmable temperature vaporizer (HT-PTV) inlet and a sulfur chemiluminescence detector (SCD) mounted in series with a flame ionization detector (FID) by use of a special mounting adapter. A simulated distillation (SimDis) software program provides an easy-to-use solution for sulfur and hydrocarbon simulated distillation. The data show that observed boiling point (BP) values agree with the ASTM D2887 consensus BP values within the allowable differences. The system also demonstrates very good repeatability for both hydrocarbon and sulfur SimDis. An example of a light cycle oil (LCO) analyzed according to D2887 is also included.



Agilent Technologies

Introduction

Sulfur and hydrocarbon simulated distillation results provide meaningful information to optimize refining processes and ensure compliance with petroleum product specifications. A previous application note [1] describes a 6890 GC based system for hydrocarbon simulated distillation by ASTM D2887 [2]. Now with the highly selective Agilent Sulfur Chemiluminescence Detector (SCD), sulfur simulated distillation is possible. This 7890A GC based simulated distillation system consists of acquiring and analyzing simultaneously the specific detector data for hydrocarbon (FID) and sulfur (SCD).

Experimental

This two-channel SimDis application uses the Agilent 7890A GC configured with a high-temperature programmable temperature vaporizer (HT-PTV) inlet, and an SCD mounted onto an FID using a special adapter. Detailed GC conditions used are listed in Table 1.

Table 1. 7890A Gas Chromatographic Conditions (1) D2887, (2) D7213

HT-PTV inlet typical temperature programs	(1) 225 to 350 °C (hold 15 min) at 200 °C/min to 225 °C at 100 °C /min (2) 50 to 420 °C (hold 15 min) at 200 °C /min to 50 °C at 100 °C /min
Split ratio	(1) 4:1 for diluted sample, 20:1 for nondiluted sample (2) 1:1
Injection volume	(1) 0.1 µL (2) 0.5 to 1 µL
Column	(1) HP-1 10 m × 530 µm × 0.88 µm (19095z-021) (2) DB-HT-SimDis 5 m × 530 µm × 0.15 µm (145-1001)
Column flow (He)	(1) 13 mL/min, constant flow mode (2) 16 mL/min, constant flow mode
FID temperatures	(1) 350 °C (2) 400 °C
H ₂ flow	40 mL/min
Air flow	400 mL/min
Make up (N ₂)	40mL/min
SCD	
Burner temperature	800 °C
Vacuum of burner	324 torr
Vacuum of reaction cell	11.6 torr
H ₂	40 SCCM
Air	8.3 SCCM
Oven programs	(1) 35 °C (hold 0.5 min) to 350 °C at 20 °C/min , hold 10 min (2) 40 to 420 °C at 20 °C/min , hold 6 min
Data acquisition rate	5 Hz typical

SimDis Software

The processes of SimDis analysis include: blank analysis for baseline subtraction, calibration for establishing the relationship between boiling point and retention time (RT), validation for verifying both the chromatographic conditions and calculations in the method, and sample analysis. The Agilent SimDis software divides these functions under separate tabs that make navigation and data processing straightforward. The software is based on four modules: Browse, Setup, SimDis, and Report. For example, the Setup module allows you to configure the files to use for BP calibration, blank selection, and QC reference. Partial integration with the GC ChemStation sequence makes automated data analysis possible.

Processing Two Signals

The software can process one or two channels of signal data (FID and SCD for example) from GC ChemStation data files. When working with dual channels, the SimDis software requires that each channel be labeled by the detector type rather than the defaults used by the GC ChemStation. Since the SCD operates off the analog input board (AIB), its signal begins with "AIB." For this reason, the post-run command macro SCDnamer.mac must be run to rename the signal file. The macro renames the AIB2B.ch channel as SCD1.ch. If the channel name is not corrected, the software will switch the FID and SCD channels during analysis, giving faulty results. The macro code to do this is shown below. It assumes the AIB is in the rear position (B).

```
!=====
! SCDNamer call this as a post run command when an SCD is
! installed
! it renames the dual channel AIB2B.ch to SCD1.ch to allow
! simdis to
! properly calibrate
!=====
NAME SCDNamer
! This macro renames the SCD files named as AIB2B.ch to
SCD1.ch
if filestat(mode,dadatapath$+dadatfile$+"\AIB2B.CH")=1
rename dadatapath$+dadatfile$+"\AIB2B.ch",dadat-
path$+dadatfile$+"\SCD1.ch"
print "File Renamed"
else
print "No AIB2B File found"
endif
RETURN
ENDMACRO
```

Results and Discussion

Calibration

A calibration mixture containing a series of known n-alkanes can be used for establishing the relationship between BP and RT. C5 to C40 is used for ASTM D2887, and Polywax 500 dissolved in toluene is used to calibrate ASTM D7213 [3]. Since both are too viscous or waxy at ambient temperature to sample with a syringe, they need to be heated manually to approximately 80 °C before injection. RT repeatability is key for consistent correlation of BP and RT. Figure 1 and Figure 2

show overlays of consecutive runs of C5 to C40 and Polywax 500, respectively. Tables 2 and 3 show repeatability for both RT and area.

Polywax 500 Sample Preparation

Place approximately 80 mg of Polywax 500 in a 2-mL vial. Add about 1.5 mL toluene followed by the addition of a suitable mixture of n-paraffins from C5 to C18 (Agilent SimDis calibration No.2). The final concentration should be approximately one part of (C5–C18) to 20 parts of toluene. Initially heat the solution to 80 °C to dissolve the Polywax 500.




Figure 1. Overlay of five consecutive runs of C5 to C40 calibration mix, vial heated to 80 °C for 3 min prior to injection. GC conditions are listed in Table 1, items (1).




Figure 2. Overlay of four consecutive runs of Polywax 500 plus C5–C18. GC conditions are listed in Table 1, items (2).

Table 2. Repeatability for C5 to C40, n = 10

	Retention Time			Area		
	Average	STDEV	RSD%	Average	STDEV	RSD%
C5	0.275	0.000	0.06	19870	126	0.64
C6	0.388	0.000	0.09	14020	83	0.60
C7	0.673	0.001	0.17	16527	108	0.65
C8	1.192	0.002	0.16	18693	81	0.43
C9	1.874	0.002	0.12	20383	107	0.53
C10	2.622	0.003	0.10	43561	280	0.64
C11	3.338	0.002	0.07	22730	158	0.69
C12	4.068	0.002	0.05	94289	714	0.76
C14	5.327	0.002	0.03	48149	393	0.82
C15	5.902	0.002	0.03	24268	199	0.82
C16	6.477	0.001	0.02	49175	408	0.83
C17	6.991	0.001	0.02	24448	201	0.82
C18	7.485	0.000	0.00	10552	84	0.80
C20	8.424	0.001	0.01	6187	53	0.86
C24	10.083	0.000	0.00	4293	17	0.40
C28	11.512	0.001	0.01	4288	45	1.06
C32	12.762	0.002	0.01	3988	66	1.66
C36	13.874	0.001	0.01	3407	66	1.94
C40	14.874	0.002	0.01	3238	69	2.14

Table3. Repeatability of Polywax 500 Plus C5 to C18, n = 10

	Retention Time			Area		
	Average	STDEV	RSD%	Average	STDEV	RSD%
C14	2.769	0.002	0.07	49126	953	1.94
C15	3.278	0.002	0.05	24337	469	1.93
C16	3.847	0.002	0.05	49304	948	1.92
C17	4.311	0.002	0.05	24597	470	1.91
C18	4.753	0.001	0.03	11374	218	1.92
C20	5.596	0.001	0.01	952	17	1.80
C22	6.424	0.001	0.01	1635	30	1.81
C26	7.904	0.001	0.01	3615	62	1.71
C32	9.783	0.001	0.01	6856	105	1.53
C36	10.858	0.001	0.01	8418	137	1.63
C40	11.823	0.001	0.01	8432	128	1.52
C44	12.690	0.002	0.01	7037	137	1.95
C48	13.480	0.001	0.01	5288	104	1.98
C52	14.208	0.001	0.01	3677	67	1.83
C60	15.512	0.001	0.01	1353	19	1.40
C70	16.931	0.002	0.01	273	5	1.92

QC Reference

A QC reference sample is the basis for quantifying total sulfur and allows the direct entry of response factors for calculation based on total area and user-entered concentrations of sulfur. In this application, a diesel sample (SDF-1X-4, AccuStandard, Inc., New Haven, CT) with a sulfur concentration of 100 µg/g

is used as the QC external reference for calibration of response factors for the SCD channel. This is needed for calculation of total sulfur in the sample. Figure 3 shows the graphic pane from the SimDis software for of the QC reference.

Reference Gas Oil Analysis

To meet the requirements of ASTM D2887, the reference gas oil (RGO) sample analysis must be performed to verify both the chromatographic performance and the calculation algo-




Figure 3. QC reference setup. GC conditions are listed in Table 1, items (1).

ritms involved in this test method. Figure 4 shows the chromatograms of RGO for both the hydrocarbon and sulfur channels. Tables 4 and 5 show the results for six runs of RGO

analysis. The data show that observed BP values agree with the ASTM D2887 consensus BP values within the allowable differences and with good repeatability.




Figure 4. Chromatograms of RGO for hydrocarbon and sulfur channels. GC conditions are listed in Table 1, items (1).

Table 4. Hydrocarbon SimDis Results for Reference Gas Oil (Six runs shown.)

ASTM D2887 Values											
OFF %	BP, °C	Allowable Difference	1	2	3	4	5	6	Average	Difference	RSD%
IBP	115	7.6	114	114	114	114	114	114	114	1	0.00
10%	176	4.1	174	174	174	174	174	174	174	2	0.00
20%	224	4.9	223	223	223	223	223	223	223	1	0.00
30%	259	4.7	258	258	258	258	258	258	258	1	0.00
40%	289	4.3	287	287	287	287	287	287	287	2	0.00
50%	312	4.3	311	311	311	311	311	311	311	1	0.00
60%	332	4.3	330	330	330	330	330	330	330	2	0.00
70%	354	4.3	352	352	351	352	352	352	352	2	0.12
80%	378	4.3	376	376	376	376	376	376	376	2	0.00
90%	407	4.3	405	405	405	405	405	405	405	2	0.00
FBP	475	11.8	471	471	471	471	471	471	471	4	0.00

Table 5. Sulfur SimDis Results for Reference Gas Oil, BP in °C

OFF%	1	2	3	4	5	6	Average	STDEV	RSD%
IBP	168	169	169	167	165	169	168	1.60	0.95
10%	265	265	265	265	265	265	265	0.00	0.00
20%	293	293	293	293	293	293	293	0.00	0.00
30%	314	314	314	314	314	314	314	0.00	0.00
40%	329	330	330	330	330	330	330	0.41	0.12
50%	344	344	344	344	344	345	344	0.41	0.12
60%	359	359	359	360	360	360	360	0.55	0.15
70%	376	376	377	377	377	377	377	0.52	0.14
80%	396	396	396	397	397	398	397	0.82	0.21
85%	408	408	408	409	409	409	409	0.55	0.13
90%	422	422	423	423	424	424	423	0.89	0.21
FBP	495	495	495	499	499	501	497	2.66	0.53

Light Cycle Oil Analysis

To illustrate repeatability, chromatographic overlays are shown in Figures 5a and 5b for an LCO sample. Tables 6 and 7 list the results for hydrocarbon and sulfur SimDis, respectively. The average total sulfur content calculated is 248 ppm with 3.5% RSD.




Figure 5a. Carbon SimDis of LCO. Five-run overlay.




Figure 5b. Sulfur SimDis of LCO. Five-run overlay.

Table 6. Carbon SimDis Results for LCO, BP in °C

OFF%	1	2	3	4	5	Average	SD	RSD%
IBP	141	140	140	140	139	140	0.71	0.51
10%	221	221	221	221	221	221	0.00	0.00
20%	233	233	233	233	234	233	0.45	0.19
30%	247	247	247	247	247	247	0.00	0.00
40%	260	260	261	260	261	260	0.55	0.21
50%	274	275	275	275	275	275	0.45	0.16
60%	291	292	292	292	292	292	0.45	0.15
70%	306	307	307	307	307	307	0.45	0.15
80%	324	324	324	324	324	324	0.00	0.00
90%	344	344	344	344	344	344	0.00	0.00
FBP	391	391	391	391	392	391	0.45	0.11

Table 7. Sulfur SimDis Results for LCO, BP in °C

OFF%	1	2	3	4	5	Average	SD	RSD%
IBP	314	314	314	314	314	314	0.00	0.00
10%	328	329	328	328	328	328	0.45	0.14
20%	329	329	329	329	329	329	0.00	0.00
30%	329	329	329	329	329	329	0.00	0.00
40%	332	332	332	332	332	332	0.00	0.00
50%	342	342	342	342	342	342	0.00	0.00
60%	345	345	345	345	345	345	0.00	0.00
70%	347	347	346	346	347	347	0.55	0.16
80%	351	351	350	350	351	351	0.55	0.16
90%	359	359	357	359	358	358	0.89	0.25
FBP	375	375	371	374	371	373	2.05	0.55
Sulfur, ppm	254	250	240	238	258	248	8.62	3.48

Conclusions

This new SimDis procedure utilizes a 7890A GC configured with the HT-PTV inlet, and an SCD mounted in series with an FID. The Agilent SimDis software is capable of processing both FID and SCD data channels, providing a solution for hydrocarbon and sulfur simulated distillation.

Sulfur simulation distillation has been demonstrated using the Agilent 355 sulfur chemiluminescence detector. With a selectivity over carbon of approximately 10^6 , reliable boiling point distributions of sulfur in petroleum fractions can be obtained.

References

1. C. Wang and R. Firor, "Simulated Distillation System for ASTM D2887," Agilent Technologies, publication 5989-2726EN
2. ASTM D2887-06 a, "Standard Test Method for Boiling Range Distribution of Petroleum Fractions by Gas Chromatography," Annual Book of Standards, Volume 05.01, ASTM, 100 Barr Harbor Drive, West Conshohocken, PA 19428 USA
3. ASTM D 7213-05, "Standard Test Method for Boiling Point Distribution of Petroleum Distillates from 100 °C to 615 °C by Gas Chromatography," Annual Book of Standards, Volume 05.04, ASTM, 100 Barr Harbor Drive, West Conshohocken, PA 19428 USA

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2008
Published in the USA
December 5, 2008
5990-3237EN



Agilent Technologies