



CDS Solutions




APPLICATIONS INFORMATION USING ADVANCED GC SAMPLE HANDLING TECHNOLOGY

Dynamic Headspace and Pyrolysis in Ink Analysis

Printing ink is a complex formulation of pigment, petroleum resins, waxes and oils. Once produced, it may be difficult to analyze since it is a combination of solid and volatile or semivolatile materials. Pyrolysis GC may be useful in the analysis of the non-volatile constituents, but the presence of volatile materials may interfere. In cases such as this, a two-step process may simplify the situation.

Figure 1 shows a pyrogram of a finished ink product with no pretreatment. The presence of a large, unresolved mass of peaks from 20 to 26 minutes makes interpretation difficult. To simplify the chromatogram, the sample was first heated to 250°C for 5 minutes in the Pyroprobe interface. This process removed the volatile hydrocarbon oil which made up a large part of the ink, leaving behind the solid resins and waxes. The volatilized hydrocarbon oil was then chromatographed, producing Figure 2.

The remaining solids were then pyrolyzed at 700°C for ten seconds, and the resulting pyrolysate chromatographed separately, as seen in Figure 3. Now there is nearly baseline resolution of the major peaks, and specific peaks may be identified and associated with individual solid constituents of the ink formulation. Many of the peaks eluting between 12 and 19 minutes, for



(Continued on back)

(Continued from front)

example, resulted from pyrolysis of the hydrocarbon resins, while the peaks eluting from 20 to 25 minutes are normal alkanes and alkenes produced from the wax.

EQUIPMENT PYROLYSIS

CDS Pyroprobe 120 with
Platinum ribbon probe
Headspace: 250°C for
5 minutes
Pyrolysis: 700°C for
10 seconds

CHROMATOGRAPHY

Column: 50m x 0.25mm SE-54
fused silica capillary
Carrier: Helium, split 70:1
Program: 50°C for 2 min, then
10°/min to 290°
Detector: FID

For more information on this
and related applications, we rec-
ommend the following readings:

Wampler, T., and Levy, E., "Pyrol-
ysis GC in the Analysis of Ink and
Paper," *LC GC*, 4, 11, (1986)
1112-1116

Levy, E., and Wampler, T., "Appli-
cations of Pyrolysis Gas Chroma-
tography/Mass Spectrometry
to Toner Material from Photo-
copiers," *J. Forensic Sci.*, 31, 1,
(1986) 258-271.

Wampler, T., and Levy, E., "Ana-
lytical Pyrolysis in the Forensic
Science Laboratory," *Crime Lab
Digest*, 12, 2, (1985) 25-28.

Additional literature may be
obtained from CDS by calling
1-800-541-6593 or in Pennsyl-
vania 215-932-3636.

ABOUT CDS

CDS Instruments, a unit of Autoclave Engineers Inc., is a worldwide leader in the manufacture of instru-
ments and equipment for research applications. With over 15 years of service to the research community,
CDS is dedicated to the development of instrumentation for use in the preparation, separation and analysis
of complex organics. We are also dedicated to providing researchers with the technical and applications
information they need to perform their work efficiently. The purpose of CDSolutions is to supply the re-
searcher with a constant stream of applications solutions for advanced GC sample handling. The article you
just read is one of many available at no charge to our customers. For a complete list of topics, please write
to our Technical Information Dept. in Oxford, Pennsylvania.



• 7000 Limestone Rd., Oxford, PA 19363, U.S.A. • Phone: (215) 932-3636 • Telex: 83-5308